

Plasmid Preparation using Corning FiltrEX™ 384 Well Filter Plates

Protocol



Life
Sciences

*Helen L. Krasnow, Linda S. Belkowski,
Ph.D. and Debra S. Hoover, Ph.D.
Corning Incorporated, Life Sciences,
45 Nagog Park, Acton, MA 01720*

*Tony West, Wellcome Trust Sanger Institute,
The Wellcome Trust Genome Campus,
Hinxton, Cambridgeshire, UK*



Introduction

This procedure describes an easily automated method for bacterial growth and lysis and plasmid purification. Bacteria are lysed under alkaline conditions and the lysates are clarified by filtration through a 0.45 µm PVDF membrane. DNA is recovered from the cleared lysate by alcohol precipitation.

Procedure

- Dispense 160 µL/well of LB or Circle Grow containing appropriate antibiotics into a 96 well plate (Corning Cat. No. 3960) Inoculate each well with a single colony of *E. coli* and incubate the plate at 37°C with vigorous shaking for 20 to 22 hours.*

*We use 96 well plates for bacterial growth and lysis followed by a 384 well filter plate for lysate clarification. Alternatively, colonies can be inoculated into 160 mL/well of the appropriate media in a 384 well block (Corning Cat. No. 3964) .

- Pellet the cells by centrifugation at 2000 x g for 2 minutes at room temperature.
- Add 25 µL/well of Solution 1 + RNase to each pellet and resuspend thoroughly.
- Add 25 µL/well of Solution 2 and vortex to mix. Incubate the plate at room temperature for 2 minutes.
- Add 35 µL/well of Solution 3 and vortex to mix.
- Transfer 80 µL of lysate to each well of a 384 well filter plate (0.45 µm PVDF, Corning Cat. No. 3531).
- Add 85 µL/well of 100% isopropanol to a 384 well collection plate (such as Corning Cat. No. 3657, 3702 or 3965).
- Place the filter plate above the collection plate and centrifuge at 2000 x g for 20 minutes at 4°C.
- Decant the isopropanol and wash the pellets with 100 µL/well 70% ethanol. Centrifuge at 2000 x g for 3 minutes at 4°C.
- Carefully decant or aspirate the alcohol from the pellets.
- Centrifuge briefly (1 minute) to assure that the pellets are located at the bottom of each well.
- Air dry and resuspend the DNA in 30 µL of the desired buffer.

Reagents

Solution 1

22.5 mM Tris pH 8.0
0.25% glucose
10 mM EDTA

1M Tris pH 8.0	22.5 mL
0.1M EDTA	100 mL
20% Glucose	12.5 mL
ddH ₂ O	865 mL
Final Volume	1000 mL

Add RNase A to Solution 1 to a final concentration of 60 µg/mL just prior to use.

Solution 2

0.2 N NaOH
0.7% SDS
0.1% TritonX-00

4N NaOH	50 mL
20% SDS	35 mL
TritonX-100	1 mL
ddH ₂ O	914 mL
Final Volume	1000 mL

Solution 3

3M Potassium acetate, pH = 5.5

Corning Incorporated Life Sciences

45 Nagog Park
Acton, MA 01720
t 800.492.1110
t 978.635.2200
f 978.635.2476

[www.corning.com/
lifesciences](http://www.corning.com/lifesciences)

Worldwide Support Offices

ASIA

Australia
t 61 2-9416-0492
f 61 2-9416-0493

China
t 86 21-6361-0826
f 86 21-6361-0827

Hong Kong
t 852-2807-2723
f 852-2807-2152

India

t 91 11 341 3440
f 91 11 341 1520

Japan

t 81 (0) 3-3586 1996/1997
f 81 (0) 3-3586 1291/1292

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-7310

Taiwan
t 886 2-2716-0338
f 886 2-2716-0339

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

LATIN AMERICA

Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589