

One-STrEP-tag for protein:protein-interaction analysis

Johansen et al., 2008

Protein:protein-interactions (PPI) govern almost all important processes in living organisms. Thus, their rapid and accurate determination and investigation is a major challenge in life sciences.

The One-STrEP-tag is a reliable tool to isolate protein complexes at high purity without losing transient binders. The advantages of the One-STrEP system in comparison to the commonly used TAP system are that only one tag is needed and the isolation of the protein complex components is carried out in one step. An One-STrEP-tagged bait protein is expressed

in the target cell and the cell lysate containing the bait and the putatively interacting preys is subjected to tag-based affinity chromatography on *Strep-Tactin*[®]. The isolated protein complexes are analyzed by SDS-PAGE and potential preys are identified by mass spectrometry (Rigaut et al., 1999). In several cases the One-STrEP-tag was already successfully applied to isolate and identify bait proteins and its interaction partners (see below for references). Johansen et al., 2008,




Figure 1: IKAP co-purifies with various cytosolic proteins.

One-STrEP-tag purification. HEK293 cells transiently transfected with either empty vector or cIKAP-strep were harvested 2 days after transfection and cytosolic extracts were prepared as described in the Materials and Methods. Purification was performed according to the manufacturer's instructions (IBA). Eluates were concentrated and run on 10% SDS-PAGE. Silver staining of the purified proteins. The indicated proteins were identified by MALDI-TOF-MS. DNPk1, DNA-dependent protein kinase; USP9X, ubiquitin-specific processing protease.

(Reproduced with permission of the Journal of Cell Science. For details see Johansen et al., 2008 in the reference list)

screened for interaction partners of the I κ B-kinase-associated protein, IKAP, which is involved in the neurodegenerative disease familial dysautonomia. "Previous studies of proteins that might associate with IKAP revealed no candidates that would fit with both its cytosolic localization and function in cell migration. Therefore we decided to use the One-STrEP-tag, since it seems to be an excellent tag for the identification of even weakly-binding interactors" (Tuula Kallunki, personal comment). The One-STrEP-tagged C-terminus of IKAP (cIKAP-strep) was used as bait and a total of 15 proteins were identified to

associate specifically with cIKAP-strep (Figure 1). Two of the potential interactors (filamin A and dynein heavy chain) had been linked to neuronal migration in previous studies, what coincides with one of the functions IKAP is involved. Further analysis of IKAP and filamin A revealed their co-localization and verified their association. This example is demonstrating the effectiveness of the One-STrEP system purifying protein complexes in combination with protein identification via MS and adjacent functional analysis.

References:

1. Junttila MR, Saarinen S, Schmidt T, Kast J, Westermarck J, 2005: *Proteomics* 5: 1199-1203. Single-step Strep-tag[®] purification for the isolation and identification of protein complexes from mammalian cells.
2. Johansen LD, Naumanen T, Knudsen A, Westerlund N, Gromova I, Junttila M, Nielsen C, Bøttzauw T, Tolkovsky A, Westermarck J, Coffey ET, Jäättelä M, Kallunki T, 2008: *J Cell Sci* 121: 854-864. IKAP localizes to membrane ruffles with filamin A and regulates actin cytoskeleton organization and cell migration. [Article](#)
3. Groth A, Corpet A, Cook AJL, Roche D, Bartek J, Lukas J, Almouzni G, 2007: *Science* 318: 1928-1931. Regulation of replication fork progression through histone supply and demand.
4. Morita E, Sandrin V, Alam SL, Eckert DM, Gygi SP, Sundquist WI, 2007: *Cell Host & Microbe* 2: 19-28. Identification of human MVB12 proteins as ESCRT-I subunits that function in HIV budding.
5. Weber M, Wehling M, Lösel R, 2008: *Am J Physiol Heart Circ Physiol* 295: H361-H365. Proteins interact with the cytosolic mineralocorticoid receptor depending on the ligand.
6. Jasencakova Z, Scharf AND, Ask K, Corpet A, Imhof A, Almouzni G, Groth A, 2010: *Mol Cell* 37:736-743. Replication stress interferes with histone recycling and predeposition marking of new histones.
7. Bekker-Jensen S, Rendtlew Danielsen J, Fugger K, Gromova I, Nerstedt A, Bartek J, Lukas J, Mailand N, 2010: *Nat Cell Biol* 12 (1): 80-86. HERC2 coordinates ubiquitin-dependent assembly of DNA repair factors on damaged chromosomes.
8. Rameix-Welti MA, Tomoiu A, Dos Santos Afonso E, Van der Werf S, Naffakh N, 2009: *J Virol* 83 (3): 1320-1331. Avian Influenza A virus polymerase association with nucleoprotein, but not polymerase assembly, is impaired in human cells during the course of infection.
9. Pegoraro G, Kubben N, Wickert U, Göhler H, Hoffmann K, Misteli T, 2009: *Nat Cell Biol* 11: 1261-1267. Aging-related chromatin defects through loss of the NURD complex.
10. Amako Y, Sarkeshik A, Hotta H, Yates J, Siddiqui A, 2009: *J Virol* 83(18): 9237-924. Role of Oxysterol Binding Protein in Hepatitis C Virus infection.
11. Gianni T, Amasio M, Campadelli-Fiume G, 2009: *J Biol Chem* 284: 17370-17378. Herpes simplex virus gD forms distinct complexes with fusion executors gB and gH/gL through the C-terminal profusion domain.
12. Rigaut G, Shevchenko A, Rutz B, Wilm M, Mann M, Séraphine B, 1999: *Nature Biotech* 17: 1030-1032. A generic protein purification method for protein complex characterization and proteome exploration.