


AGAROSE

BUFFERS

LADDERS

EQUIPMENT



Nucleic Acid Electrophoresis

APPLICATION GUIDE

REAGENTS: AGAROSE



Thermo Scientific and Fisher Scientific products deliver an end-to-end solution that can meet your most demanding electrophoresis requirements.

You can depend on our expertise in electrophoresis instruments along with ultrapure reagents that are pre-qualified for your applications. This guide is designed to help you select the right products from our best-in-class array of laboratory equipment and bioreagents.

All Fisher BioReagents® agarose are DNase- and RNase-free to ensure optimal results for your nucleic acid application.

Fisher BioReagents offers three different grades of agarose that are functionally tested and pre-qualified for specific applications.

Agarose grades used in electrophoresis of nucleic acids:

Genetic Analysis Grade—agarose that yields biologically active DNA or RNA. Testing includes enzymatic performance measurements.

Molecular Biology Grade—suitable for analytical separation of DNA or RNA.

PCR Grade—the original agarose for analytical separation of PCR amplicons (<1kb).

Agarose is a linear polysaccharide composed of alternating residues of D- and L-galactose joined by glycosidic linkages. Agarose forms gels that are both porous and resilient.

These gel properties provide a sieving matrix that allows the electrophoretic separation of charged macromolecules such as DNA or RNA according to size. Compared to polyacrylamide gel, agarose has a lower resolution but wider range of separation.

Lower grades of agarose can be contaminated with other polysaccharides, salts, and proteins. Such impurities can alter the gelling/melting temperature of agarose solutions or affect the ability to use the recovered nucleic acid sample in a post-electrophoresis application.

3-STEP Selection Process

Separation of Nucleic Acids by AGAROSE GEL ELECTROPHORESIS

1. Choose Your Reagents

- Agarose
- Buffer
- Ladders

2. Choose Your Equipment


- Power Supply
- Gel Box

3. Downstream Application Essentials

- Gel Staining
- Hybridization
- DNA Gel Extraction

Two Factors for Selecting an Agarose

1. The size of DNA or RNA fragments to be analyzed (see graph below).



2. The type of downstream application that will follow electrophoretic separation (e.g., cloning procedures directly from remelted agarose or in-gel reaction).

Agarose Selection Guide

| Type of Agarose | Low EEO | Low Melting >200bp | Low Melting <1000bp | Wide Separation Range | PCR Grade |
|---|-------------------|--------------------|---------------------|-----------------------|---------------|
| Cat. No. | BP160 | BP165 | BP1360 | BP1356 | BP2410 |
| Recovery of DNA and RNA | x | x | x | x | x |
| Southern and Northern Blots | x | | | | |
| DNA/RNA separation 50bp to 1kb | | | x | | x |
| DNA/RNA separation >1kb | x | x | | x | |
| PCR fragment analysis | x | x | x | x | x |
| In-gel reactions (ligation, transformations, PCR) | | | x | | |
| Colony lifts | x | | | | |
| Available pack sizes | 100g and 500g | 25g | 100g | 100g and 500g | 100g |
| Agarose grade | Molecular biology | Molecular biology | Genetic analysis | Genetic analysis | PCR |

REAGENTS: BUFFERS



Two buffers commonly used for DNA agarose electrophoresis are Tris-acetate with EDTA (TAE; 40mM Tris-acetate, 1mM EDTA) and Tris-borate with EDTA (TBE, 89mM Tris-borate, 2mM EDTA). Because the pH of these buffers is neutral, the phosphate backbone of DNA has a net negative charge and migrates toward the anode. TAE and TBE have different properties which makes one more suitable than the other for a specific purpose.

MOPS is a commonly used buffer system for RNA electrophoresis using formaldehyde or formamide denatured RNA. It is important to use RNase-free chemicals, water and containers when preparing the buffer solution. The typical formulation of a 10X MOPS running buffer is 0.4M MOPS (pH 7.0), 0.1M sodium acetate, and 0.01M EDTA.



The denaturing system chosen depends on the purpose of the RNA experiment and the size of the RNA fragment being separated. Formaldehyde denaturation is suitable if RNA samples are to be recovered. Formamide denaturation is suitable if the RNA needs to retain its biological activity.

| Buffer | Suggested Uses | Properties |
|--------|---|--|
| TAE | DNA recovery. Electrophoresis of large DNA (>12kb). | Low buffering capacity. Recirculation may be necessary for extended run times (>6 hr.) |
| TBE | Electrophoresis of small DNA (<1kb). Increased resolution of small DNA (<1kb). | Decreased DNA mobility. High buffering capacity – no recirculation required for extended run times. |
| MOPS | Electrophoresis of formaldehyde denatured RNA. | Buffer is low in ionic strength. Recirculation of buffer may be necessary. |

Suggested Agarose Concentrations

The optimal gel concentration depends on the size of the DNA fragments to be resolved.

| Cat. No. | Main Application | DNA Size Range in Base Pairs | Final Agarose Concentration % (W/V) | |
|----------|---|------------------------------|-------------------------------------|---------------|
| | | | 1x TAE buffer | 1x TBE buffer |
| BP1360 | Low melting temperature agarose. | 500-1,000 | 2.5 | 2.0 |
| | Certified recovery of small nucleic acid fragments. | 150-700 | 3.0 | 2.5 |
| | Outstanding resolution. | 100-450 | 3.5 | 3.0 |
| | | 70-300 | 4.0 | 3.5 |
| | | 10-100 | 4.5 | 4.0 |
| | 8-50 | 5.0 | 4.5 | |
| BP165 | Low melting temperature agarose. | 500-25,000 | 0.75 | 0.70 |
| | Broad separation range. | 300-20,000 | 1.00 | 0.85 |
| | Ideal for DNA and RNA recovery after electrophoretic separation. | 200-12,000 | 1.25 | 1.00 |
| | | 150-6,000 | 1.50 | 1.25 |
| | | 100-3,000 | 1.75 | 1.50 |
| | 50-2,000 | 2.00 | 1.75 | |
| BP1356 | Suitable for routine nucleic acid electrophoresis applications with broad separation range. | 1,000-23,000 | 0.60 | 0.50 |
| BP160 | | 800-10,000 | 0.80 | 0.70 |
| | | 400-8,000 | 1.00 | 0.85 |
| | | 300-7,000 | 1.20 | 1.00 |
| | | 200-4,000 | 1.50 | 1.25 |
| | | 100-3,000 | 2.00 | 1.75 |

NEW!

Ethanol, Molecular Biology Grade, is an ultrapure molecular biology grade ethanol that is DNase-, RNase- and Protease-free. It is used for the purification and precipitation of biomolecules such as nucleic acids and proteins.

| Cat. No. | Size |
|------------|-------|
| BP2818-100 | 100mL |
| BP2818-500 | 500mL |
| BP2818-4 | 4L |

NEW!

Water, Molecular Biology Grade, is ideal for many fundamental procedures such as PCR, electrophoresis, DNA sequencing and buffers for enzymatic analyses.

| Cat. No. | Size |
|------------|-------|
| BP2819-100 | 100mL |
| BP2819-1 | 1L |
| BP2819-4 | 4L |
| BP2819-10 | 10L |
| BP2819-20 | 20L |

Buffers for Nucleic Acid Applications

| Cat. No. | Concentration | Size |
|------------|---------------|-------|
| TBE | | |
| BP2430-1 | 1X | 1L |
| BP2430-4 | 1X | 4L |
| BP2430-20 | 1X | 20L |
| BP1396-86 | 5X | 1L* |
| BP1333-1 | 10X | 1L |
| BP1333-4 | 10X | 4L |
| BP1333-20 | 10X | 20L |
| BP1334-1 | 10X | 1L** |
| TAE | | |
| BP2434-4 | 1X | 4L |
| BP2434-20 | 1X | 20L |
| BP1335-500 | 10X | 500mL |
| BP1335-1 | 10X | 1L |
| BP1335-4 | 10X | 4L |
| BP1335-20 | 10X | 20L |
| BP1330-1 | 25X | 1L |
| BP1332-500 | 50X | 500mL |
| BP1332-1 | 50X | 1L |
| BP1332-4 | 50X | 4L |
| BP1332-20 | 50X | 20L |
| BP1331-1 | 25X | 1L** |

| Cat. No. | Description | Size |
|---------------------|---------------------|-------|
| MOPS | | |
| BP308-100 | Powder | 100g |
| BP308-500 | Powder | 500g |
| BP2900-500 | 10x Buffer Solution | 500mL |
| BP2900-1 | 10x Buffer Solution | 1L |
| WATER | | |
| BP2484-50 | Nuclease-Free | 50mL |
| BP2484-100 | Nuclease-Free | 100mL |
| BP2470-1 | DNA-Grade | 1L |
| BP561-1 | RNA-Grade | 1L |
| FORMALDEHYDE | | |
| BP531-25 | 37% by weight | 25mL |
| BP531-500 | 37% by weight | 500mL |


*Pre-weighed powder in poly bottle. Dissolve in water.

** Pre-weighed powder in foil pack. Dissolve in water.

REAGENTS: LADDERS

In addition to premixed buffers, Fisher BioReagents also offers individual buffer components for you to create your own buffer recipe

| Buffer Component | Size | Cat. No. |
|---------------------|-------|--|
| Tris Base | 500g | BP152-500 |
| | 1kg | BP152-1 |
| | 5kg | BP152-5 |
| | 10kg | BP152-10 |
| | 25kg | BP152-25 |
| Glacial Acetic Acid | 500mL | BP2401-2500, BP2401S-500 (Safe-Cote®) |
| | 2.5L | BP2401-212, BP2401SI-212 (Safe-Cote) |
| Boric Acid | 500g | BP168-500 |
| | 1kg | BP168-1 |
| EDTA Disodium Salt | 500g | BP120-500 |
| | 1kg | BP120-1 |



To achieve the most accurate qualitative and quantitative analysis via agarose gel electrophoresis, the appropriate DNA or RNA standard is required.



Fisher BioReagents provides a wide range of standards, including routine DNA ladders for quick size and quality assessment as well as *exACTGene*® DNA ladders that allow for quantitative analysis.

RiboLadders™ RNA Standards



These standards can be used to assess single-stranded RNA molecules on both native and denaturing agarose gels. These unique RNA standards are lyophilized to reduce thawing-related degradation, to prolong shelf life and to ensure consistent performance.

| Cat. No. | Application | Size Range | Number of Bands | Number of Loadings |
|-----------|-------------------------------------|------------|-----------------|--------------------|
| | Sizing unknown RNA fragments | | | |
| BP2810-50 | Small RNA fragments | 0.1 – 1kb | 8 | 50 |
| BP2811-50 | Large RNA fragments | 0.2 – 4kb | 9 | 50 |

Routine DNA Ladders



RiboLadders RNA Standards



REAGENTS: LADDERS



exACTGene® and Routine DNA Ladders

Ready-to-use (pre-mixed with the loading dye), room temperature, stable DNA ladders are available for all common electrophoresis applications.

NEW!

On Chemical Mini-site

Reagent Recommendation Tool for Nucleic Acid Electrophoresis

One click to see the ideal reagents for your specific application

www.fishersci.com/
NucleicAcidElectrophoresis

Reagent Recommendation Tool

Target Sample:
 DNA RNA





Target Fragment Size:
 Move the blue slider below to the target size.
 0 bp/b

Reagent Recommendations:
 Agarose Buffer Ladder

| Cat. No. | Application | Size Range | Number of Bands | Number of Loadings |
|------------|--|--------------|-----------------|--------------------|
| | exACTGene DNA ladders are ideal for qualitative analysis, quantitative estimation and size assessment | | | |
| BP2570-100 | PCR fragment analysis | 25-650bp | 14 | 100/10uL |
| BP2571-100 | PCR fragment analysis, small DNA digests | 25-1000bp | 12 | 100/10uL |
| BP2572-100 | Quick check of PCR or enzyme digestion results | 50-2000bp | 8 | 100/10uL |
| BP2573-100 | General purpose, small DNA fragments | 100-1000bp | 10 | 100/10uL |
| BP2574-100 | Fast run times, small DNA fragments | 100-2000bp | 11 | 100/10uL |
| BP2575-100 | Clone identification | 100-2686bp | 14 | 100/10uL |
| BP2576-100 | Large size PCR or cloning | 300-5000bp | 10 | 100/10uL |
| BP2577-100 | Small and large cloning application | 100-5000bp | 16 | 100/10uL |
| BP2578-100 | General purpose, large digested DNA | 300-10,000bp | 13 | 100/10uL |
| BP2579-100 | General purpose, wide size range | 100-10,000bp | 19 | 100/10uL |
| BP2580-100 | General purpose, extra-large fragments | 300-24,000bp | 15 | 100/10uL |
| | Routine DNA ladders are designed for qualitative analysis and size assessment | | | |
| BP2581-200 | Small fragments, quick size assessment | 50-2000bp | 11 | 200/5uL |
| BP2582-200 | Quick size assessment of broad size range | 50-10,000bp | 16 | 200/5uL |

For Lambda DNA digests or other DNA markers and ladders not containing loading dye, please visit www.fishersci.com and type BP2553-100 in the search box.

exACTGene DNA Ladders



EQUIPMENT: GEL BOXES



Submarine Electrophoresis Gel Boxes

Thermo Scientific submarine electrophoresis gel boxes are available in three sizes and sample quantities. This gel size flexibility is important in order to provide sufficient separation of closely spaced DNA fragments as well as small and large nucleic acid fragments on the same gel.

Thermo Scientific EasyCast® Systems contain:

- Compact footprint
- Thermo Scientific EasyCast gasket UVT gel tray
- Thermo Scientific SuperSafe® lid with attached power supply leads
- Optional cooling system

Thermo Scientific Millipede® Systems contain:

- Buffer chamber with 3-point leveling base
- Thermo Scientific SuperSafe lid with attached power supply leads
- UV transparent gel tray with gasket end gates
- Leveling bubble



EasyCast B1 Gel Box

| Cat. No. | Model | Cooling Option | Max. Number of Samples | Gel Size (W x L) in cm | Running Buffer Volume in mL |
|-------------|--------------------------------|----------------|--------------------------|------------------------|-----------------------------|
| 09-528-178 | Thermo Scientific EasyCast B1 | No | 34 | 9 x 11 | 600 |
| OWB1BP | Thermo Scientific EasyCast B1 | Yes | 34 | 9 x 11 | 600 |
| 09-528-110B | Thermo Scientific EasyCast B2 | No | 108 with rapid load tray | 12 x 14 | 800 |
| 09-528-118 | Thermo Scientific EasyCast B2 | Yes | 108 with rapid load tray | 12 x 14 | 800 |
| 09-528-124 | Thermo Scientific EasyCast B3 | Yes | 300 with rapid load tray | 12 x 14 | 1000 |
| OWA6 | Thermo Scientific Millipede A6 | No | 500 | 23 x 25 | 3000 |

Voltage Table

The table (below) provides recommended voltages and buffers according to DNA size and application. The distance used to determine the voltage gradients is the distance between electrodes, not the gel length. If the voltage is too high, band streaking may occur for large DNA sizes (>12kb). When the voltage is too low, the mobility of small (<1kb) DNA is reduced, and band broadening will occur due to dispersion and diffusion.

| Gel Size | Voltage | Recovery Buffer | Analytical Buffer |
|---------------|----------|-----------------|-------------------|
| <1kb | 5V/cm | TAE | TBE |
| <1kb to >12kb | 4-10V/cm | TAE | TBE |
| >12kb | 1-2V/cm | TAE | TAE |

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. or its subsidiaries.



For customer service
To fax an order
To order online

Fisher Scientific US
1-800-766-7000
1-800-926-1166
www.fishersci.com

Fisher HealthCare US
1-800-640-0640
1-800-290-0290
www.fisherhealthcare.com

Fisher Scientific Global
Export, Latin America
1-770-871-4725
1-770-871-4726
www.fishersci.com

Fisher Scientific
Canada
1-800-234-7437
1-800-463-2996
www.fishersci.ca

BN0608121 JJ/LY 5M-IW-06/12 12_5363

