# SPE Procedure – Six Steps for a Clean Extract

- 1. Sample Pre-treatment
- 2. Column Conditioning
- 3. Column Pre-equilibration
  - 4. Sample Loading
- 5. Wash Away Interferences
- 6. Elute Compounds of Interest

## 1. Sample Pre-treatment

It is important to optimize the sample for effective analyte retention. The following should be considered:

- Sample volume/analyte concentration/ matrix complexity
- · Adjust sample/matrix composition for proper dilution/ionic strength
- Sample pH for optimum retention
- Confirm that analytes are free in solution
- · Remove any unwanted particulates via filtration or centrifugation

#### 2. Column Conditioning

Prepare the sorbent for effective interaction(s) with the compounds of interest.

- Use appropriate solvent for column condition/activation
- · Prevent sorbent drying during conditioning

### 3. Column Pre-equilibration

Equilibrate with weakly eluting solvent to prepare the phase for sample addition.

- Use the same solvent as for sample pre-treatment
- Prevent sorbent drying during column equilibration

## 4. Sample Loading

Analytes are retained on the sorbent.

· Apply samples at appropriate flow rate (1mL/minute typical)

#### For Reversed-Phase Interactions

- Neutral compounds are not affected by pH
- For charged compounds, a pH at which the compound is not charged is used. Neutralize the molecule according to the following:
  - For basic compounds, the neutral molecule exists at least 2pH units below the pK<sub>a</sub> of the compound
  - For acidic compounds, the neutral molecule exists at least 2pH units above the pK<sub>a</sub> of the compound

#### For Normal-Phase Interactions

- pH is not normally an issue in normal phase interactions, as the solvents used are typically non-polar organic solvents, rather than water
- There is no need to verify the sample application pH

#### For Ion-Exchange Interactions

- pH and pK<sub>a</sub> are important considerations
- · Acidic compounds are extracted from a sample solution at least 2pH units above the pKa of the analyte
- Basic compounds are extracted from a sample solution 2 or more pH units below the pK<sub>a</sub> of the analyte
- For second (organic) wash, choose the strongest solution where no compound breakthrough occurs
- For elution step, use a solution stronger than where all the compound of interest is eluted
- · NB: when choosing these solutions allow some margin for error

### 5. Wash Away Interferences

Remove impurities bound less strongly than the compounds of interest.

- Select a strong enough wash solvent to remove interferences but weak enough to leave compounds of interest bound
- Selectively rinse away the less strongly bonded interferences
- · Wash solvent selected according to phase mechanism/analyte properties

## 6. Elute Compounds of Interest

Selectively recover the analyte(s) by disrupting the analyte-sorbent interaction.

- · Selectively elute analytes of interest using different solvents
- Smaller elution volume produces a more concentrated extract
- Select elution solvent that does not elute strongly retained impurities
- Select elution solvent according to phase mechanism/analyte properties

It is important to optimize the Wash and Elution steps in order to obtain maximum levels of recovery.

> Get your copy of our application notebook, featuring over 150 SPE applications www.thermoscientific.com/spe