Fisherbrand
QUALITY. RELIABILITY. VALUE.

Instruction Manual - Gebrauchsanweisung - Manuel d'utilisation - Manual de funcionamiento

Motorized Pipette Filler - Motorisierte Pipettierhilfe -
Pipeteur motorisé - Pipeta motorizada

15249805

Motorized Pipette Filler

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Pipeteador motorizado

Fisherbrand
QUALITY. RELIABILITY. VALUE.

INSTRUCTION MANUAL - GEBRAUCHS MANUAL DE FUNCIONAMIENTO

*Motorized Pipette Filler - Motorisierte Pipettierhilfe -
Pipeteur motorisé - Pipeta motorizada*

Table of content

1 - Introduction	5
2 - Work Safety Introductions	5
3 - Limitation of Use	5
4 - Operations Elements	5
5 - Pipetting	6
6 - Recharging the Lithium Battery	7
7 - Maintenance	8

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Fisherbrand
QUALITY. RELIABILITY. VALUE.

INSTRUCTION MANUAL - GEBRAUCHSANLEITUNG MANUAL DE FUNCIONAMIENTO

*Motorized Pipette Filler - Motorisierte Pipettierhilfe -
Pipeteur motorisé - Pipeta motorizada*

15249805

Motorized Pipette

1 - Introduction

The Fisherbrand Motorized Pipette aid is a device pipetting liquids with the use of measuring pipettes. It can work with all types of glass or plastic pipettes in the volume range from 0.1 mL to 100 mL. The Fisherbrand Motorized Pipette is a general purpose laboratory device used to transfer or measure for pipetting safely and accurately.

2 - Work Safety Introductions

Before starting the work with the Fisherbrand Motorized Pipette every user should read these operating instructions carefully. Using the device inconsistently with the operating instructions may result in damaging the device. The device should be maintained only at an authorized service center otherwise the manufacturer will be relieved from any liability under the warranty.

During the work with the Fisherbrand Motorized Pipette general safety regulations regarding risks related with laboratory work should be observed. Protective clothing, goggles and gloves should be worn.


- The Fisherbrand Motorized Pipette shall be used only for measuring liquids in conditions specified by the manufacturer, which are limited due to the chemical and mechanical resistance of the device as well as the user safety.
- The pipetting aid should not be used in an environment where explosion risk is present.
- The information and instructions provided by the manufacturers of the reagents must be observed.
- Only original spare parts and accessories, recommended by the manufacturer shall be used.
- Only the original charger, supplied by the manufacturer, shall be used for charging the lithium battery.
- In case of incorrect functioning of the pipetting aid, work shall be stopped for charging the lithium battery.
- In case of incorrect functioning of the pipetting aid, work shall be stopped. The device shall be cleaned in accordance with the operating instructions and/or sent to an authorized service center.
- In the case of mechanical damage to the casing the device shall be immediately sent for repair to an authorized service center.
- The use of excessive force during work shall be avoided.

3 - Limitation of Use

- The pipetting aid may not be used in an environment where explosion risk is present. Flammable liquid shall not be measured in particular substances with flash point below 3°C. (either, acetone.)
- The pipetting aid shall not be used for drawing acids with a concentration above 1 mol.

- The pipetting aid shall not be used for drawing solutions with the temperature above 50 °C.
- The pipetting aid may work in temperature range from +10°C. to +35°C.

4 - Operations Elements

Assembly

- Grip the lithium battery plug by the wire, confirm alignment of the plug and insert it in the socket. Push it firmly into the socket until you can hear it snap into place. Place lithium battery in compartment.
- Select lithium battery compartment cover.
 - a) Smoothly rounded
 - b) With ridges for resting the instrument inverted on the bench
- Close the lithium battery compartment.
- Recharging the lithium battery charging time is 2 hours.

Wall Support/ Invert Bench placement Mounting the wall support

The wall support is attached with snap tapes to permit easy removal for cleaning purposes. To use the snap tapes, the mounting surface must be smooth, clean and grease free. Prepare the surface for mounting by cleaning with a 50% solution of isopropanol/water. Alternatively, the wall support can be attached to the wall with two screws.

Invert Bench Placement

If the cover with ridges is mounted, instrument can be placed inverted on the bench.

Select Mode.

- Switch up: Gravity-delivery
- Switch down: blow out (with motor power)
- Adjust motor speed
Maximum motor speed is continuously adjustable from maximum (+) to minimum (-) by turning the thumbwheel. We recommend using the maximum speed setting for larger volume pipettes and a slower speed setting for lower volume pipettes.
- Pipetting Speed Control
The amount of pressure on the pipetting button provides continuously variable speeds up to the rate set by the speed

controller. Press upper pipetting button for filling and lower button for dispensing.

5 - Pipetting

- Attach the Pipette
 - Hold the pipette as near to its upper end as possible, and carefully insert into the adapter until it fits tightly.
 - Once the pipette has been securely attached, always hold the instrument in a vertical position tip down.

- Attention!
 - Be sure the pipette fits tightly into the adapter. Never use force. Thin pipettes are particularly liable to break. Avoid the risk of injury.
- Filling
 - Immerse the pipette tip into the liquid. Slowly press the upper pipetting button and fill the pipette so that the meniscus is slightly above the mark desired.

- Attention!
 - Do not over fill the pipet.
 - Adjust the volume
 - Use suitable lint free tissue to wipe the pipette tip. Slowly press the lower button and dispense liquid until the meniscuses

is adjusted exactly to the desired volume.

- Dispensing

-Slowly press the lower pipette button. The dispensing speed is dependent on the amount of pressure on the pipetting button as well as the setting of speed controller. When the mode switch is in the motor powered delivery position.

6 - Recharging the Lithium Battery

- Attention!

- Use original manufacture AC adapter only. Using a different AC adapter can damage both the instrument and the lithium battery.

- Green LED flashes : lithium battery charge is sufficient
- Red LED flashes : remaining lithium battery charge is sufficient maximum of two hours of pipetting; charging of lithium battery is recommended
- LED is not lighting : finish pipetting activity and charge lithium battery.
- Recharging procedure
 1. Temperature range for recharging: +10° C to +35°C (+50°F to 95°F)
 2. When the AC adapter is initially connect the red LED on the instrument light continuously.
 3. After 2 hours of charging the AC adaptor switches to a pulsed, long term charging mode (the green lighting)


Note:

- The charging system is design to prevent the lithium battery from over chagrining and minimizes the lazy battery effect. The maintain maximum lithium battery capacity, it is best to charge the battery only when the slowly flashing or quickly flashing LED indicates charging is needed.
- The instruments be used during recharging.

- To open the lithium battery compartment press at the textured surface and simultaneously slide cover back
- Remove the lithium battery, grip lithium battery plug by the wire and pull it gently out of the socket.
- Grip the plug of the new lithium battery by the wire and push it firmly into the socket until you can hear it is snapping in.
- Insert new battery and close lithium battery compartment.

7 - Maintenance

- Cleaning
 1. Pull out the pipette.
 2. Unscrew the adapter support.
 3. Pull out the membrane filter.
 4. Remove the adapter from the adapter support by pushing it upwards. The check valve remain in the adapter.
 5. To clean adapter with built in check valve, carefully rinse it by using an appropriate compatible solvent (e.g. water) in a wash bottle. Blow out the liquid and allow drying completely.
 6. Place the new membrane filter with the thicker end facing downwards into the adapter (with the smaller conical side up).
 7. The entire pipette adapter including membrane filter is autoclavable at 121 °C.
 8. Reassemble the instrument in reverse order and carry out leak test.

- Storage
The Fisherbrand Motorized Pipette shall be stored in a dry place. The allowable storage temperature: -20°C to +50°C.
- Note:
Do not store the pipetting aid with a filler pipette

Fisherbrand[®]
QUALITY. RELIABILITY. VALUE.

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© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.