15361901

Ultra Thin Magnetic Stirrer

Magnetrührer ultraflach

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Agitador magnético ultrafino

Table of content

Product description	5
Features	5
Technical Specifications	5
Maintenance and Safety Instructions	6
Installation and Operation	7
How to set speed ?	7
Pulse Mode	8
Stir Bar Docking	8
Warranty Statement	8

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Ultra Thin Magnetic Stirrer

PRODUCT DESCRIPTION

Congratulations! You're using an ultra thin and compact magnetic stirrer designed with modern magnetic coil technology. This is a maintenance free product with no moving parts which complies with IP protection class 65 (IP 65).

This magnetic stirrer is designed for use in laboratories, pharmacies and schools for mixing liquids. For better mixing, a pulse mode is provided in which the stir bar alternates rotating clockwise and counter clockwise every 30 seconds.

FEATURES

1. Flat surface
2. Compact and thin in size
3. Pulse function for better mixing
4. Variable speed setting from 15 to 1500 RPM
5. Indicator of ongoing speed during operation
6. Anti-skid pads
7. Water spill resistant, compliant with IP 65 standard
8. Stir bar docking location to place or attach magnetic stir bar
9. Wear-free drive with no moving parts
10. Setup plate and casing made from chemically resistant materials
11. Simple control panel with convenient and easy to use interface

TECHNICAL SPECIFICATIONS

Maximum stirring volume (H ₂ O)	800mL
Maximum stirring capability (tested in H ₂ O)	500mL at 1500 rpm with 25mm stir bar
Speed range	15 – 1500 rpm
Max. stirring bar length	25mm
Set-up plate material	Polyester
Set-up plate dimensions	Ø94mm
Dimensions (L x W x H)	185 x 115 x 13mm
Weight	0.4 – 0.5 Kg
Permissible ambient temperature	5 – 40°C
Protection class according to DIN EN 60529	IP 65
Input voltage	100 – 230 V
Frequency	50/60 Hz
Power	5W
Auto direction change/Pulse Mode	Yes
Auto direction change time interval	30 seconds
Current speed indicator	Yes

MAINTENANCE AND SAFETY INSTRUCTIONS

Read the following safety instructions before using the magnetic stirrer for your safety as well as equipment safety.


1. Do not place any steel or magnetic material on the base plate except the recommended magnetic stir bar with a beaker or flask in between. This can affect the magnetism of device.
2. Place device on an even, stable, clean, non-slippery and fire-proof surface.
3. The base plate must be clean and intact.
4. Do not try to mix more than 800mL of liquid with 25mm magnetic stir bar. It may lead to improper mixing.
5. Wear your personal protective equipment in accordance with the hazardous category of the media to be processed. Otherwise there is a risk from:
 - Splashing and evaporation of liquids.
 - Ejection of parts.
 - Release of toxic or combustible gases.
6. Do not lift or hold the device by the AC adaptor supply attached with the unit.
7. Do not use a damaged beaker, flask, stir bar or any other damaged component during operation. It may affect the efficiency of the device.
8. Reduce the speed:
 - If the medium splashes out of the vessel because speed is too high.
 - If the appliance is not running smoothly.
 - If the container moves on the base plate.
9. The base plate may heat up after continuous usage so do not touch the base plate immediately after operation.
10. Process pathogenic materials only in closed vessels under a suitable extractor hood. Do not operate the appliance in explosive atmospheres, with hazardous substances.
11. Do not use liquids which may react dangerously to the extra energy produced during operation of the device. This also applies to any extra energy produced in other ways. For example: Through light irradiation, through surrounding temperature, etc.
12. Protect the device from collision.
13. The temperature of the stirring media used in the beaker and the beaker itself should not be more than 70°C.
14. Abrasion of the dispersion equipment or the rotating accessories can get into medium you are working on.
15. A chemical reaction could occur between the PTFE stir bar and any molten or dissolved alkaline earth metals, as well as with fine partied metal powders from groups 2 and 3 on the periodic table at temperatures above 300-400°C. Only elementary fluorine, chlorine trifluoride and alkaline metals do attack PTFE, halogen hydrocarbons have a reversible swelling effect.
16. Wear protective gloves while cleaning the device.
17. Use a clean cloth and the following substances to clean the device:
 - To clean dyes → Isopropyl alcohol
 - To clean construction materials → Water containing tensile and Isopropyl alcohol
 - To clean cosmetics → Water containing tensile and Isopropyl alcohol
 - To clean foodstuffs → Water containing tensile
 - To clean fuels → Water containing tensile

INSTALLATION AND OPERATION

The magnetic stirrer is provided in the box. Open the box, then remove the packaging and gently place the device on a firm and level surface. Take care while unpacking and do not pull the AC power adaptor as it is attached with the device.

The user manual and accessories should be kept with the device.

Please keep all packaging in safe storage for at least one year for warranty purposes.

ON/OFF button

Press button to turn ON and turn OFF the unit

Pulse Mode button

Press button to rotate stirring bar in CW and CCW directions every 30 secs

Increase Speed button

Press button to increase speed value

Decrease Speed button

Press button to decrease speed value

1. Plug in the AC power adapter.
2. Fill the beaker or vessel with the medium or liquids to be processed then place the magnetic stirring bar in the beaker, flask or vessel.
3. Place the vessel or beaker on the center of the base plate.
4. Turn ON the stirrer by pressing the ON/OFF button.
5. Set the speed at which the medium or liquids are to be processed. Minimum speed is 15 RPM and maximum speed is 1500 RPM.

How to set speed ?

- Press (+) button to increase speed and press (-) button to decrease speed.
- Pressing the (+) or (-) button will increase or decrease the speed by 50 RPM on each press.
- First set the speed to any of the 5 speed values given in the control panel. Then increase or decrease the speed value to change set to the required speed.
- For Example: Setting the speed to 750 RPM.
 1. First set the speed to 500 RPM. The speed is set to 500 RPM when the LED adjacent to 500 RPM blinks and the LED below 500 RPM is constantly lit.
 2. Now to set the speed to 750 RPM press the (+) button 5 times. As each press will increase the speed by 50 RPM and 5 x 50 RPM = 250 RPM.

- Therefore $500 + 250 = 750$ RPM.
- A speed set to 750 RPM will be shown when both the LEDs adjacent to 1000 RPM and 500 RPM are blinking i.e. from the given LEDs, the LED above the required speed and the LED below the required speed will blink and the LEDs below 500 RPM will be constantly lit.

Pulse Mode

In pulse mode, the stirring direction is reversed every 30 seconds. To select pulse mode which changes the stirring direction at regular intervals, press the pulse button.

First set the required speed and then press the pulse button to start the stirring process in pulse mode. The LED adjacent to the pulse button will blink to indicate pulse mode is active.

Press the ON/OFF button again to stop the ongoing operation as well as to switch OFF the lab disc.

Stir Bar Docking

A stir bar docking location is provided as shown in the figure above. Due to magnetism in both, the stir bar and the docking place, the magnetic stir bar will stick to that place as shown in above figure. This helps users keep the stir bar with the unit and makes it easy to find for future operation.

WARRANTY STATEMENT

This product is warranted to be free from defects in material and workmanship for a period of one (1) year from the date of purchase. Your product will be duly repaired upon prompt notification in compliance with the following conditions:

- This warranty is valid if the product is used for its intended purpose and within the guidelines specified in the supplied user manual.
 - This warranty does not cover damage caused by accident, neglect, misuse, improper service, natural forces or other causes not arising from defects in original material or workmanship.
 - This warranty does not cover any other damages from the use of this product.
 - This warranty does not cover damage to paint or finish and defects or damages caused by physical and chemical abuse or normal wear and tear.
 - The warranty is invalidated by any non-factory modification, which will immediately terminate all manufacturer's liabilities for the products or damages caused by its use.
 - The customer shall be responsible for the product or use of products as well as any supervision required for safety.
-
- All items returned for the service should be sent in the original packaging or other suitable carton, padded to avoid damage.
 - Under no circumstances shall Fisher Scientific be liable for indirect, consequential or special damages of any kind.
 - Fisher Scientific reserves the right to make technical changes without prior notice.
 - For your reference, make a note of the serial number and date of purchase.
 -

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