

thermo scientific



**Thermo Scientific™**  
Nunc Serological Pipettes

Research with Confidence

Find it at [eu.fishersci.com](http://eu.fishersci.com)

 **fisher scientific**  
part of Thermo Fisher Scientific

# Thermo Scientific Nunc Serological Pipettes

Our newly expanded portfolio of serological pipettes features a full range of size and packaging options – offering the quality, accuracy and performance today’s scientists require during every stage of cell culture research.

**Cell isolation  
and preparation**



**Characterization  
and Analysis**



We are dedicated to meeting the rigorous demands of your research, at every stage of cell culture. Our new collection of serological pipettes offers the quality assurance, size variability and packaging options you need to advance your research with complete confidence.

**Growth  
and passage**



## Accurate, disposable plastic pipettes plugged and sterilized

- Sterility assurance level (SAL) of  $10^{-6}$
- 100% virgin polystyrene for maximum clarity
- Certified non-pyrogenic
- Bold clear graduations with easy-to-read scale markings
- Non-dripping tip for pipette precision
- Color-coded packaging for ease in sorting and selecting the correct size

## Packaging options that accommodate every lab!

- Individually wrapped; paper/plastic or plastic film; peel or pop through for easy opening
- Bulk packed for frequent use; reduces packaging waste
- Multiple sizes for use in a wide range of applications

### Storage



### Harvest



## Serological Portfolio Highlights

- Use new sleek shortie pipettes for easy maneuvering under the hood
- Choose multiple sizes, ranging from 1 mL to 50 mL, all with non-dripping tips and clearly marked graduations for precision
- Pick bulk packed options for higher volume use with less packaging waste
- Select individually wrapped pipettes in either paper/plastic peel or plastic film
- See the difference with new ultra clear, 100% virgin polystyrene



thermo scientific

# Thermo Scientific Nunc Serological Pipettes

## Ordering Information

Cat. No.	Description	Unit per bag	Unit per case
Individually wrapped, paper peel packaging, plugged			
13450287	1 ml pipette	100	1000
13470287	2 ml pipette	100	500
13430287	5 ml pipette	50	200
13440287	10 ml pipette	50	200
13424776	25 ml pipette	50	200
13420297	50 ml pipette	25	100
Shortie, Individually wrapped paper/plastic peel packaging, plugged			
15101539	5 ml shortie pipette	50	200
15191529	10 ml shortie pipette	50	200
Individually wrapped, plastic film packaging, plugged			
13480287	1 ml pipette	100	1000
13490287	2 ml pipette	100	500
13400297	5 ml pipette	50	200
13410297	10 ml pipette	50	200
13430297	25 ml pipette	50	200
13440297	50 ml pipette	25	100
Pipettes, Bulk Packed, Plugged			
13450297	1 ml pipette	50	1000
13460297	2 ml pipette	50	500
13470297	5 ml pipette	25	500
13480297	10 ml pipette	25	500
13490297	25 ml pipette	25	200
13400307	50 ml pipette	25	100



### Thermo Scientific S1 Pipet Filler

A lightweight, cordless pipet filler that allows longer, fatigue-free pipetting—an excellent companion to our serological pipettes.

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at **Belgium:** fishersci.be **Denmark:** fishersci.dk  
**Germany:** fishersci.de **Ireland:** fishersci.ie **Italy:** fishersci.it  
**Finland:** fishersci.fi **France:** fishersci.fr **Netherlands:** fishersci.nl  
**Norway:** fishersci.no **Portugal:** fishersci.pt **Spain:** fishersci.es  
**Sweden:** fishersci.se **Switzerland:** fishersci.ch **UK:** fishersci.co.uk

© 2022 Thermo Fisher Scientific Inc. All rights reserved.  
 Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**  
 part of Thermo Fisher Scientific