



nl new lab
su start-up
programme

Fisherbrand Favourites
2016/2017

Redeeming your voucher is **EASY!**

Place orders on your New Lab Start-Up account number



You have six months from the date of your first order to hit one of the three spend thresholds



We will email you a unique Fisherbrand voucher code



Redeem your voucher online at



eu.fishersci.com/go/nlsu

We have selected some of our top selling Fisherbrand consumable items so you can really make the most out of your Fisherbrand voucher.


*Terms and conditions apply, see back page

| | |
|-----------------------|---------|
| Autoclave Accessories | 35 |
| Cell Culture | 35 |
| Centrifugation | 36 → 38 |
| Cryogenics | 39 |
| Filtration | 40 |
| Glassware | 41 → 43 |
| Weighing Accessories | 44 |
| Liquid Handling | 46 → 48 |
| Microscopy | 49 |
| Plasticware | 49 → 53 |
| pH Testing | 53 |
| Safety | 54 → 57 |
| Thermometers & Timers | 57 → 59 |

Spend **£20,000** receive a further **£500** Fisherbrand voucher

Spend **£10,000** receive a further **£500** Fisherbrand voucher

Spend **£2,000** receive a **£500** Fisherbrand voucher



Autoclave bags

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11600312 | Autoclave bag, printed, biohazard, 40µm thickness, opening on width, polypropylene, autoclavable at high temperature 134°C (max 140°C), clear, 310mm x 660mm | 100 |
| 11553342 | Autoclave bag, printed, biohazard, 40µm thickness, opening on width, polypropylene, autoclavable at high temperature 134°C (max 140°C), clear, 410mm x 630mm | 100 |
| 11563342 | Autoclave bag, printed, biohazard, 40µm thickness, opening on width, polypropylene, autoclavable at high temperature 134°C (max 140°C), clear, 620mm x 780mm | 100 |
| 11583342 | Autoclave bag, plain unprinted, biohazard, 40µm thickness, opening on width, polypropylene, autoclavable at high temperature 134°C (max 140°C), clear, 410mm x 630mm | 100 |
| 11593342 | Autoclave bag, plain unprinted, biohazard, 40µm thickness, opening on width, polypropylene, autoclavable at high temperature 134°C (max 140°C), clear, 620mm x 780mm | 100 |



Cell Culture

Cell strainers

- Cell culture strainer, sterile, individually packed, polypropylene with nylon mesh

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11587522 | For use with 50mL conical tubes, 40µm | 50 |
| 11597522 | For use with 50mL conical tubes, 70µm | 50 |
| 11517532 | For use with 50mL conical tubes, 100µm | 50 |



Flasks, borosilicate glass, conical, narrow neck, ISO 1773

- Erlenmeyer, graduated

| Cat. No | Capacity, mL | Height, mm | Widest OD, mm | OD. Neck, mm | Pack qty |
|----------|--------------|------------|---------------|--------------|----------|
| 15499093 | 50 | 90 | 51 | 22 | 10 |
| 15409103 | 200 | 135 | 79 | 34 | 10 |
| 15419103 | 250 | 145 | 85 | 34 | 10 |
| 15429103 | 300 | 160 | 87 | 34 | 10 |
| 15439103 | 500 | 180 | 105 | 34 | 10 |
| 15449103 | 1,000 | 220 | 131 | 42 | 1 |
| 15459103 | 2,000 | 280 | 166 | 50 | 1 |
| 15469103 | 3,000 | 310 | 187 | 50 | 1 |
| 15479103 | 5,000 | 365 | 220 | 50 | 1 |



Erlenmeyer flask, plastic

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11765253 | Shaker flask, Erlenmeyer, plain bottom, vented, non-pyrogenic, non-cytotoxic, autoclavable, polycarbonate, PP closure, moulded in graduations, for suspension cultures, 250mL | 12 |




Petri dishes, plastic, round

- Inner packs of 20

| Cat. No | Ø x H, mm | Vent | Pack qty |
|----------|-----------|------|----------|
| Aseptic | | | |
| 12694785 | 90 x 14.2 | 3 | 600 |
| 12664785 | 90 x 16.2 | 3 | 600 |
| 12644785 | 90 x 16.2 | 0 | 600 |
| 12654785 | 90 x 16.2 | 1 | 600 |
| Sterile | | | |
| 12604795 | 90 x 14.2 | 3 | 600 |
| 12674785 | 90 x 16.2 | 3 | 600 |
| 12684785 | 90 x 16.2 | 1 | 600 |





Tubes, microtube


- Use between -80°C to 120°C
- Flat cap
- Withstand speeds up to 30,000xg
- DNase, RNase free
- Polypropylene
- Graduation moulded
- Autoclavable
- Non-sterile

| Cat. No | Capacity, mL | Colour | Dia. x height, mm | Max. rcf (xg) | Pack qty |
|--|--------------|---------|-------------------|---------------|----------|
| 1 Graduated microtube | | | | | |
| 11916955 | 0.6 | Natural | 8 x 30 | 30,000 | 500 |
| 11926955 | 1.5 | Natural | 11 x 40 | 26,000 | 500 |
| 11393613 | 2.0 | Natural | 11 x 40 | 25,000 | 500 |
| 2 Graduated safelock microtube | | | | | |
| 11976955 | 0.6 | Natural | - | 30,000 | 1,000 |
| 11706467 | 1.5 | Natural | - | 26,000 | 500 |
| 11966955 | 1.5 | Mix | - | 26,000 | 500 |
| 11956955 | 2.0 | Natural | - | 25,000 | 500 |
| 3 Graduated low binding microtube | | | | | |
| 11996955 | 0.6 | Natural | 10 x 30 | - | 500 |
| 11986955 | 1.5 | Natural | 13 x 40 | - | 250 |
| 11906965 | 2.0 | Natural | 13 x 40 | - | 250 |



Microcentrifuge tubes with locking snap cap

- Polypropylene (PP)
- Graduated
- Flat and pierceable cap
- Locking lid to secure the sample and avoid undesired opening
- Marking area
- Autoclavable and can withstand temperatures to -86°C
- DNase/RNase and pyrogen free
- Easy one-hand use
- Applications include:
 - Sample preparation
 - Sample storage and transportation
 - Centrifugation

Locking snap cap "hook" style



| Cat. No | Capacity, mL | Colour | Bottom shape | Max. Rcf (xg) | Diameter, mm | Length, mm | Sterile | Pack qty |
|-------------------|--------------|---------|--------------|---------------|--------------|------------|---------|----------|
| 1 15422545 | 0.5 | Natural | Conical | 30,000 | 10 | 30 | No | 500 |
| 2 15432545 | 1.5 | Natural | Conical | 26,000 | 13 | 39 | No | 500 |
| 3 15442545 | 2.0 | Natural | Round | 25,000 | 13 | 40 | No | 500 |



Tubes, centrifuge, 15mL and 50mL

- 15mL and 50mL capacity
- Available in PP and PET
- Black graduation
- White marking area
- Flat cap and plug seal cap
- Gamma irradiation sterilisation
- Bulk or rack version

Bulk

| Cap flat top | Cap plug seal | Capacity, mL | Max. rcf, xg | Material | Dia. x H, mm | Sterile | Inner pack qty | Pack qty |
|--------------|---------------|--------------|--------------|----------|--------------|---------|----------------|----------|
| 11755075 | 11765075 | 15 | 6,000 | PP | 17 x 119 | Yes | 25 | 500 |
| 11512303 | 11809650 | 50 | 9,400 | PP | 29 x 114 | Yes | 25 | 500 |
| - | 11829650 | 50 | 9,400 | PP | 29 x 114 | No | 25 | 500 |

Rack

| Cap flat top | Cap plug seal | Capacity, mL | Max. rcf, xg | Material | Dia. x H, mm | Sterile | Inner pack qty | Pack qty |
|--------------|---------------|--------------|--------------|----------|--------------|---------|----------------|----------|
| - | 11879640 | 15 | 1,800 | PET | 17 x 119 | Yes | 50 | 500 |
| 11849650 | 11889640 | 15 | 6,000 | PP | 17 x 119 | Yes | 50 | 500 |
| - | 11839650 | 50 | 1,800 | PET | 29 x 114 | Yes | 50 | 500 |
| 11819650 | 11899640 | 50 | 9,400 | PP | 29 x 117 | Yes | 50 | 500 |

Tube racks, four way

- Each rack can hold 4 x 50mL conical tubes, 12 x 15mL conical tubes, 32 x 1.5mL microtubes or 32 x 0.5mL microtubes
- Dimensions (l x w x h), mm: 174 x 95 x 52
- Autoclavable

| Cat. No | Description | Pack qty |
|----------|---|----------|
| 11750334 | Four way tube rack, natural | 5 |
| 11770045 | Four way tube rack, blue | 5 |
| 11760334 | Four way tube rack, green | 5 |
| 11780045 | Four way tube rack, pink | 5 |
| 11770334 | Four way tube rack, yellow | 5 |
| 11790045 | Four way tube rack, orange | 5 |
| 1170055 | Four way tube rack, assorted colours (blue, green, pink, yellow and orange) | 5 |



Tube racks, microcentrifuge, polypropylene, 80 well

- Polypropylene 80 well microtube racks ideal for freezer storage
- Racks hold 1.5mL to 2mL microtubes in a 5 x 16 array. Tubes sit in large wells with conical bottoms
- Dimensions (l x w x h), mm: 225 x 67 x 28
- Autoclavable

| Cat. No | Colour | Pack qty |
|----------|---------------------------------------|----------|
| 11720344 | Natural | 5 |
| 11908084 | Orange | 5 |
| 11968074 | Blue | 5 |
| 11978074 | Green | 5 |
| 11988074 | Pink | 5 |
| 11998074 | Yellow | 5 |
| 11918084 | Red | 5 |
| 11728084 | Assorted colours (one of each colour) | 5 |



Tube racks, adapt-a-rack™

- Each slot adapts to accommodate 5 to 50mL tubes sizes and holds them in a firm upright position
- Made from POM
- Racks can be linked together for additional holding capacity
- Open sides allow for label viewing
- Flat base to hold tubes upright
- Linked racks should be supported underneath when carrying
- Available in a solid or a combination of colours for easy coding identification
- Autoclavable

| Cat. No | Accommodates | Wells | Well diameter, mm | Diameter, mm (l x w x d) | Colour | Pack qty |
|----------|-----------------|-------|-------------------|--------------------------|---------------|----------|
| 15300370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | Blue/green | 2 |
| 15320370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | Blue/yellow | 2 |
| 15330370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | Pink/white | 2 |
| 15340370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | Blue/blue | 2 |
| 15350370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | White/white | 2 |
| 15360370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | Pink/pink | 2 |
| 15380370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | Yellow/yellow | 2 |
| 15390370 | 5 to 50mL tubes | 4 | 12 to 30 | 181 x 56 x 76 | Green/green | 2 |





Test tube racks


- Made from Delrin™ fibreglass-reinforced polyoxymethylene (POM)
- Tough and more highly resistant to acids, bases, solvents and heat
- Offers the highest chemical and heat resistance (-40°C to +140°C)
- Available in half and full rack sizes for five different tube diameters
- Moulded in a single, continuous piece so no assembly required
- Autoclavable

| Cat. No | Tubes | Wells | Well diameter mm | L x W x H, mm | Colour | Pack qty |
|--------------------------------|------------|-------|------------------|----------------|---------|----------|
| Full size | | | | | | |
| 15340380 | 5 to 10 mL | 72 | 13 | 104 x 202 x 59 | White | 1 |
| 15390380 | 5 to 10 mL | 72 | 16 | 127 x 250 x 70 | White | 1 |
| 15340390 | 10 to 18mL | 40 | 20 | 100 x 252 x 83 | White | 1 |
| 15370390 | 10 to 18mL | 40 | 25 | 120 x 300 x 92 | White | 1 |
| 15310400 | 15 to 50mL | 24 | 30 | 110 x 282 x 85 | White | 1 |
| 15350380 | 5 to 10 mL | 72 | 13 | 104 x 202 x 59 | Blue | 1 |
| 15300390 | 5 to 10 mL | 72 | 16 | 127 x 250 x 70 | Blue | 1 |
| 12327629 | 5 to 30mL | 40 | 20 | 100 x 252 x 83 | Blue | 1 |
| 15380390 | 5 to 30mL | 40 | 25 | 120 x 300 x 92 | Blue | 1 |
| 10257963 | 15 to 50mL | 24 | 30 | 110 x 282 x 85 | Blue | 1 |
| 15360380 | 5 to 10 mL | 72 | 13 | 104 x 202 x 59 | Green | 1 |
| 15310390 | 5 to 10 mL | 72 | 16 | 127 x 250 x 70 | Green | 1 |
| 15350390 | 10 to 18mL | 40 | 20 | 100 x 252 x 83 | Green | 1 |
| 15390390 | 10 to 18mL | 40 | 25 | 120 x 300 x 92 | Green | 1 |
| 15330400 | 15 to 50mL | 24 | 30 | 110 x 282 x 85 | Green | 1 |
| 15370380 | 5 to 10 mL | 72 | 13 | 104 x 202 x 59 | Magenta | 1 |
| 15320390 | 5 to 10 mL | 72 | 16 | 127 x 250 x 70 | Magenta | 1 |
| 15360390 | 10 to 18mL | 40 | 20 | 100 x 252 x 83 | Magenta | 1 |
| 15300400 | 10 to 18mL | 40 | 25 | 120 x 300 x 92 | Magenta | 1 |
| 15350400 | 15 to 50mL | 24 | 30 | 110 x 282 x 85 | Magenta | 1 |
| 15380380 | 5 to 10 mL | 72 | 13 | 104 x 202 x 59 | Yellow | 1 |
| 15330390 | 5 to 10 mL | 72 | 16 | 127 x 250 x 70 | Yellow | 1 |
| 15360400 | 15 to 50mL | 24 | 30 | 110 x 282 x 85 | Yellow | 1 |
| Half size | | | | | | |
| 15390400 | 5 to 10 mL | 36 | 13 | 104 x 104 x 59 | White | 1 |
| 15340410 | 5 to 10 mL | 36 | 16 | 127 x 127 x 70 | White | 1 |
| 15310420 | 10 to 18mL | 20 | 20 | 100 x 127 x 83 | White | 1 |
| 15380420 | 10 to 18mL | 16 | 25 | 120 x 122 x 92 | White | 1 |
| 15330430 | 15 to 50mL | 9 | 30 | 110 x 110 x 85 | White | 1 |
| 15300410 | 5 to 10 mL | 36 | 13 | 104 x 104 x 59 | Blue | 1 |
| 15360410 | 5 to 10 mL | 36 | 16 | 127 x 127 x 70 | Blue | 1 |
| 15340420 | 10 to 18mL | 20 | 20 | 100 x 127 x 83 | Blue | 1 |
| 15300430 | 10 to 18mL | 16 | 25 | 120 x 122 x 92 | Blue | 1 |
| 12317629 | 15 to 50mL | 9 | 30 | 110 x 110 x 85 | Blue | 1 |
| 15310410 | 5 to 10 mL | 36 | 13 | 104 x 104 x 59 | Green | 1 |
| 15380410 | 5 to 10 mL | 36 | 16 | 127 x 127 x 70 | Green | 1 |
| 15350420 | 10 to 18mL | 20 | 20 | 100 x 127 x 83 | Green | 1 |
| 15310430 | 10 to 18mL | 16 | 25 | 120 x 122 x 92 | Green | 1 |
| 15340430 | 15 to 50mL | 9 | 30 | 110 x 110 x 85 | Green | 1 |
| 15330410 | 5 to 10 mL | 36 | 13 | 104 x 104 x 59 | Magenta | 1 |
| 15390410 | 5 to 10 mL | 36 | 16 | 127 x 127 x 70 | Magenta | 1 |
| 15370420 | 10 to 18mL | 20 | 20 | 100 x 127 x 83 | Magenta | 1 |
| 15320430 | 10 to 18mL | 16 | 25 | 120 x 122 x 92 | Magenta | 1 |
| 15360430 | 15 to 50mL | 9 | 30 | 110 x 110 x 85 | Magenta | 1 |
| Polypropylene full size | | | | | | |
| 15370430 | 5 to 10mL | 72 | 13 | 104 x 202 x 59 | White | 8 |
| 15310440 | 5 to 10mL | 72 | 16 | 127 x 250 x 70 | White | 8 |
| 15380430 | 5 to 10mL | 72 | 13 | 104 x 202 x 59 | Blue | 8 |
| 15320440 | 5 to 10mL | 72 | 16 | 127 x 250 x 70 | Blue | 8 |
| 15390430 | 5 to 10mL | 72 | 13 | 104 x 202 x 59 | Green | 8 |
| 15340440 | 5 to 10mL | 72 | 16 | 127 x 250 x 70 | Green | 8 |
| 15300440 | 5 to 10mL | 72 | 13 | 104 x 202 x 59 | Magenta | 8 |
| 15360440 | 5 to 10mL | 72 | 16 | 127 x 250 x 70 | Magenta | 8 |

Cryoboxes, polypropylene, 81 place, with grid inserts

- For sample storage and transportation
- Autoclavable, temperature resistant from -90°C to +121°C
- High quality polypropylene
- Available in natural and various colours
- 81 places, 9 x 9 grid, numerically coded
- Robust hinge with safe snap-on lid
- Stackable for safe transportation
- Drain bores at the bottom for dew liquids
- Removable grid inserts
- Can also be used for larger vessels by removing the grid inserts

| Cat. No | Accommodates | Colour | Pack qty |
|--|---------------------------|---------|----------|
| Dimensions (l x w x h), 133mm x 133mm x 52mm, 81 place, 9 x 9 grid | | | |
| 11856893 | 81 x 1.2mL to 2.0mL vials | Natural | 5 |
| 11826903 | 81 x 1.2mL to 2.0mL vials | Yellow | 5 |
| 11836903 | 81 x 1.2mL to 2.0mL vials | Red | 5 |
| 11856903 | 81 x 1.2mL to 2.0mL vials | Blue | 5 |
| Dimensions (l x w x h), 133mm x 133mm x 75mm, 81 place, 9 x 9 grid | | | |
| 11866903 | 81 x 3.0mL to 5.0mL vials | Natural | 5 |
| 11876903 | 81 x 3.0mL to 5.0mL vials | Yellow | 5 |
| 11886903 | 81 x 3.0mL to 5.0mL vials | Red | 5 |
| 11896903 | 81 x 3.0mL to 5.0mL vials | Green | 5 |
| 11806913 | 81 x 3.0mL to 5.0mL vials | Blue | 5 |



Cryoboxes, polypropylene, flat-pack

- Cryoboxes are supplied flat-pack eliminating the hassle of stacking and reduces the amount of space needed for storage
- Boxes fit in standard freezer racks
- Easy to assemble
- More durable than cardboard boxes as they are not susceptible to moisture
- Significantly reduces the chance of mould
- Dividers included with the boxes
- Store down to temperatures of -80°C

| Cat. No | Accommodates | Wells | Well diameter, mm | Dimensions, mm (l x w x h) | Colour | Pack qty |
|---|-----------------|-------|-------------------|----------------------------|---------|----------|
| Corrugated polypropylene 0.2/0.5/1.5 to 2.0/5mL | | | | | | |
| 15316528 | 0.2mL tubes | 144 | 7.9 | 133 x 130 x 30 | Natural | 10 |
| 15326528 | 0.2mL tubes | 144 | 7.9 | 133 x 130 x 30 | Blue | 10 |
| 15336528 | 0.2mL tubes | 144 | 7.9 | 133 x 130 x 30 | Red | 10 |
| 15346528 | 0.2mL tubes | 144 | 7.9 | 133 x 130 x 30 | Purple | 10 |
| 15356528 | 0.5mL tubes | 81 | 12.0 | 133 x 130 x 42 | Natural | 10 |
| 15366528 | 0.5mL tubes | 81 | 12.0 | 133 x 130 x 42 | Blue | 10 |
| 15376528 | 0.5mL tubes | 81 | 12.0 | 133 x 130 x 42 | Red | 10 |
| 15386528 | 0.5mL tubes | 81 | 12.0 | 133 x 130 x 42 | Purple | 10 |
| 15226996 | 1.5/2.0mL tubes | 81 | 12.0 | 133 x 133 x 51 | Natural | 10 |
| 15396528 | 1.5/2.0mL tubes | 81 | 12.0 | 133 x 133 x 51 | Blue | 10 |
| 15306538 | 1.5/2.0mL tubes | 81 | 12.0 | 133 x 133 x 51 | Red | 10 |
| 15316538 | 1.5/2.0mL tubes | 81 | 12.0 | 133 x 133 x 51 | Purple | 10 |
| 15326538 | 5mL tubes | 25 | 22.6 | 133 x 133 x 75 | Natural | 10 |
| 15336538 | 5mL tubes | 25 | 22.6 | 133 x 133 x 75 | Blue | 10 |
| 15346538 | 5mL tubes | 25 | 22.6 | 133 x 133 x 75 | Red | 10 |
| 15356538 | 5mL tubes | 25 | 22.6 | 133 x 133 x 75 | Purple | 10 |
| Corrugated polypropylene 0.2/0.5/1.5 to 2.0/5mL | | | | | | |
| 15366538 | 15mL tubes | 36 | 19.8 | 145 x 147 x 120 | Natural | 10 |
| 15376538 | 15mL tubes | 36 | 19.8 | 145 x 147 x 120 | Blue | 10 |
| 15386538 | 15mL tubes | 36 | 19.8 | 145 x 147 x 120 | Red | 10 |
| 15396538 | 15mL tubes | 36 | 19.8 | 145 x 147 x 120 | Purple | 10 |
| 15306548 | 50mL tubes | 16 | 31.6 | 145 x 147 x 120 | Natural | 10 |
| 15316548 | 50mL tubes | 16 | 31.6 | 145 x 147 x 120 | Blue | 10 |
| 15326548 | 50mL tubes | 16 | 31.6 | 145 x 147 x 120 | Red | 10 |
| 15336528 | 50mL tubes | 16 | 31.6 | 145 x 147 x 120 | Purple | 10 |





Filter papers, qualitative, cellulose, QL and QT range

| Cat.No | Description | Whatman equivalent | Diameter, mm | Pack qty |
|----------|---|--------------------|--------------|----------|
| 11572003 | QL 100 medium retention, medium/fast flow rate | No 1 | 42.5 | 100 |
| 11546873 | QL 100 medium retention, medium/fast flow rate | No 1 | 55 | 100 |
| 11556873 | QL 100 medium retention, medium/fast flow rate | No 1 | 70 | 100 |
| 11582003 | QL 100 medium retention, medium/fast flow rate | No 1 | 85 | 100 |
| 11566873 | QL 100 medium retention, medium/fast flow rate | No 1 | 90 | 100 |
| 11576873 | QL 100 medium retention, medium/fast flow rate | No 1 | 110 | 100 |
| 11425248 | QL 100 medium retention, medium/fast flow rate | No 1 | 125 | 100 |
| 11435248 | QL 100 medium retention, medium/fast flow rate | No 1 | 150 | 100 |
| 11445248 | QL 100 medium retention, medium/fast flow rate | No 1 | 185 | 100 |
| 11592003 | QL 100 medium retention, medium/fast flow rate | No 1 | 240 | 100 |
| 11704273 | QL 115 same as QL110 already folded, cellulose | No 2V | 240 | 100 |
| 11455208 | QL 125 coarse retention and fast flow rate, cellulose | No 4 | 185 | 100 |
| 11779975 | QL 145 coarse retention and fast flow rate, folded | No 113V | 240 | 100 |
| 11794303 | QT 200 low ash, medium retention, medium flow rate | No 40 | 90 | 100 |
| 11708963 | QT 200 low ash, medium retention, medium flow rate | No 40 | 125 | 100 |
| 11718963 | QT 200 low ash, medium retention, medium flow rate | No 40 | 150 | 100 |
| 11445188 | QT 205 low ash, fine retention, slow flow rate | No 41 | 55 | 100 |
| 11752375 | QT 280 hardened low ash, coarse retention, fast flow rate | No 541 | 125 | 100 |
| 11772405 | QT 290 hardened low ash, fine retention, slow flow rate | No 542 | 125 | 100 |



Filter papers, phase separation

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11784033 | Filter paper, phase separation, 125mm diameter | 100 |

Filter papers, glass microfibre, MF range

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11724073 | Microfibre filter MF 100, fast filtration speed (60s), for scintillation counting/petroleum analysis, glass microfibre, 52g/m2 mass, 260µm thickness, 1.6µm retention, 47mm diameter | 100 |
| 11455238 | Microfibre filter MF 100, fast filtration speed (60s), for scintillation counting/petroleum analysis, glass microfibre, 52g/m2 mass, 260µm thickness, 1.6µm retention, 70mm diameter | 100 |
| 11465238 | Microfibre filter MF 100, fast filtration speed (60s), for scintillation counting/petroleum analysis, glass microfibre, 52g/m2 mass, 260µm thickness, 1.6µm retention, 90mm diameter | 100 |
| 11704083 | Microfibre filter MF 200, medium to fast filtration speed (100s), for water analysis/solids in water, glass microfibre, 52g/m2 mass, 260µm thickness, 1.2µm retention, 47mm diameter | 100 |
| 11532083 | Microfibre filter MF 200, medium to fast filtration speed (100s), for water analysis/solids in water, glass microfibre, 52g/m2 mass, 260µm thickness, 1.2µm retention, 70mm diameter | 100 |
| 11714083 | Microfibre filter MF 200, medium to fast filtration speed (100s), for water analysis/solids in water, glass microfibre, 52g/m2 mass, 260µm thickness, 1.2µm retention, 90mm diameter | 100 |
| 11724083 | Microfibre filter MF 200 medium to fast filtration speed (100s) for water analysis/solids in water, glass microfibre 52g/m2 mass 260µm thickness 1.2µm retention 110mm diameter | 100 |



Syringe filters

- Wide range of membranes and pore sizes
- Sterile and non-sterile versions
- Three different sizes 13mm, 25mm and 33mm



| Cat. No | Membrane | Diameter, mm | Pore size, µm | Sterile | Pack qty |
|----------|------------------|--------------|---------------|---------|----------|
| 15206869 | PES | 33 | 0.2 | Yes | 50 |
| 15216869 | PES | 33 | 0.45 | Yes | 50 |
| 15121499 | Nylon | 25 | 0.2 | No | 50 |
| 15131499 | Nylon | 25 | 0.45 | No | 50 |
| 15161499 | PTFE | 13 | 0.2 | No | 100 |
| 15171499 | PTFE | 13 | 0.45 | No | 100 |
| 15141499 | PTFE | 25 | 0.2 | No | 50 |
| 15151499 | PTFE | 25 | 0.45 | No | 50 |
| 15181499 | PVDF | 33 | 0.2 | Yes | 50 |
| 15191499 | PVDF | 33 | 0.45 | Yes | 50 |
| 15181489 | Hydrophilic PTFE | 25 | 0.2 | No | 50 |
| 15101499 | Hydrophilic PTFE | 25 | 0.45 | No | 50 |

Beakers, squat form

- Borosilicate glass 3.3
- Squat form with spout and graduations
- ISO 3819 DIN 12331

| Cat. No | Capacity, mL | Height, mm | Exterior diameter, mm | Pack qty |
|----------|--------------|------------|-----------------------|----------|
| 15469073 | 25 | 50 | 34 | 10 |
| 15479073 | 50 | 60 | 42 | 10 |
| 15489073 | 100 | 70 | 50 | 10 |
| 15499073 | 150 | 80 | 60 | 10 |
| 15409083 | 250 | 95 | 70 | 10 |
| 15419083 | 400 | 110 | 80 | 10 |
| 15429103 | 600 | 125 | 90 | 10 |
| 15439103 | 800 | 135 | 100 | 10 |
| 15449083 | 1,000 | 143 | 105 | 10 |
| 15459083 | 2,000 | 185 | 130 | 10 |
| 15469083 | 3,000 | 210 | 150 | 1 |
| 15479083 | 5,000 | 270 | 170 | 1 |
| 15489083 | 10,000 | 350 | 217 | 1 |



Reagent and media bottles, borosilicate clear glass, wide screw neck

- Graduated and supplied with blue polypropylene cap and pouring ring
- ISO/GL 45 screw thread standard to all sizes, conforms to ISO R1115 and DIN 168
- Chemically resistant and can be repeatedly sterilised to 140°C
- Autoclavable, do not autoclave bottles with a tightly screwed cap

| Cat. No | Capacity, mL | Height, mm (bottle only) | O.D., mm | Wall thickness, mm | Internal neck diameter, mm | Pack qty |
|----------|--------------|--------------------------|----------|--------------------|----------------------------|----------|
| 15446113 | 100 | 100 | 56 | 1.5 | 27 | 10 |
| 15456113 | 250 | 138 | 70 | 1.5 | 27 | 10 |
| 15476113 | 500 | 176 | 86 | 1.5 | 27 | 10 |
| 15486113 | 1,000 | 225 | 101 | 1.7 | 27 | 10 |
| 15406123 | 2,000 | 260 | 136 | 2.0 | 27 | 1 |
| 15416123 | 5,000 | 330 | 181 | 2.0 | 27 | 1 |
| 15436123 | 10,000 | 410 | 227 | 2.7 | 27 | 1 |

Accessories

Polypropylene screw caps, with lip seal, autoclavable to 140°C

| Cat. No | ISO thread GL | Colour | Pack qty |
|----------|---------------|--------|----------|
| 11527033 | 45 | Blue | 10 |
| 11537033 | 45 | Green | 10 |
| 11547033 | 45 | Yellow | 10 |

Polypropylene pouring rings, autoclavable to 140°C

| Cat. No | ISO thread GL | Colour | Pack qty |
|----------|---------------|--------|----------|
| 11587033 | 45 | Green | 10 |
| 11597033 | 45 | Yellow | 10 |



Burettes

Borosilicate glass, single bore, interchangeable TFE stopcock, blue ceramic graduations, Class A, BS-ISO 385: 2005

| Cat. No | Capacity, mL | Graduations, mL | Pack qty |
|----------|--------------|-----------------|----------|
| 10450893 | 25 | 0.10 | 1 |
| 10738081 | 50 | 0.10 | 1 |

Accessories

| Cat. No | Description | Pack qty |
|----------|---|----------|
| 12967730 | Burette tip for Fisherbrand burettes borosilicate glass | 5 |

Borosilicate glass, single bore, interchangeable stopcock, BS 846, blue ceramic graduations, Class B, BS-ISO 385: 2005

| Cat. No | Capacity, mL | Graduations, mL | Pack qty |
|----------|--------------|-----------------|----------|
| 10448972 | 50 | 0.10 | 1 |
| 10569771 | 50 | 0.10 | 1 |

Accessories

| Cat. No | Description | Pack qty |
|----------|---|----------|
| 12967730 | Burette tip for Fisherbrand burettes borosilicate glass | 5 |





Funnels, borosilicate glass, conical

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11572423 | Funnel, borosilicate glass, 55mm top diameter, 60mm stem | 10 |
| 11582423 | Funnel, borosilicate glass, 75mm top diameter, 80mm stem | 10 |
| 11592423 | Funnel, borosilicate glass, 100mm top diameter, 110mm stem | 10 |
| 11502433 | Funnel, borosilicate glass, 155mm top diameter, 150mm stem | 1 |
| 11512433 | Funnel, borosilicate glass, 215mm top diameter, 180mm stem | 1 |



Measuring cylinders


- Borosilicate glass measuring cylinder, spouted, with hexagonal base

| Cat.No | Description | Pack qty |
|----------|----------------------------------|----------|
| 11507702 | Measuring cylinder, 5mL x 0.1mL | 2 |
| 11517702 | Measuring cylinder, 10mL x 0.2mL | 2 |
| 11527702 | Measuring cylinder, 25mL x 0.5mL | 2 |
| 11537702 | Measuring cylinder, 50mL x 1mL | 2 |
| 11547702 | Measuring cylinder, 100mL x 1mL | 2 |
| 11557702 | Measuring cylinder, 250mL x 2mL | 2 |
| 11567702 | Measuring cylinder, 500mL x 5mL | 2 |
| 11577702 | Measuring cylinder, 1L x 10mL | 2 |



Separating funnels

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11572782 | Funnel separating, conical, borosilicate glass with glass and plastic stopper, PTFE stopcock, 100mL | 2 |
| 11582782 | Funnel separating, conical, borosilicate glass with glass and plastic stopper, PTFE stopcock, 250mL | 2 |



Test tubes and culture tubes

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 12961031 | Test tube, medium wall, rimless, soda lime glass, wall thickness 1.0mm, ISO 4142, 100mm x 12mm (len x O.D.) | 100 |
| 11972178 | Test tube, medium wall, rimless, soda lime glass, wall thickness 1.2mm, ISO 4142, 150mm x 24mm (len x O.D.) | 50 |
| 11517403 | Culture tube, rimless, borosilicate glass, disposable, 12mm x 75mm | 1,000 |
| 11527403 | Culture tube, rimless, borosilicate glass, disposable, 13mm x 100mm | 1,000 |
| 11547403 | Culture tube, rimless, borosilicate glass disposable, 16mm x 100mm | 1,000 |
| 11557403 | Culture tube, rimless, borosilicate glass, disposable, 16mm x 125mm | 1,000 |
| 11577403 | Culture tube, rimless, borosilicate glass, disposable, 16mm x 150mm | 1,000 |
| 11587403 | Culture tube, rimless, borosilicate glass, disposable, 18mm x 150mm | 500 |
| 11507413 | Culture tube, rimless, borosilicate glass, disposable, 25mm x 150mm | 500 |



Vials and caps

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11660112 | Vial, sample, with attached black poly-seal cone caps, borosilicate glass (type one class A), packed in partition trays, clear, 4mL, 15mm x 48mm | 144 |
| 11523532 | Vial, sample with attached black PTFE faced rubber lined caps, borosilicate glass (type one class A), packed in partition trays, clear, 4mL, 15mm x 48mm | 144 |
| 11571364 | Vial, 0.7mL, conical with 8mm crimp top, 40mm x 7mm, 1st hydrolytic class glass | 100 |
| 10326042 | Vial, 11mm, crimp top, wide opening, 2mL, 12mm x 32mm, 1st hydrolytic class glass | 100 |
| 11535884 | Vial, 11mm, crimp top, wide opening with label and filling lines 2mL glass 12mm x 32mm 1st hydrolytic class | 100 |
| 11545884 | Vial crimp top, 1mm wide opening, with label and filling lines, glass, amber, 12mm x 32mm, 2mL, 1st hydrolytic class | 100 |
| 11585914 | Vial, snap ring, with integrated 0.2mL glass micro insert 0.2mL, glass, 12mm x 32mm, 1st hydrolytic class | 100 |
| 11595864 | Crimp cap, centre hole, red rubber / PTFE beige, 45° shore A, 1.0mm, Agilent quality, 11mm, aluminium | 100 |
| 10457692 | Crimp cap, centre hole, with 1.0mm natural rubber/transparent TEF seal, Agilent quality, 11mm aluminium | 100 |
| 10112372 | Crimp cap, centre hole, with 3.0mm butyl/grey PTFE seal, aluminium, silver, 20mm | 100 |
| 10192172 | Snap ring cap, ND11, polyethylene, transparent, 6mm centre hole, 1.0mm thickness, with white silicone/blue PTFE cross-slit septum, durometer 55° shore A | 100 |



Volumetric flasks

- Borosilicate glass flask, volumetric class A, with polyethylene stopper, batch certificate available

| Cat.No | Description | Pack qty |
|----------|-----------------------------------|----------|
| 11576923 | Flask, volumetric, class A, 5mL | 5 |
| 11586923 | Flask volumetric, class A, 10mL | 5 |
| 11596923 | Flask, volumetric class A, 20mL | 5 |
| 11506933 | Flask, volumetric, class A, 25mL | 5 |
| 11516933 | Flask, volumetric, class A, 50mL | 5 |
| 11526933 | Flask, volumetric, class A, 100mL | 5 |
| 11536933 | Flask, volumetric, class A, 200mL | 2 |
| 11546933 | Flask, volumetric, class A, 250mL | 2 |
| 11556933 | Flask, volumetric, class A, 500mL | 2 |
| 11566933 | Flask, volumetric, class A, 1L | 2 |
| 11576933 | Flask, volumetric, class A, 2L | 1 |



Cuvettes, macro, absorption, with PTFE stopper

- Quartz windows for transmission from far ultraviolet through the infrared
- Includes: PTFE stopper

| Cat. No | Type | Material | Dimensions, mm | Capacity, μ L | Pathlength, mm | Wavelength range, nm | Internal width, mm | Pack qty |
|----------|------------------|----------|------------------|-------------------|----------------|----------------------|--------------------|----------|
| 15216023 | Macro absorption | Quartz | 12.5 x 12.5 x 45 | 3,000 | 10 | 200 to 2,500 | 10 | 4 |



Cuvette, quartz, semi-micro, UV

- Reduce reagent use with semi-micro cells
- Quartz windows for transmission from far ultraviolet through the infrared
- Includes: PTFE lid

| Cat. No | Type | Material | Dimensions, mm | Capacity, μ L | Pathlength, mm | Wavelength range, nm | Internal width, mm | Pack qty |
|----------|-----------------------|----------|------------------|-------------------|----------------|----------------------|--------------------|----------|
| 15226033 | Semi-micro absorption | Quartz | 12.5 x 12.5 x 45 | 1,000 | 10 | 200 to 2,500 | 4 | 2 |
| 15246033 | Semi-micro absorption | Quartz | 12.5 x 12.5 x 45 | 1,000 | 10 | 200 to 2,500 | 4 | 4 |



Weighing boats, polystyrene, disposable

- Available as standard or anti-static
- Diamond or square shaped
- In black or white
- 5mL to 280mL capacity

| Cat. No | Shape | Capacity, mL | Pack qty |
|-------------------|---------|--------------|----------|
| Standard black | | | |
| 12952850 | Diamond | 5 | 500 |
| 12962850 | Diamond | 30 | 500 |
| 12972850 | Diamond | 100 | 500 |
| 12982850 | Square | 10 | 500 |
| 12992850 | Square | 85 | 500 |
| 12902860 | Square | 280 | 500 |
| Anti-static black | | | |
| 12912860 | Diamond | 5 | 500 |
| 12922860 | Diamond | 30 | 500 |
| 12932860 | Diamond | 100 | 500 |
| 12942860 | Square | 10 | 500 |
| 12952860 | Square | 85 | 500 |
| 12962860 | Square | 280 | 500 |
| Standard white | | | |
| 12932840 | Diamond | 5 | 500 |
| 11573422 | Diamond | 30 | 500 |
| 12387552 | Diamond | 100 | 500 |
| 11593422 | Square | 10 | 500 |
| 11503432 | Square | 85 | 500 |
| 11513432 | Square | 280 | 500 |
| Anti-static white | | | |
| 12992840 | Diamond | 5 | 500 |
| 11533432 | Diamond | 30 | 500 |
| 11543432 | Diamond | 100 | 500 |
| 12922850 | Square | 10 | 500 |
| 11680302 | Square | 85 | 500 |
| 12608513 | Square | 280 | 500 |



Weighing funnels

| Cat.No | Description |
|----------|------------------------------|
| 11572902 | Weighing funnel, glass, 3mL |
| 11582902 | Weighing funnel, glass, 6mL |
| 11592902 | Weighing funnel, glass, 10mL |



Scoops

| Cat.No | Description |
|----------|------------------------------------|
| 11597852 | Scoop, polypropylene, white, 250mL |
| 11507862 | Scoop, polypropylene, white, 500mL |
| 11517862 | Scoop, polypropylene, white, 1L |



Pasteur pipettes, glass

- Pasteur pipettes, open point, unplugged, soda lime glass


| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11546963 | Pipette, 4 boxes of 250 pipettes, 150mm length | 1,000 |
| 11566963 | Pipette, 4 boxes of 250 pipettes, 230mm length | 1,000 |
| 11755108 | Pipette, 4 boxes of 250 pipettes, 270mm length | 1,000 |



Pipettes, transfer

- low density polyethylene
- Transparent
- Graduated or non-graduated
- Sterile options available
- Various packaging formats


| Cat. No | Description | Sterile | Length, mm | Drop volume, μ L | Drop per mL | Pack qty |
|-------------------|---|---------|------------|----------------------|-------------|----------|
| 1 13469118 | Transfer pipette PE, 1mL | No | 104 | 33 | 30 | 400 |
| 2 13499108 | Transfer pipette PE, 1mL, graduated | No | 150 | 33 | 30 | 500 |
| 3 13439118 | Transfer pipette PE, 1mL, graduated, extended tip | No | 150 | 33 | 30 | 500 |
| 13489108 | Transfer pipette PE, 1mL, graduated | Yes | 150 | 33 | 30 | 500 |
| 13419118 | Transfer pipette PE, 1mL, graduated, inner pack of 10 | Yes | 150 | 33 | 30 | 500 |
| 4 13429118 | Transfer pipette PE, 1mL, graduated, inner pack of 20 | Yes | 150 | 33 | 30 | 500 |
| 13439108 | Transfer pipette PE, 3mL, graduated | No | 155 | 40 | 25 | 500 |
| 13469108 | Transfer pipette PE, 3mL, graduated | Yes | 155 | 40 | 25 | 500 |
| 13479108 | Transfer pipette PE, 3mL, graduated, inner pack of 10 | Yes | 155 | 40 | 25 | 500 |
| 13409118 | Transfer pipette PE, 3mL, graduated, inner pack of 20 | Yes | 155 | 40 | 25 | 500 |
| 5 13459118 | Transfer pipette PE, 4mL, thin stem | No | 150 | 33 | 30 | 500 |
| 6 13459108 | Transfer pipette PE, 4mL | No | 150 | 33 | 50 | 500 |
| 7 13449108 | Transfer pipette PE, 7mL, extra long | No | 300 | 50 | 20 | 100 |
| 8 13449118 | Transfer pipette PE, 10mL, jumbo | No | 170 | 56 | 18 | 200 |



Pipettes, bulb

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11942168 | Pipette, one mark, BS 1583, class B, soda lime glass, blue ceramic graduation, 25mL | 5 |
| 11952168 | Pipette, one mark, BS 1583, class A, uncertified, soda lime glass, blue ceramic graduation, 1mL | 2 |
| 10692652 | Pipette, one mark, BS 1583, class A, uncertified, soda lime glass, blue ceramic graduation, 2mL | 2 |
| 10080251 | Pipette, one mark, BS 1583, class A, uncertified, soda lime glass, blue ceramic graduation, 5mL | 2 |
| 11982168 | Pipette, one mark, BS 1583, class A, uncertified, soda lime glass, blue ceramic graduation, 10mL | 2 |





Serological pipettes

- Serological pipette, sterile plastic with magnifier stripe

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11849181 | Serological pipette individually paper peel wrapped 1ml, plugged | 1000 |
| 11859181 | Serological pipette individually paper peel wrapped 2ml, plugged | 500 |
| 11869181 | Serological pipette individually paper peel wrapped 5ml, plugged | 200 |
| 11879181 | Serological pipette individually paper peel wrapped 10ml, plugged | 200 |
| 11839181 | Serological pipette individually paper peel wrapped 25ml, plugged | 200 |
| 11889181 | Serological pipette individually paper peel wrapped 50ml, plugged | 100 |
| 11819660 | Serological pipette individually plastic peel wrapped 1ml, plugged | 1000 |
| 11879650 | Serological pipette individually plastic peel wrapped 2ml, plugged | 500 |
| 11829660 | Serological pipette individually plastic peel wrapped 5ml, plugged | 200 |
| 11839660 | Serological pipette individually plastic peel wrapped 10ml, plugged | 200 |
| 11517752 | Serological pipette individually plastic peel wrapped 25ml, plugged | 200 |

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11879660 | Serological pipette bulk packed 1ml, plugged | 1000 |
| 11889660 | Serological pipette bulk packed 2ml, plugged | 500 |
| 11899660 | Serological pipette bulk packed 5ml, plugged | 500 |
| 11809670 | Serological pipette bulk packed 10ml, plugged | 500 |



1

Pipette fillers, universal

- The side-mounted thumbwheel is effortlessly manipulated for exact aspiration or drop-wise delivery
- The plunger is easily depressed for quick dispensing
- Colour coded barrel indicates pipette size

| Cat. No | Description |
|----------|--|
| 15209805 | Pipette filler, universal, blue, 2mL |
| 15239805 | Pipette filler, universal, green, 10mL |
| 15229805 | Pipette filler, universal, red, 25mL |
| 15219805 | Pipette controller |

1



Pipettor kit

- A set of four manual single-channel pipettors

| Cat. No | Description |
|----------|--|
| 15268638 | Pipettor kit, ELITE GLP, 4 pipettors single channel (0.2-2µL, 2-20µL, 20-200µL, 100-1000µL) plus stand |

Pipette filler, motorised

- For glass and plastic pipettes from 0.1 to 200mL
- Ergonomic handgrip
- Variable speed with two pipette buttons
- Safety valve
- NiMH battery, 2.4V, 700mAh
- Recharging indicator LED

| Cat. No | Description |
|----------|--------------------------|
| 15249805 | Motorised pipette filler |




Reservoirs

| Cat.No | Description |
|----------|--|
| 12369175 | Pipetting reservoir, non-sterile, PVC, disposable, clear, 50mL |



Pipettor tips, filter, SureOne

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11903466 | Pipettor tip, SureOne, filter, sterile, rack, micropoint tip, graduated at 2µL, polypropylene, 96 tip tray x 10, clear, 0.1 to 10µL | 960 |
| 11913466 | Pipettor tip, SureOne, filter, sterile, rack, micropoint tip, graduated at 2.5µL, extended length, polypropylene, 96 tip tray x 10 clear, 0.1 to 10µL | 960 |
| 11963466 | Pipettor tip, SureOne, filter, sterile, rack, micropoint tip, graduated at 10, 50, 100µL, polypropylene, 96 tip tray x 10, clear, 20 to 200µL | 960 |
| 11973466 | Pipettor tip, SureOne, filter, sterile, rack, micropoint tip, graduated at 100, 200, 500, 1,000µL, polypropylene, 96 tip tray x 10, clear, 100 to 1,000µL | 960 |
| 11943466 | Pipettor tip, SureOne, filter, sterile, rack, bevelled tip, polypropylene, 96 tip tray x 10 clear, 2 to 20µL | 960 |
| 11953466 | Pipettor tip, SureOne, filter, sterile, rack, bevelled tip, polypropylene, 96 tip tray x 10 clear, 10 to 100µL | 960 |



Pipette tips, universal

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11577442 | Pipettor tip, crystal, G, bulk pack, polypropylene, natural, 0.1µL to 10µL | 1,000 |
| 11587442 | Pipettor tip, yellow, E, bulk pack, polypropylene, yellow, 1µL to 200µL | 1,000 |
| 11517452 | Pipettor tip, yellow, bevelled, bulk pack, polypropylene, yellow, 1µL to 200µL | 1,000 |
| 10787524 | Pipettor tip, blue, bulk pack, polypropylene, blue, 200µL to 1,000µL | 1,000 |
| 11587452 | Pipettor tip, natural, bulk pack, polypropylene, natural, 200µL to 1,000µL | 1,000 |
| 11846943 | Pipettor tip, makro, universal, bulk pack, polypropylene, natural, 1mL to 5mL | 250 |
| 11527462 | Pipettor tip finn, bulk pack, polypropylene, natural, 1mL to 5mL | 250 |
| 11567462 | Pipettor tip, finn, bulk pack, polypropylene, natural, 2mL to 10mL | 100 |





Pipettor tips, universal, SureOne

- Micropoint tips, 0.1µL to 10µL, graduated at 2µL

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11933416 | Pipettor tip, SureOne, bulk, micropoint tip, graduated at 2µL, polypropylene, 1,000 bag, clear | 1,000 |
| 10053014 | Pipettor tip, SureOne racked sterile micropoint tip graduated at 2µL polypropylene 96 tip tray x 10 clear | 960 |

- Extended length micropoint tip, 0.1µL to 10µL, graduated at 2.5µL

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11983416 | Pipettor tip, SureOne, bulk, micropoint tip, graduated at 2.5µL, extended length, polypropylene, clear | 1,000 |

- Bevelled tips 1µL to 200µL

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11933426 | Pipettor tip, SureOne, bulk, bevelled tip, polypropylene, 1,000 bag, clear, 1 to 200µL | 1,000 |
| 11578412 | Pipettor tip, SureOne ,reload system pack. bevelled tip, polypropylene, 96 tip tray x 10, clear, 1 to 200µL | 960 |
| 10124314 | Pipettor tip, SureOne, bulk, bevelled tip, polypropylene, 1,000 bag, yellow, 1 to 200µL | 1,000 |
| 11913436 | Pipettor tip, SureOne. reload system pack, bevelled tip, polypropylene, 96 tip tray x 10, yellow, 1 to 200µL | 960 |

- Thin wall micropoint tip. 1µL to 200µL

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11963436 | Pipettor tip, SureOne, bulk, thin wall, micropoint tip, polypropylene, yellow | 1,000 |
| 11993436 | Pipettor tip, SureOne, reload system pack, thin wall, micropoint tip, polypropylene, 96 tip tray x 10, yellow | 960 |

- Bevelled tip, 1µL to 200µL, graduated at 10µL, 50µL, 100µL

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 11903446 | Pipettor tip, SureOne, bulk, yellow, bevelled tip, graduated at 10, 50, 100, 200µL | 1,000 |

- Micropoint tips, 100µL to 1,250µL, graduated at 100µL, 200µL, 500µL, 1,000µL

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11548442 | Pipettor tip, SureOne, bulk, micropoint tip, graduated at 100, 200, 500, 1,000µL, polypropylene, 1,000 bag, clear, 100 to 1250µL | 1,000 |
| 11568442 | Pipettor tip, SureOne, racked, micropoint tip, graduated at 100, 200, 500, 1,000µL, polypropylene, 96 tip tray x 10, clear, 100 to 1250µL | 960 |
| 10164694 | Pipettor tip, SureOne, racked, sterile, micropoint tip, graduated at 100, 200, 500, 1,000µL, polypropylene, 96 tip tray x 10, clear, 100 to 1250µL | 960 |
| 11588442 | Pipettor tip, SureOne, reload system pack, micropoint tip, graduated at 100, 200, 500, 1,000µL, polypropylene, 96 tip tray x 10, clear, 100 to 1250µL | 960 |
| 11963456 | Pipettor tip, SureOne, reload system pack, micropoint tip, graduated at 100, 200, 500, 1,000µL, polypropylene, 96 tip tray x 10, blue, 100 to 1250µL | 960 |



Microscope slides

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 12373118 | Ground edges, plain glass, wrapped in cellophane, 76mm x 26mm 0.8mm to 1.0mm thick | 50 |
| 12383118 | Ground edges, plain glass, wrapped in cellophane, 76mm x 26mm 1.0mm to 1.2mm thick | 50 |
| 11562203 | Ground edges, twin frosted glass, wrapped in cellophane, 76mm x 26mm 0.8mm to 1.0mm thick | 50 |
| 11572203 | Ground edges, twin frosted glass, wrapped in cellophane, 76mm x 26mm 1.0mm to 1.2mm thick | 50 |

Microscope cover slips

| Cat.No | Description | Pack qty |
|----------|-------------------------------------|----------|
| 12333128 | Cover slip, No.1 glass, 22mm x 22mm | 200 |
| 11348503 | Cover slip, No.1 glass, 24mm x 40mm | 100 |

Microscope slide holders

| Cat.No | Description |
|----------|---|
| 11741476 | Slide box for 25 slides, hinged lid, polypropylene, unlined, 8mm x 140mm x 35mm, blue |
| 11771476 | Slide box for 100 slides, hinged lid, polypropylene, unlined, 165mm x 210mm x 35mm, blue |
| 11781476 | Slide box for 100 slides, hinged lid, polypropylene, unlined, 165mm x 210mm x 35mm, red |
| 11791476 | Slide box for 100 slides, hinged lid, polypropylene, unlined, 165mm x 210mm x 35mm, white |
| 11701486 | Slide box for 100 slides, nickel-plated clasp, ABS plastic, foam lined, stackable, 165mm x 210mm x 35mm, blue |



Plasticware

Beakers

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 10792921 | Squat form, PTFE, 68mm x 54mm (h x d), 100mL | 1 |
| 10713101 | Squat form, PTFE, 97mm x 66mm (h x d), 250mL | 1 |
| 11572283 | Squat form, with spout and graduations, polypropylene, BS 5404.1, ISO 7056, 25mL | 10 |
| 11512293 | Squat form, with spout and graduations, polypropylene, 500mL | 10 |
| 11769398 | Tri-cornered, polypropylene, 250mL | 100 |
| 11759398 | Tri-cornered, polypropylene, 100mL | 100 |
| 11799398 | Tri-cornered, polypropylene, 1,000mL | 100 |



Containers, polyethylene with LDPE lid, specimen

- Designed for collection, transport and storage of liquid and dry specimens and samples.
- Moulded of autoclavable high quality HDPE, with thick walls
- Inert to most chemicals including formaldehyde, weak acids and all bases
- Resist freezing and boiling
- Easy to write on with pen or pencil
- Stackable for compact storage
- Come with LDPE lids (non-autoclavable)

| Cat. No | Style | Capacity, mL | Colour | Pack qty |
|----------|--------------|--------------|-------------|----------|
| 12029977 | Short/wide | 120 | Translucent | 300 |
| 11709388 | Tall/thin | 120 | Translucent | 100 |
| 12049977 | Short/wide | 240 | Translucent | 100 |
| 11719388 | Short/wide | 240 | White | 100 |
| 11974375 | Tall/thin | 240 | Translucent | 100 |
| 12029957 | Tall/thin | 240 | White | 100 |
| 12089977 | Multipurpose | 473 | Translucent | 100 |
| 11994375 | Multipurpose | 473 | White | 100 |
| 12089947 | Multipurpose | 1,100 | Translucent | 100 |
| 12009957 | Multipurpose | 1,900 | Translucent | 50 |
| 12009967 | Multipurpose | 1,900 | White | 50 |
| 12039957 | Multipurpose | 2,500 | Translucent | 25 |
| 12069967 | Multipurpose | 2,500 | White | 25 |
| 12099957 | Multipurpose | 5,100 | Translucent | 10 |
| 12089967 | Multipurpose | 5,100 | White | 10 |





Containers, jars, polypropylene, translucent, heavy weight

- Ideal for general laboratory applications.
- Excellent chemical resistance
- Polypropylene screw closure - leakproof
- Autoclavable with cap loosened

| Cat. No | Capacity, mL | Cap size, mm | Height, mm | O.D., mm | Pack qty |
|----------|--------------|--------------|------------|----------|----------|
| 11523242 | 30 | 33/R3 | 45 | 36 | 10 |
| 11533242 | 60 | 33/R3 | 80 | 36 | 10 |
| 11543242 | 125 | 38/R3 | 93 | 50 | 10 |
| 11553242 | 150 | 58/R3 | 68 | 60 | 10 |
| 11563242 | 250 | 58/R3 | 114 | 61 | 10 |
| 11573242 | 500 | 58/R3 | 140 | 76 | 10 |
| 11583242 | 1,000 | 70/R3 | 200 | 90 | 5 |



Bottles, wide mouth, HDPE

| Cat. No | Capacity, mL | Dimensions, mm (dia. x h) | Pack qty |
|-------------|--------------|---------------------------|----------|
| Clear | | | |
| 11775243 | 30 | 36 x 64 | 72 |
| 11785243 | 60 | 38 x 86 | 72 |
| 11795253 | 250 | 62 x 132 | 72 |
| 11775253 | 500 | 74 x 168 | 48 |
| 11745253 | 1,000 | 81 x 201 | 24 |
| Translucent | | | |
| 11917974 | 125 | 51 x 99 | 500 |
| 11947964 | 250 | 62 x 132 | 250 |
| 11907964 | 500 | 74 x 168 | 125 |
| 11987924 | 1,000 | 81 x 201 | 50 |



Bottles, wide mouth, LDPE

| Cat. No | Capacity, mL | Dimensions, mm (dia. x h) | Pack qty |
|----------|--------------|---------------------------|----------|
| 11957934 | 30 | 36 x 64 | 72 |
| 11967934 | 60 | 38 x 86 | 72 |
| 11977934 | 125 | 51 x 99 | 72 |
| 11987934 | 250 | 62 x 132 | 72 |
| 11937924 | 500 | 74 x 168 | 48 |
| 11967904 | 1,000 | 81 x 201 | 24 |



Jerricans

| Cat.No | Description |
|----------|--|
| 11597433 | Bottle, narrow neck, jerrican, polyethylene, translucent, 165mm height, 38mm neck, 1L |
| 11507443 | Bottle, narrow neck, jerrican, polyethylene, translucent, 205mm height, 45mm neck, 2,500mL |
| 11517443 | Bottle, narrow neck, jerrican, polyethylene, translucent, 282mm height, 38mm neck, 5L |
| 11527443 | Bottle, narrow neck, jerrican, polyethylene, translucent, 314mm height, 61mm neck, 10L |
| 11537443 | Bottle, narrow neck, jerrican, polyethylene, translucent, 470mm height, 61mm neck, 25L |

Carboys

- Heavy wall construction provides long lasting, rugged service
- Rectangular shape optimises bench space
- Choice of materials for general use or autoclavability


| Cat. No | Material | Capacity, mL | Spigot | Autoclavable | Pack qty |
|----------|---------------------|--------------|--------|--------------|----------|
| 13458029 | HDPE | 900 | No | No | 1 |
| 13478029 | HDPE | 2000 | No | No | 1 |
| 13438039 | HDPE | 900 | Yes | No | 1 |
| 13408039 | HDPE | 2000 | Yes | No | 1 |
| 13448029 | PP | 900 | No | Yes | 1 |
| 13468029 | PP | 2000 | No | Yes | 1 |
| 13428039 | PP | 2000 | Yes | Yes | 1 |
| 13418039 | PP | 900 | Yes | Yes | 1 |
| 13488029 | Replacement closure | | | | 6 |
| 13448039 | Replacement spigot | | | | 1 |



Cylinders, plastic, graduated, with spout

- Large rounded bases ribbed for reinforcement and stability
- Calibrated 'To Contain/To Deliver' at 20°C - meets ASTM laboratory standards
- Chemical, heat and impact resistant makes these cylinders excellent for long term use
- No meniscus to confuse readings - eliminates guesswork

| Cat. No | Capacity, mL | Pack qty |
|-----------------|--------------|----------|
| Cylinders, PPCO | | |
| 11947884 | 10 | 12 |
| 11957884 | 25 | 12 |
| 11967884 | 50 | 12 |
| 11997874 | 100 | 8 |
| 11907884 | 250 | 8 |
| 11967874 | 500 | 6 |
| 11937874 | 1,000 | 4 |
| 11927874 | 2,000 | 2 |
| Cylinders, PMP | | |
| 11907894 | 25 | 18 |
| 11917894 | 50 | 18 |
| 11977884 | 100 | 12 |
| 11917884 | 500 | 8 |
| 11977874 | 1,000 | 6 |
| 11947874 | 2,000 | 4 |



Wash bottles, LDPE

- Labelled and colour coded for most commonly used solvents
- Translucent with coloured closure
- Comes with vented 38mm closure to prevent solvent drips
- 500mL capacity

| Cat. No | Label | Closure colour | Pack qty |
|----------|-----------------|----------------|----------|
| 11755233 | Bleach | White | 6 |
| 11765233 | Deionised water | Natural | 6 |
| 11775233 | Distilled water | Natural | 6 |
| 11785233 | Ethanol | Natural | 6 |
| 11795233 | Isopropanol | Yellow | 6 |
| 11705243 | Methanol | Green | 6 |
| 11715243 | Saline | Natural | 6 |
| 11725243 | Soap | Natural | 6 |
| 11735243 | Universal | Natural | 6 |
| 11745243 | Water Solvent | Natural | 6 |
| 11765243 | Cleaning | Natural | 6 |
| 11745233 | Acetone | Red | 6 |





Bottles, LDPE, wash, safety labelled, wide neck, colour coded

- Solvent dispensing
- Leakproof dispensing of your most common solvents
- Easy to read LDPE solvent-resistant printed bottles
- Colour coded for easy identification

| Cat. No | Capacity, mL | Label | Closure colour | Pack qty |
|----------|--------------|-----------------|----------------|----------|
| 11562463 | 500 | Distilled Water | Blue | 6 |
| 11572463 | 500 | Water | Blue | 6 |
| 11582463 | 1,000 | Acetone | Red | 6 |
| 11592463 | 500 | Acetone | Red | 6 |
| 11502473 | 500 | Methanol | Green | 6 |
| 11512473 | 500 | Isopropanol | Yellow | 6 |
| 11522473 | 500 | Toluene | Red | 6 |
| 11532473 | 500 | Deionised water | Blue | 6 |



Bottles, LDPE, wash, round, wide neck, venting, solvent labelled

- Labelled and colour coded for some of the most commonly used solvents
- Unique patented two piece venting assembly prevents build up of solvents within the bottle, eliminating solvent drips from the spout due to internal evaporation
- LDPE bottle with polypropylene closures and venting assembly

NB: Only use with solvents compatible with LDPE refer to Table 7 Chemical Compatibility Chart page 14 to 15

| Cat. No | Capacity, mL | Label | Closure colour | Pack qty |
|----------|--------------|-------------|----------------|----------|
| 11567153 | 500 | Acetone | Red | 3 |
| 11577153 | 500 | Methanol | Green | 3 |
| 11587153 | 500 | Isopropanol | Yellow | 3 |
| 11507163 | 500 | DIY/custom | Natural | 3 |



Bottles, LDPE, wash, wide neck, venting, labelled, Right-to-Know

- Pre-printed with chemical name and formula, NFPA section 704 four colour diamond, health hazard, CAS number and suggested protective clothing and equipment
- Wide neck for easy refilling
- Vented to prevent pressure build up, eliminating solvent drips
- Valve allows solvent to be dispensed in both an upright or inverted position
- Leakproof
- Translucent with colour coded polypropylene closures (except sodium hypochlorite bottle, which is white LDPE for light protection between 230nm to 450nm)

| Cat. No | Capacity, mL | Label | Closure colour | Pack qty |
|----------|--------------|------------------------------|----------------|----------|
| 12368616 | 500 | Ethanol | Natural | 3 |
| 12359195 | 500 | Assorted pack | Assorted | 6 |
| 12369195 | 500 | Ethanol | Natural | 6 |
| 12309205 | 500 | Isopropanol | Yellow | 6 |
| 12319205 | 500 | Methanol | Green | 6 |
| 11532463 | 250 | Sodium hypochlorite (bleach) | Yellow | 6 |
| 11542463 | 500 | Distilled water | Natural | 6 |
| 11938485 | 1,000 | Acetone | Red | 4 |

Funnels, polypropylene

| Cat. No | Top diameter, mm | Stem length, mm | Stem outer diameter, mm | Capacity, mL | Pack qty |
|---------------------|------------------|-----------------|-------------------------|--------------|----------|
| Funnels, analytical | | | | | |
| 11947914 | 34 | 52 | 6 | 12 | 36 |
| 11957914 | 48 | 50 | 7 | 23 | 36 |
| 11967914 | 55 | 61 | 7 | 37 | 36 |
| 11977914 | 66 | 65 | 7 | 50 | 36 |
| 11987914 | 77 | 80 | 7 | 100 | 36 |
| 11977904 | 104 | 99 | 9 | 150 | 24 |
| 11987904 | 158 | 151 | 14 | 225 | 24 |
| Funnels, powder | | | | | |
| 11997914 | 65 | 22 | 16 | 50 | 36 |
| 11907924 | 79 | 29 | 16 | 100 | 36 |
| 11997904 | 104 | 33 | 21 | 225 | 24 |
| 11907914 | 147 | 30 | 27 | 750 | 24 |
| Funnels, utility | | | | | |
| 11997934 | 68 | 21 | 14 | 40 | 72 |
| 11907944 | 87 | 22 | 13 | 100 | 72 |
| 11917944 | 109 | 37 | 18 | 140 | 72 |
| 11927944 | 127 | 46 | 19 | 210 | 72 |
| 11917914 | 160 | 66 | 48 | 410 | 24 |
| 11997884 | 203 | 46 | 24 | 500 | 12 |



pH

pH Indicator paper

| Cat.No | Description | Pack qty |
|----------|---|----------------|
| 10642751 | pH indicator paper sticks, non-bleeding, supplied with colour comparison chart 0pH to 14pH, 85mm x 6mm | 100 sticks |
| 10018080 | pH indicator paper sticks, non-bleeding, supplied with colour comparison chart 0pH to 6pH, 85mm x 6mm | 100 sticks |
| 10353641 | pH indicator paper sticks, non-bleeding, supplied with colour comparison chart 3.6pH to 6.1pH, 85mm x 6mm | 100 sticks |
| 10333501 | pH indicator paper sticks, non-bleeding, supplied with colour comparison chart 4.5pH to 10pH, 85mm x 6mm | 100 sticks |
| 10271751 | pH indicator paper sticks, non-bleeding, supplied with colour comparison chart 7.0pH to 14pH, 85mm x 6mm | 100 sticks |
| 10482341 | pH indicator paper reel, 1pH to 14pH, 5m x 7mm | 4 reels of 5m |
| 11532413 | pH indicator paper litmus, 2 to 12pH, 50mm x 10mm | 10 x 20 leaves |
| 11542413 | pH indicator paper litmus, colour change red-blue, 5pH to 8pH, 50mm x 10mm | 10 x 20 leaves |
| 11552413 | pH indicator paper litmus, colour change blue-red, 5pH to 9pH, 60mm x 10mm | 10 x 20 leaves |
| 11592413 | Indicator paper starch / iodide, colour change white-blue, oxidising agents, 60mm x 10mm | 10 x 20 leaves |



Gloves, latex, disposable, PPE Category III

- Latex examination gloves provide a protective barrier and contain soothing Aloe vera to help prevent skin from drying, cracking and chapping
- Textured to provide excellent grip
- Powder free
- EN 374:2003 (AQL 1.5)

| Cat. No | Size | Pack qty |
|----------|------|----------|
| 15153269 | S | 100 |
| 15173269 | M | 100 |
| 15163269 | L | 100 |
| 15183269 | XL | 100 |

Gloves, nitrile, disposable, PPE Category III

- Gloves feature a textured surface to provide ideal grip for dry and wet working conditions
- Aloe vera coating on the inner surface protects hands from drying, chapping and cracking, making this glove ideal for prolonged use
- Beaded cuff
- Mint green for easy identification

| Cat. No | Size | Pack qty |
|----------|------|----------|
| 15193269 | S | 100 |
| 15103279 | M | 100 |
| 15123279 | L | 100 |
| 15113279 | XL | 100 |

Gloves, nitrile, disposable, powder-free

- Disposable examination gloves, non-sterile, powder free purple nitrile, 250mm length

| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11502723 | Gloves, 250mm length, extra extra small XXS | 100 |
| 11512723 | Gloves, 250mm length, extra small XS | 100 |
| 11522723 | Gloves, 250mm length, small S | 100 |
| 10223882 | Gloves, 250mm length, medium M | 100 |
| 11542723 | Gloves, 250mm length, large L | 100 |
| 11552723 | Gloves, 250mm length, extra large XL | 100 |

- Extra protection, disposable examination gloves, non-sterile, powder free purple nitrile, 300mm length

| Cat.No | Description | Pack qty |
|----------|--|----------|
| 10090812 | Gloves, EP, 300mm length, extra small XS | 50 |
| 10611783 | Gloves, EP, 300mm length, small S | 50 |
| 10000822 | Gloves, EP, 300mm length, medium M | 50 |
| 10769824 | Gloves, EP, 300mm length, large L | 50 |
| 10580043 | Gloves, EP, 300mm length, extra large XL | 50 |

Glove box holder

| Cat.No | Description |
|----------|---|
| 11895843 | Glove box holder, epoxy-coated steel wire, wall mount screws included, 3 box capacity |



Laboratory coat, unisex

- Unisex laboratory coat with concealed stud front and single rear vent. Three patch pockets.

65/35 pre-shrunk polyester/cotton, white

| Cat.No | Description |
|----------|--|
| 11558102 | Laboratory coat, unisex, 92cm, S, 100cm long |
| 11568102 | Laboratory coat, unisex, 100cm, M, 100cm long |
| 11578102 | Laboratory coat, unisex, 108cm, L, 100cm long |
| 11588102 | Laboratory coat, unisex, 116cm, XL, 100cm long |
| 11598102 | Laboratory coat, unisex, 124cm, XXL, 100cm long |
| 11508112 | Laboratory coat, unisex, 132cm, XXXL, 100cm long |



Disposable Polypropylene Lab Coats

Latex-free spunbond polypropylene construction lab wear

- Elastic cuffs hold sleeves in place
- Five snaps for a secure closure
- Standard collar
- No pockets
- Fabric weight: 50g/m²

| Cat. No | Colour | Size | Pack qty |
|----------|--------|----------|----------|
| 15431922 | White | Small | 25 |
| 15441922 | White | Medium | 25 |
| 15451922 | White | Large | 25 |
| 15461922 | White | X-Large | 25 |
| 15471922 | White | XX-Large | 25 |
| 15481922 | White | 3X-Large | 25 |



Disposable Polypropylene Lab Coats with Pockets

Single-use lab coats provide economical, lightweight and breathable protection

- Latex-free spunbond polypropylene construction
- Knit cuffs hold sleeves in place
- Three pockets and a five snap closure
- Standard collar
- Fabric weight: 50g/m²

| Cat. No | Colour | Size | Pack qty |
|----------|--------|----------|----------|
| 15491922 | White | Small | 30 |
| 15401932 | White | Medium | 30 |
| 15411932 | White | Large | 30 |
| 15421932 | White | X-Large | 30 |
| 15431932 | White | XX-Large | 30 |
| 15441932 | White | 3X-Large | 30 |





Disposable Polypropylene Facemasks

Provides advanced protection

- Disposable, single use
- Choice of ear loops or tie closure
- Nonwoven polypropylene outer layer
- Thermally bonded nonwoven polypropylene inner layer
- Latex-free
- Choice of blue or white
- Universal size/fit

| Cat. No | Description | Colour | Thickness | Pack qty |
|-------------------|---------------------------------------|--------|-----------|----------|
| 1 15451932 | Disposable face mask with ear loops | White | 2-ply | 500 |
| 15461932 | Disposable face mask with ear loops | Blue | 3-ply | 500 |
| 15471932 | Disposable face mask with ear loops | White | 3-ply | 500 |
| 2 15481932 | Disposable face mask with tie closure | Blue | 3-ply | 500 |

Disposable Polypropylene Bouffant Caps

Single-use cap provides basic protection

- Disposable
- 14gsm polypropylene
- Diameter: 53.3cm or 60.9cm
- Latex-free
- Choice of blue or white

| Cat. No | Colour | Size, cm | Pack qty |
|-------------------|--------|----------|----------|
| 1 15491932 | Blue | 53.3 | 1,000 |
| 2 15401942 | White | 53.3 | 1,000 |
| 15411942 | Blue | 60.9 | 1,000 |
| 15421942 | White | 60.9 | 1,000 |

Disposable Pleated Bouffant Caps

Bouffant's elastic headband is encapsulated in soft poly thread for great comfort and basic protection

- Disposable
- Pleated
- Polypropylene
- Choice of white or blue

| Cat. No | Colour | Size | Pack qty |
|-------------------|--------|---------|----------|
| 1 15431942 | White | Small | 1,000 |
| 15441942 | Blue | Small | 1,000 |
| 15451942 | White | Large | 1,000 |
| 2 15461942 | Blue | Large | 1,000 |
| 15471942 | White | X-Large | 1,000 |

Disposable Beard Cover

Polypropylene facial hair cover for basic protection

- Disposable
- 10gsm spunbonded polypropylene
- Single-use
- Latex-free
- White

| Cat. No | Colour | Size | Pack qty |
|----------|--------|-----------|----------|
| 15471952 | White | Universal | 500 |

Disposable Shoe Covers

Spunbond polypropylene shoe covers protect against dirt, grime and certain dry particulates

- Single-use
- Blue
- Polypropylene
- Latex-free
- With or without anti-skid sole

| Cat. No | Description | Colour | Size | Pack qty |
|----------|------------------------------|--------|-----------|----------|
| 15481942 | Shoe cover without anti-skid | Blue | Universal | 300 |
| 15491942 | Shoe cover without anti-skid | Blue | X-Large | 300 |
| 15401952 | Shoe cover with anti-skid | Blue | Universal | 300 |
| 15411952 | Shoe cover with anti-skid | Blue | X-Large | 300 |



Disposable High-Tack Shoe Covers

High-Tack, disposable shoe covers provide superior protection

- Single-use
- Polyurethane
- Blue
- Latex-free
- High-Tack sole

| Cat. No | Description | Colour | Size | Pack qty |
|----------|-----------------------|--------|-----------|----------|
| 15421952 | Shoe cover, High-Tack | Blue | Universal | 300 |
| 15431952 | Shoe cover, High-Tack | Blue | X-Large | 300 |



Disposable Polyethylene Shoe Covers

Single use, disposable covers provide convenient protection

- Polyethylene
- Blue
- Latex-free

| Cat. No | Description | Colour | Size | Pack qty |
|----------|--------------------------|--------|-----------|----------|
| 15451952 | Shoe cover, polyethylene | Blue | Universal | 300 |
| 15461952 | Shoe cover, polyethylene | Blue | X-Large | 300 |



Thermometers and Timers

Traceable™ thermometer, lollipop, stainless steel

- Featherweight, wide-range lab partner covers 99% of routine temperature measurement needs
- Shock/waterproof unit delivers reading updates every second
- One-key operation recalls MIN/MAX temperature results, °C or °F switchable
- Supplied with battery, probe, probe cover with pocket clip, Traceable™ Certificate
- Temperature range -50 to 300°C, resolution 0.1° from -20 to 200°C; or 1°C elsewhere, accuracy ±1°C

| Cat.No | Description |
|----------|--|
| 11785853 | Thermometer, lollipop, stainless steel, Traceable™ |



Thermometers

| Cat.No | Description |
|----------|---|
| 11561993 | Thermometer, general purpose/stirring, 76 mm immersion, red spirit filled, amber graduations, white backed, -10°C to 110°C x 1°C, 305mm |
| 11587273 | Thermometer, general purpose/stirring, total immersion, red spirit filled, amber graduations, white backed, -20°C to 110°C x 1°C, 305mm |





Traceable™ thermometer/clock/calendar

- Clock has stylish two-tone colour and slim silhouette
- Features: 12/24-hour clock, date/day, temperature
- Supplied with battery, wall-mount, Traceable™ Certificate
- Temperature range -5 to 50 °C, resolution 1 min, accuracy 0.01%

| Cat.No | Description |
|----------|---|
| 11795833 | Thermometer/clock/calendar, 235 x 254mm, Traceable™ |



Traceable™ thermometer/clock/humidity monitor

- Simultaneously displays time of day, temperature and humidity
- Features: MIN/MAX memory, memory clear, °C/°F switchable, 12/24-hour clock
- At the touch of a button, memory recalls the highest and lowest temperature and humidity readings
- Internal sensors make it ideal for use in hoods, storerooms, cleanrooms, incubators, drying chambers and environmental cabinets
- Supplied with stand, wall mount, double-back tape, batteries, Traceable™ Certificate
- Temperature range 0 to 50 °C, resolution 1 min, accuracy 0.01%
- Humidity range 20-90% RH, resolution 1%, accuracy ±5% RH midrange

| Cat.No | Description |
|----------|--|
| 11725843 | Thermometer/clock/humidity monitor, 108 x 57 x 6mm, Traceable™ |



Traceable™ digital humidity/temperature meter

- Continuously and simultaneously displays RH and temperature with dual-recall memories
- Ideal for routine measurements, monitoring round-the-clock, quality control needs, and critical experimental requirements
- Minimum and maximum memories can be displayed or cleared at any time with the touch of a key
- Fulfills CLIA requirements, °C/°F switchable
- Supplied with wall mount, bench, stand, batteries, Traceable™ Certificate
- Temperature range 0 to 50°C, resolution 1°C, accuracy ±1°C
- Humidity range 25-95% RH, resolution 1%, accuracy ±2% RH midrange
- Features min/max memory

| Cat.No | Description |
|----------|--|
| 11765843 | Digital humidity/temperature meter, Traceable™ |



Traceable™ refrigerator/freezer thermometer

- Triple display simultaneously shows MIN/MAX and current temperatures, °C/°F switchable
- Accurately monitor temperatures in freezers, water baths, heating blocks, incubators and refrigerators with the enclosed temperature-buffered sensor
- Bottle insulates sensor from rapid temperature changes when refrigerator door is opened
- Minimum/maximum monitors high and low readings overnight, on weekends, or for any time period
- Alarm is programmable in 1° increments, sounds continuously when temperature rises above/falls below set points and can be turned off manually
- Visual and audible signals continue even if temperature returns to non-alarm range
- High-impact, chemical-resistant ABS plastic case: 70mm x 108mm x 19mm, weight 113g
- Supplied with stand, 3m micro cable, wall mounting, Velcro™, magnetic strips, Traceable™ Certificate
- Temperature range -50 to 70°C, resolution 1°, accuracy ±1°C
- Probe length 19mm, diameter 5mm

| Cat.No | Description |
|----------|--|
| 11873460 | Refrigerator/freezer thermometer, Traceable™ |

Refrigerator/freezer plus™ thermometer, Traceable™

- Temperature-buffered bottle insulates sensor from rapid temperature changes
- Triple display simultaneously shows MIN/MAX and current temperatures, °C/°F switchable
- Monitors high and low readings overnight, on weekends, or for any time period - a significant advantage over current read-only temperatures
- Alarm is programmable in 1° increments and alerts user when temperature rises above or falls below a set point
- Supplied with battery, Velcro™, magnetic strips, 3m cable, Traceable™ Certificate
- Temperature range -50 to 70°C, resolution 0.1°, accuracy ±0.5°C
- Probe length 19mm, diameter 5mm



| Cat.No | Description |
|----------|--|
| 11705853 | Refrigerator/freezer thermometer, Traceable™ |

Four-channel alarm timer, Traceable™

- Features: count up/down, time in/out, memories, alarm, 12/24 hour clock
- Large high display allows viewing from a distance
- Remarkable memory returns display to previously programmed countdown time at the touch of a button
- Extra-loud, high decibel alarms for 1 minute or can be silenced manually
- Timing capacity: 99 hours, 59 minutes, 59 seconds
- Resolution: 1 second
- Accuracy: 0.01%
- Supplied with clip, stand, magnet, battery, Traceable™ Certificate



| Cat.No | Description |
|----------|--------------------------------------|
| 11745863 | Four-channel alarm timer, Traceable™ |

Multicoloured timer, Traceable™

- Three bright, assorted coloured, interchangeable frames supplied with each timer
- Counts down and alarms on four separate channels
- Counts up during alarming indicating elapsed time since alarming
- Counts up in stopwatch mode
- Automatic Bounceback™ memory for repetitive timing
- Times up to 99hr, 59min, 59s
- Settable to 1s with 1s resolution
- Time of day clock (12/24 hour)
- Easel bench stand and magnetic back
- Accuracy: 0.01%



| Cat.No | Description |
|----------|---------------------------------|
| 11775863 | Multicoloured timer, Traceable™ |

You may also be interested in...

Spatula balance




| Cat.No | Description |
|----------|----------------------------|
| 11744416 | Spatula balance Traceable™ |

Foil dispenser



| Cat.No | Description |
|----------|--------------------------|
| 11957994 | Foil dispenser, mini bin |

Storage set, MagLab™



| Cat.No | Description | Pack qty |
|----------|---|----------|
| 11942581 | Storage set, Maglab™, complete lab set of 4 MagLab™ magnetic holders, ABS, assorted colours: blue, green, purple, orange, 522mm x 393mm x 102mm (w x d x h) | 4 |

Transport box, polycarbonate



- Ideal for transporting biological and clinical samples, as well as instruments and products that you want to keep clean and dry under testing environmental conditions
- Polycarbonate construction
- Available in three colours
- Silicone water-tight seal and three lid clasps provide secure closure
- Carry handle folds neatly into body space allowing the units to be stacked
- Lid opens a full 180°, allowing total access to contents and making it easier to clean
- Accommodates 13mm and 16mm tubes in a 72 place Delrin rack or similar sized rack
- Use the separators provided to create your own compartment sizes, or use empty to transport larger products and equipment
- Autoclavable

| Cat. No | Dimensions, (w x d x h), mm | Colour | Pack qty |
|----------|-----------------------------|---------------------------|----------|
| 15165564 | 380 x 196 x 160 | Clear with blue handle | 1 |
| 15175564 | 380 x 196 x 160 | Red with red handle | 1 |
| 15155564 | 380 x 196 x 160 | Yellow with yellow handle | 1 |

Storage tray, polypropylene, sample



- Constructed from corrugated polypropylene
- Moisture resistant, significantly reducing the chance of mold
- Lightweight and sturdy to hold up to 24 full specimen containers
- Stackable with other trays
- Removable dividers allow outer shell to be used for storage or transport of laboratory products
- Supplied flat-packed
- Pack includes: 10 trays and 10 dividers
- Disposable

| Cat. No | Wells | Well diameter, mm | Dimensions, mm (l x w x h) | Colour | Pack qty |
|----------|-------|-------------------|----------------------------|--------|----------|
| 15376548 | 24 | 58 | 406 x 260 x 89 | White | 10 |



Green credentials and services

Please talk to your local representative about our following services:

-  **Pipettor tip packaging recycling service**
-  **Waste packaging collection service**

Here are some of
the ways we are
working to
protect the
environment...



Recycle 22 tonnes of glass every month ●

Collect 40,000 empty bottles from customers every month, these are sent to a recycling company ●

Recycle plastic bottles – 1 tonne of HDPE ●
bottles per month are sent for reprocessing

Send 20 tonnes of cardboard ●
per month to mill for re-pulping

Green-fill™ totally biodegradable packaging ●

Monitor CO² emissions per vehicle on a monthly basis ●

ISO 14001:2004 Environment Accreditation ●

WEEE Compliant ●



Notes

A series of horizontal dotted lines for taking notes.



Dotted lines for form input.

Lab Start-Up Terms & Conditions

1) The only terms and conditions relating to the Fisher Scientific New Lab Start-Up Programme (the Programme) are those set out below. In relation to any question affecting eligibility or entitlement under the Scheme, the discretion of the Fisher Scientific authorised representative shall be final.
2) To be eligible for the Programme, individuals must either be:
a) Setting up a new laboratory; or
b) Moving to a new laboratory; or
c) A first-time grant holder; or
d) Renovating your lab
3) Eligible individuals must place an initial New Lab Start-Up order of at least £2000 with Fisher Scientific to be entitled to claim the promotions as detailed in individual product offers.
4) New Lab Start-Up offers are only valid up until 6 months after the first order has been placed.
5) When redeeming Fisherbrand Vouchers the value of free product must be calculated using the current Fisher Scientific list prices for the year. Voucher cannot be used in conjunction with any other Fisherbrand New Lab Start-Up offer. Voucher cannot be used in conjunction with any other Fisher Scientific promotion. Voucher only applies to Fisherbrand consumable products. One Voucher only per single order. Voucher is available to New Lab Start-Up customers only and all free product must be claimed or redeemed after the sixth month period of the account and before four weeks after account deactivation.

6) When placing a New Lab Start-Up order, please include your unique account number on all New Lab Start-Up paperwork. Please also include Product Claim Forms where relevant.
7) Prices are valid until 31st March 2017 and may be subject to change without prior notice. The stocking order must be placed within 6 months of the account being opened.
8) While Fisher Scientific will take all reasonable steps to avoid alterations to current offers, no guarantee of this is given as offers are subject to change without notice.
9) In the event of cancellation of a New Lab Start-Up order, Fisher Scientific shall be entitled at its option to either:
a) Reclaim all free goods supplied under the Programme; or
b) Issue an invoice for the value of all free goods supplied under the Programme at Fisher Scientific list prices.
10) All sales will be subject to our standard Terms & Conditions with payment due within 30 days of the invoice date.
11) The Programme is available to direct customers of Fisher Scientific and may not be used in conjunction with any other Fisher Scientific promotional offer.



nl new lab
su start-up
programme

Set Your Lab in Motion
2016/2017

© 2016 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

In the United Kingdom:

For customer service, call 01509 555500

To fax an order, use 01509 555111

To order online: www.fisher.co.uk



**Fisher
Scientific**

A Thermo Fisher Scientific Brand