

# Automated DNA isolation

## Automated DNA extraction using InnuPure® systems

### innuPREP/ INSTANT

x – Recommended Kit  
(x) – Recommended with limitations

#### Sample type

	innuPREP Blood DNA Mini Kit-IP	innuPREP Blood DNA Midi Kit-IP	innuPREP DNA Kit-IP	innuPREP DNA Easy Kit-IP	innuPREP Swab DNA Kit-IP	innuPREP Swab DNA Easy Kit-IP	innuPREP Stool DNA Kit-IP	innuPREP Forensic Kit-IP	innuPREP Forensic Easy Kit-IP	innuPREP Plant DNA Kit-IP	innuPREP Bacteria DNA Kit-IP	innuPREP Mycobacteria DNA Kit-IP	innuPREP Virus DNA Kit-IP	INSTANT Bacteria DNA Kit-IP	INSTANT Virus DNA Kit-IP
Bacterial cells (gram+ & gram-)											x			x	
Blood	x	x						(x)	(x)						
Bronchoalveolar lavage (Mycobacteria)												x			
Cell culture supernatants (Virus)													x		x
Cell-free bodily fluids (Virus)													x		
Cerebrospinal fluid (Virus)															
Chewing gum								x	x						
Cigarette butts								x	x						
Dental floss								x	x						
Eucaryotic cells															
Finger nails								x	x						
Forensic materials								x	x						
Hairs, barb hairs, hair roots								x	x						
Liquor (Virus)															
Mycobacteria												x			
Mycobacterial cell pellets												x			
Paraffin embedded tissue			x	x											
Plant material										x					
Plasma (Virus)													x		x
Rodent tails			x	x	x	x									
Serum (Virus)													x		x
Sperm stains								x	x						
Sputum (Mycobacteria)												x			
Stamps and envelopes															
Stool							x								
Swabs			x	x	x	x		x	x						
Tissue samples			x	x	x	x		(x)	(x)		(x)			(x)	
Virus (from various sources)													x		x

Automated DNA extraction using KingFisher® systems

	KFmL						KF96		KFFLX		
	innuPREP/ INSTANT						innuPREP		innuPREP		
	innuPREP DNA Kit I-KFmL	innuPREP Bacteria DNA Kit-KFmL	innuPREP Virus DNA Kit-KFmL	innuPREP Virus DNA/RNA Kit-KFmL	INSTANT Virus DNA Kit-KFmL	INSTANT Virus DNA/RNA Kit-KFmL	innuPREP Stool DNA Kit-KF96	innuPREP Tissue DNA Kit-KF96 & FLX	innuPREP Blood DNA Kit-KFFLX	innuPREP Blood DNA Midi Kit-KFFLX	innuPREP Issue DNA Kit-KF96 & FLX
	x	x							x	x	
			x		x						
			x		x						
			x		x						
			x		x						
	x							x			x
			x		x						
	x							x			x
			x		x						
							x				
	x		x		x						
	x		x		x			x	x		x
			x		x						