

fisherbrand



Fisherbrand Laboratory Guide

This brochure features over 450 specially selected Fisherbrand products, taken from the wider Fisherbrand portfolio of well over 11,000 items, that have all been chosen to make your laboratory life easier! So, we are pleased to offer not only a selection of brand new items (highlighted below), but also an exclusive collection of popular products from our existing catalogue including many devoted to important application areas such as cryogenics, sealing and labelling, weighing and microbiology to name but a few.

For further information on the full Fisherbrand range visit eu.fishersci.com/go/fisherbrand

Featured New Products




Workstation Lab Mat
Pg. 8



Reservoir Chillers
Pg. 11



100-Place Premium Plus Slide Box
Pg. 34



Slide Folder Clasps
Pg. 37



Clinical 4-Way Tube Rack
Pg. 16



PCR Cooler, 96 Well
Pg. 31



Dual Purpose Scoops
Pg. 49



Box Top Dispenser
Pg. 58

Contents



Benchtop Essentials

Centrifuges, Stirrers, Lab Mats & Notebooks

4

Liquid Handling

Pipettors, Reservoirs & Chillers, Pipettor Racks & Stands

10

Racks

Tube, Microtube, Freezer & Storage Racks

16

Histology and Microscopy

Slide Boxes & Mailers, Staining Trays

34

Sample and Blood Collection

Counters, Collection Trays & Storage, Transport Boxes

39

Cryogenics

Boxes, Vial Canes & Sleeves, Storage Rack Systems

43

Weighing

Papers, Dishes, Boats, Scoops & Spatulas

47

Microbiology

Lead Rings, Cell Spreaders & Inoculating Picks

51

Organisation

Bins, Dispensers, Containers & Systems

52

Sealing and Labeling

Parafilm & Foil Dispensers, Heat Sealer, Label Station

58



Centrifuge/Stirrer

Multi-Functional Lab Solution
Saves Valuable Benchtop Space
Perfect for Small Research & Education Labs



Holds 5mL Tubes

High-Speed Mini Centrifuge



Magnetic Stirrer

High-Speed Centrifuge is easy to use and designed to deliver maximum speed (12,500rpm/9,783x g) and flexibility with rapid acceleration and fast controlled breaking

Ideal for most protocols requiring fast spins including PCR, micro-filter cell separation and HPLC. Versatility to handle 0.2, 0.5, 1.5, 2.0 and 5.0mL tube sizes.

Intuitive One-Dial Programming

- Set speed and time with ease
- Adjustable speed from 500 to 12,500rpm
- Speed setting for rpm or RCF
- Adjustable time from 30 seconds to 60 minutes

Safe to Use

Safety features help ensure it will not operate when the lid is open, a rotor is jammed, or the unit is out of balance

High-Performance Magnetic Stirrer is designed using high quality rare earth magnets creating a robust magnetic coupling that significantly reduce the chance of spin-outs.

Perfect for work with sensitive samples and protocols up to 3L requiring aggressive mixing for full dilution.

Easy to Use

- Set stir mode, speed and time with easy one-dial programming:
- Select Mode: clockwise, counter-clockwise and oscillate
 - Adjustable speed from 50 to 2,500rpm
 - Continuous and adjustable time settings from 30 seconds to 60 minutes

Durable Accessories

- Re-usable, easy to clean silicone mats included to keep vessel stable and prevent spills and two rare-earth stir bars (40x12mm) with chemical resistant PTFE coating included

Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
15909350	Grey	9.2 x 7.3 x 4.7	23.5 x 18.5 x 12	1 each
Replacement parts				
16380886		High Speed Centrifuge Rotor 12x1.5/2mL tube		1 each
16390886		High Speed Centrifuge Rotor 6 x 5mL tube		1 each
16300896		Stirrer Rotor		1 each
16320896		Silicone Lid Mat		1 each
16310896		Stir Bar 40x12mm for Centrifuge / Stirrer		2
16370886		Tube Adapters, 0.2mL, 0.5mL		6 each

Package Includes		
Centrifuge/Stirrer	Stirrer Rotor	12 Tube Adapters 0.2mL
5mL Tube Rotor	2 Stir Bars	12 Tube Adapters 0.5mL
1.5mL Tube Rotor	2 Silicone Mats for Stirrer Lid	
1 low voltage, double insulated power adapter with 4 interchangeable plugs		

Specifications: Centrifuge	
Speed Range	500 – 999rpm (increments of 1); 1,000 – 5,000rpm (increments of 5); 5,000 – 12,500rpm (increments of 25)
RCF	16 – 9,783 x g (max RCF calculated)
Run Time Ranges	30 seconds – 60 minutes (5 second increments)

Specifications: Magnetic Stirrer	
Speed Range	50 – 999rpm (increments of 1); 1,000 – 2,500rpm (increments of 5)
Run Time Ranges	Stirrer Mode: 30 seconds to 60 minutes or continuous (5 second increments)
	Oscillating Stirrer Mode: 1 second – 60 minutes or continuous
	Oscillating Period: 30 seconds – 15minutes
Capacity	3L Plus. Diameter of vessel not to exceed 130mm

Patent Information
Community Design No. 003515337-0001-0002
Euro. Pat. App. 3 246 088
US Design Patent No. D814,045
Utility Patent Pending



High-Speed Mini-Centrifuge

Maximum speed with minimal space



This compact centrifuge is ideal for protocols calling for fast spins with a higher RCF requirement. Made of durable ABS with UV-resistant coating.

User Friendly

Intuitive push-button programming and bright LCD screen make it easy to use.

Quiet Operation

Brush-less motor enables quiet operation and low vibration.

Safety Features


Lid safety feature prevents opening while rotor is turning and clear polycarbonate allows user to see rotor movement for added visual safety. Also includes automatic shut-off feature for sample imbalance and tilt.



Cat. No.	Power Supply	L x W x H in	L x W x H cm	Pack qty
12972041	Mini Centrifuge 100-240VAC, 50/60 Hz	8 x 6.7 x 4.5	20.3 x 17.1 x 11.4	1 each
12992041	Universal power cord set			
12902051	Rotor 12 x 1.5/2mL tube (with cover and knob)			
12912051	Rotor 4 x 8 strip 0.2mL tube (with knob)			
12922051	Rotor knob			
12932051	Rotor cover (for 12 x 1.5/2mL tube rotor only)			
12982041	Power adapter			

Package includes:

- High-Speed Mini-Centrifuge
- 12 place 1.5/2.0mL tube rotor
- 32 place 0.2mL tube rotor (singles, or 4 x 8 Strip)
- 12 tube adapters 0.2mL
- 12 tube adapters 0.5mL
- Universal double insulated,
- Power adapter 110-240V 50/60 Hz, with 4 interchangeable cord types



Tool-free rotor changes

Specifications:

Speed: 500-12,500rpm
 Timer Range: 1 to 30 minutes
 RCF: 15-9,783 x g
 Five-year warranty



Micro & Conical Centrifuge Tubes

Pure black for light sensitive samples

Microcentrifuge Tubes:

1.5mL (a)

- Attached flat caps for easy one-handed use
- Frosted writing area
- Syringe pierceable cap

Conical Centrifuge Tubes:

15 and 50mL (b-c)

- Sterile with lot-to-lot traceability
- Writing area and flat polyethylene lids for easy labeling



Certified DNase, RNase and Pyrogen free

Autoclavable

Microcentrifuge Tubes

	Cat. No.	Colour	Vol. mL	Outer Dimensions mm	Closed Height mm	Graduations mL	Max RCF Rating	Packaging	Pack qty
a	15386548	Black	1.5	13.4	40.6	0.1, 0.5, 1, 1.5	18,000 x g	Bulk	500pk

Conical Centrifuge Tubes

	Cat. No.	Colour	Vol. mL	Outer Dimensions mm	Closed Height mm	Graduations mL	Max RCF Rating	Packaging	Pack qty
b	15396548	Black	15	17	120	0.5	9,400 x g	Bulk wrap sterile	500pk
c	15306558	Black	50	29	114	5.0	9,400 x g	Bulk wrap sterile	500pk



Mini-Centrifuge

Mini but Mighty - ideal for small runs

The Mini Centrifuge is ideal for microfilter cell separations and HPLC samples, with adapters to accommodate all major microtube and PCR tube styles (1.5/2.0, 0.5, 0.2mL) for customised use

Easy to Use

Tool-free SnapSpin rotors and quick-spin options make the centrifuges fast and easy to use.

Ergonomic Design

The lid was functionally designed to fit in the palm of the hand for one handed lock and unlock allowing for quick tube extraction. Larger lid size accommodates longer and skirted tubes and tall tube lids.

The housing is made from ABS plastic and lid is polycarbonate plastic enabling the unit to be easily cleaned.

Safety Measures

Fitted with a braking function, the unit stops quickly when powered off or when lid is opened.

A weighted rubberised base eliminates potential movement and dampens vibration and noise.



Accessories and Storage Case Included



Specifications:


Maximum Speed: Approx. 6,000rpm (fully loaded)
 Maximum RCF: 2000 x g
 Power: 100–240 VAC 50/60 Hz
 0.6A, Approx. 7.3 W

Five-year warranty



Package includes:

- Mini-Centrifuge
- 6 place 1.5/2.0mL tube rotor
- 16 place 0.2mL tube rotor (singles, or 2 x 8 strip)
- 6 tube adapters 0.2mL
- 6 tube adapters 0.5mL
- 4 O-rings
- Storage case for rotors and adapters
- Tube rack
- Low voltage, double insulated power adapter with 4 interchangeable plugs



SnapSpin Rotors

change out with the push of a button.
 No tools are needed!

Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
15358266	Grey	6 x 5 x 4.3	15.3 x 12.8 x 10.9	1 each

Replacement Parts

12803342	Replacement 'O' Ring 10 pack
12883332	Rotor SnapSpin 6 x 1/1.5mL or 1.5/2mL tubes
12893332	Rotor SnapSpin 2 x 8 Strip 0.2mL



Specially Designed Inserts
seat the centrifuge, all its accessories,
and includes a tube rack

12 Volt Car Adapter
for on-the-go testing

Portable Centrifuge Kit

Spin Samples On-the-Go

The Portable Centrifuge Kit offers a convenient way to spin down and prepare samples outside of a conventional laboratory environment

Immediate Field Processing

Process samples in un-powered remote areas before deterioration can take place.

Allows for immediate on-site testing for faster results.

Made of durable ABS plastic for easy clean in field environments.

Easy to Use

Centrifuge is equipped with tool-less SnapSpin rotors for easy change out.

Simple operation with lock/unlock activation

Portable Convenience

All items in kit are housed in a handled transport box for portability.

The handled box has an air-tight seal for optimal sample protection.



Cat. No.	Description	Pack qty
15817944	Portable Centrifugation Kit, 100-240VAC, 50/60Hz Universal Plug + 12V Power Adapter	1 each
15827944	12-Volt Mini Centrifuge Power Adapter*	1 each

**Adapter can be used with Mini Centrifuge units manufactured Jan 1, 2016 - present*

Package includes:

- Mini Centrifuge
 - 6 Place 1.5/2.0mL SnapSpin tube rotor
 - 16 Place 0.2mL SnapSpin tube rotor (singles or 2x8 strip)
 - 6 Tube adapters (0.2mL)
 - 6 Tube adapters (0.5mL)
 - 4 O-rings
 - Storage case for rotors and adapters
 - 16 Place tube rack
 - Low voltage, double insulated power adapter with 4 interchangeable plugs



- Transport Box
 - 16 Place tube rack insert
 - Transport box dividers
 - Mini Centrifuge cradle
- Portable 12-volt centrifuge adapter
- Five-year warranty



Reusable Lab Mats

Sustainable, Safe Solutions

Reusable lab mats are an essential lab solution to help keep benchtops organised and clean. They help facilitate workflow processes and are ideal for every type of laboratory

Save Money & Eliminate Waste

Since the lab mats are reusable, they reduce cost and waste making it a more environmentally friendly option over the alternative paper mats.

Adds Stability & Safety

These mats offer a large working surface that is anti-skid which keeps delicate glassware and expensive equipment in place.

The mat's surface is also noise dampening, reducing sound from the vibration of common benchtop equipment.

Made from durable food grade silicone, they are chemically resistance and have a raised edge to contain chemical spills.

Easy to Clean & Store


When the mats are done being used, they can be easily cleaned using any standard disinfectant by simply wiping it, or putting it in a dishwasher and it's ready for repeated use.

Workstation Lab Mat

Safely facilitate workflow processes

Gridded Surface

for tube staging, preventing movement of objects



Reversible,

Side One: combination gridded and smooth surface
Side Two: smooth surface only

- Protect counter tops from hot items up to 200°C (392°F)
- Lip around the edge can hold approx. 250mL of liquid
- Small footprint, 12 x 14.9 x 0.1 in (30.5 x 37.9 x 0.2cm)
- **Autoclavable** at 121° for 15 minutes

Item No.	Colour	L x W x H in	L x W x H cm	Pack qty
16361648	Blue	12 x 14.9 x 0.1	30.5 x 37.9 x 0.2	1 each

Silicone Lab Mat

Large working space


The lab mat can be used with items such as centrifuges, vortex mixers, hot plates, stirrers, pipettors and so much more

Mat features designated magnetic area to keep stir bars in place.

- Protect countertops from hot items up to 200°C (392°F)
- **Autoclavable** at 121° for 15 minutes

Large Working Surface

accommodates benchtop equipment, weigh boats, stir bars, scales, etc.



Raised Edges

enable liquid retention of approximately 300mL of fluid



Reversible

Blue and Grey


Cat. No.	Side 1 Colour	Side 2 Colour	L x W x H in	L x W x H cm	Pack qty
15837944	Blue	Blue/Grey	14 x 23.5 x .1	35 x 60 x .2	1 each



Keep Information Safe
with a durable hardbound, water-resistant cover

Easy Identification
with a foil stamped title

Marker Ribbon
to keep your place



Lined: college-ruled pages for easy report writing

Grid: 6 x 6mm grid size

fisherbrand

Laboratory Notebooks

Keep information safe and well organised with these laboratory notebooks. Record research data, experimental procedures, design and development ideas, and more

- Quality 26 Smyth-sewn binding holds pages securely and lies flat when opened
- Made of acid-free paper, pages are numbered and include space for title, date, project number, book number, date and witness signature
- User data, notebook instructions, table of contents, periodic table, and conversion tables are included in every book

Cat. No.	Colour	Pages	Page Format.	L x W in	L x W cm	Pack qty
15801658	Black	100	Grid	8.6 x 11.1	21.8 x 28.1	6cs
15841658	Blue					
15881658	Red					
16216870	Green	200	Grid	8.6 x 11.1	21.8 x 28.1	6cs
15811658	Black					
15851658	Blue					
15891658	Red	100	Lined	8.6 x 11.1	21.8 x 28.1	6cs
16226870	Green					
15821658	Black					
15861658	Blue	200	Lined	8.6 x 11.1	21.8 x 28.1	6cs
15801668	Red					
16236870	Green					
15831658	Black	100	Lined	8.6 x 11.1	21.8 x 28.1	6cs
15871658	Blue					
15811668	Red					
16246870	Green	200	Lined	8.6 x 11.1	21.8 x 28.1	6cs
	Black					
	Blue					

Fillers and Bulbs

Rubber Bulbs (Natural Rubber)

Cat. No.	Colour	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm	Pack qty
11967984		12	2	51	1.3	34	0.2	4	1 each
11977984		15	2.2	56	1.3	34	0.2	4	1 each
11987984	Red ■	30	2.3	66	1.8	46	0.2	6	1 each
11997984		60	3.1	80	2.1	54	0.2	6	1 each
11314055		150	4.1	103	2.9	73	0.2	6	1 each

- Contoured shape and heavy-wall neck increase grip for easier use with pipe



Safety Bulb Pipettor Filler (Natural Rubber)

Cat. No.	Colour	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm	Pack qty
11779165	Red ■	15	4.2	106	1.6	41	0.2–0.3	4.5–6.5	3pk

- Features a white ABS attachment containing a silicone conical adapter that seals around pipettors without twisting
- Pipettor does not contact hard plastic, which reduces the risk of cracking and chipping
- Elastic band for hanging on reagent bottle or wrist
- Filler fits standard size pipettors

Reservoir Chillers


Modern Design & Enhanced Features

Ice-free, non-toxic way to extend working period and preserve integrity for temperature sensitive reagents

- Reagents are kept at less than 5°C (41°F) for approximately 3 hours
- Engineered custom fit design for 50/55mL or 100mL reservoirs to ensure a maximum contact for longer uniform cooling
- Eliminates the risk of contamination from messy ice bath

Stackable

Saves storage and work space



Symmetrical cooler design

Allows reagent reservoirs to be mounted in either direction

50/55mL



100mL



Reservoirs Chillers

Item No.	Colour	Volume mL	L x W x H in	L x W x H cm	Pack qty
16330412	Purple ■	50/55	5.6 x 2.9 x 1.75	14.4 x 7.4 x 4.5	2 ea/pk
16340412	Blue ■				
16350412	Assorted* ■■	100	6 x 3.5 x 1.75	15 x 9 x 4.3	2 ea/pk
16360412	Purple ■				
16370412	Blue ■				
16380412	Assorted* ■■				

*Assorted colours (Purple and Blue)

Multi-Channel Reagent Reservoirs

Designed to streamline pipetting procedures

- Superior durability with enhanced walls and crush-resistant “turrets” at the corners
- Easy handling with tabs to prevent stacked reservoirs from sticking together
- Sterile versions gamma irradiated and individually sealed in sani-tear HDPE packets in an ISO-certified facility. Traceability certificate available upon request

Cat. No.	Volume mL	Material	Sterility	L x W x H in	L x W x H cm	Pack qty
16371076	55mL	Clear PVC	Non-Sterile	4 x 1.1 x 1.2 in	10.2 x 2.8 x 3.0cm	100pk
16361076	55mL	White Polystyrene	Non-Sterile			
16351076	55mL	White Polystyrene	Sterile			
16344862	100mL	White Polystyrene	Non-Sterile	4.3 x 1.9 x 1.3 in	10.9 x 4.8 x 3.2cm	
16334862	100mL	White Polystyrene	Sterile			



100mL

Easy pour spouts



Tapered "V" Design
Allows for full sample recovery

Dual Solution Reservoir

Innovative shape guides and protects pipettor tips

Multi-functional reagent reservoir has a single compartment 50mL basin on one side and 12 slotted 5mL basins on the other side

- Angled wells guide pipettor tips toward a centre position to help avoid cross-contamination and waste
- Individual 5mL basin wells numerical imprinted to identify samples
- Cover included to prevent contamination

Cat. No.	Volume mL	Sterility	L x W x H in	L x W x H cm	Pack qty
11978084	5/50	Non-sterile	4.9 x 2 x 1.5	12.5 x 5 x 3.7	25pk



Autoclavable



Universal Pipettor Stand

Pipettors never looked so good on your bench top

Stylish pipettor stand holds and stores up to 8 single channel, 4 multi channel, or a combination of both pipettors

Versatile

Patent pending design enables this stand to securely hold virtually every known brand of pipettor in the marketplace

Modular upper holder inserts accommodate the fit of different pipettor widths. Inserts come in two separate colours allowing for colour coding configurations

Durable Utility

Stable weighted base allows for easy one handed removal and replacement of pipettors

Rotating head makes pipettor selection fast and easy

Moving the stand while keeping pipettors in place is simple with a quick grip handle allowing a user to grab and go





Made of Durable ABS

Easy to wipe clean

Vertical Storage

Helps prevent contamination



Inserts

on the upper holders can be added or removed to accommodate the fit of different pipettor widths

US Patent Pending

Cat. No.	Upper Holder Colour	Quantity	Lower Holder Colour	Quantity	Base Width in	Top W x H in	Base Width cm	Top W x H cm	Pack qty
15650785	Blue/Green	8 each/colour	Blue	8 each	6.2	9 x 16.1	15.8	22.9 x 40.8	1 each

Universal Pipettor Wall-Mount

Saves valuable bench top space

The pipettor holder comes with two wall mounting brackets that each hold two pipettors

Secure

Upper and lower segment inserts easily slide in place to securely hold virtually every known brand of pipettor

Segments unique design includes traction grips that centre pipettors allowing for easy one handed removal and replacement


Adaptable

Inserts come in two separate colours allowing for colour coding and quick identification or add a pop of colour to your workspace



Vertical Storage

Helps prevent contamination



US Patent Pending

Cat. No.	Lower Holder	Upper Holder	L x H in	L x H cm	Pack qty
15847944	2 Blue	2 Blue	6.9 x 1.8	17.5 x 4.5	2pk
	2 Green	2 Green			
	2 Clear	2 Purple			
		2 Yellow			



Secure Clamping
for tabletops, benches and shelves

Clamp-ette™ Pipettor Holders

Secure, easy clamping

- Triple format fits most major pipettor brands
- Made of chemical-resistant polypropylene

Autoclavable

Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
11821582	Blue	6.2 x 4.7 x 3.7	15.7 x 12 x 9.5	1 each

Flip & Grip™ Pipettor Holder

Giving labs a choice

Unique pipettor holder that gives you the option of using on the bench top, or hanging under shelf to help maximise your working space


- Holds 12 single-channel, or a combo of 3 single and 3 multi-channel most common brands of pipettors
- Weighted base and rotating rack make for easy pipettor selection
- Rubber grips prevent scratching pipettor



Cat. No.	Colour	Bench Top Style in	W x H cm	Under Shelf Style in	W x H cm	Pack qty
15291588	Blue ■					
15380360	Lime ■	13.1 x 7.6	33.2 x 19.3	7.4 x 7.6	18.7 x 19.3	1 each
15390360	Magenta ■					



Pipettors can be held in the rubber grips and/or the plunger hooks



Mount base undershelf or counter with included hardware



US Design Patent No. D754872
US Design Patent No. D763464



3-Place Stand Acrylic Clear

4-Place Stand Acrylic Clear

6-Place Stand Acrylic Blue

6-Place Stand ABS Blue

Pipettor Stands

Convenient and safe storage

Pipettor stands are designed to hold most major pipettor brands

- Acrylic stands are transparent/clear
- ABS stands are solid opaque colour

Cat. No.	Colour	Material	Pipettor Capacity	Slot Upper in	mm	Slot Lower in	mm	Pack qty
10760774	Clear <input type="checkbox"/>	Acrylic	3					
10122952	Clear <input type="checkbox"/>	Acrylic	4					
11394045	Clear <input type="checkbox"/>	Acrylic	6	1	25	0.8	20	1 each
11384045	Blue ■	Acrylic	6					
10563893	Blue ■	ABS	6					

Cargo Pipettor Rack

Meets the needs of any lab

Customisable for any situation

Racks securely lock together horizontally or vertically to desired configuration or can be used individually

Can be used on top of benchtop or conveniently hung on any metallic surface with strong magnets that are designed into the racks

- Angled shelves store 0.1 to 50mL pipettors
- Protective lid encloses the unit to prevent dust build-up and provides additional shelf space


Made of Durable ABS
Easy to wipe clean



Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
11958004	Blue	11.9 x 3.9 x 5	30.1 x 10 x 12.8	3pk

Magnetic

Can be hung on any metallic surface with strong magnets that are designed into the racks




Pipettor Racks

Protects and organises glass and plastic pipettors

Four angled-shelved compartments keep your 0.1 to 50mL pipettors or cans within easy reach while preventing damage and clutter

- Heavy-duty acrylic version provides crystal-clear viewing from any angle
- ABS versions have cut-out viewing window and come with two removable attachments that hold pipettor bulbs, fillers, and pumps
- Magnetised ABS version has strong magnets on both sides and mounts securely to steel surfaces



Cat. No.	Colour	Material	L x W x H in	L x W x H cm	Pack qty
11703926	Clear	Acrylic			
11793849	Blue	ABS			
11791455	Red	ABS	11.4 x 3.5 x 15.8	29.1 x 8.8 x 40.1	1 each
11757325	Blue	ABS (mounting magnets)			



Pipettor Attachments

ABS version includes 2 removable pipettor attachments

Clinical 4-Way Tube Rack

Multiple combinations for any situation

The only tube rack with a six way connection to create various configurations and designed to hold round, flat & conical blood collection and culture tubes

Versatile

- Use as standalone unit or infinite units can be linked to create configurations to fit number of wells and tube sizes needed
- Accommodates all brands of 13 and 17mm round, flat and conical bottom blood collection & culture tubes
- Assorted colour pack available for colour coding options

Robust Functional Design

- Uniformly spaced wells help eliminate pipetting errors which enables more accurate results
- Moulded alphanumeric identification for quick location and retrieval of tubes
- Made from polypropylene that is chemically-resistant and easy to clean
- Racks can be used in water baths and placed in freezer


Hold Multiple Tube Sizes without multiple racks!



Easy to Configure

Racks slide together to create many combinations

Autoclavable



Blood Tube Brands (include & not limited to)

- Becton Dickinson
- Greiner
- Covidien
- Thermo Scientific™
- Applied Biosystems
- Sarstedt

Recommended Tube Dimensions:

- 13 x 75mm
- 13 x 100mm
- 16 x 100mm
- 16 x 125mm
- 17 x 120mm

Tube Volume Guide:

- 3.5mL
- 7.0mL
- 10mL
- 15mL



Tube Bottoms Nest


in well pockets keeping them from going through the rack and onto the bench top

Holds 15mL Conical Tubes!

Patent Pending

Item No.	Colour	Item No.	Colour
16350402	Blue	16380402	Natural
16360402	Purple	16340402	Assorted*
16370402	Green		

*Assorted colours (Blue, Purple, Green, Red, Natural)



Individual Rack Specifications

Side	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
A	12	13	9.1 x 2.8 x 3.8	23.2 x 7.0 x 9.5	5 ea
B	18	13			
C	12	17			
D	18	17			



fisherbrand

4-Way Tube Rack

Link together to meet your application needs

Versatile

- Use as standalone unit or infinite units can be linked to create configurations to fit number of wells and tube sizes needed

Robust Functional Design

- Moulded alphanumeric identification for quick location and retrieval of tubes
- Made from polypropylene that is chemically-resistant and easy to clean
- Racks can be used in water baths and placed in freezer

Easy to Configure

Racks slide together end to end or side to side to create many combinations

Autoclavable

Cat. No.	Colour	Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11750334	Natural <input type="checkbox"/>	11770334	Yellow <input type="checkbox"/>	50	4	30	4.3 x 4.3 x 4.5	10.8 x 10.8 x 11.3	5pk
11770045	Blue <input type="checkbox"/>	11790045	Orange <input type="checkbox"/>	15	12	17.5			
11760334	Green <input type="checkbox"/>	11700055	Assorted* <input type="checkbox"/>	0.5	32	8.2			
11780045	Pink <input type="checkbox"/>			1.5/2.0	32	11.6			

**Assorted Colours (Blue, Green, Pink, Yellow, Orange)*

fisherbrand

Cube Tube Rack

Interlocking multi-tube racks

Create custom configurations by linking cubes together

- Assorted colours make colour coding easy
- Made from durable polypropylene that can be easily cleaned and stand up to ice bath and freezer conditions



Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11784726	<input type="checkbox"/>	50	4	30	4.3 x 4.3 x 4.5	10.8 x 10.8 x 11.3	5pk
		15	10	17			
		3-7	12	12.5			
		1.5/2.0	16	11.5			

Autoclavable

**Assorted colours (Blue, Green, Pink, Yellow, Orange)*

Tube Racks 5mL/15mL/17mm and 10-17mm

Exclusive design, optimal functionality


Unique design allows for flexibility while working with a multitude of tube diameters and heights

Customisable

- Rack has two hinged shelves that can be pivoted in either an upper or lower position to securely accommodate the size of tube being used. Lower position is ideal height when working with shorter microcentrifuge tubes

Improved Functionality

- Open design allows for better filling visibility and tube graduations
- Rack wells have a unique teardrop design that enables better retention of tubes over traditional designs




Autoclavable

Embossed Index Numbers for quick sample ID and retrieval

Cat. No.	Colour	No. of Included Segments	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
15569115	Green/Blue	2	5-15	20	17	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	1 each
15728177	Purple/Natural	2	5-15	20	13	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	1 each
		2	5-10	20	17			

Alphanumeric well identification



Easy tube removal
Tapered well design suspends tube for easy removal

Autoclavable

Rota-Rack™ Duo Tube Rack

Rotate to customise workflow processing

Transfer solutions between 0.2mL PCR tubes and 0.5 and 1.5mL microcentrifuge tubes

- Purple section holds 48 x 0.5mL and 48 x 1.5/2.0mL microfuge tubes on the reverse side
- Blue section holds 96 x 0.2mL tubes and 48 x 1.5/2.0mL tubes on the reverse side

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11304095	Blue/	0.2	96	5	11.6 x 4.6 x 1.3	29.6 x 11.8 x 3.4	1 each
	Purple	0.5	48	8			
		1.5/2.0	96	11			

Rota-Rack™ Tube Racks

Modular, fun rotating rack

Each module is interlocked and adjustable with just a twist allowing you to work with multiple sized tubes in the same rack

- Suitable for variety of culture, centrifuge, micro-centrifuge, round, flat or conical bottomed tubes
- Two versions to accommodate small and large tube sizes

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11384085	Yellow	50	2	29	8.8 x 4.3 x 4.3	22.4 x 11 x 11	1 each
	Green	10-15	8	20			
	Blue	15-20	6	16			
	Purple	10-13	10	13.5			
11394085	Yellow	15	6	16.4	8.2 x 3.9 x 3.9	21 x 10 x 10	1 each
	Green	1.5/2.0	9	11.2			
	Blue	0.5	12	8			
	Purple	0.2	32	6.5			

Alphanumeric well identification



Autoclavable



Adapt-a-Rack™ Adaptable Multi-Tube Rack

Flexible & linkable to your needs

Made of POM material with flexible fingers that adapt to accommodate 12-30mm tube diameters and hold them firmly in place

- Racks can be linked together for additional capacity
- Open sides allow for better label and fill level viewing
- Available in solid or a combination of colours for easy coding identification

Fits
25mL
Conical
Tubes



Autoclavable

Cat. No.	Rack No.1 Colour		Rack No.2 Colour	
15300370	Blue/Green	■ ■	Green/Blue	■ ■
15320370	Blue/Yellow	■ ■	Yellow/Blue	■ ■
15330370	Pink/White	■ □	White/Pink	□ ■
15340370	Blue/Blue	■ ■	Blue/Blue	■ ■
15350370	White/White	□ □	White/White	□ □
15360370	Pink/Pink	■ ■	Pink/Pink	■ ■
15380370	Yellow/Yellow	■ ■	Yellow/Yellow	■ ■
15390370	Green/Green	■ ■	Green/Green	■ ■

Tubes mL	Wells	Well Dia. mm	L x W x H		Pack qty
			in	cm	
7-50	4	12-30	7.1 x 2.2 x 3	18.1 x 5.6 x 7.6	2pk



Chain Tube Rack

Thinking outside the traditional tube rack "box"

Unique & Versatile Design

- Individual racks link together to create a chain that you can configure in multiple shapes and sizes improving workflow process
- Each well has the flexibility to hold 5mL tubes on one side or rotate it 180° to hold 15mL or 17mm tubes (round and conical bottom).
- Open design allows for better visibility of graduations when filling
- Assorted colours to colour code tubes or testing protocols

US Patent pending. European Patent App. 3 248 684

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H mm	Pack qty
15660785	■ □ ■ ■*	5-15	16	17.2	3.7 x 1	94 x 27	16pk

*Assorted colours (Green, Clear, Yellow and Blue)

One-Well Connecting Tube Racks

Holds 50mL tubes

Form a line or use as a single, stand-alone rack

- Open design enables easy identification of samples
- Moulded circular base ring keeps conical or round-bottom tubes in place

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11727285	■ ■ ■ ■ ■*	50	1	29	3.2 x 2.8 x 2.6	8.3 x 7 x 6.7	5pk

*Assorted colours (Blue, Green, Yellow, Orange, Pink)




Two Tier Microtube Rack 20-Well

Interlock multiple racks to increase tube storage

Autoclavable

Improved Experience

- Tiered rack allows for improved staging and visibility for easier, more accurate sample filling and retrieval
- Made of easy to clean polypropylene
- Easily remove tubes using finger slot on front of the rack
- Embossed index numbers for quick sample ID and retrieval



Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11784016	Blue	1.0	20	10	8.3 x 3.5 x 2	21.1 x 9 x 5	6pk
		1.5					
		2.0					

Universal Tube Rack

Ideal for holding a wide range of tubes

Improved Experience

- Holds up to 12 similar size tubes
- Alphanumeric identification for quick location and retrieval
- Side handles for easy transport
- Made of polypropylene



Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11987994	White/Purple	15-50	12	16-30	7.6 x 5.5 x 2	19.4 x 14 x 5.1	1 each

Transformable Cube Tube Rack

Lidded tube rack & storage

Autoclavable

Quick change over with two snap-in inserts that hold 15 and 50mL tubes

- Rack keyed and numeric wells for tube identification
- Polypropylene construction capable of handling extreme storage temperatures down to -80°C when used with appropriate temperature viable tubes

Compact Storage

- Clear lid included for viewable, stackable storage
- Open space available under insert for extra storage



Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11700634	Purple	15	16	18	4.9 x 4.9 x 5	12.5 x 12.5 x 12.9	5pk
		50	9	31			

Snap-Together Conical Tube Rack

Space Saver - shipped flat and easy-to-assemble

Conveniently holds 15 and 50mL tubes simultaneously

- ABS racks snap together securely and are designed for stable use in waterbaths
- Alphanumeric identification for quick location and retrieval



Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
10753257	Orange	15	30	17	8.1 x 6.7 x 2.4	17.1 x 17 x 6.1	1 each
		50	20	28			

Pop-Up Tube Racks

Collapse for compact storage and "pop-up" for quick use

Designed to save space!

- Easy to open, holds 15 and 50mL tubes
- Easy to close, stackable for compact storage
- Alphanumeric well identification for quick location and retrieval of samples

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm
11710634	Blue	15, 50	21, 12	17, 30
12269309	Pink			
12693466	Purple	15	45	17
12603476	Green	50	18	30

L x W x H				Pack qty
(Expanded) in	(Collapsed) in	(Expanded) cm	(Collapsed) cm	
10 x 5.4 x 2.8	10 x 5.4 x 0.8	25.5 x 13.7 x 7.2	25.5 x 13.7 x 2.1	2pk



Tube Racks

Tough and highly resistant to acids, bases, solvents and heat

Durable racks are moulded in a single, continuous piece so no assembly required and no detachable pieces

Intuitive Design

- Unique tear-drop shaped openings offer better retention of tubes
- Moulded in coordinates for tube orientation
- Flat labeling area for easier identification
- Stable design for use with water baths

Available in two materials:

Fibreglass-Reinforced Polyoxymethylene (POM)

- Available in half and full rack sizes for 5 different tube diameters
- Offers the highest chemical and heat resistance (-40 to +121°C)
- **Autoclavable**

Polypropylene

- Available in full rack sizes for 2 different tube diameters
- Offers good chemical and heat resistance (0 to 110°C)





Full-Size POM

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H	Pack qty
15340380	White					
15350380	Blue					
15360380	Green	5–10	72	13	20.2 x 10.4 x 5.9	1 each
15370380	Magenta					
15380380	Yellow					
15390380	White					
15300390	Blue					
15310390	Green	5–10	72	16	25 x 12.7 x 7	1 each
15320390	Magenta					
15330390	Yellow					
15340390	White					
12327629	Blue					
15350390	Green	10–18	40	20	25.2 x 10 x 8.3	1 each
15360390	Magenta					
15370390	White					
15380390	Blue					
15390390	Green	10–18	40	25	30 x 12 x 9.2	1 each
15300400	Magenta					
15310400	White					
10257963	Blue					
15330400	Green	15–50	24	30	28.2 x 11 x 8.5	1 each
15350400	Magenta					
15360400	Yellow					

Half-Size POM

Cat. No.	Colour	Tubes	Wells	Well Dia. mm	L x W x H cm	Pack qty
15390400	White					
15300410	Blue					
15310410	Green	5–10	36	13	10.4 x 10.4 x 5.9	1 each
15330410	Magenta					
15519821	Yellow					
15340410	White					
15360410	Blue					
15380410	Green	5–10	36	16	12.7 x 12.7 x 7	1 each
15390410	Magenta					
15529821	Yellow					
15310420	White					
15340420	Blue	10–18	20	20	10 x 12.7 x 8.3	1 each
15350420	Green					
15370420	Magenta					
15380420	White					
15300430	Blue	10–18	16	25	12 x 12.2 x 9.2	1 each
15310430	Green					
15320430	Magenta					
15330430	White					
12317629	Blue	15–50	9	30	11 x 11 x 8.5	1 each
15340430	Green					
15360430	Magenta					

Full-Size Polypropylene

Cat. No.	Colour	Tubes	Wells	Well Dia. mm	L x W x H cm	Pack qty
15370430	White					
15380430	Blue					
15390430	Green	5–10	72	13	20.2 x 10.4 x 5.9	1 each
15300440	Magenta					
15310440	White					
15320440	Blue	5–10	72	16	25 x 12.7 x 7	1 each
15340440	Green					
15360440	Magenta					

Reversible Tube Racks 96-Well

Easy-to-grasp tabs for sample handling and storage

96 wells for 0.5mL tubes on one side and 96 wells for 1.5/2.0mL tubes on the reverse side

Optimum Function & Identification

- Wells arranged in standard 8 x 12 array
- Moulded-in gridlines separate 0.5mL wells making them easier to see
- Raised rings around the 1.5mL/2.0mL wells make tubes easier to remove
- Alphanumeric identification for quick location and retrieval of tubes

Secure lid

- Translucent lid snaps onto rack tabs for secure, visible storage
- Spacious enough to fit most screw-top micro-centrifuge tubes
- Stackable, saving valuable space

Cat. No.	Colour	Cat. No.	Colour
11728174	Assorted*	11968014	Green
11738174	Natural	11978014	Yellow
11958014	Blue		

*Assorted colours (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
1.5/2	96	11	9.6 x 4.7 x 2	24.4 x 12 x 5.1	5pk
0.5		7.6			



Autoclavable

Microtube Racks 80-Well

Popular fraction-collector-plate style rack


Compact 5x16 array makes filling and transfer easy

- Frosted labeling area on bottom of rack
- Embossed index numbers for quick sample ID and retrieval
- Made of polypropylene

Cat. No.	Colour
11720344	Natural
11968074	Blue
11978074	Green
11988074	Pink
11998074	Yellow
11908084	Orange
11728084	Assorted*
11918084	Red

*Assorted colours (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
1.5/2.0	80	11	8.9 x 2.6 x 1	22.5 x 6.7 x 2.8	5pk




Autoclavable

Microtube Racks

Handles for portability

Durable, Simple to Use

- Made from reinforced nylon offering superior chemical and heat-resistance over most other plastic based racks
- Handles provide easy and safe transportation
- Stable design can be used in water baths
- Alphanumeric identification for quick location and retrieval of tubes



Autoclavable


Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
15377997	Blue ■	1.5/2.0	96	11.1	8.3 x 5.2 x 2.1	21 x 13.1 x 5.4	4pk

Mega Rack™

High-capacity storage for your tubes

Fully-assembled, polypropylene mega racks are great for high-throughput labs requiring large sample storage

- Sturdy, material for easy and secure transport of glass tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Place different test groupings into one rack eliminating the need for multiple racks
- Not recommended for autoclaving.




Double Rack

Single Rack

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11720634	Blue ■	5-7	216	13	13.3 x 9.3 x 3.6	33.7 x 23.5 x 9.1	1 each
14801442	White □	5-7	216	13	13.3 x 9.3 x 3.6	33.7 x 23.5 x 9.1	1 each
11730634	Blue ■	10	120	16	13.3 x 9.3 x 3.6	33.7 x 26.7 x 9.1	1 each
14811442	White □	10	120	16	13.3 x 9.3 x 3.6	33.7 x 26.7 x 9.1	1 each

Double Rack

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
15509811	White □	5-7	432	13	26.5 x 9.3 x 3.1	67.3 x 23.5 x 9	1 each
15519811	Blue ■	5-7	432	13	26.5 x 9.3 x 3.1	67.3 x 23.5 x 9	1 each
15529811	White □	10	240	16	26.5 x 10.5 x 3.5	67.3 x 26.7 x 9	1 each
15539811	Blue ■	10	240	16	26.5 x 10.5 x 3.5	67.3 x 26.7 x 9	1 each



Single Rack

One-Well Friction-Fit Tube Rack

- Sturdy base made from polystyrene
- Open design enables easy identification of samples
- Holds conical and round bottom tubes

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11354055	Blue ■	50	1	29	3.6 X 1.2	9.2 x 3	5pk



Expanding Storage Box

Single rack, double function

Streamline and reduce sample handling time with this convenient working and storage rack, eliminating the need for several racks when working with samples for cold temperature storage

Efficient Handling

- Expand rack to opening position instantly creating more working space. Enough room for sample manipulation with no interference from open microfuge lids
- Collapse rack into closed position when it's ready for freezer storage
- Compact design will fit into standard 53mm freezer racking system
- Alphanumeric identification for quick location and retrieval of tubes

Safe and Stackable Storage

- Clear lids locks the rack together for safe handling and transportation
- Low profile, easily stackable



Durable Construction

Sides made from polycarbonate, tube holders made of polypropylene

Autoclavable



Clear Lid

Locks rack together for safe handling and transportation

Expanding Storage Box

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H (Closed) in	L x W x H (Open) in	L x W x H (Open) cm	L x W x H (Closed) cm	Pack qty
15162855	White	1.5/2.0	64	12.3	12.2 x 5 x 1.3	5.2 x 5.2 x 2	31 x 12.8 x 3.4	13.3 x 13.3 x 5.3	1 each

Expanding Storage Box

- Holds 64 x 1.5/2.0mL cryo and/or microtubes
- Expanding wells have a centre depression to stabilise conical tubes and a cryogenic vial lock for one handed cap removal





TracRack™ Convenient Rack System

Comprised of a storage box and eight removable tube segments allowing you to easily take what you need - an entire segment or individual sample - and store the rest!

- Printed index letters on each segment make for quick sample ID and retrieval

Robust Storage

- Moisture and mould resistant
- 2 in. (5.3cm) box made from thin film polypropylene material fits in standard freezer racks
- Store down to -80°C. Must use tubes capable of being stored at the required temperature of the user

Available in two sizes:

TracRack™ Tube Rack

- Eight polypropylene segments hold 1.5–2.0mL cryo and microcentrifuge tubes (64 tubes)
- Wells feature a centre depression stabilising conical tubes and a cryogenic vial lock for one-handed cap removal

TracRack™ MicroTube Rack

- Eight reversible polypropylene segments that hold 64 x 0.5mL tubes, 96 x 0.2mL micro-tubes, or a combination of both



TracRack™ Tube Rack

Cat. No.	Colour	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
15346548	Purple	■	1.5/2.0	64	12.3	5.1 x 5.1 x 1.4	12.9 x 12.9 x 3.5	1 each

TracRack™ Microtube Rack

Cat. No.	Colour	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
15356548	Lt. Purple	■	0.2 0.5	96 64	5.1 7.8	5.1 x 5.1 x 1.4	12.9 x 12.9 x 3.5	1 each

Stylish & Identifiable

- Made of translucent, vibrantly coloured polypropylene allowing for visibility inside box and colour coding identification
- Tube capacity is indicated on the outside the box within a distinctive wave design
- Quick visual identification with uniform wells and alphanumeric location identification along with a frosted labeling area

Robust

- Stackable, saving valuable space in the lab
- Store down to -80°C with temperature viable tubes

Autoclavable

fisherbrand

50-Well Microtube Storage Boxes

Ideal for sample or case study storage

- For smaller storage requirements of 1.5/2.0mL microtubes or cryogenic vials.
- Easy-open hinged lid with snap clasp

Cat. No.	Colour
11704316	Natural <input type="checkbox"/>
10707992	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

*Assorted colours (Blue, Green, Purple, Yellow, Orange)

Tubes mL	Wells	Well Dia. mm	L x W x H in	cm	Pack qty
1.5/2.0	50	12.3	5.6 x 3.6 x 2.2	14.2 x 9.2 x 5.6	5pk



fisherbrand

81-Well Microtube Storage Boxes

Compact footprint for microtube and cryogenic vial storage

- Store 1.5/2.0mL microtubes or cryogenic vials in polypropylene storage boxes
- Easy open, friction-fit lid
- Black, light blocking version available for light sensitive samples

Cat. No.	Colour	Cat. No.	Colour
10500203	Natural <input type="checkbox"/>	10325002	Yellow <input type="checkbox"/>
10273222	Blue <input type="checkbox"/>	10385042	Orange <input type="checkbox"/>
10050182	Green <input type="checkbox"/>	15558319	Black <input type="checkbox"/>
10070182	Pink <input type="checkbox"/>	10243272	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>

*Assorted colours (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	cm	Pack qty
1.5/2.0	81	12.3	5.1 x 5.1 x 1.9	13 x 13 x 4.7	5pk



fisherbrand

100-Well Microtube Storage Boxes

Secure storage for microtube and cryogenic vials

- For larger storage requirements of 1.5/2.0mL microtubes or cryogenic vials
- Black, light blocking version available for light sensitive samples
- Easy-open hinged lid with snap clasp

Cat. No.	Colour	Cat. No.	Colour
11790334	Natural <input type="checkbox"/>	15579811	Green <input type="checkbox"/>
11918014	Black <input type="checkbox"/>	15589811	Yellow <input type="checkbox"/>
11700344	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	15459123	Orange <input type="checkbox"/>
15569811	Blue <input type="checkbox"/>		

*Assorted colours (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	cm	Pack qty
1.5/2.0	100	12.3	5.6 x 6 x 2.2	14.1 x 15.1 x 5.7	5pk



Low-Temp PCR Rack: 96-Well

Slide, set, work & store PCR samples

Multi-purpose rack converts from compact workstation to lidded storage solution accommodating individual and strip 0.2mL tubes

Versatile & Robust

- Polycarbonate lid serves as cover and elevated rack base - just slide open and pivot lid under rack for compact workstation
- Lid can be used in right or left handed orientation and is detachable
- Stackable saving valuable bench-top and freezer space
- Multiple colours offered for coding options
- Store down to -80°C with appropriate temperature viable tubes

US Patent Pending, Community Design No: 004101756-0001-0002



Cat. No.	Colour	Cat. No.	Colour
15758177	Assorted*	15778177	Green
15768177	Blue	15788177	Purple

*Assorted colours (Purple, Blue, Green, Orange, Yellow)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
0.2	96	5.2	4.9 x 3.6 x 1.5	12.5 x 9.2 x 3.7	5pk

PCR Tube Racks 96-Well

Perfect for pre- or post-PCR sampling

- Made of polypropylene
- Accommodates individual tubes or strips
- Removable hinged lid lies flat
- Alphanumeric identification for quick location and retrieval of tubes

Autoclavable

Cat. No.	Colour	Cat. No.	Colour
11710344	Assorted*	11948074	Blue
11527593	Natural	11958074	Green

*Assorted colours (Blue, Green, Pink, Yellow, Orange)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
0.2	96	5.2	5 x 3.9 x 1.2	12.7 x 9.9 x 3.0	5pk



PCR Tube Racks Reversible

Removable hinged lid fits both sides

- Made of polypropylene
- PCR side of the rack holds 8 or 12 x 0.2mL tube strips or individual tubes
- Reverse side holds 0.5mL tubes
- Both sides of the rack have 12 wells for 1.5mL tubes

Autoclavable

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11787265		0.2	168	5.5	8.1 x 5.2 x 2.1	20.6 x 13.1 x 5.3	5pk
		0.5	40	7.6			
		1.5/2.0	12	11.4			

*Assorted colours (Blue, Green, Pink, Yellow, Orange)



PCR COOLER, 96-WELL

Ice-Free Cooling System For PCR Preparation and Storage

Cools industry standard format PCR tubes, strips and single format as well as most 96-well plates, in both standard and low profile well shape


Consistent & Reliable

- Modern design cools faster in freezers and stays cold longer on the bench top
- Maintains the entire 96-well plate at a temperature of 0°C for over an hour
- Maintains sample temperature of less than 5°C for up to three hours (with lid off)
- When the temperature ceiling of 6°C has been exceeded, users are signaled by a contrasting colour change of the cooler or a single well

96 well 8x12 configuration



Changes Colour when the Temperature Changes!



Two Colour Contrasting Options Available

purple to pink or dark blue to light blue allowing easy monitoring of samples



Ergonomic Textured Finger Grips

provide a secure hold, reducing the chance of dropping the PCR Cooler, and makes handling and transportation easy, especially when using gloves

Vented Lid

enables faster cooling of samples while in freezer

Lid doubles as a stand

When used as a stand the space between the cooler and bench top enables the cooler to stay cold longer as there is no conductive heat transfer

Alphanumeric Well Identification

for quick location and retrieval of tubes

Item No.	Warm Colour	Cold Colour	Wells	Tube Sizes mm	L x W x H in	L x W x H cm	Pack qty
16300412	Pink ■	Purple ■	96	0.2, 0.5 Assay Plates	4.3 x 5.8 x 1.9	11 x 14.8 x 4.8	2 ea/pk
16310412	Light Blue ■	Dark Blue ■					
16320412	Assorted*						

*Assorted colours (Blue, Pink)

Refrigerator Tube Racks

Contoured design for safe handling

Fill lower section with water or crushed ice to chill sample tubes quickly

- Made from polystyrene
- Top inset tray secures tubes in place
- Stackable, saves bench space

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11304085	White	<input type="checkbox"/> 15 & 50	15	29.2	10.9 x 6.4 x 5.3	27.6 x 16.4 x 13.4	1 each



Cool Containers

Durable containers for chilling temperature-sensitive samples

Heavily insulated for ultra-cold conditions, these Cool Containers are ideal for ice, ice-salt slurries and more

- Made from single-moulded polyurethane
- Close-fit lid keeps evaporation to a minimum and features comfortable recessed handles
- Lid nests under bucket
- For low-temperature use to -196°C at 1 bar
- Stackable
- Pour spout and large handles

Cat. No.	Style	Capacity	L x W x H	Colour	Pack qty
11314085	Bucket	4 L	181 (H) x 234 (D)	Black onyx ■	1 each
11324085	Pan	1 L	267 x 195.5 x 114	Emerald ■	
11334085	Pan	4 L	381 x 254 x 140	Arctic blue ■	
11344085	Pan	9 L	521 x 292 x 163	Ruby ■	



Lid fits below
bucket to
save space



Foam Tube Racks

High storage, minimal footprint

Two-piece racks ideal for long-term studies or tube storage

- Comprised of lightweight expanded polystyrene (EPS) foam
- Stack securely with nesting features on lid and base
- Extra spacing between wells allows for easy gripping with fingertips
- Store down to temperatures of -80°C when used with temperature viable tubes



Cat. No.	Colour	Tubes mL	Wells	Well Dia mm	L x W x H in	L x W x H cm	Pack qty
11774016	White □	1.5/2	100	11.5	13.2 x 3.7 x 2.9	33.6 x 9.5 x 7.3	10pk
11948084	White □	1.5/2	50	11.5	8.3 x 4.3 x 2.7	21 x 11 x 7.1	10pk


Floating Microtube Rack

For use on bench and in ice and water baths

- Made of polypropylene

Detachable legs create multi-functional usage

- Attach legs and rack stands upright for use on bench top
- Detach legs and samples can be quickly removed by pushing down onto any surface



Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
13436959	White □	1.5/2.0	16	10.9	Square 4	10.2	4pk
13456959			20		Round 3.8	9.8 dia.	

Floating Tube Racks

Reusable and flexible for a variety of applications

- Ideal for floating tubes in water baths, tubs, or beakers
- Brightly-coloured foam racks made of high-density polyethylene that resists water absorption
- Racks include detachable carrying handles—easily lift samples out of baths



Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty
11704736	Green ■	50	4	29	6.4 x 5.3 x 2.5	16.5 x 13.7 x 6.5	5pk
11714736	Blue ■	15	8	16	6.7 x 2.2 x 2.5	17.2 x 5.7 x 6.5	
11724736	Yellow ■	1.5/2.0	24	9.3	6 x 3.2 x 2.7	15 x 8.3 x 7	
11736565	Blue ■	0.2/0.5/1.5/2.0	18	9.3	3.2	8.2 dia.	

"X" Design
to hold a variety of tube sizes



Sliding Lock
provides secure storage with easily identified lock and unlock symbols.

Slide Rail Design
ensures a secure fit allowing easy placement and removal of slides

Centrally Positioned Numbering
enables simple slide organisation with easy viewing of location numbers

Patent Pending

fisherbrand

100-Place Premium Plus Slide Box

Modern design for maximum functionality

Optimal Organisation

- Organise by "project" and "date" moulded-in on lid
- Inventory card index corresponds to slide numbers and is removable for trouble-free writing on flat surface.
- Vibrant translucent colours are easy to identify and enable coding options

Functional Design

- Flexible storage, store vertically or horizontally to efficiently use space
- Moulded feet allow secure horizontal stacking of multiple units
- Unique side storage can accommodate a pen, desiccant or other small items needed for working
- Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement

Safe & Secure

- EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- Stylish wave surface texture is easy to hold for safe sample transport

100-Place Premium Plus Slide Box

Cat. No.	Colour	Cat. No.	Colour
15810086	Blue	15840086	Red
15820086	Orange	15850086	Purple
15830086	Green	15860086	Assorted*

*Assorted colours (Purple, Blue, Green, Red, Orange)

Lining	Hinge	L x W x H in	L x W x H cm	Pack qty
EVA Foam	Sliding Latch	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2	5pk

Inventory Card

Cat. No.	L x H in	L x W cm	Pack qty
15870086	8.5 x 11	21.6 x 28	10pk (Polybag with Label)

Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (26 x 76mm and 1 x 3 inch)

Organise & Identify

- Lining and slots are numbered to correspond to the slide inventory sheet
- Organise by "Project" and "Date" with pre-moulded labeling area
- Available in multiple colours for project coding

Durable & Safe

- ABS plastic boxes available with either cork or foam lining
- Designed to securely store slides in grooved slots
- Moulded feet secure stacking for easy transport and storage surfaces and prevents misplacement



25-Place Slide Box

Easily access slides with friction-fit cover



Cat. No.	Colour	Lining	Hinge	L x W x H in	L x W x H cm	Pack qty
11354135	White <input type="checkbox"/>	Cork	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5	1 each
11314135	Blue <input type="checkbox"/>					
11324135	Green <input type="checkbox"/>					
11334135	Red <input type="checkbox"/>					
11711476	Blue <input type="checkbox"/>	Foam	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5	1 each
11324135	Green <input type="checkbox"/>					
11721476	Red <input type="checkbox"/>					
11731476	White <input type="checkbox"/>					

50-Place Slide Box

Thumb-latch lock mechanism



Cat. No.	Colour	Lining	Hinge	L x W x H in	L x W x H cm	Pack qty
11731486	Blue <input type="checkbox"/>	Cork	Thumb-Latch	8.2 x 3.4 x 1.4	20.9 x 8.6 x 3.5	1 each
11771486	Red <input type="checkbox"/>					
11781486	White <input type="checkbox"/>					

100-Place Slide Box

Closes securely with a rust-resistant nickel-plated clasp and pin hinges



Cat. No.	Colour	Lining	Hinge	L x W x H in	L x W x H cm	Pack qty
11314145	Blue <input type="checkbox"/>	Cork	Nickel Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4	1 each
11374145	Black <input type="checkbox"/>					
11324145	Green <input type="checkbox"/>					
11334145	Red <input type="checkbox"/>					
11344145	Yellow <input type="checkbox"/>					
11354145	White <input type="checkbox"/>					
11364145	Grey <input type="checkbox"/>					
11701486	Blue <input type="checkbox"/>	Foam	Nickel Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4	1 each
10428241	Green <input type="checkbox"/>					
11711486	Red <input type="checkbox"/>					
11394135	Yellow <input type="checkbox"/>					
11721486	White <input type="checkbox"/>					

100-Place Tall Slide Box

Additional space for larger 3 x 2 in. slides.
Stackable with 100-place Premium Slide Boxes



Cat. No.	Colour	Lining	Hinge	L x W x H in	L x W x H cm	Pack qty
15186823	Blue <input type="checkbox"/>	Cork	Nickel Clasp	8.2 x 6.9 x 2.4	20.8 x 17.5 x 6	1 each

Standard Microscope Slide Boxes

Constructed of durable polypropylene

Improved Identification

- Moulded slot numbers correspond directly to the slide inventory card
- Organise with "Project" and "Date" pre-moulded labeling area

Secure & Stackable

- Hinged lids are held securely closed with a latch lock
- Moulded feet secure stacking for easy transport and storage

25-Place Slide Box

Cat. No.	Colour	Lining	Hinge	L x W x H in	L x W x H cm	Pack qty
11741476	Blue	<input type="checkbox"/>	None	5.6 x 3.6 x 1.2	14.1 x 9.2 x 3	1 each
11751476	Red	<input type="checkbox"/>	Latch			
11761476	White	<input type="checkbox"/>	Lock			

100-Place Slide Box

Cat. No.	Colour	Lining	Hinge	L x W x H in	L x W x H cm	Pack qty
11771476	Blue	<input type="checkbox"/>	None	8.3 x 6.7 x 1.5	21 x 16.9 x 3.7	1 each
11781476	Red	<input type="checkbox"/>	Latch			
11791476	White	<input type="checkbox"/>	Lock			



Holds standard slides
(26 x 76mm or 1 x 3 inches)

Cardboard Slide Tray

Cardboard slide trays hold standard slides
(26 x 76mm and 1 x 3 inches)

Cardboard slide folders offer a convenient solution for safely working, transporting and short-term storing of slides

Organise & Protect

- Index on the front cover corresponds with the number grid inside for identification
- Recessed slide compartments and inside covers for specimen protection
- Slides cannot be accessed when cover is closed reducing contamination

Sturdy Construction

- Compartment holder is made from a double layer of densely pressed cardboard. Cover is single layer cardboard
- Front protection flaps are attached with durable tape



Cat. No.	Colour	Capacity	L x W x H in	L x W x H cm	Pack qty
11703217	Black	<input type="checkbox"/>	8.2 x 13.4 x 0.4	20.8 x 34.1 x 0.9	72cs
11790465	Blue	20 slides			
12942017	Green	<input type="checkbox"/>	8.2 x 12.9 x 0.4	28.1 x 32.8 x 1	40cs
11713217	Black	30 slides			

Slide Folder Clasps

Keep slide folder flaps securely closed

- Ideal for histology, microbiology and pathology as well as multiple user research environments
- Clasps can be labeled with a permanent marker or affixed label

Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
15813977	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	2.2 x 1.7 x 0.6	5.6 x 4.4 x 1.6	100cs
15823977	Clear <input type="checkbox"/>			


*Assorted colours (Clear, Grey, Blue, Green, Red)



Slide Holders

Hold slides during evaluation, cleaning, and drying

- Made from chemical and stain-resistant polyethylene
- Solid, strong construction with rubber feet for benchtop stability
- Angled slide holder slots allow easy access to individual slides and viewing of slide labels



Cat. No.	Colour	Capacity	L x W x H in	L x W x H cm	Pack qty
11937984	Natural <input type="checkbox"/>	40 Slides	11.7 x 4.1 x 0.9	29.8 x 10.5 x 2.4	1 each
11947984	Natural <input type="checkbox"/>	78 Slides	11.8 x 7 x 0.8	29.9 x 17.8 x 2.1	1 each

Slide Staining Tray Set

Reduce possible cross-contamination and messy work areas

Rugged staining tray is an economical alternative to metal versions

- Made from a polypropylene and polyethylene blend for stability
- Stain, rinse, and dry slides on a single working tray
- Fits up to eight slides and deepwell holds up to 38mL
- Compact size and recessed handles for easy transport from work area to sink - pack of 4 base trays and 1 lid



Cat. No.	Colour	Capacity	L x W x H in	L x W x H cm	Pack qty
11968004	Black <input type="checkbox"/>	8 Slides	12 x 5 x 1.2	30.5 x 12.7 x 3	1pk

Slide Mailers

Safe, secure transport of microscope slides

Cardboard Slide Mailers

Cat. No.	Colour	Capacity	L x W x H in	L x W x H mm	Pack qty
12323107	Natural <input type="checkbox"/>	1-place	3.7 x 1.7 x 0.2	95 x 42 x 5	25pk
15873967	Natural <input checked="" type="checkbox"/>	4-place	4.7 x 5.7 x 0.3	120 x 140 x 6	20pk

- Heavy cardboard mailers for safe transport
- Cover that opens for viewing and thumb grooves for easy slide removal
- Securely holds slides with or without cover glass
- Flat surface area is ideal for labeling or writing necessary information
- Each will hold slides 26 x 76mm (1 x 3 in)



Plastic 2-Place Slide Mailers

Cat. No.	Colour	Capacity	L x W x H in	L x W x H mm	Pack qty
11745088	Assorted* <input checked="" type="checkbox"/>	2-place	3.7 x 2.7 x 0.2	93 x 68 x 5	25pk
11755088	Natural <input type="checkbox"/>	1-place	3.3 x 1.6 x 0.2	83 x 41 x 5	50pk

*Assorted colours (Blue, Purple, Red, Orange, Yellow)

- Ideal for secure shipping and storing of slides
- Polypropylene protectors hold two standard (26 x 76mm and 1 x 3 in) slides securely
- Large, easy-to-secure snaps, provide safe storage for valuable slides
- Variety of colours brings some fun to the lab and enables colour coding




Plastic 5-Place Slide Mailers

Cat. No.	Colour	Capacity	L x W x H in	L x W x H mm	Pack qty
11714316	Natural <input type="checkbox"/>	5-place	3.2 x 0.7 x 1.2	81 x 19 x 31	25pk
11719885	Natural <input type="checkbox"/>	5-place	1.7 x 0.9 x 3.5	43 x 24 x 88	25pk
16310736	Assorted* <input checked="" type="checkbox"/>	5-place	1.7 x 0.9 x 3.5	43 x 24 x 88	25pk

*Assorted colours (Blue, Purple, Red, Orange, Yellow)

- Ideal for secure shipping, storing, or staining slides
- Polypropylene containers hold five standard slides (26 x 76mm)
- Inside wall features slotted channels to keep slides separated and to make removal easier
- Lid closes securely to keep slides safe





Complete system for easy counting

fisherbrand

Colony Counter

Rapid accurate colony counter

Combines a multi-function electronic counter with a Sharpie™ pen for marking to prevent missing or double-counting colonies

- Light, ergonomic body is balanced for easy manipulation
- Select between an audible beep or an LED flash to acknowledge count
- Store up to 32 separate counts and display the total of all counts
- Count up/down verification
- Quick and easy to change between different coloured markers to check different colony counts in the same view
- Lock counter off to allow use of the marker without registering more counts

Cat. No.	Description	Pack qty
15253305	Electronic colony counter package	1



Weighted stand for bench top use



Thumb button allows you to use the ecount™ as a digital tally counter

Package includes:

- Colony counter
- Black Sharpie™ pen for highly visible marks on both plastic and glass
- Long-life button battery
- Hand tally adapter counts without marking
- Weighted stand for convenient bench top use



fisherbrand

Mini-Magnifier

Perfect accessory for the Colony Counter

- Swiveling, moulded lens. 1.75x magnification
- Enlarges objects for easier viewing

Cat. No.	Description	L x W x H	Pack qty
12816355	Mini magnifier	161 x 135 x 207mm	1 each


fisherbrand

Benchtop Counters

Mechanical benchtop counters are perfect for blood cells (label included), bacterial colonies or other repetitive counting

- Each key records up to 999
- Total window keeps track of total number of strokes made on all other keys; bell rings when total reaches 100
- Blood cell key labels

Cat. No.	No. of Keys	L x W x H in	L x W x H cm	Pack qty
11334095	Five key	9.8 x 3 x 2.2	25 x 7.6 x 5.6	1 each
11344095	Eight key	12.5 x 3 x 2.2	32 x 7.6 x 5.6	





Ratcheting Handle
folds out of the way for easy access,
even over tall sharp containers



Supply Drawer

Side Pockets

fisherbrand

Blood Collection Tray

Designed with phlebotomists in mind

- Keeping a clean and organised tray has never been easier. Wash the ABS tray as needed, and replace the insert periodically for a fresh start.
- Multiple compartments keep a variety of supplies neatly at hand, and organised
- Built in side pockets offer additional storage space
- Handle includes soft rubber grip to reduce hand and finger strain
- Each tray includes either a 13 or 16mm tube rack (page 21) and two disposable polystyrene inserts

Droplet Configurations

Cat. No.	Tray Insert	Tube Rack mm	L x W x H in	L x W x H cm	Handle Height in	Handle Height cm	Pack qty
11740634	11760634	13	16.8 x 10.4 x 6.7	42.9 x 26.5 x 17.2	11	28	1 each
11750634	11770634	16					

Droplet Inserts

Cat. No.	Description	L x W in	L x W x H cm	Pack qty
11760634	Tray insert for a 13 tube rack	14.2 x 8.2	36 x 21	4pk
11770634	Tray insert for a 16 tube rack			

Tray Inserts

come in two configurations and include an extra-deep compartment for a sharps container.



Disposable Insert A
for 13mm tube rack



Disposable Insert B
for 16mm tube rack

US Patent No. D544703



Tray Lid

Cover for the Droplet Blood Collection Tray

Protect the contents of the Droplet Blood Collection Tray with custom designed, secure lid

- Notched design provides access to tray handle to maintain portability
- Made of easy to clean polystyrene
- Can be secured with cable ties for extra security

Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
11780634	White <input type="checkbox"/>	16.5 x 9.4 x 8	42 x 24 x 20.5	1 each

Sharps Container

Safe & Portable

Provides safe disposal of syringes, needles and other small sharps

- Capacity 1L with a clear lid to view fill level
- Slide access lid for temporary and final closures for safety and security
- Built-in needle removal port, preventing direct contact with sharps



Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
15290009	Yellow ■	4.5 x 4.5 x 7.5	11.5 x 11.5 x 19.2	1 each



Droplet™ Sample Storage Tray

Sample organisation for the busy lab

Hard wearing sample storage trays that can be used for specimen containers with the insert provided, or for general lab products without the insert

- Lightweight and sturdy to hold up to 24 full specimen containers (up to 58mm in diameter). Specimen containers not included.
- Removable insert to allow outer shell to be used for storage or transportation of laboratory products
- Stackable with other trays to support as many as four high
- Integrated label holder and cards to enable easily changed identification

Cat. No.	Wells	L x W x H in	L x W x H cm	Pack qty
15524149	24	18.8 x 13 x 3.3	48 x 33.1 x 8.4	1 each


Corrugated Sample Storage Tray

Disposable sample organisation and storage

Corrugated storage trays can be used with the included dividers for specimen containers or remove them for general lab product storage or transport

- Made from corrugated polypropylene
- Moisture resistant, significantly reducing the chance of mould

Cat. No.	Colour	Wells	Well Dia mm	L x W x H in	L x W x H cm	Pack qty
15376548	White □	24	58	16 x 10.2 x 3.5	40.6 x 26 x 8.9	10pk



Trays are supplied "flat-packed"
Eliminating the hassle of stacking and reduces the amount of space needed for storage

Lightweight and Sturdy
to hold up to 24 full specimen containers



Removable Inserts
allow outer shell to be used for storage or transport of laboratory products



Carry-Handle Folds
neatly into body space allowing the units to be stacked

Moulded Slots
to hold tubes while working



Silicone Water-tight Seal
and three lid clasps provide secure closure

fisherbrand

Transport Box Specimen or Sample Transport Box

Ideal for transporting biological and clinical samples, as well as instruments and products that you want to keep clean and dry under testing environmental conditions

The transparent box design helps comply with OSHA 29 CFR 1910.1030 standards without having to re-label containers, and enabling the user to visibly check sample integrity.


- Polycarbonate construction and available in three colours for user/owner identification
- Will accommodate 13 and 16mm tubes in a 72-place tube rack (page 21), or similar sized/brand of rack
- Use the separators provided to create your own compartment sizes
- Take out separators to transport larger products and equipment
- **Autoclavable** 120°C for 15 min., 16 psi (1.1 bar)




A silicone water-tight seal and three clasps provide secure closure




Load



Lock



Transport



Store




Cat. No.	Colour		L x W x H in	L x W x H cm	Pack qty
15165564	Clear with Blue Handle	■			
15175564	Red with Red Handle	■	15 x 7.7 x 6.3	38 x 19.6 x 16	1 each
15155564	Yellow with Yellow Handle	■			

Arctic Squares™ Polycarbonate Cryostorage Boxes

Safely store cryogenic vials

Our polycarbonate boxes are designed with unique features to help orient, identify and access vials with ease

- Liquid nitrogen vent and drainage holes
- Magnet on top of forceps enables easy storing on metal surfaces
- Stackable, saves valuable space in the lab
- Store down to temperatures of -196°C with appropriate temperature viable tubes



US Patent No. D529622; EP Registered Design No.: 00551429-0001, 00551429-0002, 00551429-0003


	Cat. No.	Colour	Tube Array	Tubes mL	Wells	Well Dia. mm	Lid	L x W x H in	L x W x H cm	Forceps	Pack qty
a	11938084	Ruby Red	5 x 5	2.0	25	12.5	Slide On	3 x 3 x 2	7.6 x 7.6 x 5.3		8pk
b	11998004	Assorted*	9 x 9	2.0	81	12.5	Hinged	5.2 x 5.2 x 2	13.3 x 13.3 x 5.3	✓	4pk
c	11978004	Arctic Blue	10 x 10	2.0	100	12.5	Hinged	5.2 x 5.2 x 2	13.3 x 13.3 x 5.3	✓	4pk

*Assorted colours (Ruby Red, Amethyst Purple, Arctic Blue, Emerald Green)

Freezer Boxes

Reusable, economical and moisture resistant plastic

- Easy to assemble boxes made from polypropylene material
- Moisture resistant, more durable than cardboard boxes
- Dividers included with the boxes
- Corrugated polypropylene: store down to temperatures of -80°C with appropriate temperature viable tubes
- Thin film: can be used in mechanical freezers & vessels down to -196°C with appropriate temperature viable tubes. Vapour phase only.



Freezer boxes are supplied “flat-packed”
Eliminating the hassle of stacking and reduces the amount of space needed for storage

Corrugated Polypropylene 0.2/0.5/1.5–2.0/5mL/Vial


Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty	
15316528	White	<input type="checkbox"/>	0.2	144	7.9	5.2 x 5.1 x 1.2	13.3 x 13 x 3	10pk
15326528	Blue	<input checked="" type="checkbox"/>	0.2	144	7.9	5.2 x 5.1 x 1.2	13.3 x 13 x 3	10pk
15356528	White	<input type="checkbox"/>	0.5	81	12	5.2 x 5.1 x 1.6	13.3 x 13 x 4.2	10pk
15366528	Blue	<input checked="" type="checkbox"/>	0.5	81	12	5.2 x 5.1 x 1.6	13.3 x 13 x 4.2	10pk
15226996	White	<input type="checkbox"/>	1.5/2.0	81	12	5.2 x 5.1 x 2	13.3 x 13 x 5.1	10pk
15396528	Blue	<input checked="" type="checkbox"/>	1.5/2.0	81	12	5.2 x 5.1 x 2	13.3 x 13 x 5.1	10pk
15326538	White	<input type="checkbox"/>	5	25	22.6	5.2 x 5.1 x 2.9	13.3 x 13 x 7.5	10pk
15336538	Blue	<input checked="" type="checkbox"/>	5	25	22.6	5.2 x 5.1 x 2.9	13.3 x 13 x 7.5	10pk
15413790	White	<input type="checkbox"/>	Vial	100	17	7.4 x 7.8 x 1.6	18.8 x 19.8 x 4.1	10pk

Corrugated Polypropylene 15mL and 50mL

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty	
15366538	White	<input type="checkbox"/>	15	36	19.8	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12	10pk
15376538	Blue	<input checked="" type="checkbox"/>	15	36	19.8	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12	10pk
15306548	White	<input type="checkbox"/>	50	16	31.6	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12	10pk
15316548	Blue	<input checked="" type="checkbox"/>	50	16	31.6	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12	10pk

Thin Polypropylene film 1.5–2mL

Cat. No.	Colour	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	Pack qty	
15610785	White	<input type="checkbox"/>	1.5/2.0	81	12.1	5 x 4.9 x 1.9	12.5 x 12.4 x 4.9	10pk
15620785	Blue	<input checked="" type="checkbox"/>	1.5/2.0	81	12.1	5 x 4.9 x 1.9	12.5 x 12.4 x 4.9	10pk
15630785	Lt. Green	<input checked="" type="checkbox"/>	1.5/2.0	81	12.1	5 x 4.9 x 1.9	12.5 x 12.4 x 4.9	10pk
15640785	Purple	<input checked="" type="checkbox"/>	1.5/2.0	81	12.1	5 x 4.9 x 1.9	12.5 x 12.4 x 4.9	10pk



Corrugated Polypropylene 50mL

Thin Polypropylene film 1.5–2mL

Corrugated Polypropylene 1.5–2.0/5mL

Freezer Microscope Slide Box

Store in ultra-low temperatures with confidence

Protect your slides in low temperature ranges down to -80°C

- Durable polycarbonate and stainless steel clasp/friction fit lid
- Safe storage and transport of 26 x 76mm slides
- Uniform grooved slide slots separate slides preventing slide contact
- Stay organised with a moulded labeling area with "Project" and "Date"
- Low temperature-resistant numbered foam lining and identification sheet
- Moulded feet enable secure stacking for easy transport and storage



Cat. No.	Colour	Capacity	Hinge	L x W x H in	L x W x H cm	Pack qty
15366548	Blue ■	25 Slides	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5	1 each
11730594	Purple ■	100 Slides	Stainless Steel	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4	1 each

Cryogenic Vial Canes

Organise cryogenic samples in dewars

Aluminium canes hold vials for freezer storage or immersion in LN₂ dewars

- Use in liquid and vapour phase; use approved vials only
- Canes bend to fit all major cryogenic vial brands

Cat. No.	Tubes mL	No. of Vials	Vial Size mL	Length in	Length cm	Pack qty
11958084	15 & 50	5	Up to 2.0	1.1	2.9	12pk
11968084	15 & 50	6	Up to 5.0	1.2	3	12pk



PVC Cryogenic Cane Sleeves

For extra security and sample safety

Clear plastic sleeves protect vials in aluminium canes during handling and storage

- Safe for cryogenic freezing
- More durable than cardboard sleeves
- Transparency allows sample identification without removing sleeve

Cat. No.	Diameter in	Length in	Diameter cm	Length cm	Pack qty
11734673	0.6	10.8	1.5	27.5	100pk



Universal Freezer Rack

Innovative horizontal drawer system for upright freezers

- Constructed of corrosion-resistant stainless steel with four drawers to hold 16 boxes
- Drawers feature a unique stop, that when fully extended pivot down allowing for easy box access
- Steel drawer handles provide a user-friendly grip
- Freezer rack also includes a pull-out wire handle for easy removal off the freezer shelf
- Compatible with most standard 5.3cm freezer storage boxes



Unique drawer stop
that pivots down allowing for easy box access


Steel drawer handles
have a space to insert a label, and provide a user-friendly grip

Cat. No.	Config.	Box Height		L x W x H	L x W x H	Pack qty
		in	cm	in	cm	
11786013	4 x 4	2	5.3	5.4 x 22.4 x 9.3	13.9 x 57.1 x 23.8	1 each

Cube Freezer Rack – Upright

For 15 and 50mL tube box storage

- Constructed of corrosion-resistant stainless steel with 6 rack spaces
- Provides safe, convenient storage of your samples
- Handles at each end accommodate freezers with right or left handed doors
- Compatible with 15.2cm, or smaller boxes



Handles for easy transfer
Robust mounting brackets assure handles stay securely in place for easy use

Cat. No.	Config.	Opening L x W x H		L x W x H	L x W x H	Pack qty
		in	cm	in	cm	
12643476	2 x 3	6 x 6.1 x 5.2	15.2 x 15.7 x 13.4	18 x 6.2 x 10.6	45.9 x 15.9 x 27.1	1 each



fisherbrand

Paper Boats

Green solution to sample weighing

Made from FSC™ Certified (FSC™ Mix Credit) paper (Eco-IN) and when disposed of in the standard waste stream they will naturally biodegrade (Eco-OUT)

- A cost-effective solution to reduce plastics in the lab
- Made of grease-proof paper that meets FDA 21 CFR 176.170 (components of paper in contact with aqueous and fatty foods) and with 21 CFR 176.180 (contact with dry foods)
- Control contents like a polystyrene weigh boat because they are easier to pour
- An alternative to regular weighing papers because they are easier to handle

Cat. No.	Colour	Size	L x W x H in	L x W x H mm	Pack qty
15422214	Brown <input checked="" type="checkbox"/>	Small	13.8 x 13.8 x 5.9	35 x 35 x 15	500pk
15408113	White <input type="checkbox"/>				
15433214	Brown <input checked="" type="checkbox"/>	Medium	17.7 x 17.7 x 5.9	45 x 45 x 15	500pk
15442214	White <input type="checkbox"/>				
15893967	Brown <input checked="" type="checkbox"/>	Large	3.2 x 3.2 x 1	95 x 95 x 25	250pk
15853967	White <input type="checkbox"/>				

Functional in the Lab

Available in a small and medium in both brown for light coloured samples and white for dark samples



fisherbrand

Disposable Spatulas


Designed to fit perfectly with standard 0.2, 0.5, 1.5/2.0, 15 and 50mL tubes

- Made of virgin, food-grade polypropylene
- Ideal for dispensing or sample collection of powders or soft solids
- Durable to work with liquid nitrogen frozen samples
- Hollow core minimises heat transfer—good for hot and cold liquids
- Antistatic, sterile and eco-friendly versions available

Always there and always clean!

Autoclavable

Lab tools that scoop, stir and scrape.



Cat. No.	Description	Colour	Suggested Tube mL	L x W mm	Spoon (Flare) mm	Scoop/Pick L x W mm	Taper	Pack qty
a 15356558	Micro Antistatic	Opaque <input type="checkbox"/>	0.2 and 0.5	149 x 4	—	11 (V-Tip)	15°	300pk
b 15326134	Eco-Friendly	Opaque <input type="checkbox"/>	15 and 1.5	140 x 7	25 x 14	30 x 3		300pk
c 15336558	Standard	Blue <input checked="" type="checkbox"/>						300pk
d 15346558	Standard	Opaque <input type="checkbox"/>	15 and 1.5	210 x 7	25 x 14	30 x 3		300pk
15376558	Std. EtO Sterile	Opaque <input type="checkbox"/>						100pk
e 15366558	Macro	Green <input checked="" type="checkbox"/>	50	310 x 10	42 x 23.5	30	35°	150pk

Weighing Boats

Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss

- Ideal for weighing, dispensing, or storing applications
- Select white, black, or blue for best colour contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Antistatic versions available making the boats nonreactive with most substances and suitable for weighing static-affected samples
- Wide flat bottoms ensure stability and prevent tipping
- Rounded corners allow simple flexing for controlled pouring of liquid, powdered or granular samples
- Flat bottoms and sloping sides allow easy stacking for convenient storage
- Can be used with microwave moisture analysers
- Made from FDA/food-grade compliant material



Standard



Hexagonal



Pour-boat



Square



Diamond

Shape	Regular Material Cat. No.	Antistatic Material Cat. No.	Colour	Size	Volume mL	L x W x H in	L x W x H cm	Pack qty
Standard	11593422	12922850	White <input type="checkbox"/>	Small	10	1.8 x 1.8 x .3	4.6 x 4.6 x .8	500pk
	15728187	15768167	Blue <input checked="" type="checkbox"/>					
	12982850	12942860	Black <input checked="" type="checkbox"/>					
	11503432	11680302	White <input type="checkbox"/>	Medium	85	3.4 x 3.4 x .9	8.5 x 8.5 x 2.4	
	15738187	15778167	Blue <input checked="" type="checkbox"/>					
	12992850	12952860	Black <input checked="" type="checkbox"/>					
	11513432	12608513	White <input type="checkbox"/>	Large	280	5.5 x 5.5 x .9	14 x 14 x 2.2	
	15748187	15788167	Blue <input checked="" type="checkbox"/>					
	12902860	12962860	Black <input checked="" type="checkbox"/>					
Hexagonal	—	15778187	White <input type="checkbox"/>	Small	6	1.5 x 1 x .4	3.8 x 2.5 x .9	500pk
	—	15788187		Medium	46	2.5 x 2 x .7	6.4 x 5.1 x 1.9	500pk
	—	15798187		Large	194	4 x 3 x .9	10.2 x 7.6 x 2.2	500pk
	—	15708197		X-Large	364	5 x 3.5 x 1.2	12.7 x 8.9 x 3.1	500pk
Pour-boat	15798177	15798167	White <input type="checkbox"/>	Small	14	1.7 x 2.3 x .5	4.3 x 5.8 x 1.3	500pk
	15708187	15708177		Medium	137	3.3 x 5.2 x 1	8.3 x 13.2 x 2.6	500pk
	15718187	15718177		Large	259	4.3 x 7.2 x 1	10.8 x 18.3 x 2.6	500pk
Square	15758187	—	White <input type="checkbox"/>	Medium	51	3.5 x 3.5 x .3	9 x 9 x .7	500pk
Diamond	12932840	12992840	White <input type="checkbox"/>	Small	5	2.1 x 1.3 x .2	5.3 x 3.3 x .6	500pk
	12952850	12912860	Black <input checked="" type="checkbox"/>					
	15768187	15718197	White <input type="checkbox"/>	Medium	25	2.6 x 1.9 x .6	6.7 x 4.7 x 1.4	
	15738197	15728197	Black <input checked="" type="checkbox"/>					
	11573422	11533432	White <input type="checkbox"/>	Large	30	3.1 x 2.2 x .6	7.8 x 5.6 x 1.4	
	12962850	12922860	Black <input checked="" type="checkbox"/>					
	12387552	11543432	White <input type="checkbox"/>	X-Large	100	4.7 x 3.5 x .7	11.9 x 9 x 1.9	
	12972850	12932860	Black <input checked="" type="checkbox"/>					

Weighing Boats & Dishes

For general-purpose weighing, dispensing, storage, food processing, evaporating, or drying

Aluminium Weighing Boats

- Sturdy yet lightweight boats feature a smooth, flat bottom for sample handling and a flanged, crimped wall for rigidity
- Finger tab handle can be etched for identification prior to weighing



Aluminium Weighing Boats



Aluminium Weighing Dish

Aluminium Weighing Dish


- Designed for moisture determination balances with slightly tapered walls for stacking
- Use as dust covers for open vessels



Cat. No.	Size	Dia. x H in	Dia. x H mm	Pack qty
11998084	Small	1.7 x 0.5	43 x 12	
11908094	Medium	2.2 x 0.6	57 x 14	100pk
11918094	Large	2.8 x 0.6	70 x 16	
11988084	Large	0.4 x 0.3	10 x 8	50pk

Aluminium Weighing Dish

- Ideal for microbalances, organic micro elementary and thermal analysis, and pharmaceutical use.
- Made of Ultra High Vacuum (UHV) Foil - "A" wettable and completely free from oils normally found on other grades of aluminium foil
- Foil meets ASTM B-479, specification 3.1.4 and 10.3.1
- "Cahn style" accessory pans are compatible with most commercially available analysers and microbalances



Cat. No.	Dia. x H in	Dia. x H mm	Pack qty
15803977	1 x 0.1	25 x 3.6	20pk

Dual Purpose Scoops

Scoop + Weigh Boat

Unique economical solution that combines the function of a scoop and a weigh boat into one, streamlining the process. Just scoop to retrieve, place on scale to weigh!

- Made from anti-static polystyrene that is waterproof and reduces the chance of materials sticking to the scoop
- Ships and stores flat, easy to assemble for use
- Locking tab feature enables a scoop shape to help retain materials
- Coated to stand up to powders, gels, and semi-solids
- Electronic beam sterilisation procedures compliant with ISO11137-3:2017
- Certificate of Sterility available upon request



Cat. No.	Colour	Size	L x W in	cm	Sterile	Pack qty
15833977	White	<input type="checkbox"/> Large	5 x 3.5	12.7 x 9	No	100pk
15843977	White	<input type="checkbox"/> Small	5 x 2.5	12.7 x 6.3	No	100pk
16286860	Black	<input checked="" type="checkbox"/> Large	5 x 3.5	12.7 x 9	Yes	100pk
16296860	Black	<input checked="" type="checkbox"/> Small	5 x 2.5	12.7 x 6.3	Yes	100pk



fisherbrand

Heat Sealer with Cutter

Quickly and efficiently seal polyethylene or polypropylene bags, sheets or flat tubing from 0.01 to 0.8mm thick


Sliding blade removes excess material 19mm from the seal

- Can seal a bag up to 12 in (30.5cm) long
- Heating element measures 31cm long and is only hot when sealing
- Includes one replacement part kit
- Additional replacement part kits can be ordered separately

Cat. No.	Power Supply	D x W x H cm	Pack qty
15350360	220/230 VAC, U.K	8.6 x 44.6 x 16.6	1 each
15360360	220/230 VAC, Cont. Europe/Israel		1 each
15350407	Heat Sealer Replacement Parts Kit		1 each

Replacement Parts Kit contains:

- Blades (2), Heating element (2), PTFE cloths (2)



fisherbrand

Safe Slicer

Laboratory Safety Razor Blade

Say no to dangerous sharps

- At just 0.27mm the US Patent Pending product is as thin as a single edge industrial razor blade, but much safer
- Cleanly slice directly on light box without scratching the surface
- Adjust the size required by simply cutting with scissors
- Made from food grade PET

Cat. No.	L x W x H		Pack qty
	in	mm	
15316558	0.8 x 1.5 x 0.01	19 x 38 x 0.3	100pk
15326558	1 x 2 x 0.01	25 x 50 x 0.3	

fisherbrand

Paper Disposable Funnel

A clean transfer every time

- Supplied flat-pack to minimise storage
- With a minimal amount of finger pressure the funnel takes its working shape
- Maximum flexibility – cut the stem to any size or shape to meet with your needs
- Fits a wide variety of openings and can be used with virtually any liquid* or powder
- Disposable

*Not recommended for hot liquids



Cat. No.	Capacity mL	Top Dia. in	Top Dia. cm	Bottom Dia. in	Bottom Dia. cm	Height in	Height cm	Pack qty
15346568	186	2.4	6	0.3	0.6	5	12.8	25pk

fisherbrand

Lead Rings

Stabilise vessels and labware or immerse in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customised fit

Cat. No.	Shape	Colour	Flask Sizes mL	ID cm	OD cm	Weight g	Pack qty
11726725	C	Blue ■	125–500	4.2	7.5	215	1 each
10195573	O			4.8	7.7	210	
11736725	C	Blue ■	250–1000	5.1	9	465	1 each
11797884	O			5.1	9	490	
11746725	C	Blue ■	500–2000	5.4	9.8	595	1 each
11707894	O			5.7	9.8	610	
11756725	C	Blue ■	1000–4000	6.6	11.3	890	1 each
11717894	O			7	11.5	935	

*Weights may vary by +/- 3%



fisherbrand

Cell Spreaders

Delta™

- Made from ABS
- Flexible rippled handle and domed spreading bar to help reduce agar tears
- Gamma sterilised
- Disposable



Domed Headpiece
promotes smooth spreading

Stainless Steel

- Sterilisable
- Reusable
- Flame and **Autoclavable**



Cat. No.	Shape	L x W x H in	L x W x H cm	Packaging	Sterile	Pack qty
11364035	Delta™	6.8 x 1.5 x .2	17.2 x 3.7 x .4	1 each/bag	✓	500pk
13569240	Stainless Steel	8 x 1.3 x .2	20.3 x 3.3 x .4	—		1 each

Inoculating Picks

The 'smart' way to isolate individual colonies

Available in three lengths and two needle sizes to suit a wide range of application needs

- Naturally biodegradable, made from 100% U.S.A. grown white birch wood
- Manufactured in an FDA-registered factory and complies with FDA standards for food contact
- Available in semi or needle point styles for delicate operations



Autoclavable

Cat. No.	Size	Point	Tip	Suggested Tube mL	Length in	Length cm	Pack qty
15386558	Mini	Semi Point		Fits 1.5 tube (closed)	1.3	3.3	2000pk
15396558	Mini	Needle Point					
15306568	Standard	Semi Point	Single Tip		2.5	6.4	1600pk
15316568	Standard	Needle Point					
15326568	Large	Semi Point	Single Tip		6	15.2	500pk
15336568	Large	Needle Point					



Mini Sized SmartPicks™ fit neatly into a 1.5mL microfuge tube allowing the lid to be closed for full sample transfer/recovery

Drop it in & Close it up!



Workstation Dispensing Bins

Crystal-clear bins let you see your small items inventory

- See inventory levels of pipettor tips, tubes, fittings, finger cots, and other bulk small items at a glance
- Flip-top covers on generous openings make filling and dispensing easy and keep products dust-free
- Four sizes available to cater to a wide variety of storage needs
- Includes non-skid rubber feet

Cat. No.	Colour	Capacity	L x W x H cm	Pack qty
11764834	Clear <input type="checkbox"/>	1	18 x 17.5 x 23.5	1 each
11774834		2	13.9 x 16 x 36.8	
11907994		3	21 x 15.5 x 21.2	
11784834		4	27.8 x 15.8 x 21.7	

Petri Dish Dispenser

- Sturdy acrylic holds up to 30 x 100mm dishes in two easy-access partitions
- Bench top or wall mountable (screws provided)



Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
11384145	Clear <input type="checkbox"/>	8.7 x 4.2 x 9.3	22.2 x 10.7 x 23.6	1 each

Lab Wipe Holders

Holds lab wipes such as Kimwipes™ and other major brands

Clear acrylic holders fit boxes of small or large lab wipes.

- Mount on a shelf or wall
- The smaller version attaches with double-stick tape and the larger version with screws
- Mounting materials included

Cat. No.	Size	L x W x H in	L x W x H cm	Pack qty
12377619	Small	4.6 x 4.8 x 3.3	1.8 x 12.2 x 8.4	1 each
12387619	Large	12.5 x 5 x 3.8	32 x 12.7 x 9.7	



Clear Acrylic Glove Box Holders

Keep gloves organised & ready for immediate use

Free up bench-top space by mounting holders on wall for easy dispensing of gloves

- Universal fit for most major brands of glove boxes
- Mounting screws included
- Includes a polyester foam insert to securely hold a range of glove box sizes—snug fit ensures easy removal of gloves from the box

Cat. No.	Capacity	L x W x H in	L x W x H cm	Pack qty
11897102	1	10.3 x 5.9 x 4.5	26.3 x 14.6 x 11.4	1 each
11807112	2	11.4 x 10.3 x 4.5	29.2 x 26 x 11.4	1 each
11817112	3	11.6 x 15.2 x 4.3	29.4 x 38.8 x 11.6	1 each



Epoxy Coated Glove Box Holders

Keep gloves organised & ready for immediate use

- Epoxy-coated steel wire racks are an economical solution
- Universal fit for most major brands of glove boxes
- Mounting screws included

Cat. No.	Capacity	L x W x H in	L x W x H cm	Pack qty
11885843	1	5.5 x 8.2 x 4.3	14 x 20.9 x 10.8	1 each
11895843	3	9.8 x 17.9 x 4.3	25 x 45.4 x 10.8	1 each

Modular Glove Box Holders

Available infused with Microban™ for anti-microbial protection

The anti-microbial version (blue) is ideal for use in clean rooms, hospitals, and microbiological work areas. The Microban additive prevents growth of moulds and bacteria

- High-impact, durable resin is resistant to harsh cleaning solutions
- Modular design allows single boxes or horizontal stacking for multiple box holders
- Universal fit for most major brands of glove boxes
- Mounting screws included

Cat. No.	Description	Colour	L x W x H in	L x W x H cm	Pack qty
15380440	Standard	White	10 x 6.1 x 3.8	25.4 x 15.6 x 9.7	3pk
11710644	Anti-microbial	Blue			3pk



US Patent No. 6669157

Microban™ antimicrobial protection has been built-in to our anti-microbial glove boxes during manufacturing. When microbes, such as bacteria, mould and mildew that can cause stains, odours and product deterioration come in contact with the glove box surface, Microban protection penetrates the cell wall of the microbe and disrupts key cell functions so that the microbe cannot function, grow, or reproduce. Microban protection is an intrinsic part inside and at the surface and will not wash off or wear away for the useful lifetime of the glove box.



Adaptable
Store different size boxes all in one holder

Flexible Glove Box Holder

Glove box holder ideal for gloves, face masks, shoe covers...

The adjustable holding points provide true flexibility to hold up to four different boxed protective products and other lab essentials (size dependent)

- Ideal for holding different boxed sizes of gloves, mop caps, ear plugs, face masks, shoe covers
- Adjustable universal fit
- Easy to clean, with a low profile
- Wall or door-mountable
- Mounting hardware included

Cat. No.	Colour	L x W x D in	L x W x D cm	Pack qty
15279568	White	22.3 x 5.3 x 1.8	56.5 x 13.3 x 4.6	1 each



Workstation Storage Bin

The essential combination accessory holder

Organise your most commonly used laboratory disposables in one central location

- Acrylic unit has room to store gloves, Kimwipes™, and Parafilm™
- Use bins to store tubes, tips, etc.

Cat. No.	Colour	L x W x H cm	Pack qty
11927994	Clear	26.5 x 13.7 x 44	1 each

Parafilm™ is a registered trademark of Bemis Inc.

Drawer Organiser

Keeps supplies organised and easily accessible

Polystyrene organiser perfect for storing and organising small laboratory supplies

- Protect small and fragile items
- Fit in most standard drawers or store on benches and shelves



Cat. No.	Colour	L x W x H in	L x W x H cm	Pack qty
11747444	White	17.9 x 12.1 x 2.5	45.5 x 30.7 x 6.3	1 each

Compartments	L x W in	L x W cm
2 Small Compartments	17 x 3.1	43 x 8
1 Large Compartment	17 x 4.3	43 x 11

Storage Racks for 80-Well Racks

Strong, clear acrylic construction

These racks help you keep large numbers of microtubes organised on your benchtop

- Holds up to eight 80-well microtube racks filled with either 1.5 or 2.0mL microtubes
- Plenty of room to slide the racks in and out



Cat. No.	Description	Colour	L x W x H in	L x W x H cm	Pack qty
11708214	Vertical	Clear	6.2 x 9.3 x 9.2	15.7 x 23.6 x 23.4	1 each
11710055	Horizontal	Clear	12.2 x 9.3 x 4.8	31 x 23.6 x 12.2	

Benchtop Dispensing Bins

Work great with gloved hands, bright colours

Ergonomically designed with large, angled openings and oversized lip

- Lid remains open when tilted back
- Nonskid rubber feet and pre-drilled holes for wall mounting
- Green bin includes a removable divider
- Mounting screws included



Cat. No.	Colour	Compartments	L x W x H in	L x W x H cm	Pack qty
11304065	Neon Green	2	7.1 x 6.1 x 6.5	18 x 15.5 x 16.5	1 each
11314065	Neon Blue	1	6.1 x 6.1 x 6.7	15.5 x 15.5 x 17	

Earplug Dispenser

Translucent acrylic dispensers provide easy access to safety accessories

- Holds approximately 200 pairs of foam earplugs
- Set on counter or wall mount (screws included)

Cat. No.	Colour	L x W x D in	L x W x D cm	Pack qty
11507473	Green	8.9 x 6.4 x 6.1	22.7 x 16.2 x 15.6	1 each



Soft Cover Dispenser

Store and dispense with ease

Translucent acrylic dispensers provide easy access to safety accessories

- Hinged lid protects contents from dust and debris
- Set on counter or wall mount (screws included)
- Great for holding other personal items, too!

Cat. No.	Colour	L x W x D in	L x W x D cm	Pack qty
15599801	Natural	8.9 x 6.3 x 6	22.7 x 16.2 x 15.6	1 each



3-in-1 Dispenser

Three separate lift-up lids make loading and dispensing supplies easy

Three-compartment dispenser can be used to organise safety glasses, earplugs, soft covers (including bouffant caps, beard covers, shoe covers) or any small items that need to be kept within reach

The unit includes screws for wall mounting or can be set on a counter

- Middle compartment has reach-in hole for dispensing
- Durable acrylic is see-through so product levels are always visible

Cat. No.	Colour	L x W x D in	L x W x D cm	Pack qty
11927984	Grey	16 x 16 x 8	40.6 x 40.6 x 20.3	1 each



Safety Glasses Dispensers

Ergonomically angled for easy dispensing

Translucent acrylic dispensers provide easy access to safety accessories

- Hinged lid protects glasses from dust and debris
- Set on counter or wall mount (screws included)
- Large dispenser holds approx. 20 glasses, small approx. 10 glasses
- Great for holding other personal items

Cat. No.	Colour	L x W x D in	L x W x D cm	Pack qty
11803470	Blue	8.9 x 6.4 x 6.1	22.7 x 16.2 x 15.6	1 each
11893460	Blue	8 x 8 x 16	20.3 x 20.3 x 40.6	1 each

fisherbrand

Handy Container

Sturdy storage tub with lid and divider

Our handy stackable storage container is designed to help organise a laboratory

- The removable dividers can create as big or as small a compartment as needed. Perfect for keeping tubes, tips, pipettors, gloves, and other small items organised and out of the way.
- The lid snaps in place to secure the contents and for safe transport.



Constructed of high-impact polypropylene

Cat. No.	L x W x H in	L x W x H cm	Pack qty
11938014	13 x 7.9 x 4.5	33 x 20 x 11.5	5pk



fisherbrand

MagLab

Magnetic lab storage series

Easily organise common lab products off your bench top with the MagLab storage system

Vibrant & Functional

- Made of robust ABS plastic in bright translucent colours
- Secure direct to metal with built-in powerful magnets or non-metal surfaces with included self-adhesive mounting plate
- Available individually or as a pack of four

MagRack holds four 50mL tubes and nine 15mL tubes

MagWipe holds a standard (110 x 210mm) Kimwipes™ box


MagPette holds two pipettors (fits most popular styles)

MagBox holds miscellaneous accessories

Cat. No.	Description	Colour	L x W x H in	L x W x H cm	Pack qty
11942581	MagLab Series	Assorted ■ ■ ■ ■ *	—	—	4pk
11982581	MagRack	Teal ■	5.1 x 4.4 x 2.5	13.1 x 11.2 x 6.4	1 each
11992581	MagWipe	Green ■	5.1 x 3.5 x 3.8	13 x 9.1 x 9.8	
11932581	MagPette	Purple ■	5.1 x 4.4 x 2.9	13.1 x 11.2 x 7.6	
11922591	MagBox	Orange ■	5.1 x 3 x 4	13 x 7.8 x 10.2	

Kimwipes™ is a registered trademark of Kimberly-Clark Corp.

*Assorted colours: Teal, Green, Purple and Orange



Box Top Dispenser

Easy application with no mess

Easy application, simply place the dispenser on top of the Parafilm™ M box and pull the sealing film through, there is no need to remove rolls from the packaging

- Accommodates two 50mm, or one 100mm Parafilm™ M rolls
- Built in safety blade slides and cuts with no danger to fingertips and eliminates the need for scissors
- Supplied with two safety blades, additional replacement blades are available for purchase
- Made from ABS plastic making the unit easy to clean



Item No.		L x W x H in	L x W x H cm	Pack qty
16390412	Natural <input type="checkbox"/>	4.9 x 6.8 x 3.9	12.3 x 17.2 x 9.9	1 each
16390402	Replacement blades for Box Top Dispenser - Black <input checked="" type="checkbox"/>			3 ea/pk

Parafilm™ is a registered trademark of Bemis Company Inc.



ABS dispenser for Parafilm™ Sealing Film

Unique design for maximum safety

Store, dispense, and cut Parafilm™, tape, labels, and Tough-Tags™

- ABS plastic dispenser designed with a safety razor that slides and cuts with no danger to fingers, supplied with two easily exchanged cutting blades
- Dispenser holds one or two 50mm rolls, or one 100mm roll

Cat. No.		L x W x H in	L x W x H cm	Pack qty
11350040	Blue <input checked="" type="checkbox"/>			
11865993	Green <input checked="" type="checkbox"/>	4.7 x 6.1 x 6.7	12 x 15.6 x 17.1	1 each
11875993	Red <input checked="" type="checkbox"/>			
11928014	Replacement blades for ABS dispenser			5pk

Parafilm™ is a registered trademark of Bemis Company, Inc.

Acrylic Dispenser for Parafilm™ Sealing Film

Heavy-duty construction for increased durability

Store, dispense, and cut Parafilm™, tape, labels, and Tough-Tags™

- Strong, crystal-clear acrylic dispenser makes your film supply visible and features a serrated cutting blade
- Dispenser holds one or two 50mm (2-inch) rolls, or one 100mm (4-inch) roll

Cat. No.	L x W x H in	L x W x H cm	Pack qty
11937994	4.7 x 6.6 x 6.9	11.9 x 16.8 x 17.5	1 each

Parafilm™ is a registered trademark of Bemis Company, Inc.



Label Station

Organise your labels for easy accessibility

Portable label station holds up to six standard size boxes of Tough-Tags™, Cryo-Babies™, and Tough-Spots™

Just drop the box in the label station and feed the paper through the slot. Labels shown not included, order separately

- Acrylic station includes a writing surface, a cutting edge for easy dispensing, and non-skid rubber feet
- Counter-sunk holes provided for wall mounting
- Translucent blue



Cat. No.	Colour	L x W x H in	L x W x H mm	Pack qty
11947994	Blue	11 x 4.7 x 3.5	28 x 12 x 8.9	1 each

Unique safety lid



Mini Bin Foil Dispenser

Safely cut foil to any size

The mini bin foil dispenser is a great alternative to wasteful giant-sized foil rolls and expensive pre-cut foil squares

- Compact to fit in limited work space areas and easy to transport
- Safety guard feature to prevent injuries
- Foam strip prevents foil retracting back into unit
- For use with 102mm or 153mm length foil rolls (sold separately)
- Made from ABS


Cat. No.	L x W x H in	L x W x H cm	Pack qty
11957994	7.5 x 3.9 x 4.4	19 x 10 x 11.2	1 each

Foil Rolls

Designed for easy use and minimal waste

- Pieces of foil are ideal for covering small surfaces such as beakers, flasks, tubes and Petri dishes
- Ideal for autoclaving, freezing, incubating, protecting labware, lining, sealing, and general storage
- Two sizes available. designed to fit the Mini Bin Foil Dispenser

Cat. No.	Roll Width cm	Roll Length m	Pack qty
11967994	10.2	152.4	1 each
11977994	15.3	152.4	





Find out more at eu.fishersci.com/go/fisherbrand

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303 Sweden: +46 31 352 32 00
Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500 International: +44 (0)1509 555 939

© 2020 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**
part of Thermo Fisher Scientific