



Disposable Safety Apparel

Disposable Polypropylene Lab Coats

Latex-free spunbond polypropylene construction lab wear

- Elastic cuffs hold sleeves in place
- Five snaps for a secure closure
- Standard collar
- No pockets
- Fabric weight: 50g/m²

Cat. No	Colour	Size	Pack qty
15431922	White	Small	25
15441922	White	Medium	25
15451922	White	Large	25
15461922	White	X-Large	25
15471922	White	XX-Large	25
15481922	White	3X-Large	25



Disposable Polypropylene Lab Coats with Pockets

Single-use lab coats provide economical, lightweight and breathable protection


- Latex-free spunbond polypropylene construction
- Knit cuffs hold sleeves in place
- Three pockets and a five snap closure
- Standard collar
- Fabric weight: 50g/m²

Cat. No	Colour	Size	Pack qty
15491922	White	Small	30
15401932	White	Medium	30
15411932	White	Large	30
15421932	White	X-Large	30
15431932	White	XX-Large	30
15441932	White	3X-Large	30



Lab coat sizing guide

Size	Body Length	Chest Circumference	Sleeve Length
Small	89cm (35in)	119cm (47in)	52cm (20.5in)
Medium	97cm (38in)	125cm (49in)	52cm (20.5in)
Large	99cm (39in)	130cm (51in)	52cm (20.5in)
X-Large	104cm (41in)	135cm (53in)	52cm (20.5in)
XX-Large	109cm (43in)	140cm (55in)	53cm (21in)
3X-Large	114cm (45in)	145cm (57in)	53cm (21in)



A: Body Length
B: Chest Circumference
C: Sleeve Length

Material type quick reference

Material	recommended when a customer needs:	Breathability	Protection Rating
Spunbond Polypropylene	<ul style="list-style-type: none"> • Economical, lightweight, breathable protection • Lower cost material 	Highest level of breathability	Basic protection



Disposable Polypropylene Facemasks

Provides advanced protection

- Disposable, single use
- Choice of ear loops or tie closure
- Nonwoven polypropylene outer layer
- Thermally bonded nonwoven polypropylene inner layer
- Latex-free
- Choice of blue or white
- Universal size/fit

	Cat. No	Description	Colour	Thickness	Pack qty
1	15451932	Disposable face mask with ear loops	White	2-ply	500
	15461932	Disposable face mask with ear loops	Blue	3-ply	500
	15471932	Disposable face mask with ear loops	White	3-ply	500
2	15481932	Disposable face mask with tie closure	Blue	3-ply	500



Disposable Polypropylene Bouffant Caps

Single-use cap provides basic protection

- Disposable
- 14gsm polypropylene
- Diameter: 53.3cm or 60.9cm
- Latex-free
- Choice of blue or white

	Cat. No	Colour	Size, cm	Pack qty
1	15491932	Blue	53.3	1,000
2	15401942	White	53.3	1,000
	15411942	Blue	60.9	1,000
	15421942	White	60.9	1,000



Disposable Pleated Bouffant Caps

Bouffant's elastic headband is encapsulated in soft poly thread for great comfort and basic protection

- Disposable
- Pleated
- Polypropylene
- Choice of white or blue

	Cat. No	Colour	Size	Pack qty
1	15431942	White	Small	1,000
	15441942	Blue	Small	1,000
	15451942	White	Large	1,000
2	15461942	Blue	Large	1,000
	15471942	White	X-Large	1,000





Disposable Shoe Covers

Spunbond polypropylene shoe covers protect against dirt, grime and certain dry particulates

- Single-use
- Blue
- Polypropylene
- Latex-free
- With or without anti-skid sole

Cat. No	Description	Colour	Size	Pack qty
15481942	Shoe cover without anti-skid	Blue	Universal	300
15491942	Shoe cover without anti-skid	Blue	X-Large	300
15401952	Shoe cover with anti-skid	Blue	Universal	300
15411952	Shoe cover with anti-skid	Blue	X-Large	300



Disposable High-Tack Shoe Covers

High-Tack, disposable shoe covers provide superior protection

- Single-use
- Polyurethane
- Blue
- Latex-free
- High-Tack sole

Cat. No	Description	Colour	Size	Pack qty
15421952	Shoe cover, High-Tack	Blue	Universal	300
15431952	Shoe cover, High-Tack	Blue	X-Large	300



Disposable Anti-Static Shoe Covers

Conductive carbon strip helps eliminate static buildup

- Single-use
- Elastic ankle
- 42gms polypropylene
- With carbon dissipative strip sewn in
- Latex-free
- Blue

Cat. No	Description	Colour	Size	Pack qty
15441952	Shoe cover, with conductive carbon strip	Blue	Universal	300





Disposable Polyethylene Shoe Covers

Single use, disposable covers provide convenient protection

- Polyethylene
- Blue
- Latex-free

Cat. No	Description	Colour	Size	Pack qty
15451952	Shoe cover, polyethylene	Blue	Universal	300
15461952	Shoe cover, polyethylene	Blue	X-Large	300



Disposable Beard Cover

Polypropylene facial hair cover for basic protection

- Disposable
- 10gsm spunbonded polypropylene
- Single-use
- Latex-free
- White

Cat. No	Colour	Size	Pack qty
15471952	White	Universal	500



Disposable Polyethylene Sleeve Covers

Disposable general purpose sleeves

- Polyethylene
- Length: 45.7cm
- Latex-free
- Blue
- Thickness: 6 mil
- Elastic cuffs

Cat. No	Colour	Length, cm	Size	Pack qty
15481952	Blue	45.7	Universal	200

