

# WHEATON® COMPLETEPAK

READY TO USE SINGLE SOURCE PRIMARY PACKAGING SOLUTION

The WHEATON COMPLETEPAK product line is a range of tubular glass serum vials, rubber stoppers and seals supplied ready-to-use in one complete kit, certified to European Pharmacopoeia and US Pharmacopoeia, and available off the shelf. The components are also available to be purchased separately.

18 months  
shelf life  
from date of  
processing



Enhancing the Efficiency of Drug Discovery and Drug Compounding

## Benefits

### Easy

- Supplied ready-to-use
- One source for total primary packaging solutions

### Efficient

- Reduces lead times
- Saves processing time
- Allows focus on core competencies

### Cost and space saving

- Eliminates capital investments in equipment required for cleaning and preparing product
- Frees up valuable clean room space

## Target Applications:

- Pharma pilot batch fills
- Ophthalmic packaging
- Small scale production of Biotech Reagents including Controls, Standards & Antibodies
- Freeze dry/lyophilisation drug products
- Compound Pharmaceuticals
- Drug discovery & research
- Validation batch runs

Cat. No	Description	Pack qty	List Price, €
WHEATON COMPLETEPAK Sterile Kit			
<b>16297460</b>	2mL serum vials 13mm ultra pure serum stopper and blue flip seal	210	<b>1303,00</b>
<b>16227300</b>	2mL serum vials 13mm OmniFlex serum stopper and blue flip seal	210	<b>1427,00</b>
<b>16237300</b>	2mL serum vials 13mm OmniFlex Igloo stopper and blue flip seal	210	<b>1466,00</b>
<b>16247300</b>	2mL serum vials 13mm ultra pure serum stopper and red flip seal	210	<b>1303,00</b>
<b>16257300</b>	2mL serum vials 13mm OmniFlex serum stopper and red flip seal	210	<b>1427,00</b>
<b>16267300</b>	2mL serum vials 13mm OmniFlex Igloo stopper and red flip seal	210	<b>1466,00</b>
<b>16277300</b>	5mL serum vials 20mm ultra pure serum stopper and blue flip seal	200	<b>1368,00</b>
<b>16287300</b>	5mL serum vials 20mm OmniFlex serum stopper and blue flip seal	200	<b>1492,00</b>
<b>16297300</b>	5mL serum vials 20mm OmniFlex Igloo stopper and blue flip seal	200	<b>1531,00</b>
<b>16207310</b>	5mL serum vials 20mm ultra pure serum stopper and red flip seal	200	<b>1368,00</b>
<b>16217310</b>	5mL serum vials 20mm OmniFlex serum stopper and red flip seal	200	<b>1492,00</b>
<b>16227310</b>	5mL serum vials 20mm OmniFlex Igloo stopper and red flip seal	200	<b>1531,00</b>
<b>16237310</b>	10mL serum vials 20mm ultra pure serum stopper and blue flip seal	158	<b>1433,00</b>
<b>16247310</b>	10mL serum vials 20mm OmniFlex serum stopper and blue flip seal	158	<b>1557,00</b>
<b>16257310</b>	10mL serum vials 20mm OmniFlex Igloo stopper and blue flip seal	158	<b>1596,00</b>
<b>16267310</b>	10mL serum vials 20mm ultra pure serum stopper and red flip seal	158	<b>1433,00</b>
<b>16277310</b>	10mL serum vials 20mm OmniFlex serum stopper and red flip seal	158	<b>1557,00</b>
<b>16287310</b>	10mL serum vials 20mm OmniFlex Igloo stopper and red flip seal	158	<b>1596,00</b>
WHEATON COMPLETEPAK Sterile Components			
<b>16297280</b>	2mL Type 1 clear serum vials	210	<b>513,00</b>
<b>16207290</b>	5mL Type 1 clear serum vials	200	<b>553,00</b>
<b>16217290</b>	10mL Type 1 clear serum vials	158	<b>714,00</b>
<b>16227290</b>	13mm OmniFlex 3G serum stoppers	230	<b>540,00</b>
<b>16237290</b>	13mm ultra pure serum stoppers	230	<b>429,00</b>
<b>16247290</b>	13mm OmniFlex 3G Igloo Lyo stoppers	230	<b>577,00</b>
<b>16257290</b>	20mm OmniFlex 3G serum stoppers	220	<b>577,00</b>
<b>16267290</b>	20mm ultra pure serum stoppers	220	<b>365,00</b>
<b>16277290</b>	20mm OmniFlex 3G Igloo Lyo stoppers	220	<b>625,00</b>
<b>16287290</b>	13mm flip off red seal	230	<b>397,00</b>
<b>16297290</b>	13mm flip off blue seal	230	<b>397,00</b>
<b>16207300</b>	20mm flip off red seal	220	<b>421,00</b>
<b>16217300</b>	20mm flip off blue seal	220	<b>421,00</b>

Find it at [fishersci.be](http://fishersci.be)

# WHEATON® E-Z Ex-Traction Vials

20% off

The WHEATON biopharmaceutical E-Z Ex-Traction Vials combine the best attributes of high recovery vials with a unique patented bottom design, providing optimum benefits to users of automated or manual handling and storage systems in compound or biological sample applications. Precision manufactured from low extractable borosilicate glass, E-Z Ex-Traction Vials meet ASTM Type 1 Class A and USP Type 1 standards.

- Quality engineered for automated sample handling systems
- Made from low extractable borosilicate glass
- Conical well ensures maximum recovery
- Excellent for chemical compound or biological sample storage



Cat. No	Description	Pack qty	Promo Price, €
10426453	Type I clear borosilicate glass 2mL E-Z Ex-Traction™ vial, without 13-425 screw cap	250	497,79
10373843	Type I clear borosilicate glass 4mL E-Z Ex-Traction™ vial, without 13-425 screw cap	250	508,66
10598013	Type I clear borosilicate glass 8mL E-Z Ex-Traction™ vial without 15-425 screw cap	250	527,13
10292743	Type 1 clear borosilicate glass 10mL E-Z Ex-Traction vial, without 22-400 cap	100	208,76
10661844	Type I clear borosilicate glass 20mL E-Z Ex-Traction™ vial, without 24-400 screw cap	100	242,58
10757443	Black phenolic screw cap with PTFE faced (14B) styrene-butadiene rubber liner, size 13-425	250	155,76
10008752	White polypropylene screw cap with bonded PTFE faced silicone liner, size 13-425	250	155,76
10393363	Black phenolic screw cap with PTFE faced (14B) styrene-butadiene rubber liner, size 15-425	250	171,22
10190343	White polypropylene screw cap with bonded PTFE faced silicone liner, size 15-425	250	171,22
10303563	Black phenolic screw cap with PTFE faced (14B) styrene-butadiene rubber liner, size 22-400	100	133,19
10506103	White polypropylene screw cap with bonded PTFE faced silicone liner, size 22-400	100	133,19
10742973	Black phenolic screw cap with PTFE faced (14B) styrene-butadiene rubber liner, size 24-400	100	144,50
10028082	White polypropylene screw cap with bonded PTFE faced silicone liner, size 24-400	100	135,10

## WHEATON® THE LABWARE YOU TRUST

### PRE-BAR CODED AND READY-TO-USE

It is critical that bar codes are read correctly every time, remain adhered to the container and are durable enough to withstand severe environments. We can customise a bar code solution to meet any specific application utilising a variety of symbologies and materials to provide unique identification solutions for a range of applications.

#### Applications

- Compound storage
- Biobanking
- Diagnostic kitting
- Clinical trials
- Petrochemistry
- Agrochemicals
- Flavours and Fragrances

#### Benefits

- Increased workflows
- Reduce errors
- Save time and money

Allowing you and your team to focus on your core competencies

Contact us to find out how bar coding vials and laboratory glassware can help with traceability and tracking through your workflow process



### Prices valid until 31st December 2019

Enter Promotion Code 13705\* to access the promotional prices

\*Discount will apply to qualifying orders received by Fisher Scientific between 1st June 2019 to 31st December 2019. Promotion code must be presented at the time of order. Discount applies to the list price in effect at the time the order is received. Promotion cannot be combined with other discounts or promotions. Offer void where prohibited, licensed, or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply.

#### Contact us today:

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20  
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478  
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00  
 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303  
 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

