Thermo Scientific Flammable Storage Lab Refrigerators, Freezers and Combination Units

Designed to be used for the storage of volatiles, compounds and solvents where a spark-free internal cabinet environment is required.



5.6 cu. ft. lab flammable storage refrigerator



50 cu. ft. flammable storage refrigerator

Thermo Scientific flammable storage laboratory refrigerators, freezers are available from undercounter to full height freestanding.

Available in North America and select Latin America countries, consult your sales representative.

- Sparkfree interior to reduce the risk of internal explosion
- · All models are manual defrost

- Tough white exterior (except model no. 3575-10)
- CFC- and HCFC-free refrigeration system and insulation
- Designed to meet NFPA 45 and 99 along with OSHA article 1910.307

Warranty: 13 months parts and labor **Certifications:** CSA and or UL





10 cu. ft. flammable storage combination refrigerator and freezer





Thermo Scientific Flammable Storage Lab Refrigerators (Temperature Range 1° to 12°C)

Cat. No.	Capacity cu. ft.	Defrost	Voltage	Shelves (interior)	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
22-299-776	5.6	manual	120V / 60Hz	4 adjustable	yes	16 x 20 x 27.5**	28.75 x 23.5 x 33.5	110
11-675-953*	21	manual	120V / 60Hz	4 adjustable	yes	16.5 x 27 x 55	30.5 x 32 x 74.5	230
11-675-955	21	manual	120V / 60Hz	4 fixed	yes	16.5 x 27 x 55	31.75 x 32 x 74.5	230
NC0163773	50	manual	115V / 60Hz	8 adjustable	yes	26.5 x 45.75 x 59	34.25 x 52 x 77.25	400

^{*} Corrosion-resistant model, features stainless steel interior door liner and stainless steel shelves

Thermo Scientific Flammable Storage Lab Freezers (Temperature Range -12° to -20°C)

Cat. No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock	Interior Dimensions D x W x H (inches)	Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
3556-4	5	manual	115V / 60Hz	2 fixed, 3 door	no	15 x 18 x 28.5	31.4 x 21 x 33.5	84
11-675-957	20.6	manual	115V / 60Hz	4 fixed	yes	20 x 27 x 57	36 x 32 x 73	200

Thermo Scientific Flammable Storage Combination Lab Refrigerator/Freezers (Temperature Range: 1° to 12°C [ref), -12° to -20°C [frz])

Cat. No.	Capacity cu. ft.	Defrost	Voltage	Shelves	Key Lock		Exterior Dimensions D x W x H (inches)	Shipping Weight (lbs.)
22-037-263	10	manual	120V / 60Hz	Ref: 4 door, 3 adjustable Frz: 1 fixed	no	13.5 x 19.5 x 38.5 (ref.) 15 x 17.5 x 11.5 (frz.)	30 x 23.5 x 60	150





© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.





^{** 10&}quot; (for bottom 8"), 16" (for top 19.5")