



Thermo Scientific
GP Series Laboratory
Refrigerators and Freezers

your samples
our obsession



Introducing the Thermo Scientific GP Series Refrigerators and Freezers

Reliable Sample Protection and Storage

Ideal for the demanding requirements of daily lab applications, our range of general purpose (GP) refrigerators and freezers include:

▶ Thermo Scientific GP Series Refrigerators and Freezers

▶ Thermo Scientific GP Series Chromatography Refrigerators

All models of our Thermo Scientific GP Series refrigerators and freezers deliver the performance, security and quality features you have come to depend on from us.

- Wide range of capacities – from 12 cu. ft. to 72 cu. ft.
- Exceptional security features
- Outstanding temperature stability
- Intuitive digital controls
- Quality, robust construction



Thermo Scientific
GP Series Freezer



Thermo Scientific
GP Series Refrigerator



Thermo Scientific
GP Series Chromatography Refrigerators

Our complete line of GP Series refrigerators and freezers provide the reliability and security essential for the protection of sensitive lab samples in the +4° to -30°C range.



Thermo Scientific GP Series Refrigerators and Freezers

Reliable cold storage for
general-purpose applications

Our GP Series refrigerators and freezers offer the dependability and safety required for your daily lab applications.

- Range of sizes to suit your lab space needs
- Choice of painted or stainless steel interior finish (refrigerators only)
- Available with or without chart recorders



12 cu. ft. lab refrigerator with glass door, cat. no. 11670301



23 cu. ft. lab freezer with solid stainless steel door, cat. no. 11670337



27 cu. ft. lab refrigerator with solid stainless steel door, cat. no. 11670317



Performance

- Refrigerator cabinet temperature range: 1° to 12°C (factory setpoint 4°C)
- Freezer cabinet temperature range: -12° to -30°C (factory setpoint -30°C)
- Digital temperature controller with intuitive interface
- Excellent temperature stability of $\pm 3^{\circ}\text{C}$
- High-density, fluorine free insulation
- Automatic defrost

Security

- Keyed on/off switch
- Optional chart recorder
- High/Low Temperature, Visual/Audio alarms
- Choice of solid or glass lockable doors (refrigerators only)
- Easily accessible contacts for backup monitoring
- Door ajar alarm (except sliding glass doors)
- Control panel with key lock

Quality

- Painted, high-quality steel exterior
- Choice of painted or stainless steel interior (refrigerators only)
- Casters for easy mobility with leveling for better stability
- Noise levels – quiet refrigeration system
- 1" rear access port
- Certifications: UL listed



49 cu. ft. lab refrigerator with solid stainless steel doors and chart recorder, cat. no. 13991895



72 cu. ft. lab refrigerator with glass doors, cat. no. 11670335

Thermo Scientific GP Series Refrigerators and Freezers

SPECIFICATIONS & ORDERING INFORMATION

	Cat. No.	Capacity	Voltage	Plug Type	Interior Finish	Door[s]
Thermo Scientific GP Series Freezers	11670337	24 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670338	24 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670339	34 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670340	34 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670341	48 cu. ft.	208-230V 60Hz	L14-20	Stainless Steel	Solid stainless-steel door, double
	11670342	48 cu. ft.	208-230V 60Hz	L14-20	Stainless Steel	Solid stainless-steel door, double
Thermo Scientific GP Series Refrigerators	11670301	12 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	11670302	12 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	11670307	23 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, single
	11670308	23 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, single
	11670309	23 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670310	23 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670303	23 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	11670304	23 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	11670305	23 cu. ft.	115V 60Hz	5-15P	White	Solid stainless-steel door, single
	11670306	23 cu. ft.	115V 60Hz	5-15P	White	Solid stainless-steel door, single
	11670315	27 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, single
	11670316	27 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, single
	11670317	27 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670318	27 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, single
	11670311	27 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	11670312	27 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	11670313	27 cu. ft.	115V 60Hz	5-15P	White	Solid stainless-steel door, single
	11670314	27 cu. ft.	115V 60Hz	5-15P	White	Solid stainless-steel door, single
	11670319	38 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	11670320	38 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	11670323	45 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Sliding glass door, double
	13991893	45 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Sliding glass door, double
	11670321	45 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	11670322	45 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	11670329	49 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, double
	11670330	49 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, double
	11670331	49 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, double
	11670332	49 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Solid stainless-steel door, double
	11670325	49 cu. ft.	115V 60Hz	5-15P	White	Glass door, double
	13991894	49 cu. ft.	115V 60Hz	5-15P	White	Glass door, double
	13991896	49 cu. ft.	115V 60Hz	5-15P	White	Solid stainless-steel door, double
	13991895	49 cu. ft.	115V 60Hz	5-15P	White	Solid stainless-steel door, double
	11670335	72 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, triple
	11670336	72 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, triple
	11670333	72 cu. ft.	115V 60Hz	5-15P	White	Glass door, triple
	11670334	72 cu. ft.	115V 60Hz	5-15P	White	Glass door, triple

Add 3" to exterior height if using provided standard casters

	No. of Shelves	Interior Dimensions D x W x H	Exterior Dimensions* D x W x H	Chart Recorder	Shipping Weight
	3	28 x 22 x 60"	32 x 26 x 78.5" (Add 1 3/4" to Depth for handle)	No	392 lb.
	3	28 x 22 x 60"	32 x 26 x 78.5" (Add 1 3/4" to Depth for handle)	Yes	392 lb.
	3	28 x 31 x 60"	32 x 35 x 78.5" (Add 1 3/4" to Depth for handle)	No	429 lb.
	3	28 x 31 x 60"	32 x 35 x 78.5" (Add 1 3/4" to Depth for handle)	Yes	429 lb.
	6	28 x 48 x 60"	32 x 52 x 78.5" (Add 1 3/4" to Depth for handle)	No	630 lb.
	6	28 x 48 x 60"	32 x 52 x 78.5" (Add 1 3/4" to Depth for handle)	Yes	630 lb.
	4	20.75 x 21 x 45.75"	24 x 24 x 72"	No	230 lb.
	4	20.75 x 21 x 45.75"	24 x 24 x 72"	Yes	230 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78" (Add 1 3/4" to Depth for handle)	No	300 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78" (Add 1 3/4" to Depth for handle)	Yes	300 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78" (Add 1 3/4" to Depth for handle)	No	300 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78" (Add 1 3/4" to Depth for handle)	Yes	300 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78"	No	300 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78"	Yes	300 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78" (Add 1 3/4" to Depth for handle)	No	300 lb.
	4	28.5 x 24 x 61.75"	32 x 27 x 78" (Add 1 3/4" to Depth for handle)	Yes	300 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	No	386 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	Yes	386 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	No	386 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	Yes	386 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78"	No	386 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78"	Yes	386 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	No	386 lb.
	4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	Yes	386 lb.
	8	27 x 40.5 x 61.75"	32 x 43.5 x 78"	No	524 lb.
	8	27 x 40.5 x 61.75"	32 x 43.5 x 78"	Yes	524 lb.
	8	27 x 48.63 x 61.75"	32 x 52 x 78"	No	561 lb.
	8	27 x 48.63 x 61.75"	32 x 52 x 78"	Yes	561 lb.
	8	27 x 48.63 x 61.75"	32 x 52 x 78"	No	561 lb.
	8	27 x 48.63 x 61.75"	32 x 52 x 78"	Yes	561 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	No	597 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	Yes	597 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	No	597 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	Yes	597 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78"	No	597 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78"	Yes	597 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	No	597 lb.
	8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	Yes	597 lb.
	12	28.5 x 72 x 61.75"	32 x 75 x 78" (Add 1 3/4" to Depth for handle)	No	730 lb.
	12	28.5 x 72 x 61.75"	32 x 75 x 78" (Add 1 3/4" to Depth for handle)	Yes	730 lb.
	12	28.5 x 72 x 61.75"	32 x 75 x 78"	No	730 lb.
	12	28.5 x 72 x 61.75"	32 x 75 x 78"	Yes	730 lb.

Thermo Scientific GP Series Chromatography Refrigerators

Secure storage for chromatography applications

Our Thermo Scientific GP Series chromatography refrigerators combine performance and security required for everyday chromatography applications.

- Capacities from 27 cu. ft. to 72 cu. ft.
- Range of sizes: from single-door to large-capacity, three-door uprights
- Choice of painted or stainless interior



Side access ports

45 cu. ft. chromatography refrigerator with sliding glass doors, cat. no. 11670349



Performance

- Refrigerator cabinet temperature range: 1° to 12°C (factory setpoint 4°C)
- Excellent temperature stability of $\pm 3^{\circ}\text{C}$
- Bright digital display
- Automatic defrost

Security

- Optional chart recorder
- High/low temperature alarms
- Keyed on/off switch
- Door ajar alarm (except sliding glass door models)

Quality

- Painted, high-quality steel exterior
- Choice of painted or stainless steel interior
- Two 2-inch side access ports
- Energy-efficient, double pane, argon-filled glass doors
- Interior duplex outlet
- Certifications: UL listed



Backwall plugs



27 cu. ft. chromatography refrigerator
with glass door and stainless steel interior,
cat. no. 11670345

Thermo Scientific GP Series Chromatography Refrigerators

SPECIFICATIONS & ORDERING INFORMATION

	Cat. No.	Capacity	Voltage	Plug Type	Interior Finish	Door[s]
Thermo Scientific GP Series Chromatography Refrigerators	11670345	27 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, single
	11670346	27 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, single
	11670351	45 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Sliding glass door, double
	11670352	45 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Sliding glass door, double
	11670355	49 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, double
	11670356	49 cu. ft.	115V 60Hz	5-15P	Stainless Steel	Glass door, double
	11670359	72 cu. ft.	115V 60Hz	5-20P	Stainless Steel	Glass door, triple
	11670360	72 cu. ft.	115V 60Hz	5-20P	Stainless Steel	Glass door, triple
	11670343	27 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	11670344	27 cu. ft.	115V 60Hz	5-15P	White	Glass door, single
	13991891	38 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	13991892	38 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	11670349	45 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	11670350	45 cu. ft.	115V 60Hz	5-15P	White	Sliding glass door, double
	11670353	49 cu. ft.	115V 60Hz	5-15P	White	Glass door, double
	11670354	49 cu. ft.	115V 60Hz	5-15P	White	Glass door, double
	11670357	72 cu. ft.	115V 60Hz	5-20P	White	Glass door, triple
	11670358	72 cu. ft.	115V 60Hz	5-20P	White	Glass door, triple

Add 3" to exterior height if using provided standard casters

▶ Services designed to help you achieve your goals!

Responsive and knowledgeable laboratory service is only a phone call or email away, providing you with peace of mind:

- Increased equipment uptime, ensuring continuity of operations
- Broad network of factory trained and certified technicians
- Convenient and dependable access to our expertise and original parts
- A fixed cost solution with service maintenance options

No. of Shelves	Interior Dimensions D x W x H	Exterior Dimensions* D x W x H	Chart Recorder	Shipping Weight
4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	No	300 lb.
4	28.5 x 27 x 61.75"	32 x 30 x 78" (Add 1 3/4" to Depth for handle)	Yes	300 lb.
8	27 x 48.63 x 61.75"	32 x 52 x 78"	No	561 lb.
8	27 x 48.63 x 61.75"	32 x 52 x 78"	Yes	561 lb.
8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	No	597 lb.
8	28.5 x 48.63 x 61.75"	32 x 52 x 78" (Add 1 3/4" to Depth for handle)	Yes	597 lb.
12	28.5 x 72 x 61.75"	32 x 75 x 78" (Add 1 3/4" to Depth for handle)	No	730 lb.
12	28.5 x 72 x 61.75"	32 x 75 x 78" (Add 1 3/4" to Depth for handle)	Yes	730 lb.
4	28.5 x 27 x 61.75"	32 x 30 x 78"	No	300 lb.
4	28.5 x 27 x 61.75"	32 x 30 x 78"	Yes	300 lb.
8	27 x 40.5 x 61.75"	32 x 43.5 x 78"	No	524 lb.
8	27 x 40.5 x 61.75"	32 x 43.5 x 78"	Yes	524 lb.
8	27 x 48.63 x 61.75"	32 x 52 x 78"	No	561 lb.
8	27 x 48.63 x 61.75"	32 x 52 x 78"	Yes	561 lb.
8	28.5 x 48.63 x 61.75"	32 x 52 x 78"	No	597 lb.
8	28.5 x 48.63 x 61.75"	32 x 52 x 78"	Yes	597 lb.
12	28.5 x 72 x 61.75"	32 x 75 x 78"	No	730 lb.
12	28.5 x 72 x 61.75"	32 x 75 x 78"	Yes	730 lb.



Thermo Scientific Labware

Thermo Scientific plastic and glass labware offer a complete range of tubes, bottles, flasks and vials for all of your cold storage applications. Our labware products are designed and manufactured for tough applications under rigorous laboratory conditions. For more information, visit www.fishersci.com.



Warranty: 2 years parts and labor

© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



**For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com**

**In Canada: For customer service, call 1-800-234-7437
To fax an order: 1-800-463-2996
To order online: www.fishersci.ca**