

your samples our obsession

Temperature variation can have an impact on the viability of samples and the efficacy of reagents. These variations can shorten product shelf-life, degrade enzymes and reduce vaccine effectiveness. Therefore, it's important to choose a refrigerator that has been designed for the demanding standards of your laboratory.

Our Revco® high-performance lab refrigerators with glass doors are available in a wide range of cabinet sizes to fit virtually any application need. Each refrigerator comes with an advanced microprocessor temperature control, high-quality cabinet construction and advanced alarm systems. Positive, forced-air circulation maintains temperature uniformity, protecting your valuable samples and inventory.

- An advanced microprocessor control allows for precise regulation of internal conditions
- Directional airflow for uniform temperature throughout the cabinet
- Key-operated alarms and setpoint security
- Seven storage capacities ranging from undercounter to a spacious, three-door model
- Environmentally-friendly, mercury-free, LED interior lighting
- Easy-to-read digital displays
- Multi-pane, glass doors
- Standard 1-inch access ports
- Optional temperature recording charts to confirm storage parameters
- Compatible with optional, remote wireless monitoring



designed for sample protection



Powerful industrial compressors:
provide fast recovery from
door openings

Optional chart recorder: inkless,
seven day, graphic recorder

Positive forced-air circulation:
maintains temperature uniformity


Heavy-walled, CFC-free insulation:
ensure temperature stability and
minimize energy use

Adjustable shelves

**Robust, spring-loaded,
self-closing doors with
90° stay-open feature**

Set of four 2" casters:
front locking, rear non-locking

Microprocessor-based control:
for simplified operation and improved reliability



SPECIFICATIONS AND ORDERING INFORMATION

Temperature Range: +1°C to +8°C, Factory Pre-Set to +4°C (Automatic Defrost)

Cat. No.	Cu. Ft. (Liters)	Voltage (Hz)	Amps/ Breaker (Plug)	Cabinet	Door	Adjustable Shelves	Interior Dimensions* H x D x W inches (mm)	Exterior Dimensions* H x D x W inches (mm)	Shipping Weight lbs. (kg)
RGL404A RGL404D RGL404V RGL404W	4.9 (133)	115 (60) 208/230 (60) 230 (50) 230 (50)	8.0/15 (P1) 6.0/15 (P4) 6.0/15 (P7) 6.0/15 (P8)	Undercounter	Single, Glass	3	20.0 x 20.5 x 20.0 (508 x 521 x 508)	33.4 x 26.0 x 24.0 (848 x 660 x 610)	200 (91)
RGL1204A RGL1204D RGL1204V RGL1204W	11.5 (326)	115 (60) 208/230 (60) 230 (50) 230 (50)	9.4/15 (P1) 6.0/15 (P4) 6.0/15 (P7) 6.0/15 (P8)	Space Saver	Single, Glass	4	52.4 x 21.75 x 20.0 (1331 x 553 x 508)	73.6 x 29.4 x 24.0 (1869 x 748 x 610)	280 (127)
RGL2304A RGL2304D RGL2304V RGL2304W	23.3 (659)	115 (60) 208/230 (60) 230 (50) 230 (50)	9.4/15 (P1) 6.0/15 (P4) 6.0/15 (P7) 6.0/15 (P8)	Upright	Single, Glass	4	58.0 x 29.0 x 24.0 (1473 x 737 x 610)	79.2 x 37.2 x 28.0 (2012 x 945 x 711)	370 (168)
RGL3004A RGL3004D RGL3004V RGL3004W	29.2 (826)	115 (60) 208/230 (60) 230 (50) 230 (50)	9.4/15 (P1) 6.0/15 (P4) 6.0/15 (P7) 6.0/15 (P8)	Upright	Single, Glass	4	58.0 x 29.0 x 30.0 (1473 x 737 x 762)	79.2 x 37.2 x 34.0 (2012 x 945 x 864)	435 (197)
REL4504A REL4504D REL4504V REL4504W	45.8 (1297)	115 (60) 208/230 (60) 230 (50) 230 (50)	15.0/20 (P2) 8.5/15 (P4) 8.5/15 (P7) 8.5/15 (P8)	Upright	Double, Sliding Glass	8	58.0 x 26.0 x 52.5 (1473 x 660 x 1334)	79.1 x 34.5 x 56.5 (2009 x 876 x 1435)	670 (304)
RGL5004A RGL5004D RGL5004V RGL5004W	51.1 (1447)	115 (60) 208/230 (60) 230 (50) 230 (50)	15.0/20 (P2) 8.5/15 (P4) 8.5/15 (P7) 8.5/15 (P8)	Upright	Double, Glass	8	58.0 x 29.0 x 52.5 (1473 x 737 x 1334)	79.2 x 37.2 x 56.5 (2012 x 945 x 1435)	585 (265)
RGL7504A RGL7504D RGL7504V RGL7504W	78.8 (2231)	115 (60) 208/230 (60) 230 (50) 230 (50)	16.0/20 (P2) 10.5/15 (P4) 10.5/15 (P7) 10.5/15 (P8)	Upright	Triple, Glass	12	58.0 x 29.0 x 81.0 (1473 x 737 x 2057)	79.2 x 37.2 x 85.0 (2012 x 945 x 2159)	805 (365)

* Published volumes are nominal; useable cu. ft. (L) capacities may be offset by shelves, drawers, blower coils or cabinet protrusions. For specific cabinet dimensions and allowances, or for verification of clearance tolerances, contact your sales representative. Exterior dimensions do not reflect standard 2" (5cm) casters shipped with unit.

OPTIONS AND ACCESSORIES FOR GLASS DOOR REFRIGERATORS

For Use With Cat. No.:	RGL404	RGL1204	RGL2304	RGL3004	REL4504	RGL5004	RGL7504
Shelves							
Epoxy-Coated, Full Size	12370653	12370653	13990163	12370607	13990163	13990163	13990163
Epoxy-Coated, Half Size	N/A	N/A	13990164	12370609	13990164	13990164	13990164
Door Options and Accessories*							
LH Swing, Glass	N/A	12370627	12370629	12370631	N/A	N/A	N/A
LH Swing, Heated Glass	N/A	12370640	12370642	12370644	N/A	N/A	N/A
RH Swing, Heated Glass	N/A	13990601	12370643	12370645	N/A	N/A	N/A
LH/LH Swing, Glass	N/A	N/A	N/A	N/A	N/A	12370637	N/A
LH/LH Swing, Heated Glass	N/A	N/A	N/A	N/A	N/A	12370646	N/A
RH/RH Swing, Glass	N/A	N/A	N/A	N/A	N/A	12370639	N/A
RH/RH Swing, Heated Glass	N/A	N/A	N/A	N/A	N/A	12370647	N/A
LH/RH Swing, Heated Glass	N/A	N/A	N/A	N/A	N/A	12370648	N/A
Roll-Out Drawers							
Complete Set	12370654	12370655	13990167	12370616	12370615	13990170	12370613
Single Drawers	N/A	N/A	12370614	50868161	12370620	12370614	12370614
Casters (2" casters standard)							
3" Casters (76 mm), Set of 4**	50868420	50868420	50868420	50868420	13990162	50868420	50868420
Finishes (adds 1" or 25 mm to unit height)							
Stainless Steel Interior	12370649	12370657	12370600	12370617	12370618	12370601	12370602
Temperature Chart Recorders							
Built-in, Factory Installed	N/A	13990217	13990217	13990217	13990217	13990217	13990217
Free-standing, Field Installed***	13990220	13990220	13990220	13990220	13990220	13990220	13990220
Replacement Chart Paper, -40 to +25°C Box of 50	1399016A	1399016A	1399016A	1399016A	1399016A	1399016A	1399016A
Surge Suppressor							
Specify Voltage with Order	13990603	13990603	13990603	13990603	13990603	13990603	13990603
Deluxe Electronic Remote Alarm							
Requires Dedicated Phone Line	12370604	12370604	12370604	12370604	12370604	12370604	12370604
Remote Alarm Module							
Specify Voltage with Order	1398632	1398632	1398632	1398632	1398632	1398632	1398632
Alarm Delay Module							
Specify Voltage with Order	1399291	1399291	1399291	1399291	1399291	1399291	1399291
Digital Temperature Monitor							
Specify Voltage with Order	1399286	1399286	1399286	1399286	1399286	1399286	1399286
4-20mA Transmitter (factory installed)							
	13990604	13990604	13990604	13990604	13990604	13990604	13990604
+5°C Setpoint							
Must be Specified with Order	13990605	13990605	13990605	13990605	13990605	13990605	13990605

* Standard door configuration is; single door RH; double door LH/RH; triple door LH/RH/RH. LH or RH designates hinge on left or right side. Designations read from left to right facing unit.

** 2" casters are standard on all units except RGL404 where casters are shipped uninstalled.

*** Must be installed by a qualified professional

© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific