

LabServ Centrifuge Tubes




Priced to save you money without sacrificing quality.

Ranging from \$109.75 – \$142.75 for a case of 500.

Features:

- Made from polypropylene
- Available in rack or sleeve
- Flat screw top
- Sterile
- Capacity of 15mL or 50mL
- DNase/RNase
- RCF rating of 12,000 x g



LabServ Centrifuge Tubes

Learn more from your sales representative or by visiting www.fishersci.com today.

LabServ Polypropylene Centrifuge Tubes	Capacity	Material	Sterile	Closure Type	Graduations	Packaging Type	Quantity
14-955-159	15mL	Polypropylene	Yes	Flat Screw Top	1mL	Rack	500
14-955-160	15mL	Polypropylene	Yes	Flat Screw Top	1mL	Sleeve	500
14-955-161	50mL	Polypropylene	Yes	Flat Screw Top	2.5mL	Rack	500
14-955-162	50mL	Polypropylene	Yes	Flat Screw Top	2.5mL	Sleeve	500

© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

In the United States:

For customer service, call 1-800-766-7000
To fax an order, use 1-800-926-1166
To order online: www.fishersci.com

In Canada:

For customer service, call 1-800-234-7437
To fax an order, use 1-800-463-2996
To order online: www.fishersci.ca



**Fisher
Scientific**

A Thermo Fisher Scientific Brand