

*fisher*brand



PHARMA & BIOPHARMA

Chromatography Consumables for Pharma and Bio-Pharma Selections

Not every chromatographer has the need for the same vial & closure or the same sample handling container. The here shown selection represent the most used and helpful autosampler products for this segment of the business.

- Certified KITs
- Screw, Crimp and Snap Ring Vials & Closures
- Plastic and Microsampling Vials
- Well Plates and Sealmats
- Vial Rack

The selection has been made by experts in order to meet the special requirements of this segment and the connected challenges of the “to be analyzed” samples and molecules. It might not be complete but represents 90% of the core products.

LC/MS-GC/MS certified KITs & HPLC/GC certified KITs



13429748

HPLC/GC Certified Vial Kit: 1.5ml Short Thread Vial, clear glass, 1st hydrolytic class, label; UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1.0mm

13439748

HPLC/GC Certified Vial Kit: 1.5ml Short Thread Vial, clear glass, 1st hydrolytic class, label; UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1.0mm

15562320

LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, clear glass, wide opening, with overwind-barrier; Ultra High Performance Seal: PP Short Thread Cap, blue, centre hole; Silicone darkblue-translucent/PTFE natural, 35° shore A, 1,0mm

15572320

LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, amber glass, wide opening, with overwind-barrier; Ultra High Performance Seal: PP Short Thread Cap, blue, centre hole; Silicone darkblue-translucent/PTFE natural, 35° shore A, 1,0mm

Short Thread Vials ND9

fisherbrand



11575884

1.5ml Short Thread Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, label and filling lines



10080952

1.5ml Short Thread Vial, 32 x 11.6mm, amber glass, 1st hydrolytic class, wide opening, label and filling lines



15166425

1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, with overwind-barrier

Microsampling Vials ND9



12951011

Short Thread Vial with integrated Micro-Insert, 32 x 11.6mm, clear glass, 1st hydrolytic class, Base bonded



15592310

Short Thread Vial with integrated Micro-Insert, 32 x 11.6mm, amber glass, 1st hydrolytic class, Base bonded



11515894

1.1ml Microliter Short Thread Vial ND9, 32 x 11.6mm, clear glass, 1st hydrolytic class

Plastic Vials ND9



11727597

1.5ml PP Short Thread Vial, 32 x 11.6mm, transparent, with filling lines



11824940

0.7ml PP Short Thread Micro-Vial, 32 x 11.6mm, transparent



11707597


0.7ml PP Short Thread Micro-Vial, 32 x 11.6mm, transparent

PHARMA & BIOPHARMA

fisherbrand

PP Short Thread Seals ND9

fisherbrand



11863021

9mm Combination Seal: PP Short Thread Cap, blue, with centre hole; RedRubber / PTFE beige, 45° shore A, 1,0mm



11787567

UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1.0mm



11797567

9mm Combination Seal: PP Short Thread Cap, blue, centre hole; Silicone white/PTFE blue, 55° shore A, 1.0mm, slit

Snap Ring Vials and Seals ND11



11535894

1.5ml Snap Ring Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, label and filling lines



11545894

1.5ml Snap Ring Vial, 32 x 11.6mm, amber glass, 1st hydrolytic class, wide opening, label and filling lines

PE Snap Ring Seals ND11



15522320

11mm Combination Seal: PE-Snap Ring Cap, blue, with centre hole; RedRubber/PTFE beige, 45° shore A, 1.0mm



11814910

UltraClean Closure: 11mm PE Snap Ring Cap, blue, centre hole; Silicone white/PTFE red, 45° shore A, 1.3mm



11834910

11mm Combination Seal: PE Snap Ring Cap, blue, centre hole; Silicone white/PTFE blue, 55° shore A, 1.0mm, cross-slitted

**16233151**

Deep Well Microplate, PP, 96 Positions, certified, height 14.7mm, V-shape, 7mm dia., 220µl total volume (non coated, non sterile)

16213151

Deep Well Microplate, PP, 96 Positions, certified, height 41.6mm, U-shape, 7mm dia., 1000µl total volume (non coated, non sterile)

16223151

Square Well Microplate, PP, 96 Positions, certified, height 44.4mm, V-shape, 7mm dia., 2000µl total volume (non coated, non sterile)

**16263151**

Sealmat, MicroMat CLR, clear, Silicone, for 96 position Deep Well Microplate, round well - flat base, 7mm diameter

16243151

Sealmat, blue, Silicone/PTFE, for 96 position Deep Well Microplate, round well, flat base, 7mm diameter (non steril)


16253151

Sealmat, MicroMat CLR, clear, Silicone, for 96 position Square Well Microplate (non sterile)

15359939

Sealmat, blue, Silicone/PTFE, for 96 position Square Well Microplate

PP Storage Box and Vial Rack

**12692495**

PP Storage Box for 1,5ml vials or 2ml shell vials, blue, with cover (130 x 130 x 45mm), 81 cavities with alphanumeric coding of all 4 margins as well as the cavities at the bottom

11767746

PP Vial-Rack, (200 x 105 x 17mm), for 1.5ml vials, 50 cavities, blue, stackable

fisherbrand



PHARMA & BIOPHARMA

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Czech Republic: +420 466 798230
Denmark: +45 70 27 99 20 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Lithuania: +370 5 244 4442 Netherlands: +31 (0)20 487 70 00 Norway: +47 22 95 59 59
Portugal: +351 21 425 33 50 Spain: +34 902 239 303 Sweden: +46 31 352 32 00
Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

Other Central & Eastern Europe: +49 2304 932 715 Middle East: +44 1509 555939
Africa (Anglophone): +44 1509 555939 Afrique (Francophone): +33 3 8867 5327

© 2019 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**
part of Thermo Fisher Scientific