


Bata de aislamiento desechable

- Tela en láminas
- Polipropileno: 28 g/m²
- Polietileno: 15 g/m²
- Puños tejidos
- Azul
- PPE de categoría I
- No estéril
- No quirúrgica
- Abierta por detrás
- Con tiras para el cuello y la cintura
- Talla: universal



N.º cat.	Descripción	Cantidad
16804471	Bata de aislamiento desechable	2000

Dimensiones	Especificaciones	Tolerancia
Universal		
A. Ancho, parte superior	70cm	±2
B. Longitud, mangas	62cm	±2
C. Longitud, bata	118cm	±2
D. Ancho, abertura para brazos	28cm	±2
E. Longitud, abertura de los puños tejidos	5,5cm	±1.3
F. Longitud, del cuello hasta la parte inferior	109cm	±2
G. Longitud, tira para la cintura	175cm	±8.5
H. Longitud, tira para el cuello	32cm	±2.5
Peso	42gsm	±3

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
 Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
 Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
 Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
 Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**
 part of Thermo Fisher Scientific