


Blouse de protection à usage unique

- Tissu laminé
- Polypropylène : 28 g/m²
- Polyéthylène : 15 g/m²
- Poignets tricotés
- Bleu
- EPI de catégorie I
- Non stérile
- Non chirurgical
- Fente à l'arrière
- Avec attache au cou et à la taille
- Taille : Universelle



Ref.	Désignation	Qté
16804471	Blouse de protection à usage unique	2000

Dimensions	Spécification	Tolérance
Universelle		
A. Largeur, haut	70 cm	±2
B. Longueur, manche	62 cm	±2
C. Longueur, blouse	118 cm	±2
D. Largeur, emmanchure	28 cm	±2
E. Longueur, ouverture des poignets tricotés	5,5 cm	±1,3
F. Longueur, de l'encolure jusqu'au bas	109 cm	±2
G. Longueur, attache à la taille	175 cm	±8,5
H. Longueur, attache au cou	32 cm	±2,5
Poids	42gsm	±3

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
 Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
 Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
 Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
 Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**
 part of Thermo Fisher Scientific