


Cubremangas desechable de polietileno

Cubremangas desechable fabricado en polietileno de alta densidad. Protege de las salpicaduras en los sectores industrial, de fabricación y de asistencia sanitaria.

- Tela de PEAD
- Grosor: 18 µm
- Peso: 3 g
- Azul
- PPE de categoría I
- No estéril
- Muñeca y codo elásticos



N.º cat.	Descripción	Cantidad
16844471	Cubremangas con muñeca y codo elásticos	2000

Especificaciones técnicas

Dimensiones	Tolerancia
A. 40 cm de largo	±2
B. 22 cm de ancho	±2

Item Specification	Método de prueba (300mm/min)	Valor
Resistencia de tensión a la ruptura (N/5cm)	MD ISO 13934-1	>19
	CD ISO 13934-1	>12
Tasa de elongación (%)	MD ISO 13934-1	>50
	CD ISO 13934-1	>50

Envasado	
Paquete / Cantidad	100 unidades por bolsa 20 unidades por envase de cartón
Dimensiones del envase de cartón	40cm x 25cm x 3cm

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**
part of Thermo Fisher Scientific