

# Cubrezapatos SPP con suela antideslizante

Cubre-zapatos desechables de polipropileno, con suela antideslizante. Para usar sobre los zapatos en ambientes no peligrosos para protegerlos de la suciedad. Ideal para áreas no críticas o salas de entrada previas a la zona limpia.

- Tejido SPP
- Largo trecho: 44cm (A)
- Muy alto: 16cm (B)
- Peso: 38gsm
- Azul/blanco
- PPE de categoría I
- No estéril



N.º cat.	Descripción	Cantidad
16814471	Cubrezapatos SPP con suela antideslizante	300

Artículo	Especificaciones	Método de prueba	Valor
Peso (gsm)		ISO 3801	38±0.2
Propiedades de resistencia de tensión a la ruptura (N/5cm)		MD ISO 13934-1	>55
		CD ISO 13934-1	>45
Tasa de elongación (%)		MD ISO 13934-1	>30
		CD ISO 13934-1	>30

Envasado	
Paquete / Cantidad	100 fundas de zapatos por bolsa 3 unidades por envase de cartón
Dimensiones del envase de cartón	46 x 25 x 20cm

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**  
part of Thermo Fisher Scientific