

# Couvre-chaussure avec semelle antidérapante

Couvre-chaussure jetable en polypropylène, avec semelle antidérapante. À porter par-dessus vos chaussures dans des environnements non dangereux pour les protéger de la saleté et de la crasse. Idéal pour les zones non critiques ou les pièces d'entrée avant l'habillage.

- Tissu SPP
- Longue durée : 44 cm (A)
- Hauteur du talon : 16 cm (B)
- Poids : 38 g/m<sup>2</sup>
- Bleu/blanc
- EPI de catégorie I
- Non stérile



Ref.	Désignation	Qté
16814471	Couvre-chaussures avec semelle antidérapante	300

Spécification de l'article	Méthode de test	Valeur
Poids (g/m <sup>2</sup> )	ISO 3801	38+0.2
Propriétés de résistance à la tension (N/5 cm)	MD ISO 13934-1	>55
	CD ISO 13934-1	>45
Taux d'allongement (%)	MD ISO 13934-1	>30
	CD ISO 13934-1	>30

Conditionnement	
Emballage / quantité	100 couvre-chaussures par sac 3 sacs par cartons
Dimensions de la boîte	46 x 25 x 20cm

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**  
part of Thermo Fisher Scientific