# **KIMTECH**



# Why KIMTECH™?

Kimberly-Clark Professional<sup>TM</sup> partners with you to create exceptional workplaces, making them healthier, safer and more productive. When you work in challenging environments with no margin for error, count on Kimtech<sup>TM</sup> products to protect your critical processes and the people who perform them. We understand the risks and costs associated with contamination, Kimtech<sup>TM</sup> products are designed for precision tasks and exacting performance. We partner with you to deliver seamless, worry-free performance that minimise risk and maximise performance.



#### **Protect the Scientist**

Laboratory Safety is your #1 Priority. Kimtech™ gloves are specifically designed to protect Scientists from Chemical Splash and Micro-Organism Hazards.

- **EN ISO 374-1 Chemical splash protection certification**
- ▶ EN 374-4: 2003 Resistance to degradation by chemicals
- ➤ EN ISO 374-5:2016 Micro Organism and Virus protection, the highest level of micro-organism protection possible for laboratory gloves
- Non detectable level of chemical accelerators to mitigate risk of skin irritation



#### **Protect the Science**

Protection of procedures and assays from contamination is vital. Contamination costs time, money and credibility.

- Manufactured to reduce levels of residuals, outperforming competitor laboratory gloves on cleanliness
- > Certificates of Analysis for each production lot demonstrates low residue levels1
- Non detectable levels of silicone and amide



#### **Assured Compliance**

Kimberly-Clark Professional™ makes it easy to demonstrate compliance with Laboratory gloves.

- ▶ PPE Category III Certified for protecting from Chemical Splash and Micro-Organism
- ▶ Publish Certificates of Analysis (CoA) with every production lot to provide you with assurance of barrier quality and cleanliness of the gloves¹
- Our website www.kimtech.eu is a rich source of on-line documents, including Technical Data Sheets, Declaration of Conformity and Certificates of Analysis

## Kimtech™ Purple Nitrile™ Gloves



Industry-leading gloves offering unrivalled protection, cleanliness and quality.

Nitrile construction results in products that are stronger and leaner than latex gloves, and feature better protection against a wider range of chemicals, including cytotoxic drugs.

Gloves are anti-static tested to protect the wearer and equipment, and ambidextrous.

Textured fingertips enhance grip and tactile sensitivity for safer and more efficient processes.

Beaded cuffs add strength to the gloves, reducing the risk of tearing and increasing their durability, while also reducing roll down for easier donning and doffing.

Contain no natural rubber latex, silicone or powder, reducing the risks of skin irritation for the wearer.

Available in two lengths.



The high quality nitrile material delivers seamless protection whenever and wherever it is needed. The synthetic nitrile polymer material is designed for fit and reliability, with textured fingertips for improved grip and excellent water tightness (AQL 0.65) that results in a low risk of pinholes. The gloves are ambidextrous and incorporate a beaded cuff for added strength and ease of donning, so the wearer can simply grab and go without any fear of ripping the material. Our nitrile safety gloves are also anti-static tested to ensure no disruption to samples or equipment and are latex, silicone and powder-free.

Purple Nitrile<sup>™</sup> gloves keep hands comfortable and protected while ensuring that research applications can be carried out contamination-free. The gloves are designated as PPE Cat III according to (EU) Regulation 2016/425, and are ideal for use in higher-risk applications as well as being food contact approved.

Assured Compliance: PPE Cat III according to Regulation (EU) 2016/425, EN ISO 374-1:2016 Type B (JKT) Chemical Splash protection, EN 374-4:2003 Resistance to degradation by chemicals, EN ISO 374-5:2016 Micro Organism and Virus protection, Food contact approved

Description	XS	S	M	L	XL	Gloves / Box	Total / Case
Kimtech™ Purple Nitrile™ Gloves, 24/25 cm	90625	90626	90627	90628	90629	100 / 10	1,000
Kimtech™ Purple Nitrile™ Xtra™ Gloves, 30 cm	97610	97611	97612	97613	97614	50 / 10	500









#### Kimtech™ Sterling™ Nitrile Gloves

- Manufactured using proprietary Sterling<sup>™</sup> technology that combines security and comfort. The efficient, environmentally-friendly construction minimises waste without compromising safety.
- ➤ Nitrile construction results in products that are stronger and leaner than latex gloves, and feature certified protection against a wide range of contaminants while also being food contact approved.
- Gloves are anti-static tested to protect the wearer and equipment, ambidextrous and grey in colour.
- Textured fingertips enhance grip and tactile sensitivity for safer and more efficient processes.
- ▶ Beaded cuffs add strength to the gloves, reducing the risk of tearing and increasing their durability, while also reducing roll down for easier donning and doffing.
- ➤ Contain no natural rubber latex, silicone or powder, reducing the risks of skin irritation for the wearer.
- Available in two lengths.



Assured Compliance: PPE Cat III according to Regulation (EU) 2016/425, EN ISO 374-1:2016 Type C (K) Chemical Splash protection, EN 374-4:2003 Resistance to degradation by chemicals, EN ISO 374-5:2016 Micro Organism and Virus protection, Food contact approved

Description	XS	S	М	L	XL	Gloves / Box	Total / Case
Kimtech™ Sterling™ Nitrile Gloves, 24/25 cm	99210	99211	99212	99213	99214	150 / 10 (XL= 140/10)	1,500 (XL=1400)
Kimtech™ Sterling™ Nitrile-Xtra™ Gloves, 30 cm	98341	98342	98343	98344	98345	100 / 10 (XL = 90/10)	1,000 (XL=900)

# K










#### Certificate of Analysis (CoA) with every production lot

(for Kimtech™ Purple Nitrile / Purple Nitrile™ Xtra™ / Sterling™ Nitrile / Sterling™ Nitrile-Xtra™ Gloves only)









## Kimtech™ Green Nitrile Gloves

- > Textured finish enhances grip and tactile sensitivity.
- **)** Beaded cuffs reducing the risk of tearing and increasing durability.

Kimtech<sup>TM</sup> Green Nitrile Ambidextrous Gloves offer superb protection with an environmentally-friendly glove thickness of just 0.06 mm. The comfortable design is powder- and latex-free and guards against a number of common contaminants.



Description	XS	S	М	L	XL	Gloves / Box	Total / Case
Kimtech™ Green Nitrile Gloves, 24/25 cm	99850	99851	99852	99853	99854	250 / 6 (XL = 225/6)	1,500 (XL = 1,350)







## Kimtech™ Comfort Nitrile Gloves

- Innovative nitrile technology for gloves engineered for comfort and reduced hand fatigue.
- > Beaded cuffs reducing the risk of tearing and increasing durability.

Latex free Kimtech™ Comfort Nitrile Ambidextrous Gloves are designed to maximise comfort, quality and cleanliness. The ambidextrous gloves are suitable for low risk applications and protective against chemical splash, viruses and micro-organisms.



Description	XS	S	М	L	XL	Gloves / Box	Total / Case
Kimtech™ Comfort Nitrile Gloves, 24/25 cm	47672	47673	47674	47675	47676	150 / 10	1,500









## Kimtech™ Blue Nitrile Gloves

- Smooth bisque finish to enhance grip and tactile sensitivity in wet or dry conditions.
- Beaded cuffs reducing the risk of tearing and increasing durability.

Kimtech™ Blue Nitrile Ambidextrous Gloves provide excellent tactile sensitivity and superb grip in wet and dry environments. The latex, silicone and powder free gloves are engineered for a comfortable, skin friendly fit and extended use.



Description	S	М	L	XL	Gloves / Box	Total / Case
Kimtech™ Blue Nitrile Gloves, 24/25 cm	97982	97983	97984	97985	100 / 10 (XL= 90 /10	1,000 (XL=900)







## Kimtech™ PFE Latex Gloves

- > Textured palm and fingertips enhance grip and tactile sensitivity.
- **>** Powder free rubber latex to reduce the risks of skin irritation.

Kimtech<sup>TM</sup> PFE Latex Ambidextrous Gloves provide comfortable, skin-friendly contamination protection. The non-sterile PFE latex gloves also feature a textured finish and beaded cuff to add strength and simplify donning.



Description	XS	S	М	L	XL	Gloves / Box	Total / Case
Kimtech™ PFE Latex Gloves, 24/25 cm	E110	E220	E330	E440	E550	100 / 10	1,000
Kimtech™ PFE-Xtra Latex Gloves, 30 cm		50501	50502	50503	50504	50 / 10	500







## Kimtech™ Satin Plus Latex Gloves

- > Smooth finish for excellent tactile sensivity.
- **>** Powder free rubber latex to reduce the risks of skin irritation.

Kimtech<sup>TM</sup> Satin Plus Latex Ambidextrous Gloves feature a smooth finish that is ideal for double donning and delivers sensitivity and performance in many applications. The non-sterile latex gloves are also rigorously tested and feature a beaded cuff.



Description	xs	S	М	L	XL	Gloves / Box	Total / Case
Kimtech™ Satin Plus Latex Gloves, 24/25 cm	SP2110E	SP2220E	SP2330E	SP2440E	SP2550E	100 / 10 (XL = 90/10)	1,000 (XL=900)







## **All-round protection for Scientific Applications**

Applications									
Molecular Biology Biochemistry	Chemistry	Analytical Chemistry	Virology	Genomics	Proteomics	Forensic Science	Non-Sterile Pharmaceutical Production		
Main Technique	Primary Risk		Secondary Risk		Primary/ Alter	native Choi	ce		
Electrophoresis	Chemicals		Contamination		-10		dis		
Polymerase chain reaction	Contamination		Biohazard						
DNA extraction	Chemicals		Contamination						
Chemical Synthesis & Derivatization	Chemicals		Contamination						
Pathology	Biohazard		Contamination			,	<u> </u>		
Toxicology	Chemicals		Biohazard		JW-0-	,	***		
Cell Culture	Contamination		Sensitive Instru	ments	The same of the sa				
Microscopy	Contamination		Sensitive Instru	ments					
Spectroscopy	Contamination		Sensitive Instru	ments					
Chromatography	Chemicals		Contamination						

# Kimtech<sup>™</sup> Laboratory Gloves Selection Guide

	Product Code	Description	Sizing	Chemical Protection EN ISO 374-1	Biohazard Protection EN ISO 374-5	AQL	Glove Length		Chemical
	97610- 97614	Purple Nitrile™ Xtra™ Gloves	XS - XL	Type B (JKT)	Virus	0.65	30cm	+	py & cations
	90625 -90629	Purple Nitrile™ Gloves	XS - XL	Type B (JKT)	Virus	0.65	24-25cm		HIGHEST Chemotherapy & High Risk Applications
	98341 -98345	Sterling™ Nitrile Xtra™ Gloves	XS - XL	Type C (K)	Virus	0.65	30cm		Ξ
*	99210 -99214	Sterling™ Nitrile Gloves	XS - XL	Type C (K)	Virus	0.65	24-25cm	ÆL	ılash
	99850 -99854	Green Nitrile Gloves	XS - XL	Type C (K)	Virus	0.65	24-25cm	CHEMICAL PROTECTION LEVEL	Moderate chemical splash
	47672 -47676	Comfort Nitrile Gloves	XS - XL	Type C (K)	Virus	0.65	24-25cm	EMICAL PRO	Moderate
*	97982 - 97985	Blue Nitrile Gloves	S-XL	Type C (K)	Virus	1.5	24-25 cm	CHI	
	50501 -50504	PFE-Xtra Latex Gloves	S - XL	Type C (K)	Virus	1.5	30cm		<b>T</b> it splash
	E110 -E550	PFE Latex Gloves	XS - XL	Type C (K)	Virus	1.5	24-25cm		<b>LOWEST</b> Low risk - Light splash
	SP2110E -SP2550E	Satin Plus Latex Gloves	XS - XL	Type C (K)	Virus	1.5	24-25cm		Ľ



## Kimtech™ A7 P+ Laboratory Coat

- High quality disposable laboratory coats with triple stitched seams and breathable, multi-layered material.
- Enhanced all-round protective features are included, such as high mandarin collars, thumb-loops and extra-length arms to fit lab gloves over coats.
- The coats are individually packed making them suitable for ISO 7 (Grade C) areas.

With a cloth-like feel and a reduced risk of heat stress, Kimtech™ A7 P+ lab coats ensure contamination control during extensive periods of use. Made from non-woven SMS fabric with triple stitched seams, they are silicone and latex free.



Description	S	М	L	XL	XXL
Kimtech™ A7 P+ Laboratory Coat	97700	97710	97720	97730	97740









## Kimtech™ Science Wipes

- For delicate equipment, surfaces, optics and materials, our low-lint and high absorbency white wipes have you covered.
- Light duty, white-coloured wipes made from 100% pure fibre.
- > Provided in a hygienic, user-friendly pop-up box.

#### **Precision Wipes**

Low-lint, high absorbency and chemically inert wipes for delicate tasks and instrumentation in a range of applications. The 1-ply white-coloured wipes are provided in a convenient pop-up box.

Product Code	Sheets / Box	Ply	Colour	Size W x L	Total / Case
7551	196 x 15	1	White	30.5 x 30.5 cm	2,940
7552	280 x 30	1	White	11.5 x 21 cm	8,400



#### **Delicate Task Wipes**

Low-lint, high absorbency and chemically inert wipes for delicate tasks and instrumentation in a range of applications. The 2-ply white-coloured whites are provided in a convenient pop-up box.

Product Code	Sheets / Box	Ply	Colour	Size W x L	Total / Case
7557	100 x 24	2	White	20 x 21 cm	2,400
7558	200 x 15	2	White	20 x 20.5 cm	3,000



## WypAll® L Wipers



This range of wipers excels when a single-use wiper for mopping up spills, cleaning hands and wiping surfaces is needed.

The unique AIRFLEX<sup>™</sup> Fabric structure, with large air pockets, absorbs more, is stronger and thicker than a conventional tissue wiper structure and requires fewer sheets for each job. Ideal for cleaning large surface spills, the WypAll® L Wiper range is both task and cost efficient.

Product Description	Code	Ply	Case Count	Sheet size L x W in cm	Dispenser Code
WypAll® L20 EXTRA PLUS Wipers - Large Roll Blue	7301	2	1x500	38 x 33	6146, 6154, 6155
WypAll® L20 EXTRA PLUS Wipers - Large Roll Blue	7300	2	1x500	38 x 23.5	6146, 6154, 6155
WypAll® L20 EXTRA Wipers - Centerfeed Roll White	7303	2	6x300	42.5 x 18.5	7087
WypAll® L30 Ultra + Wipers - Large Roll Blue	7426	3	1x750	38 x 33	6146, 6154, 6155
WypAll® L40 Folded Wipers	7471	1	18x56	30.4 x 31.7	7969

## WypAll® X Cloths

The reusable WypAll® X Cloths are absorbent and solvent resistant and strong enough to tackle tough tasks, whilst also designed to be cost and time efficient. Using the patented HYDROKNIT™ Material, WypAll® X Cloths keep their shape and can be used again and again. In addition, the unique PowerPockets technology helps trap dirt and grime quicker and more effectively, leaving surfaces cleaner and reducing waste.

Product Description	Code	Ply	Case Count	Sheet size L x W in cm	Dispenser Code
WypAll® X60 Cloths - Large Roll White	6036	1	1x750	38 x 42	6146, 6154, 6155
WypAll® X70 Cloths - Quarter Fold White	8387	1	12x76	30.5 x 31.8	7969
WypAll® X70 Cloths - Large Roll White	8384	1	1x500	38 x 42	6146, 6154, 6155
WypAll® X80 Cloths - Large Roll White	8377	1	1x475	34 x 31.5	6146, 6154, 6155



#### **Assured Compliance**

European Personal Protective Equipment Regulation (EU) 2016/425 – PPE Cat III Complex Design: For use in applications where effects are irreversible or mortal risk is present.

#### Hand Protection - Chemical Splash Protection

Old:



EN374-1:2003 Low Chemical Protection Protection against Chemical Splash. Ask Manufacturer for breakthrough times against of chemicals being used.

Test: Permeation EN374- 3:2003

EN 374-3 :20

XYZ

EN374-3:2003 Chemical Protection Protection against Chemical Permeation. List of 12 chemicals, each assigned a letter a to l, breakthrough time of >30 minutes on at least 3 chemicals. Test: Permeation EN374-3:2003

New:

EN ISO 374-1:2016 / Type C



Type C: At least Level 1 performance (more than 10 minutes) against at least 1 chemical on the list. Test: Permeation EN 16523-1:2015 EN ISO 374-1:2016 / Type B



Type B: At least Level 2 performance (more than 30 minutes) against at least 3 chemicals on the list. Test: Permeation EN 16523-1:2015 EN ISO 374-1:2016 / Type A



UVWXYZ

Type A: At least Level 2 performance (more than 30 minutes) against at least 6 chemicals on the list. Test: Permeation EN 16523-1:2015

#### Existing List of 12 Chemicals

- A Methanol
- **B** Acetone
- C Acetonitrile
- **D** Dichloromethane
- E Carbon disulphide
- F Toluene

- **G** Diethylamine
- H Tetrahydrofuran
- I Ethyl acetate
- J n-heptane
- K Sodium hydroxide 40%
- L Sulphuric acid 96%

#### Additional 6 Chemicals

- M Nitric acid 65%
- N Acetic acid 99%
- O Ammonium hydroxide 25%
- P Hydrogen peroxide 40%
- S Hydrofluoric acid 40%
- T Formaldehyde 37%

#### Hand Protection - Micro-organisms hazard protection

EN 374-2 :2003



> Level 2

EN374-2:2003 Micro-organism hazard protection against micro-organisms, Level 3 - AQL 0.65 / Level 2 - AQL 1.5 Test: water-tightness, ISO2859 sampling by lot.

EN ISO 374-5:2016



VIRUS

For protection against bacteria and fungi a penetration test is required using EN 374-2:2014. For protection against viruses, compliance to ISO 16604:2004 (method B) standard is necessary.

#### **Chemical Permeation Classifications**

Classification	Unclassified	Level 1	Level 2	Level 3	Level 4	Level 5	Level 6	
Time	<10	10-30	30-60	60-120	120-240	240 - 480	>480	
Usage	Not recommended	Splash protection		Medium p	rotection	High protection		




## **Chemical Permeation Data**

- Visit www.kimtech.eu, from the top menu select resources and then chemical permeation
- Select up to 4 gloves to compare. Select either breakthrough times in minutes, or classification



> Submit and download your PDF






## **General Product Selection Guide**

	Product Code	Size /Format	Total Case content	Material	Characteristics	ISO or Class	Regulatory	Chemical Splash Protection	Micro-organisms Protection
Sterile Cleanroom Ra	ange								
Kimtech™ Sterile Cleanroom									
G3 sterile Sterling™ Nitrile Gloves	11821 to 11828	6.0 to 10.0	300 pairs	Nitrile	slick - double donning	ISO 5 Class A/B	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus
G3 sterile White Nitrile Gloves	HC61160 to HC61110	6.0 to 10.0	200 pairs	Nitrile	slick - double donning	ISO 5 Class A/B	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus
G3 sterile Latex Gloves	HC1360S to HC1310S	6.0 to 10.0	200 pairs	Latex	slick - double donning	ISO 5 Class A/B	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus
G5 sterile Latex Gloves	HC1160S to HC1110S	6.0 to 10.0	200 pairs	Latex	slick - double donning	ISO 5 Class A/B	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus
Kimtech™ Sterile Cleanroom	Apparel								
A5 Sterile Cleanroom Coverall	88800 to 88806	S to 4XL	25 garments	SMS	vacuum packed	ISO 5 Class A/B	PPE Cat III	Type 5 Type 6	
A5 Sterile Cleanroom Hood with ties	25797	Universal	100 hoods	SMS	vacuum packed	ISO 5 Class A/B	PPE Cat III	Type 6 (PB)	
A5 Sterile Cleanroom Hood & Mask	36072	Universal	75 hoods	SMS	vacuum packed	ISO 5 Class A/B	PPE Cat III	Type 6 (PB)	
A5 Sterile Overboots with vinyl sole	12922 / 88808	S/M and univ.	100 pairs	SMS	vacuum packed	ISO 5 Class A/B	PPE Cat I		
A5 Sterile Overboots with wrap-around vinyl foot	31683 / 31696	S/M and univ.	100 pairs	SMS	vacuum packed	ISO 5 Class A/B	PPE Cat I		
A5 Sterile sleeves	36077	Length 45 cm	100 pairs	SMS	vacuum packed	ISO 5 Class A/B	PPE Cat I		
Kimtech™ Sterile Face Masks	;								
M3 Certified Sterile Pouch-Style Face Mask	62483	Width 25 cm	200 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with headbands	ISO 5 Class A/B			
M3 Certified Sterile Pleat-Style Face Mask	62467	Width 18 cm	200 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with ties	ISO 5 Class A/B			
M3 Certified Sterile Pleat-Style Face Mask	62470	Width 18 cm	200 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with earloops	ISO 5 Class A/B			
M3 Certified Sterile Pleat-Style Face Mask	62494	Width 23 cm	200 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with ties	ISO 5 Class A/B			
Non-Sterile Cleanroo	m Range								
Kimtech™ Cleanroom Gloves									
	99234 to					ISO 3 or higher		EN ISO 374-1	
G3 Sterling™ Nitrile Gloves	99239 HC61010 to	XS to XL	1500 gloves	Nitrile	tacky	Class C/D	PPE Cat III	Type C (K) EN ISO 374-1	EN ISO 374-5 Virus
G3 White Nitrile Gloves	HC61014 62990 to	XS to XL	1000 gloves	Nitrile	tacky slick -	Class C/D	PPE Cat III	Type C (K) EN ISO 374-1	EN ISO 374-5 Virus
G3 NxT™ Nitrile Gloves	62995 HC225 to	XS to XL	1000 gloves	Nitrile	double donning	Class C/D	PPE Cat III	Type C (K) EN ISO 374-1	EN ISO 374-5 Virus
G3 Latex Gloves	HC555	S to XL	1000 gloves	Latex	smooth	ISO 3 or higher Class C/D	PPE Cat III	Type C (K)	EN ISO 374-5 Virus
G5 Sterling™ Nitrile Gloves	98184 to 98189	XS to XL	1500 gloves	Nitrile	tacky	ISO 3 or higher Class C/D	PPE Cat III	EN ISO 374-1 Type C (K) EN ISO 374-1	EN ISO 374-5 Virus
G5 Latex Gloves	HC2211 to HC5511	S to XL	1000 gloves	Latex	smooth	ISO 3 or higher Class C/D	PPE Cat III	Type C (K)	EN ISO 374-5 Virus
Kimtech™ Apparel									
A8 Particles Protection Coverall	47691 to 47697	S to 4XL	25 garments	SMS <sup>(4)</sup>	individually packed	ISO 6 or higher Class C/D	PPE Cat III	Type 5 Type 6	
A8 overshoes	39370 to 39372	S / Univ/ XL	300 overshoes	SMS <sup>(4)</sup>	fully coated	ISO 6 or higher Class C/D	PPE Cat I		
Kimtech™ Cleanroom Face M	asks								
M3 Certified Pouch-Style Face Mask	62484	Width 25 cm	300 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with headbands	ISO 3 or higher Class C/D			
M3 Certified Pleated-Style Face Mask	62451	Width 23 cm	500 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with earloops	ISO 3 or higher Class C/D			
M3 Certified Pleated-Style Face Mask	62452	Width 23 cm	500 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with ties	ISO 3 or higher Class C/D			
M3 Certified Pleated-Style Face Mask	62465	Width 18 cm	500 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with earloops	ISO 3 or higher Class C/D			
M3 Certified Pleated-Style Face Mask	62466	Width 18 cm	500 masks	PE <sup>[1]</sup> /[PE/ PET] <sup>[2]</sup> /PP <sup>[3]</sup>	with ties	ISO 3 or higher Class C/D			
M6 Certified Pleated-Style Face Mask	62477	Width 18 cm	500 masks	PP <sup>[3]</sup> /[PE/ PET] <sup>[2]</sup>	with earloops	ISO 3 or higher Class C/D			
Kimtech™ Cleanroom Face Ve	eils								
Kimtech™ Cleanroom Face Ve M7 Face Veils	eils 62751	33 x 38 cm	150 veils	Low lint PE <sup>[1]</sup>	with snaps	ISO 6 or higher			

	Product Code	Size / Format	Total Case content	Material	Characteristics	ISO or Class	Regulatory	Chemical Splash Protection	Micro-organisms Protection	Food Contact
Laboratory Range										
Kimtech™ Laboratory Gloves										
Purple™ Nitrile Gloves	90625 to 90629	XS to XL	10 x 100 <sup>(5)</sup> gloves	Nitrile	textured fingertips	ISO 6 or higher Class C/D	PPE Cat III	EN ISO 374-1 Type B (JKT)	EN ISO 374-5 Virus	Approved
Purple™ Xtra™ Nitrile Gloves	97610 to 97614	XS to XL	10 x 50 gloves	Nitrile	textured fingertips	ISO 6 or higher Class C/D	PPE Cat III	EN ISO 374-1 Type B (JKT)	EN ISO 374-5 Virus	Approved
Sterling™ Nitrile Gloves	99210 to 99214	XS to XL	10 x 150 <sup>(6)</sup> gloves	Nitrile	textured fingertips	ISO 6 or higher Class C/D	PPE Cat III	EN ISO 374-1 Type C(K)	EN ISO 374-5 Virus	Approved
Sterling™ Xtra™ Nitrile Gloves	98341 to 98345	XS to XL	10 x 100 gloves	Nitrile	textured fingertips	ISO 6 or higher Class C/D	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus	Approved
Green Nitrile Gloves	99850 to 99854	XS to XL	10 x 250 <sup>(7)</sup> gloves	Nitrile	textured fingertips	ISO 8 or higher	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus	
Comfort Nitrile Gloves	47672 to 47676	XS to XL	10 x 150 gloves	Nitrile	smooth	ISO 8 or higher	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus	Approved
Blue Nitrile Gloves	97982 to 97985	S to XL	10 x 100 <sup>(5)</sup> gloves	Nitrile	smooth	ISO 8 or higher	PPE Cat III	EN ISO 374-1 Type C(K)	EN ISO 374-5 Virus	
Satin Plus Latex Gloves	SP2110E to SP2550E	XS to XL	10 x 100 <sup>(5)</sup> gloves	Latex	smooth	ISO 8 or higher	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus	
PFE Latex Gloves	E110 to E550	XS to XL	10 x 100 <sup>(5)</sup> gloves	Latex	textured fingertips	ISO 8 or higher	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus	
PFE Xtra Latex Gloves	50501M to 50504M	S to XL	10 x 50 gloves	Latex	textured fingertips	ISO 8 or higher	PPE Cat III	EN ISO 374-1 Type C (K)	EN ISO 374-5 Virus	
Kimtech™ Lab Coats										
A7 Lab Coat	96700 to 96740	S to 2XL	15 lab coat	SMS <sup>(4)</sup>	individually packed	ISO 8 or higher	PPE Cat III	Type 6 (PB)		
A7P+ Lab Coat	97700 to 97740	S to 2XL	15 lab coat	Laminated SMS <sup>(4)</sup>	individually packed	ISO 8 or higher	PPE Cat III	Type 6 (PB)	EN 14126:2003	
Kimberly-Clark™ Face Mask										
Procedure Masks Pleated-Style (Blue)	47080	Width 18 cm	500 masks	PP <sup>(3)</sup> / wCellulose/ PU <sup>(8)</sup>	with earloops	ISO 8 or higher	MD Class I			
Procedure Masks Pleated-Style (White)	47090	Width 18 cm	500 masks	PP <sup>(3)</sup> / Cellulose/ PU <sup>(8)</sup>	with earloops	ISO 8 or higher	MD Class I			
Fluidshield Masks Pleated-Style (Orange)	47147	Width 18 cm	100 masks	PP <sup>(3)</sup> / Polyester/ PE <sup>(1)</sup> / Cellulose/ PU <sup>(8)</sup>	with wraparound visor and earloops	ISO 8 or higher	MD Class I			
Procedure Masks Pleated-Style (Blue)	48100	Width 18 cm	300 masks	PP <sup>(3)</sup> / Cellulose / Polyester	with ties	ISO 8 or higher	MD Class I			
Kimtech™ Science Wipes										
Precision Wipes (White)	7551	30 x 30 cm	15 boxes x 196 sheets	Paper	1 ply	ISO 8 or higher				
Precision Wipes (White)	7552	11,5 x 21 cm	30 boxes x 280 sheets	Paper	1 ply	ISO 8 or higher				
Delicate Task Wipes (White)	7557	20 x 21 cm	24 boxes x 100 sheets	Paper	2 ply	ISO 8 or higher				
Delicate Task Wipes (White)	7558	20 x 21 cm	15 boxes x 200 sheets	Paper	2 ply	ISO 8 or higher				



<sup>\*</sup>For more details on the products please refer to the corresponding technical data sheet

<sup>[1]</sup> PE: polyethylene

 $<sup>^{\</sup>text{\tiny{[2]}}}\,\text{PE/PET:}$  polyethylene/ polyethylene terephtalate

<sup>[3]</sup> PP :polypropylene

<sup>[4]</sup> SMS : spunbond - meltblown- spunbond

 $<sup>^{[5]}</sup>$  10 x 90 for size XL

<sup>&</sup>lt;sup>(6)</sup> 10 x 140 for size XL

 $<sup>^{\</sup>text{\tiny{[7]}}}$  10 x 225 for size XL

<sup>[8]</sup> PU : polyurethane

# KIMTECH

# BE YOUR BEST

> WE'VE GOT YOU COVERED

The Kimtech™ Brand offers you additional resources to help you select the right product for the right application.

For more information visit www.kimtech.eu

➤ For technical assistance contact us at kimtech.support@kcc.com



It is the employer's responsibility to assess the risk of the task to be undertaken and determine the correct choice of personal protective equipment for the task. The manufacturer, Kimberly-Clark, does not accept any responsibility for the incorrect choice or misuse of the personal protective equipment shown in this brochure. All care has been taken to ensure that the information contained herein is as accurate as possible at the time of publication, however errors may occur and legislation concerning personal protective equipment is under constant review and may change in the lifetime of this brochure. Accordingly, the specification for the products may be subject to change. We would advise you to contact your local sales representative if you have any queries concerning the products shown or the suitability of such products for a particular task. Always dispose of used protective equipment in a safe and appropriate manner in accordance with European, National and Local environmental regulations.

 $^{\text{M}}$  Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates.  $^{\text{CKWW}}$ . The colour Purple and Purple Nitrile  $^{\text{M}}$ , the colour Grey and Sterling  $^{\text{M}}$  are trademarks of Owens & Minor, Inc. or its affiliates and used under license.

**& Kimberly-Clark**PROFESSIONAL™