

# Nitrile Disposable Gloves

The Fisherbrand™ range of nitrile single-use gloves are designed to provide maximum protection and comfort. They are the ideal gloves for everyday protection in the lab with new packaging designed to reduce waste and environmental impact.

## Nitrile Indigo Disposable Gloves, PPE Category III

Choose a glove with excellent protection for a wide range of applications and general use. Fisherbrand Indigo gloves are powder-free and latex-free, with textured fingertips.

- Powder-free, latex-free
- Ambidextrous
- Non-sterile
- 200 gloves per standard box, reduces waste and environmental impact
- Thickness: 0.06mm to 0.09mm (palm), 0.08mm to 0.14mm (finger)
- Colour: indigo
- Conforms to EN 420:2003, EN ISO 374-1/Type C, EN ISO 374-5, EN ISO 374-4



### Ordering Information

Cat. No.	Size	Pack qty
17112192	XS	200
17192182	S	200
17182182	M	200
17172182	L	200
17102192	XL	200

## Nitrile Disposable Gloves, Comfort, PPE Category III

Choose a glove with enhanced comfort and tactile sensitivity. Fisherbrand comfort nitrile gloves are powder-free and latex-free with a beaded cuff and textured fingers.

- Powder-free, latex-free
- Exceptional fit, feel and durability
- Ambidextrous with textured fingers and beaded cuff
- Non-sterile
- 200 gloves per standard box (size XL 180 gloves per box) reduces waste and environmental impact
- Thickness: 0.07mm (palm), 0.09mm (finger)
- Colour: blue
- Length 241.3mm
- Conforms to EN ISO 374/Type B
- Certified for food contact



### Ordering Information

Cat. No.	Size	Pack qty
15622367	XS	200
15632367	S	200
15642367	M	200
15652367	L	200
15662367	XL	180

## Nitrile Disposable Gloves, Extended Cuff, PPE Category III

Choose a glove with exceptional durability and added protection. The Fisherbrand extended cuff nitrile gloves are extra-thick and fully textured for a superior grip.

- Powder-free, latex-free
- Fully textured finish for a superior wet or dry grip
- Ambidextrous and non-sterile with beaded cuff
- Thickness: 0.15mm (palm), 0.18mm (finger)
- Colour: teal green
- Length: 305mm
- Conforms to EN ISO 374/Type B
- Certified for food contact

### Ordering Information

Cat. No.	Size	Pack qty
15612367	XS	50
15662357	S	50
15672357	M	50
15682357	L	50
15692357	XL	50



## Nitrile Disposable Gloves with Aloe Vera, PPE Category III

These gloves feature a textured surface to provide ideal grip for dry and wet working conditions. They are green for easy identification, and have a soothing Aloe vera coating on the inner surface making them ideal for prolonged use.

- Premium nitrile gloves offering the comfort and feel of latex
- Comfortable, strong and durable
- Aloe vera coating on the inner surface to protect hands from drying, chapping and cracking
- 100 ambidextrous gloves per standard box reduces waste and environmental impact
- Colour: mint green
- PPE complies with Category III, 1.5 AQL, EN ISO 374/Type B
- Thickness 0.06mm

### Ordering Information

Cat. No.	Size	Pack qty
15193269	S	100
15103279	M	100
11772779	L	100
15113279	XL	100



Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2021 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**  
part of Thermo Fisher Scientific