



**PHARMA &  
BIOPHARMA**

## Essential Vial and Closure Solutions

Most used and helpful auto-sampler products  
for Pharma and BioPharma

# Chromatography Consumable for Pharma and BioPharma Selections

Not every chromatographer has the need for the same vial & closure or the same sample handling container. The here shown selection represent the most used and helpful autosampler products for this segment of the business.

- Certified KITs
- Screw, Crimp and Snap Ring Vials & Closures
- Plastic and Microsampling Vials
- Well Plates and Sealmats
- Vial Rack

The selection has been made by experts in order to meet the special requirements of this segment and the connected challenges of the “to be analyzed” samples and molecules. It might not be complete but represents 90% of the core products.

## LC/MS-GC/MS certified KITs & HPLC/GC certified KITs



**16475359**

HPLC/GC Certified Vial Kit: 1.5ml Short Thread Vial, clear glass, 1st hydrolytic class, label; UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1.0mm

**13439748**

HPLC/GC Certified Vial Kit: 1.5ml Short Thread Vial, clear glass, 1st hydrolytic class, label; UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1.0mm

**15562320**

LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, clear glass, wide opening, with overwind-barrier; Ultra High Performance Seal: PP Short Thread Cap, blue, centre hole; Silicone darkblue-translucent/PTFE natural, 35° shore A, 1,0mm

**15572320**

LC/MS and GC/MS Certified Vial Kit: 1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, amber glass, wide opening, with overwind-barrier; Ultra High Performance Seal: PP Short Thread Cap, blue, centre hole; Silicone darkblue-translucent/PTFE natural, 35° shore A, 1,0mm

## Short Thread Vials ND9



**11575884**

1.5ml Short Thread Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, label and filling lines



**10080952**

1.5ml Short Thread Vial, 32 x 11.6mm, amber glass, 1st hydrolytic class, wide opening, label and filling lines



**16388357**

1.5ml Short Thread SureStop™ Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, with overwind-barrier

## Microsampling Vials ND9



**12951011**

Short Thread Vial with integrated Micro-Insert, 32 x 11.6mm, clear glass, 1st hydrolytic class, Base bonded



**15592310**

Short Thread Vial with integrated Micro-Insert, 32 x 11.6mm, amber glass, 1st hydrolytic class, Base bonded



**11515894**

1.1ml Microliter Short Thread Vial ND9, 32 x 11.6mm, clear glass, 1st hydrolytic class

## Plastic Vials ND9



**11727597**

1.5ml PP Short Thread Vial, 32 x 11.6mm, transparent, with filling lines



**11824940**

0.7ml PP Short Thread Micro-Vial, 32 x 11.6mm, transparent



**11707597**

0.7ml PP Short Thread Micro-Vial, 32 x 11.6mm, transparent

## PP Short Thread Seals ND9



**11863021**

9mm Combination Seal: PP Short Thread Cap, blue, with centre hole; RedRubber / PTFE beige, 45° shore A, 1,0mm



**11787567**

UltraClean Closure: 9mm PP Short Thread Cap, blue, centre hole; Silicone white/PTFE red, 55° shore A, 1.0mm



**11797567**

9mm Combination Seal: PP Short Thread Cap, blue, centre hole; Silicone white/PTFE blue, 55° shore A, 1.0mm, slit

## Snap Ring Vials and Seals ND11



**11535894**

1.5ml Snap Ring Vial, 32 x 11.6mm, clear glass, 1st hydrolytic class, wide opening, label and filling lines



**11545894**

1.5ml Snap Ring Vial, 32 x 11.6mm, amber glass, 1st hydrolytic class, wide opening, label and filling lines

## PE Snap Ring Seals ND11



**16308377**

11mm Combination Seal: PE-Snap Ring Cap, blue, with centre hole; RedRubber/PTFE beige, 45° shore A, 1.0mm



**11814910**

UltraClean Closure: 11mm PE Snap Ring Cap, blue, centre hole; Silicone white/PTFE red, 45° shore A, 1.3mm



**11834910**

11mm Combination Seal: PE Snap Ring Cap, blue, centre hole; Silicone white/PTFE blue, 55° shore A, 1.0mm, cross-slitted

## Well Plates and Sealmats



**16243151**


Sealmat, blue, Silicone/PTFE, for 96 position Deep Well Microplate, round well, flat base, 7mm diameter (non sterile)



**16223151**

Square Well Microplate, PP, 96 Positions, certified, height 44.4mm, V-shape, 7mm dia., 2000µl total volume (non coated, non sterile)

## PP Storage Box and Vial Rack



**12692495**

PP Storage Box for 1,5ml vials or 2ml shell vials, blue, with cover (130 x 130 x 45mm), 81 cavities with alphanumeric coding of all 4 margins as well as the cavities at the bottom



**11767746**

PP Vial-Rack, (200 x 105 x 17mm), for 1.5ml vials, 50 cavities, blue, stackable

Find out more at [eu.fishersci.com/go/fisherbrandchromatography-vials-closures](https://eu.fishersci.com/go/fisherbrandchromatography-vials-closures)



Distributed by Fisher Scientific. Contact us today:

**Austria:** +43(0)800-20 88 40 **Belgium:** +32 (0)56 260 260 **Denmark:** +45 70 27 99 20  
**Germany:** +49 (0)2304 9325 **Ireland:** +353 (0)1 885 5854 **Italy:** +39 02 950 59 478  
**Finland:** +358 (0)9 8027 6280 **France:** +33 (0)3 88 67 14 14 **Netherlands:** +31 (0)20 487 70 00  
**Norway:** +47 22 95 59 59 **Portugal:** +351 21 425 33 50 **Spain:** +34 902 239 303  
**Sweden:** +46 31 352 32 00 **Switzerland:** +41 (0)56 618 41 11 **UK:** +44 (0)1509 555 500

© 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representatives for details.