



SPOTLIGHT

Unlock your Life Science Research & Testing solutions

MICROBIOLOGY, CELL CULTURE, GENOMICS, PROTEOMICS,
DRUG DISCOVERY, TOXICOLOGY, BIOPRODUCTION



INDUSTRY



Genomics p.6 — Microbiology p.44 — PPE p.56

Genomics p.6 — Toxicology p.27 — PPE p.56

Microbiology p.44 — PPE p.56

Genomics p.6 — Microbiology p.44 — PPE p.56

Microbiology p.44 — PPE p.56

Genomics p.6 — Toxicology p.27 — Microbiology p.44 — PPE p.56

Genomics p.6 — Proteomics p.17 — PPE p.56

Genomics p.6 — Microbiology p.44 — PPE p.56

Cell culture p.30 — PPE p.56

Drug discovery p.14 — Proteomics p.17 — Cell culture p.30 — PPE p.56

Bioproduction p.40 — PPE p.56

Proteomics p.17 — Cell culture p.30 — Microbiology p.44 — PPE p.56

SPOTLIGHT

Index

GENOMICS

6-13

Agarose gels.....	6
Assay tanks.....	12
Balances	7
Baths.....	11
Blotters.....	12
Containers	6
DNA amplification kits.....	9
Electrophoresis gel tanks.....	11
Microtubes	10
Non indicating FTA cards.....	6
PCR beads.....	9
PCR kits.....	8
PCR plates and sealing film.....	10
Pipettors.....	10
Power supplies	10
Purification kits	9
Sample preparation systems.....	6
Shakers	13
Transilluminators	13
UV sterilisation cabinets	8,13
Vacuum concentrators	7
Water purification systems.....	8

Filters.....	17
Filtration systems.....	18
Imaging systems	23
Magnetic beads.....	19
Microplates	24
Power supplies	21
Speciality application particles	25
Spectrophotometers	26
Syringe filters.....	18
Transfer stacks.....	22
Transilluminators	23
Water purification systems.....	17
Western blotting accessories.....	22

TOXICOLOGY

27-29

Burettes.....	29
Funnel stands.....	29
Pipettors and tips	27
Samplers	28
Serological pipettes.....	28
Stirrers and stirring hotplates	29
Transfer pipettes	28

DRUG DISCOVERY

14-16

Baths.....	16
Bottles	14-15
Microplates and readers	14
Pipettors, tips and kits	15-16

PROTEOMICS

17-26

Antibodies.....	26
Cell assay kits.....	25
Cell stains.....	25
Centrifugal concentrators.....	18
Chemiluminescent substrates	22
Chillers	22
Columns	17
Concentrators	21
Electrophoresis units	21

CELL CULTURE

30-38

Aspiration systems.....	33
Baths.....	32,34
Cell counters.....	36
Cell culture dishes and flasks	32
Cell culture media and reagents...	32-33
Centrifuge tubes.....	31
Centrifuges.....	31
Cryogenic storage systems	37
Detection kits and reagents	36
Filter units.....	30
Freezer racks	38
Freezers, ultra low temperature...	38
Imaging systems	36, 37
Incubators	34
Media bottles	31
Microscopes and microscopy reagents.....	35

SPOTLIGHT

Index

Pipette tips.....	35	Cutlery	45
Serological pipettes.....	33	Cuvettes	52
Shakers	34	Dry baths.....	49
Sonic dismembrators.....	30	Incubators	46-47
Spatulas	30	Labels.....	52
Streptavidin conjugates	36	Membranes and filters.....	50
Thermometers	37	Microcentrifuges.....	44
Trays	31	Mixers, vortex	44,50
Wash bottles.....	30	Petri dishes.....	48
		Plate readers.....	54
		Racks.....	48
BIOPRODUCTION	40-43	Refrigerators and freezers.....	54,58
Biological safety cabinets.....	41	Rotator	45-46
Bottles	40	Shakers	47,51
Carboys	40	Spectrophotometers	52
Centrifuges.....	41	Tape dispensers	52
Chillers	42	Thermometers	53
Circulators	42-43	Transfer membranes.....	50
Incubators	41	Trays	48
Jerrycans.....	40	Tubes	44-45
pH meters.....	41-42		
Pipettors.....	43	PERSONAL	
Pumps, solvent.....	40	PROTECTIVE	
Tubes	43	EQUIPMENT	56-61
Vials and vial storage	43		
		Coats	56
MICROBIOLOGY	44-54	Coveralls.....	56
Autoclaves.....	54	Face masks.....	57
Autoclave pads	54	Funnels, safety waste.....	60
Bag holders	52	Gloves	58-59
Bags	54	Goggles.....	56
Balances	45	Laminar flow cabinets.....	60
Baths.....	53	Respirators.....	57-58
Blenders and diluters.....	46	Safety storage cabinets.....	61
Bottles, centrifuge	47	Shoe covers	60
Cell spreaders.....	48	Spectacles	56
Cell strainers.....	45	Surface protectors	60
Centrifuges.....	51		
Containers	44		
Control organisms	51		
Cooling racks and boxes	49		
Culture media.....	49		



FASTPREP-24™ 5G SAMPLE PREPARATION SYSTEM

- The most advanced sample preparation system available!
- Delivering the most DNA, RNA and proteins from the most resistant samples in less than 40 seconds
- Intuitive software, user-friendly touch screen
- Highest speed and most power
- Shortest processing
- Widest range of adapters for sample flexibility
- Highest yield of DNA, RNA and proteins
- Process any tissue or cellular sample type with highest quality results
- The kits include the lysing matrix, binding solutions, washing buffer, elution buffer, spin filters and catch tubes



FEATURED PRODUCT

Cat No.	Description	Tubes	Pack qty	List price
15260488	FastPrep-24 5G instrument		each	94210,56
11482370	QuickPrep™	2 x 24mL	each	7313,77
11472380	BigPrep™	2 x 50mL	each	12408,31
10429774	TeenPrep™	12 x 15mL	each	12408,31
11490576	HiPrep™	2 x 48mL	each	12408,31
11471525	CoolBigPrep™	2 x 50mL	each	15302,95
11461525	CoolTeenPrep™	12 x 15mL	each	15302,95
11482360	CoolPrep™ Cryogenic	24 x 2mL	each	15302,95

Cat No.	Description	Extraction type	Sample type	Pack qty	List price
---------	-------------	-----------------	-------------	----------	------------

Lysing matrix

10591758	Lysing matrix A	DNA/Protein	Animals, bacteria, yeast, fungi, plants	100pcs	4350,00
11462410	Lysing matrix B	DNA/Protein	Bacteria, spores	100pcs	7231,80
11422420	Lysing matrix D	DNA/RNA/Protein	Animal or plant tissues	100pcs	3056,94
11462420	Lysing matrix E	DNA/RNA	Soil, sediments, water, faeces	100pcs	3039,60
11401555	Lysing matrix S	RNA/Protein	Tough Tissues Skeletal Muscle Lung Heart Bone Seeds Spores	100pcs	3039,60

Kits (extracted DNA or RNA is suitable for PCR / RT-PCR analysis and other downstream applications)

11442400	FastDNA™ SPIN Kit	DNA	Plants, animals, bacteria, yeast, algae, fungi	100prep	4139,32
11492400	FastDNA™ SPIN Kit	DNA	Plants, animals, bacteria, fungi	50prep	9029,04
11481535	FastDNA™ SPIN Kit	DNA	Fresh or frozen human and animal	1pcs	3725,04
11482380	FastRNA™ Pro Blue Kit	RNA	Gram positive and gram negative bacteria	50prep	5622,00
11492380	FastRNA™ Pro Red Kit	RNA	Yeast, fungi	50prep	5627,00
11402390	FastRNA™ Pro Green Kit	RNA	Plant, animal tissue or cultured cells	50prep	2374,56
12703057	FastRNA™ Pro Soil-Direct Kit	RNA	Soil	50prep	4911,22

thermo scientific

CERTIFIED UNIVERSAL CONTAINERS



- Manufactured under controlled cleanroom conditions, this clear robust polypropylene container has a wealth of features which make it a laboratory essential
- With a working volume of 25mL and convenient self-standing conical base, it's perfect for a wide range of general laboratory uses from sample handling to centrifugation applications
- Free from RNase, DNase, human DNA and pyrogens
- ½ turn Quickstart cap for ease of handling
- Excellent leak free performance – tested in accordance with BS EN 14254 Annexe D

SPECIAL PRICE

Cat No.	Label	Volume, mL	Pack qty	List price	Promo price
15663934	Plain	30	400pcs	1312,11	1075,93
15673934	Yes, graduated & marker spot	30	400pcs	1367,16	1121,07



Distributor
GE Healthcare

WHATMAN™ NON-INDICATING FTA CLASSIC CARDS



- Collect, stabilise, process, transport, and archive blood samples
- Improves process efficiency by enabling direct amplification for STR analysis
- Ideal for coloured samples, like blood
- Immediate stabilisation and protection of nucleic acids due to FTA technology
- Different card options to support manual or automated processing for fast and efficient sampling
- Long term archiving of samples at room temperature due to FTA card technology

SPECIAL PRICE

Cat No.	No. of Samples	Sample Size, µL	Volume (Metric) Sample, µL	Pack qty	List price	Promo price
11310854	4	125	500	25pcs	1187,98	1129,00
10223352	4	125	500	100pcs	4329,33	4116,00



Distributor
GE Healthcare

ILLUSTRA™ SEPHADEX™ G-50 AND G-100 DNA GRADE AGAROSE GELS



- Used for purification of DNA from small molecules. GE Healthcare illustra™ Sephadex™ G-50 and G-100 DNA grade agarose gels purify DNA by gel filtration
- Tested to ensure reproducibly high recovery of DNA
- Suitable for researchers who prefer to prepare their own columns for purifying nucleic acids

FEATURED PRODUCT

Cat No.	Particle Size	pH Range	Bed Volume	Pack qty	List price
11584925	50 to 150µm	2 to 10	9 to 11mL/g in distilled water	100gr	2839,00



BALANCES, PRECISION, EXPLORER™ SERIES

- Ohaus™ Explorer™ has advanced applications to simplify even the most complex laboratory measurements. Whether it's determining the difference between initial and residual weights or calculating the density of solids and liquids, the Explorer™ eliminates the need for time consuming manual calculations and data logging. The high resolution display and innovative user interface make balance setup and application use effortless.
- High resolution multi-lingual colour VGA touch display with icons for simple operation
- Programmable multi-functional touch free sensors for cleanliness, speed and safety
- Multiple application modes plus library for storage of personal application settings



SPECIAL PRICE

Cat No.	Capacity, g	Readability, mg	Pan diam., mm	Pack qty	List price	Promo price
15436101	120/220	0,01/0,1	225	1unit	63767,00	54074,40
15456101	220	0,01	225	1unit	100956,00	85610,40
15416101	120/200	0,01/0,1	225	1unit	60949,00	51684,80
15446101	225	0,01	225	1unit	98138,00	83220,80



PURIFICATION SYSTEMS, KINGFISHER™

- Prepare samples efficiently from a variety of materials with Thermo Scientific™ KingFisher™ Flex purification systems, which offer highly versatile, automated magnetic-particle processing for DNA/RNA, protein or cell purification from virtually any source. Using revolutionary magnetic particle separation technology, this system provides excellent reproducibility and quality



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
11849520	KingFisher Flex with 96 PCR head	1pcs	454045,00
11839520	KingFisher Flex with 96 KF head	1pcs	454045,00
11829520	KingFisher Flex with 96 deep well head	1pcs	436408,01
11819520	KingFisher Flex with 24 deep well head	1pcs	454045,00
15337317	KingFisher Duo Prime	1pcs	210097,15



DNA SPEEDVAC™ CONCENTRATOR KITS

- Save time and money with compact Thermo Scientific™ Savant™ DNA SpeedVac™ Concentrator Kits.
- These complete systems provide everything you need to quickly dry and concentrate small sample volumes in one fully integrated package, sold under one part number.
- Consistent performance: designed to ensure consistent drying times and reproducible results
- Simplicity: single program supports all DNA/RNA applications
- User-friendly interface: with temperature/time settings
- Three drying rates: ambient, 43°C, 65°C
- Dual digital timers: independent control of heat/run times



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
10068669	DNA SpeedVac™ Concentrator Kits	1pcs	57412,00	38925,34



CENTRIVAP DNA VACUUM CONCENTRATOR

- Designed to speed evaporation of solvents from DNA samples and other very small samples
- Compact benchtop cabinet houses the concentrator, DNA Rotor, diaphragm vacuum pump, two glass traps and all connecting tubing to begin processing
- Quick-Start™ one button start up saves time. Pressing just one button starts the rotor, the heater, the timers and the vacuum pump. Three separate Quick-Start buttons store one user-set program each
- Quick-Stop™ system. Once run time has expired or the RUN/STOP button is pressed, within seconds the vacuum pump stops, the valve opens to bleed air into the chamber, and the rotor stops
- Compliance : UL 3101-1/61010-1



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12347223	DNA concentrator 230V EU	1pcs	88995,00

Grant

UV CABINET/PCR WORKSTATION

- DNA/RNA UV-cleaner box UVC/T-M-AR is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.
- Ozone free high density UV decontamination
- Long living UV lamps (9000 hours average)
- Automatic switch off of UV-lamps when the protective screen is opened
- Bactericidal flow-type recirculator providing permanent decontamination inside UV - cleaner box during operation
- Shockproof glass walls



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
11441048	UV cabinet/PCR workstation, s. steel work surface, glass sides & front, dual UV protection lamp and UV recirculator	1pcs	22964,23

thermo scientific

WATER PURIFICATION SYSTEMS SMART2PURE™

- Tap to Type 1 water purification system
- 3 or 6L/hr capacities
- Easy to replace cartridges
- Variable flow rate
- Tilting display
- Type II water produced in tank
- Includes a further RO membrane with integrated pre-treatment, ultrapure polisher cartridge and sterile 0,2µm filter for free!



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
11337021	Smart2Pure™ 3L/hr, with UV oxidation, ultrafiltration membrane and 6L integrated tank	1pcs	48139,14
13237449	Smart2Pure™ 6L/hr, with UV oxidation, ultrafiltration membrane and 6L integrated tank	1pcs	57859,00

thermo scientific

DREAMTAQ™ POLYMERASE

- Robust amplification with minimal optimisation
- High yields of PCR products
- Higher sensitivity compared to conventional Taq DNA polymerase
- Amplification of long targets up to 6 kb from Genomics DNA and up to 20 kb from viral DNA
- Incorporates modified nucleotides, but does not incorporate dUTP



invitrogen

PLATINUM™ SUPERFI™ DNA POLYMERASE

- Exceptional >100x Taq fidelity
- High specificity and increased yields with Platinum hot-start technology
- Robust amplification of difficult-to-amplify targets including those of sub-optimal purity or with 65% GC content
- Convenient workflow with room temperature reaction setup and direct gel loading with green buffer formats



SPECIAL PRICE

Cat No.	Format	Pack qty	List price	Promo price
15525033	Platinum SuperFi DNA Polymerase	100unit	1820,00	1310,40
15535033	Platinum SuperFi DNA Polymerase	500unit	8180,00	5889,60
15545033	Platinum SuperFi DNA Polymerase	2500unit	34760,00	25027,20
15585033	Platinum SuperFi Green DNA Polymerase	100unit	1820,00	1310,40
15595033	Platinum SuperFi Green DNA Polymerase	500unit	8180,00	5889,60
15505043	Platinum SuperFi Green DNA Polymerase	2500unit	34760,00	25027,20
15525043	Platinum SuperFi PCR Master Mix	100rxn	2365,00	1702,80
15535043	Platinum SuperFi PCR Master Mix	500rxn	10640,00	7660,80
15545043	Platinum SuperFi PCR Master Mix	2500rxn	45180,00	32529,60
15565043	Platinum SuperFi Green PCR Master Mix	100rxn	2365,00	1702,80
15575043	Platinum SuperFi Green PCR Master Mix	500rxn	10640,00	7660,80
15585043	Platinum SuperFi Green PCR Master Mix	2500rxn	46920,00	33782,40
15555033	SuperFi Buffer	5mL	612,00	440,64
15565033	SuperFi Green Buffer	5mL	614,00	442,08

SPECIAL PRICE

Cat No.	Format	HotStart	Green	Pack qty	List price	Promo price
10559120	Standalone			200unit	345,00	265,65
10160790	Standalone			500unit	654,00	503,58
11883813	Standalone			5 x 500unit	2700,00	2079,00
15282086	Standalone			5000unit	4585,00	3530,45
10679320	Standalone			20 x 500unit	8580,00	6606,60
10689320	Standalone		x	200unit	344,00	264,88
10776411	Standalone		x	500unit	654,00	503,58
11803823	Standalone		x	5 x 500unit	2700,00	2079,00
10262450	Standalone		x	20 x 500unit	8570,00	6598,90
11896833	Master Mix			200rxn	836,00	643,72
11806843	Master Mix			1000rxn	3330,00	2564,10
11813933	Master Mix		x	200rxn	836,00	643,72
11816843	Master Mix		x	1000rxn	3320,00	2556,40
15699364	Standalone	x		200unit	1020,00	785,40
15609374	Standalone	x		500units	2300,00	1771,00
15619374	Standalone			2500units	9640,00	7422,80
15636154	Standalone	x		4 x 2500unit	36060,00	27766,20
15629374	Standalone	x	x	200unit	1020,00	785,40
15639374	Standalone	x	x	500unit	2300,00	1771,00
15649374	Standalone	x	x	2500unit	9640,00	7422,80
15659374	Standalone	x	x	4 x 2500unit	36060,00	27766,20
15669374	Master Mix	x		200rxn	1482,00	1141,14
15679374	Master Mix	x		1000rxn	6600,00	5082,00
15689374	Master Mix	x	x	200rxn	1482,00	1141,14
15699374	Master Mix	x	x	1000rxn	6600,00	5082,00



Distributor
GE Healthcare

ILLUSTRATM GENOMIPHI V2 DNA AMPLIFICATION KIT

- Starting material can be purified DNA or nonpurified cell lysates
- Simple protocol provides typical yields of 4 to 7µg DNA and average product lengths >10kb
- Less hands-on time compared to traditional isolation methods
- Representative amplification of the whole genome
- Automation-friendly protocol with great reproducibility
- High-quality DNA for array CGH, high-throughput genotyping, DNA archival, PCR, restriction enzyme digestion, hybridisation, and cloning
- Contains all of the components necessary for Genomics DNA preparation by isothermal strand displacement amplification



FEATURED PRODUCT

Cat No.	Includes	Pack qty	List price
10309494	5 x 0,9mL Sample Buffer, 5 x 0,9mL Reaction Buffer, 5 x 100µL Enzyme mix, 20µL Control DNA (lambda)	500rxn	20727,00



Distributor
GE Healthcare

ILLUSTRATM PURETAQ READY-TO-GO™ PCR BEADS

- PuReTaq DNA polymerase and high-purity reagents provide robust performance with very low levels of contaminating DNA
- Manufacturing process locks in performance for excellent reproducibility with less risk of pipetting errors and contamination
- Long-term ambient-temperature stability simplifies shipping and removes any need for freezer storage
- Verified for real-time PCR
- Ready-To-Go means simply add template DNA and primers and cycle



FEATURED PRODUCT

Cat No.	Includes	Format	Pack qty	List price
10678095	Room temperature stable beads containing ~ 2,5 units of recombinant PuReTaq DNA Polymerase, dATP, dCTP, dGTP, dTTP, stabilisers, BSA and reaction Buffer. When a bead is reconstituted to a 25µL final volume, the concentration of each dNTP	0,2mL Tube with hinged cap	1pcs	1284,00



Distributor
GE Healthcare

ILLUSTRATM READY-TO-GO™ RT-PCR BEADS

- Reagents are optimised for full-length cDNA synthesis (>7,5 kb) and optimal sensitivity from PCR
- Provided in either thin walled 0,5 ml or 0,2 ml tubes/plates compatible with most thermocyclers
- Supplied with both pd(N)6 and oligo(dT) primers for flexibility
- The Ready-To-Go manufacturing process locks in performance for excellent reproducibility with less risk of pipetting errors and contamination
- Long-term ambient-temperature stability simplifies shipping and removes any need for freezer storage



FEATURED PRODUCT

Cat No.	Includes	Format	Pack qty	List price
10187044	Room temperature stable beads containing Buffer, dATP, dCTP, dGTP, dTTP, approximately 2 units of Taq DNA Polymerase, M-MLV reverse Transcriptase, RNAGuard (porcine) and stabilisers (including RNase/DNase free BSA), when a bead is reconstituted to 50	0,2mL Tube with hinged cap	1pcs	3271,00



Distributor
GE Healthcare

ILLUSTRATM GFX™ PCR DNA AND GEL BAND PURIFICATION KIT

- Typical recoveries range from 60 to 80% for DNA fragments from agarose gel to as high as 95% for PCR products from solution
- Rapid and easy-to-use purification method with less than 10 minutes hands-on time
- Flexible 10 to 50µL elution volume for different DNA concentration needs
- Improved column design with increased maximum volume holding of up to 800µL liquid volume (950µL without a cap)
- DNA of exceptional purity (to 99,5%) is ready for direct use in sequencing, PCR, labelling, restriction enzyme digestion, and cloning
- Visual colour indicator added to the capture buffer to ensure optimal pH for maximum DNA binding and recovery
- Kit combines a versatile chaotropic buffer with a glass fibre matrix supported in a spin-column for the purification of DNA from both solution and agarose gel



FEATURED PRODUCT

Cat No.	Includes	Formulation	Pack qty	List price
10218134	GFX™ columns, collection tubes, colour-coded bottles of capture Buffer, Wash Buffer, two Elution Buffers	Tris HCl and sterile Water	250rxn	4098,00



BRAND

PCR COMBO PACK



- PCR-96 plates
- Including PCR-sealing film

FEATURED PRODUCT

Cat No.	Skirted	Profile	Pack qty	List price
15594409	Non-skirted	Low profile plates	1set of 50	1036,84
15504419	Non-skirted	Standard profile plates	1set of 50	1036,84
15514419	Semi-skirted	Elevated rim, LP plates	1set of 50	1036,84
15524419	Semi-skirted	Low profile plates	1set of 50	1036,84
15534419	Semi-skirted	Standard profile plates	1set of 50	1036,84



PIPETMAN M ELECTRONIC PIPETTERS



- New, guaranteed performance in standard and repetitive pipetting
- Extended volume range, now available in 20 models including 8 single-channel (from 0,5µL to 10mL) and 12 multi-channel (from 0,5µL to 1200µL)
- Complete pipetting operation with only two buttons
- Intuitive interface with five pipetting modes: Pipette, Repetitive, Mix, Reverse, and new Custom
- Ideal for many applications: qPCR, ELISA, NGS, cell culture, and more
- The electronic pipettor with new features to help you achieve exceptional performance, reproducible results and data traceability, quickly and easily
- With a.c. adapter

FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
8-channel			
15560500	P8x10M, 0,5 to 10µL	1pcs	7851,02
15580500	P8x20M, 1 to 20µL	1pcs	7851,02
15590490	P8x100M, 5 to 100µL	1pcs	7851,02
15500510	P8x200M, 5 to 200µL	1pcs	7851,02
15520510	P8x300M, 10 to 300µL	1pcs	7851,02
15510500	P8x1200M, 50 to 1200µL	1pcs	7851,02
12-channel			
15570500	P12x10M, 0,5 to 10µL	1pcs	8761,96
15590500	P12x20M, 1 to 20µL	1pcs	8761,96
15500500	P12x100M, 5 to 100µL	1pcs	8761,96
15510510	P12x200M, 5 to 200µL	1pcs	8761,96
15530510	P12x300M, 10 to 300µL	1pcs	8761,96
15520500	P12x1200M, 50 to 1200µL	1pcs	8761,96



MICROPLATE SEALING FILM AND TAPES



- For applications ranging from PCR and Real-Time PCR to ELISA and cell culture, Axygen films are an easy and cost-effective way to seal plates and enhance automation. Used to seal multi-well microplates. Sealing film temperature ranges from -40°C to +104°C
- Provides exceptional seal on multi-well plates with varied temperatures and reagents
- With compression mats, films and tapes seal during thermocycling
- Ideal for light-sensitive samples and its uniformly applied adhesive provides optimum well sealing
- Pierceable, PCR-AS-200

GREAT SAVINGS

Cat No.	Sterile	Temperature Range, °C	Thickness, µm	Pack qty	List price	Promo price
12577947	No	-80 to 104	35	100pcs	922,31	727,00
11360215	No	-80 to 104	50	100pcs	799,66	631,00
11304325	No	-80 to 104	70	100pcs	2487,01	1960,00



MICROTUBES, SNAPLOCK, MAXYCLEAR™



- Certified DNase, RNase and endotoxin free
- Polypropylene
- Colour: clear

GREAT SAVINGS

Cat No.	Volume, mL	Max. rcf, xg	Sterile	Pack qty	List price	Promo price
15227179	5	14000	Yes	250pcs	454,78	336,00
15237179	5	14000	No	250pcs	500,25	370,00



POWER SUPPLY, MIDI 300V/4



- Stackable design
- In-built safety systems
- Dual voltage compatibility
- Current, mA: 1 to 700
- Power (max.), W: 150
- Dimensions, mm: 190 x 305 x 95

FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12643546	Midi 300V/4	1pcs	3411,80



SUB-GEL TEACHING GEL TANKS

- Two new teaching tanks:
 - SUB-GEL 4X Mini featuring four 70mm x 70mm gel trays, maximum samples 64
 - SUB-GEL 6X Mini featuring six 70mm x 70mm gel trays, maximum samples 96
- Saves bench space, time and cost; teachers are able to demonstrate horizontal gel electrophoresis to many students using just a single gel tank and power supply
- Features all of the benefits of the wider SUB-GEL range including
 - High quality injection moulded construction
 - Leak free 'Plug-and-Go' casting dams
 - UV transparent gel trays
 - Height adjustable combs
 - Economic low gel and buffer volumes
 - Cassette-style electrodes – difficult to break, but inexpensive and easy to change
- Each tank is available with or without power supply
- Complete classroom workstation available, features four SUBGEL 6X Mini units together with power supply, to demonstrate electrophoresis for up to 24 students at a time



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15341357	SUB-GEL 4X Mini, Teaching Gel Tank, includes four 70mm x 70mm UV trays, eight 8 sample combs, eight casting dams, loading guides and cables	1pcs	3247,00
15351357	SUB-GEL 4X Mini, Teaching Gel Tank, includes four 70mm x 70mm UV trays, eight 8 sample combs, eight casting dams, loading guides and cables plus power supply. (Cat. No 12643546)	1pcs	6493,00
15321357	SUB-GEL 6X Mini, Teaching Gel Tank, includes six 70mm x 70mm UV trays, twelve 8 sample combs, twelve casting dams, loading guides and cables	1pcs	4329,00
15331357	SUB-GEL 6X Mini, Teaching Gel Tank, includes six 70mm x 70mm UV trays, twelve 8 sample combs, twelve casting dams, loading guides and cables plus power supply (Cat. No 12643546)	1pcs	7575,00
15361357	SUB-GEL 6X Mini classroom workstation for 24 students, includes four SUB-GEL 6X Mini units (Cat. No 15321357) and power supply (Cat. No 12613546)	1pcs	22182,00



BATHS, CIRCULATORS, REFRIGERATED/HEATED

- For laboratory and industrial applications
- The robust refrigeration system ensures that your samples and applications can be cooled quickly
- Stainless steel tanks for durability and ease of cleaning
- Fine tuned refrigeration metering for excellent temperature stability
- Bath drains are located on the front of the unit
- Capacity, L: 5,4 to 6,5
- Temperature range, °C: -20 to +100
- Stability, °C: ±0,05
- Heater, kW: 1,2
- Cooling capacity, W (at 20°C): 250



GREAT SAVINGS

Cat No.	Description	Overall dimensions (H x W x L), mm	Pack qty	List price	Promo price
11366641	Refrigerated/heated bath circulator Isotemp 4100 R20	622,3 x 203,2 x 416,6	1pcs	19135,00	15013,32



MINI HORIZONTAL GEL SYSTEM

- Designed for quick separation of DNA or nucleic acids in agarose mini-gels
- The Fisher Scientific Mini Horizontal Gel System is a complete, compact mini gel electrophoresis system for easy nucleic acid separation. The horizontal system comes with an attached power supply, two casting stands, six gel trays and four combs — everything necessary for DNA or nucleic acid horizontal electrophoresis
- Quick and easy gel casting - no spacers or clamps required
- Compact design - 19 x 12,7cm footprint
- Simple on-off switch
- Membrane keypad with selectable parameters including 3 preset voltages and a 99 minute timer
- Vented safety lid helps dissipate heat generated during electrophoresis
- Power supply connects to gel tank without cables or leads



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15530340	Mini Horizontal Gel System	1pcs	4020,00	3536,19



SUB-GEL ELECTROPHORESIS UNITS, HORIZONTAL

- The Fisherbrand SUB-GEL range features Mini, Midi, Midi-Plus and Maxi units. The relatively compact size of each unit results in economical buffer and gel consumption without compromising resolution and separation speed. All units are provided with loading guides and dams



FEATURED PRODUCT

SUB-GEL	Mini	Midi	Midi-Plus	Maxi
UV trays, mm	70x 70 100 x 70	70 x 70 100 x 70	70 x 15 100 x 150 150 x 150	200 x 100 200 x 200
Combs (1mm)	2 x 8 well sample	2 x 16 well sample	2 x 20 well sample	2 x 20 well sample
Cat. No.	11863303	11853303	11833293	11843303
List price	3525,12	3469,00	3922,00	6225,00



HORIZONTAL GEL TANK

- Includes: 1 x gel casting tray with gates, 2 x 1,0mm 20 sample combs, power supply connectors and loading strips
- Buffer recirculation ports are included into the Midi Plus Wide unit



FEATURED PRODUCT

Cat No.	Description	Max. number of samples	Gel dimensions, (L x W), mm	Buffer volume, mL	Pack qty	List price
11553352	Mini Plus Wide	80	102 x 144	500	1pcs	3875,00
11563382	Midi Plus Wide	200	140 x 230	800	1pcs	5989,00



RUNVIEW™ MINI AND MIDI SYSTEMS

- Three types of runVIEW™ packages are available:
 - runVIEW™ CONVERTER package includes runVIEW™ adjustable viewing platform and bluVIEW Mini or Midi emission filter lid (allows standard SUB-GEL Mini or SUB-GEL Midi tanks to be converted to real-time electrophoresis)
 - runVIEW™ STANDARD package includes runVIEW™ adjustable viewing platform and runVIEW™ SUB-GEL Mini or Midi tank
 - runVIEW™ COMPLETE package includes STANDARD package in Mini or Midi format with power supply
- The runVIEW™ portable adjustable blue light illuminator viewing platform can be used with SUB-GEL Mini and SUB-GEL Midi electrophoresis tanks
- bluVIEW emission filter lid with built-in extractor fan enables condensation free viewing of blue light fluorescently stained 70mm or 100mm wide agarose gels
- Blue light is completely safe and results in improved cloning efficiency compared to UV
- Typical applications - Education - ideal for quick checks of low and medium numbers of samples following PCR and cloning
- runVIEW™ is an innovative new system designed for real-time size fractionation and recovery of nucleic acids. runVIEW™ can maximise the efficiency of DNA recovery from EtBr and SYBR stained gels by minimising the number of steps involved in postelectrophoretic purification



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15391357	c. CONVERTER package includes RunVIEW™ viewing platform and bluVIEW lid for use with customers own SUB-GEL Mini unit	1pcs	6696,00
15301367	RunVIEW™ CONVERTER package includes RunVIEW™ viewing platform and bluVIEW lid for use with customers own SUB-GEL Midi unit	1pcs	6769,00
15321367	RunVIEW™ STANDARD package includes RunVIEW™ viewing platform, plus runVIEW™ SUB-GEL Mini System	1pcs	9763,00
15331367	RunVIEW™ STANDARD package includes RunVIEW™ viewing platform, plus runVIEW™ SUB-GEL Midi System	1pcs	10007,00
15341367	RunVIEW™ COMPLETE includes RunVIEW™ viewing platform, plus RunVIEW™ SUB-GEL Mini System plus Mini 300V Plus power supply	1pcs	12564,00
15351367	RunVIEW™ COMPLETE includes RunVIEW™ viewing platform, plus RunVIEW™ SUB-GEL Midi System plus Mini 300V Plus power supply	1pcs	12806,00



COMET ASSAY TANKS

- The Comet Assay or Single Cell Gel Electrophoresis (SCGE) was first introduced to quantify double stranded DNA breakages in single cells exposed to irradiation. It has since been adapted to analyse specific DNA lesions and repair processes to analyse the potential for DNA damage of new drugs and environmental conditions.
- For high throughput Single Cell Gel Electrophoresis (SCGE)
- Unique patent pending design employs two carriers, each capable of accommodating 25 slides in a vertical laminar orientation
- Slide carriers facilitate batch processing of multiple slides simultaneously, eliminating the need for manual handling of individual slides
- Ebony acrylic construction of tank and staining dishes minimises exposure to light and reduces background DNA damage
- Enhanced cooling is provided by a ceramic cooling base with a sliding drawer designed to accommodate pre-frozen cool packs in order to maintain optimal buffer temperature
- Typical applications:
 - Quantifying DNA damage and repair in drug development



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15381347	High throughput Comet system for 50 slides includes two 25 slide carriers, four staining dishes, tank with ceramic cooling platform and cool pack, lid and power cables	1pcs	16112,00
15301357	High throughput Comet system for 50 slides includes two 25 slide carriers, four staining dishes, tank with ceramic cooling platform and cool pack, lid and power cables plus power supply	1pcs	21076,00

Accessories

Cat No.	Description	Pack qty	List price
15579901	Vertical slide carrier for 25 slides	1pcs	1137,00
15589901	Ebony acrylic stain dish	1pcs	632,00
15599901	Ebony acrylic stain dish	4pcs	1794,00
15569901	Positive electrode	1pcs	812,00
15559901	Negative electrode	1pcs	812,00



SEMI DRY BLOTTER

- Rapid transfer times
- Western, Southern and Northern blots
- Economic transfers due to very low buffer volumes
- Screw down lid
- Gels from 0,25 up to 10mm thick can be blotted
- Uniform heat dispersion
- Long life electrodes



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12357297	Semi dry blotter, Maxi, 100mm x 100mm	1pcs	4987,00
12367297	Semi dry blotter, Maxi, 200mm x 200mm	1pcs	10507,00



SAFEVIEW BLUE LIGHT TRANSILLUMINATOR

- Offers the user a safe way to view and document their samples
- This light source also has the added advantage in that it does not damage DNA or RNA that you would normally associate with UV light
- Supplied as a standalone unit and can be used with the microDOC, as part of a fully integrated gel documentation system
- No DNA damage to samples
- Safer for user - No UV light
- High purity LED light
- Strong metal enclosure with stainless steel filter frame
- Fast start up



FEATURED PRODUCT

Cat No.	Description	Filter size, mm	Pack qty	List price
15502506	Safeview BLUE Light Transilluminator	154 x 154	1pcs	3846,00



UV STERILISATION CABINETS

- Inactivates nucleic acids in 5 to 30 minutes
- Doors fitted with safety switches
- Complete with four powerful, timer controlled UV bulbs
- Efficient decontamination of the complete work surface
- Suitable for work with 32P
- White light provides excellent visibility when working within the cabinets
- For use with trays of size 68 x 54cm
- Fisherbrand UV sterilisation cabinets are a low cost alternative to a cleanroom. The cabinet incorporates safety features to prevent user exposure to UV light. The UV lights themselves are timer controlled and there are safety switches on the cabinet doors which turn off the UV lights when opened. The cabinet also acts as an efficient shield of beta radioactive emissions and can therefore be safely used with isotopes such as 32P
- The new mini UV cabinet is the perfect solution for laboratories constrained by space and budget. It has all the benefits of the standard UV sterilisation cabinet, however its reduced height and footprint allows it to fit snugly and conveniently into those small areas within the laboratory that would otherwise remain unused



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15562496	UV sterilisation cabinet with timer, four UV lights and white light, no tray	1pcs	14260,00
15592496	UV Sterilisation cabinet with timer, four UV lights and white light, includes tray	1pcs	14882,00
15572496	Mini UV sterilisation cabinet with timer, UV light and white light, no tray	1pcs	12469,00
15582496	Mini UV sterilisation cabinet with timer, UV light and white light, includes tray	1pcs	13092,00



GENOMICS PACKAGE: MAXQ SHAKERS

- Each shaker has package configurations optimised for 250mL or 1L flasks to maximise your research
- Incubated shaker packages are available in temperatures up to 80°C
- Open air shaker packages are available in a variety of sizes or capacities to fit your lab's needs
- **Thermo Scientific™ MaxQ™ Digital Shakers**
 - A digital shaker designed to provide reliable results required for bacteria and yeast cultures
 - Durable triple eccentric drive backed by 10 year parts warranty
 - Maximise capacity of 250mL and 2L flasks in a compact footprint
 - Digital display with simultaneous viewing of temperature, speed and time with audible and visual alarms
- **Thermo Scientific MaxQ Digital High Performance Shakers**
 - A digital high performance shaker designed for large or heavy loads to maximise your production
 - Increased weight capacity of 32,2kg for heavy loading
 - Doubled capacity of 250mL and 2L flasks for high throughput shaking
 - Options available for continuous or timed runs for added flexibility
 - Each package includes shaker, platform and clamps in one part number for easy ordering
 - Select a package ideally suited to protect your valuable samples by delivering outstanding performance, accuracy and reproducible results



FEATURED PRODUCT

Cat No.	Application Package	Shaker Model	Shaker Description	Accessories Description	Pack qty	List price
15689045	Small Open Air Package A	MaxQ 2000	Digital Shaker	457 x 457mm universal platform	1pcs	16323,35
				250mL flask clamps	16pcs	
15699045	Small Open Air Package B			457 x 457mm universal platform	1pcs	16323,35
				1L flask clamps	9pcs	
15629055	Mid-Sized Open Air Package A	MaxQ HP416	High Performance Digital Shaker	250/300mL dedicated platform with 35 pre-installed clamps	1pcs	30820,76
15639055	Mid-Sized Open Air Package B			1L dedicated platform with 12 pre-installed clamps	1pcs	30820,76
15669055	Large Open Air Package A	MaxQ HP430	High Performance Digital Shaker	250/300mL dedicated platform with 70 pre-installed clamps	1pcs	41161,63
15679045	Large Open Air Package			1L dedicated platform with 24 pre-installed clamps	1pcs	41161,63
15679055	Small Incubated Package A	MaxQ 4450	Digital Shaker	330 x 279mm universal platform	1pcs	33839,52
15689055	Small Incubated Package B			1L flask clamps	4pcs	
				330 x 279mm universal platform	1pcs	33839,52
				250mL flask clamps	9pcs	
15649055	Large Incubated Package A	MaxQ HP420	High Performance Digital Shaker	250/300mL dedicated platform with 25 pre-installed clamps	1pcs	54209,29
15659055	Large Incubated Package B			1L dedicated platform with 8 pre-installed clamps	1pcs	54209,29



NUNC™ 2.0 EDGE 96 WELL PLATES

- Prevent edge effect with a surrounding reservoir that serves as an evaporation barrier
- Achieve full-plate cell viability and consistency by reducing the rate of overall plate evaporation to <2
- Optimise cell adherence with Nunclon™ Delta cell culture treated surface
- Enabling the use of all 96 wells allows researchers to plan experiments more efficiently, resulting in potentially less cost and related product waste
- Efficient, economical and simple solution to prevent edge effect for full-plate results during cell-based assays



GREAT SAVINGS

Cat No.	Type	Sterile	Lid	Pack qty	List price	Promo price
15533115	Nunclon delta	Yes	With lid	50pcs	967,29	605,00
15543115	Nunclon delta	Yes	With lid	160pcs	2798,57	1758,00
15553115	Nunclon delta	Yes	Without lid	160pcs	2601,35	1635,00
15563115	Nunclon delta	Yes	Without lid	50pcs	829,30	520,00
15573115	Non-treated	Yes	With lid	50pcs	771,02	482,00
15583115	Non-treated	Yes	With lid	160pcs	2319,62	1455,00
15593115	Non-treated	Yes	Without lid	160pcs	2159,97	1361,00
15513125	Non-treated	Yes	Without lid	50pcs	674,29	425,00
15503125	Non-treated	No	Without lid	160pcs	2122,40	1332,00
15523125	Non-treated	No	With lid	160pcs	2084,84	1314,00



MICROPLATE READERS, MULTISKAN™ GO + FREE DROP PLATE

- Very fast plate measurements and a full spectrum of a sample in less than 10sec
- Incubation up to 45°C and variable shaking speeds to optimise sensitive applications such as enzyme kinetics and cellular assays
- A unique power save function for reduced energy consumption
- Applications:
 - Immunoassays
 - Protein assays, enzymes assay
 - Growth curve and hormone assays
 - Endotoxins, cytotoxicity
 - Food diagnostics
 - HIV assays
 - Hybridisation assays
 - Mini sequencing assays
 - Cell proliferation, cell adhesion
 - Signal transduction



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15244687	Multiskan GO without cuvette + FREE µDrop plate	1set	117992,70
15224687	Multiskan GO with cuvette + FREE µDrop plate	1set	140086,80



HDPE NATURAL 250ML SQUARE BOTTLES

- Closure type: screw cap with wad
- Closure colour: red
- Sterile



GREAT SAVINGS

Cat No.	Capacity, mL	Dimensions, diameter x height, mm	Pack qty	List price	Promo price
11793054	250	37 x 109	1set of 210	1843,74	1430,42



HDPE NATURAL 250ML ROUND BOTTLES

- Closure type: screw cap with wad
- Closure colour: red
- Sterile



GREAT SAVINGS

Cat No.	Capacity, mL	Dimensions, diameter x height, mm	Pack qty	List price	Promo price
11730555	250	37 x 125	1set of 210	2170,00	871,63



NALGENE™ STERILE NARROW-MOUTH HDPE PACKAGING BOTTLES WITH CLOSURE

- Shrink-wrapped trays are packed in double-polybagged cartons
- Durable construction
- Narrow mouth for easy pouring
- Natural HDPE bottles with white PP closures
- 10⁻⁶ SAL
- Sterile
- Leakproof



UP TO 15% OFF

Cat No.	Capacity, mL	Dimensions, diameter x height, mm	Pack qty	List price	Promo price
10304103	30	34 x 61	864pcs	6385,00	5427,25
10621294	60	39 x 86	540pcs	4780,00	4063,00
10190173	125	50 x 104	240pcs	1803,10	1555,69
10474743	250	61 x 133	180pcs	2685,00	2394,48
10262903	500	73 x 171	120pcs	1906,41	1644,82
10465123	1000	92 x 217	24pcs	640,00	544,00



BOTTLES, WASH, GRAM STAIN, VENTING, LDPE

- Pre-labelled with Gram stain solutions, hazard warning labels and instructions for use
- With patented DripLok™ vapour venting technology
- The set of five wash bottles are printed with:
 - Crystal violet solution
 - Lugol's iodine solution
 - Gram iodine solution
 - Gram colour/de-colourising solution
 - Safranin solution



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
12831464	Set of five 250mL Gram stain venting wash bottles	5pcs	249,94	225,49



BOTTLES, WASH, DIY, WIDE NECK, LDPE

- Multi-lingual wide neck wash bottles
- Indelibly printed in multiple languages
- Choice of three water or five solvent labels with new Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information
- Colour coded one-piece PP closures with precision jet-tip aperture



GREAT SAVINGS

Cat No.	Description	Capacity, mL	Cap type	Pack qty	List price	Promo price
12540066	Custom/DIY with white cap	250	Vapour-venting	5pcs	250,62	90,00
12590056	Custom/DIY with white cap	500	Vapour-venting	5pcs	301,19	90,00
15307047	Custom/DIY with mixed colour caps	500	Vapour-venting	1pcs	320,21	115,00
12510066	Custom/DIY with white cap	250	Non-venting	5pcs	218,95	90,00
12520066	Custom/DIY with white cap	500	Non-venting	5pcs	237,64	92,00
12530066	Custom/DIY with white cap	1000	Non-venting	5pcs	275,68	140,13

Accessories

10339694	Indelible black marker pen for custom/DIY wash bottles	2pcs	150,15	75,00
----------	--	------	--------	-------



BOTTLES, WASH, LDPE, MULTI-LINGUAL, VENTING

- Multi-lingual wide neck wash bottles
- Indelibly printed in multiple languages
- Choice of three water or five solvent labels with new Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information
- Colour coded one-piece PP closures with precision jet-tip aperture



GREAT SAVINGS

Cat No.	Description	Capacity, mL	Closure colour	Language	Pack qty	List price	Promo price
11373514	Acetone	250	Red	Region 1	5pcs	251,85	90,00
12126230	Ethanol	250	Orange	Region 1	5pcs	251,06	90,00
12116230	I.M.S	250	Yellow	Region 1	5pcs	251,06	125,00
12196220	Isopropanol	250	Blue	Region 1	5pcs	251,06	125,00
11383514	Methanol	250	Green	Region 1	5pcs	251,06	125,00
11329643	Acetone	500	Red	Region 1	5pcs	267,29	90,00
11349643	Ethanol	500	Orange	Region 1	5pcs	267,29	90,00
11393514	I.M.S	500	Yellow	Region 1	5pcs	266,73	130,00
11303524	Isopropanol	500	Blue	Region 1	5pcs	267,29	90,00
12550456	Methanol	500	Green	Region 1	5pcs	267,29	130,00



E1-CLIPTIP™ EQUALIZER ADJUSTABLE TIP SPACING MULTICHANNEL ELECTRONIC PIPETTES

- NEW Program multiple pipettes from your PC with the My Pipette™ Creator app
- ClipTip™ interlocking tip interface ensures complete and secure seal on every channel with lowest tip attachment and ejection forces
- Electronic tip ejection and index finger pipetting action enables comfortable and safe pipetting
- Customised and flexible user interface offers optimal and consistent operation for every application



10% OFF

Cat No.	Description	Pack qty	List price	Promo price
15623046	E1-ClipTip™ Equalizer 384 8-channels 0,5-12,5µL	1pcs	13212,83	12170,96
15633046	E1-ClipTip™ Equalizer 384 12-channels 0,5-12,5µL	1pcs	14222,14	13100,68
15643046	E1-ClipTip™ Equalizer 384 8-channels 1-30µL	1pcs	13212,83	12170,96
15653046	E1-ClipTip™ Equalizer 384 12-channels 1-30µL	1pcs	14222,14	13100,68
15663046	E1-ClipTip™ Equalizer 8-channels 2-125µL	1pcs	13212,83	12170,96
15673046	E1-ClipTip™ Equalizer 384 8-channels 2-125µL	1pcs	13212,83	12170,96
15683046	E1-ClipTip™ Equalizer 384 12-channels 2-125µL	1pcs	14222,14	13100,68
15693046	E1-ClipTip™ Equalizer 8-channels 10-300µL	1pcs	13212,83	12170,96
15603056	E1-ClipTip™ Equalizer 6-channels 50-1250µL	1pcs	12175,99	11215,87
15613056	E1-ClipTip™ Equalizer 8-channels 50-1250µL	1pcs	13212,83	12170,96

Introducing My Pipette™ Creator

Revolutionise the way you pipette with a new era of technology.

www.thermofisher.com/mypipette

thermo**scientific**

FI-CLIPTIP™ GLP PIPETTOR KITS

- ClipTip is a unique, innovative interlocking technology ensuring a complete secure seal with minimal tip attachment and ejection force
- Upgrade to exceptional performance with comfort, redefined
- Volume adjustment is virtually effortless and silent, enabling also one-handed volume changes for convenience
- A textured adjustment knob and optional grip strips enable a firmer grip that makes a huge difference in your work
- Lowered pipetting forces now reduce Repetitive Stress Injury (RSI) even further without compromising results
- For your convenience, we have packaged a variety of F1-ClipTip pipettes in a range of volumes, compatible ClipTip pipette tips and a pipette stand to support good laboratory pipetting in your lab



15% OFF

Cat No.	Description	Pack qty	List price	Promo price
15581945	F1-ClipTip GLP Kit 1 includes:			
	<ul style="list-style-type: none"> • F1-ClipTip 1-10µL single channel • F1-ClipTip 10-100µL single channel • F1-ClipTip 100-1000µL single channel • F-Series stand • ClipTip 20 rack • ClipTip 200 rack • ClipTip 1000 rack • 25mL reservoir • 25mL divided reservoir • 100mL divided reservoir 	1pcs	5202,74	4422,30
15591945	F1-ClipTip GLP Kit 2 includes:			
	<ul style="list-style-type: none"> • F1-ClipTip 2-20µL single channel • F1-ClipTip 20-200µL single channel • F1-ClipTip 100-1000µL single channel • F1-ClipTip 30-30µL 8-channel • F-Series stand • ClipTip 20 rack • ClipTip 200 rack • ClipTip 300 rack • ClipTip 1000 rack • 25mL reservoir • 25mL divided reservoir • 100mL divided reservoir 	1pcs	8320,58	7072,49

thermo**scientific**

BATHS, CIRCULATORS, REFRIGERATED/HEATED, ARCTIC SERIES

- Featuring superior cooling power, expansive temperature ranges, powerful force/suction pumps, and sophisticated digital control to ensure accuracy and reproducibility. Choose from three different controllers that allow you the flexibility to create the perfect fit for your application.
- Energy-saving feature limits power consumption: tight control of cooling capacity dependent on application needs
- Small refrigeration circulators are optimised for small footprint and low noise
- Robust stainless-steel baths
- Thermostat can be indexed 90° on all sides of the bath for easy viewing
- Three-side ventilation concept: up to two sides can be blocked with full refrigeration performance and no degradation of equipment



FEATURED PRODUCT

Cat No.	Dimensions, (L x W x H), mm	Temperature range, °C	Cooling capacity, W	Capacity, L	Pack qty	List price
11690901	483 x 273 x 748,7	-25 to 150	500	12	1pcs	39920,00

thermo**scientific**

CLIPTIP PIPETTOR TIPS

- Feature clips that lock tips in place and form a complete seal
- Each time ClipTip pipette tips are on, they are sealed and will not loosen regardless of application pressure
- Clipped on the F1-ClipTip pipette with a light touch, eliminating the frustration of banging tips on your pipette
- Utilises low retention technology to enable maximum sample recovery
- Certified free of RNase, DNase, DNA, ATP and endotoxin contamination
- Each tip insert has ink-jetted lot and volume information
- Convenient colour key makes it easy to match the colour of the F1-ClipTip pipette plunger cap with the corresponding ClipTip pipette tip rack insert



FEATURED PRODUCT

Cat No.	Filtered	Volume, µL	Pack qty	List price
ClipTip 200				
13226259	-	200	960pcs	386,92
13236259	-	200	960pcs	498,67
13286269	Filtered	200	960pcs	621,70
13206279	Filtered	200	960pcs	600,10
ClipTip 300				
13256259	-	300	960pcs	539,99
13276259	-	300	960pcs	498,67
13296269	Filtered	300	960pcs	621,70
13216279	Filtered	300	960pcs	600,10
ClipTip 384				
15216293	-	12,5	3840pcs	3249,35
15226293	Filtered	12,5	3840pcs	4385,68
15246293	-	30	3840pcs	3249,35
15256293	Filtered	30	3840pcs	4385,68
15276293	-	125	3840pcs	3249,35
15286293	Filtered	125	3840pcs	4385,68
ClipTip 1250				
15236303	-	1250	768pcs	590,70
15266303	Filtered	1250	768pcs	681,80
15256303	-	1250	768pcs	538,25
15276303	Filtered	1250	768pcs	662,08

thermo**scientific**

BATHS, CIRCULATORS, ULTRA LOW TEMPERATURE, GLACIER SERIES

- With heated tank top to avoid ice build up
- Effective cooling capacity at ultra low temperatures allows you to reach your specific application temperature requirement
- Clamped work area cover
- Insulated supply and return ports eliminate ice build up and process temperature variation
- 10 ramp programs
- Colour screen with soft touch keys



FEATURED PRODUCT

Cat No.	Dimensions, (L x W x H), mm	Temperature range, °C	Cooling capacity, W	Capacity, L	Pack qty	List price
11317100	419 x 554 x 851	-50 to 200	1000	12	1pcs	63680,00

thermo scientific

SMART2PURE™ WATER PURIFICATION SYSTEMS

- Easy-to-operate:
 - Dispensing is easy and features variable speed to control flow
 - The display can be tilted for optimal reading
 - Optional UV bulb and UF ultrafilter are available to customise ultrapure water
- Smart integrated 6L reservoir for the Smart2Pure 3 and 6:
 - Constructed with pigment-free materials
 - Conical bottom outlet for complete draining and efficient cleaning and disinfection



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
11327021	Smart2Pure 3L/hr, with UV oxidation and 6L integrated tank	1pcs	42786,17
11367021	Smart2Pure 6L/hr, with UV oxidation and 6L integrated tank	1pcs	52506,03


Distributor
GE Healthcare

PD MINITRAP™ G-25 SAMPLE PREPARATION COLUMNS

- Reliable and reproducible cleanup of proteins and other large biomolecules (>5000 Mr).
- Range of applications such as desalting, buffer exchange, and removal of low-molecular weight compounds.
- High desalting capacity.
- Flexible protocols with ability to use gravity flow or centrifugation.



FEATURED PRODUCT

Cat No.	Description	Sample preparation	Pack qty	List price
11570654	PD MiniTrap™ G-25 Sample Preparation Columns	0,1 to 0,5mL for gravity and 0,2 to 0,5mL with centrifugation,	50pcs	2723,00
11580654	PD MidiTrap™ G-25 Sample Preparation Columns	0,5 to 1,0mL for gravity and 0,75 to 1,0mL with centrifugation,	50pcs	2723,00
11783329	PD SpinTrap™ G-25 Desalting Columns	100 to 180µL without sample dilution,	50pcs	2225,00



Distributor
GE Healthcare

DISPOSABLE PD-10 DESALTING COLUMNS

- Packed with Sephadex™ G-25 medium
- Designed for ease of use and equilibration
- Prepare samples in range of 1,0 to 2,5mL by gravity flow or 1,75 to 2,5mL by centrifugation
- Desalt multiple samples in parallel with desalting capacity >90%
- Sample recovery in range of 70 to >95%
- Includes adapters for centrifugation



FEATURED PRODUCT

Cat No.	Particle Size	Diameter Inner Column, cm	Pack qty	List price
11768488	>50µm (Dry), 38 to 235µm (wet)	1,5	30pcs	2696,00



Distributor
GE Healthcare

DESALTING COLUMNS

- HiPrep 26/10 Desalting is a convenient, easy-to-use 53mL column prepacked with Sephadex G-25 Fine for desalting and buffer exchange
- Provide reliable and reproducible desalting and buffer exchange with sample sizes of up to 15mL per column
- HiTrap™ Desalting Prepacked Columns are prepacked, ready to use columns for group separation between high and low molecular weight substances



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
10460505	HiPrep™ 26/10 Desalting Columns	1pcs	4755,00
15147956	HiTrap™ Desalting Prepacked Columns	1 x 5mL	516,00
11525015	HiTrap™ Desalting Prepacked Columns	5 x 5mL	2042,00



Distributor
GE Healthcare

MINI-UNIPREP™ SYRINGELESS FILTERS

- Provides a faster and easier way to remove particulates from UHPLC samples
- Compatible with most major autosamplers
- Fewer consumables required - reduce costs by up to 40%



SPECIAL PRICE

Cat No.	Pore size, µm	Membrane	Pack qty	List price	Promo price
11305124	0,45	Cellulose	100pcs	1478,23	1786,00



VIVASPIN™ TURBO 15 CENTRIFUGAL CONCENTRATORS

- Vivaspin Turbo 15 centrifugal concentrators offer the optimal solution to any concentration or buffer exchange application. In a single spin, 15mL solutions can be concentrated up to 150x in 10–30 minutes with macromolecular recoveries in excess of 95%, making it the fastest ultrafiltration device in the market
- The new sleek internal profile keeps the filtration pressure high (i.e. speed) throughout the concentration. This superior design allows optimal crossflow conditions over the large twin vertical PES membranes that minimises protein polarisation and membrane fouling
- Made from highly inert and stable materials, Vivaspin Turbo 15 has the highest chemical resistance and provides ultra-low protein binding. The unique angular dead stop pocket (patent pending) has a pipette friendly design and allows easy and complete recoveries of the concentrate
- Vivaspin Turbo 15 can handle up to 15mL sample volume in swing bucket rotors and 11mL in fixed angle rotors accepting 50mL centrifuge tubes



SPECIAL PRICE

Cat No.	MWCO	Pack qty	List price	Promo price
12663465	10000	12pcs	1314,76	1256,85
12673465	10000	48pcs	4789,50	4578,53



VIVASPIN™ TURBO 4 CENTRIFUGAL CONCENTRATORS

- Vivaspin™ Turbo centrifugal concentrators offer the optimal solution to any concentration or buffer exchange application with their broad range of MWCOs prior to further sample analysis. Vivaspin Turbo 4 is suitable for initial sample volumes of 2 to 4mL
- Active membrane area: 3,2 cm²
- Hold up volume of membrane: <10µL
- Dead stop volume in swing out: 40µL
- Dead stop volume in fixed angle (25°): 30µL
- Membrane: polyethersulfone
- Recommended Spin Speed (xg)

Rotor	Swing Bucket	Fixed Angle
	4000	7500



SPECIAL PRICE

Cat No.	MWCO	Pack qty	List price	Promo price
15340364	10000	1pcs	1943,97	1858,34
15498801	10000	100pcs	5625,31	5377,52



CLARISTEP™ FILTRATION SYSTEM

- No syringe required – for less waste and no hand strain
- Up to eight samples are processed simultaneously – saving time
- Hold-up volume < 30µL – for maximum sample recovery



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15539714	Sartorius Claristep™ filtration station and tray	1pcs	3131,64	2835,00
15549714	Sartorius Claristep™ tray only	1pcs	730,03	670,95
15564899	Sartorius Claristep™ filters, RC 0,2µm non-sterile	96pcs	2019,10	1828,58
15559714	Sartorius Claristep™ filters, RC 0,2µm non-sterile	480pcs	8686,83	7867,13
15544899	Sartorius Claristep™ filters, RC 0,45µm non-sterile	96pcs	2019,10	1828,58
15554899	Sartorius Claristep™ filters, RC 0,45µm non-sterile	480pcs	8686,83	7867,13



MILLEX SYRINGE FILTERS, NON STERILE

- Millex™ Syringe filters set the standard for reliable small volume filtration (1mL to 200mL)



SPECIAL PRICE

Cat No.	Diameter, mm	Membrane material	Pore size, µm	Colour	Pack qty	List price	Promo price
10747984	33	Hydrophilic Millipore Express™ PLUS PES	0,22	Green	50pcs	1802,00	1711,90
10509582	33	PE Housing - PES	0,22	Green	50pcs	1802,00	1711,90
10040752	33	Hydrophilic Durapore™ PVDF	0,22	Yellow	50pcs	1720,00	1634,00
10142642	33	Hydrophilic Durapore™ PVDF	0,45	Yellow	50pcs	1731,00	1644,45
10611585	33	Nylon	0,2	Purple	50pcs	1751,00	1663,45
10769503	33	Nylon	0,45	Purple	50pcs	1751,00	1663,45

invitrogen

**DYNABEADS™,
MAGNETIC-BEADS
LIGAND-COUPLED**

- Very low non-specific binding with high signal-to-noise ratio
- No columns, centrifugations, or time-consuming pre-clearing required
- High reproducibility and high throughput compatible with KingFisher™ instruments
- IP in less than 40 minutes
- High target protein yield with low antibody consumption



20% OFF

Cat No.		Pack qty	List price	Promo price
10334693	Dynabeads™ Protein A for Immunoprecipitation	1mL	1612,00	1289,60
10746713	Dynabeads™ Protein A for Immunoprecipitation	5mL	6340,00	5072,00
10348814	Dynabeads™ Protein A for Immunoprecipitation	50mL	34930,00	27944,00
10442945	Dynabeads™ Protein A Immunoprecipitation Kit	40rxn	3000,00	3040,00
10446293	Dynabeads™ Protein G for Immunoprecipitation	1mL	1594,00	1275,20
10765583	Dynabeads™ Protein G for Immunoprecipitation	5mL	6310,00	5048,00
10607605	Dynabeads™ Protein G for Immunoprecipitation	50mL	34930,00	27944,00
10278723	Dynabeads™ Protein G Immunoprecipitation Kit	40rxn	3810,00	3048,00

- Isolate protein in less than 40 minutes
- Binds most mouse IgG subclasses

Cat No.		Pack qty	List price	Promo price
10515883	Dynabeads™ M-280 Sheep Anti-Mouse IgG	2mL	3505,00	2804,00
10077512	Dynabeads™ M-280 Sheep Anti-Mouse IgG	10mL	10260,00	8208,00
10716653	Dynabeads™ M-280 Sheep Anti-Rabbit IgG	2mL	3515,00	2812,00
10557013	Dynabeads™ M-280 Sheep Anti-Rabbit IgG	10mL	10310,00	8248,00

- Direct and fast isolation of any biotinylated molecule
- Low-charged and neutral beads, optimal for binding of DNA fragments, proteins, peptides, and antibodies
- Flexible protocols with gentle and efficient liquid-phase reaction kinetics
- Biomagnetic protocols are easily adapted to automated platforms

Cat No.		Pack qty	List price	Promo price
10292593	Dynabeads™ M-280 Streptavidin	2mL	4395,00	3516,00
10465723	Dynabeads™ M-280 Streptavidin	10mL	16050,00	12840,00
10202983	Dynabeads™ M-280 Streptavidin	100mL	91700,00	73360,00



Distributor
GE Healthcare

**HITRAP™ MABSELECT™
COLUMNS**

- HiTrap MabSelect SuRe columns are prepacked with MabSelect SuRe, a BioProcess medium for purification of monoclonal antibodies from large sample volumes. The novel alkali-tolerant recombinant Protein A ligand allow cleaning-in-place (CIP).



FEATURED PRODUCT

Cat No.	Description	Volume, mL	Pack qty	List price
15137956	HiTrap Mabselect SuRe	1	1pcs	2458,00
10064394	HiTrap Mabselect SuRe	1	5pcs	9756,00
11758408	HiTrap Mabselect SuRe	5	1pcs	9161,00
11574865	HiTrap Mabselect SuRe	5	5pcs	38344,00



Distributor
GE Healthcare

**HITRAP™ TALON CRUDE
COLUMNS**

- Columns prepacked with cobalt-based chromatography media designed for purification of histidine-tagged recombinant proteins by immobilised metal affinity chromatography (IMAC).



FEATURED PRODUCT

Cat No.	Description	Volume, mL	Pack qty	List price
15147946	HiTrap TALON crude	1	1pcs	401,00
11618893	HiTrap TALON crude	1	5pcs	1513,00
12034118	HiTrap TALON crude	5	5pcs	5181,00



Distributor
GE Healthcare

HISTRAP™ COLUMNS

- HisTrap columns are Ni Sepharose affinity media for capture and purification of histidine-tagged proteins by immobilised metal ion affinity chromatography (IMAC)



FEATURED PRODUCT

Cat No.	Description	Media	Pack qty	List price
15187946	HisTrap excel	Ni Sepharose excel	1x1mL	487,00
12632966	HisTrap excel	Ni Sepharose excel	5x1mL	1862,00
12683466	HisTrap excel	Ni Sepharose excel	5x5mL	6277,00
10126414	HisTrap FF	Ni Sepharose 6 Fast Flow	5x1mL	1384,00
10571680	HisTrap FF	Ni Sepharose 6 Fast Flow	5x5mL	4802,00
15197946	HisTrap FF Crude	Ni Sepharose 6 Fast Flow	1x1mL	401,00
10685355	HisTrap FF Crude	Ni Sepharose 6 Fast Flow	5x1mL	1513,00
11723219	HisTrap FF Crude	Ni Sepharose 6 Fast Flow	5x5mL	5162,00
15167956	HisTrap HP	Ni Sepharose High Performance	1x1mL	458,00
10431065	HisTrap HP	Ni Sepharose High Performance	5x1mL	1721,00
11763209	HisTrap HP	Ni Sepharose High Performance	1x5mL	1531,00
11773209	HisTrap HP	Ni Sepharose High Performance	5x5mL	5899,00



Distributor
GE Healthcare

SUPEROSE™ 6 INCREASE COLUMNS

- 'Increase' size exclusion (SEC) chromatography columns from GE Healthcare are the latest generation of SEC columns prepacked with agarose-based media offering higher resolution and shorter run times than their predecessors. They comprise smaller, more rigid beads with a narrower particle size distribution and a higher selectivity. These media are available in three different prepacked column formats for sample volumes between 4µL and 500µL, suitable for both preparative and analytical purposes
- Superdex 75 Increase is designed for small-scale preparative purification and analysis of recombinant tagged or other proteins with molecular weights from Mr 3000 up to 70 000. Compared with its predecessor Superdex 75, new Superdex 75 Increase offers
- Runtime reduced to one third for faster separations and results
- Up to 50% higher resolution
- Improved lot-to-lot consistency
- Superdex 200 Increase columns are excellent for detection and separation of antibody monomers, dimers, and aggregates in monoclonal antibody (MAb) preparations.
- Superose 6 Increase columns are designed for purification of large proteins and protein complexes



FEATURED PRODUCT

Cat No.		Pack qty	List price
15303234	Superose™ 6 Increase 3,2/300	1pcs	19871,00
15383224	Superose™ 6 Increase 5/150 GL	1pcs	18850,00
15377653	Superose™ 6 Increase 10/300 GL	1pcs	18850,00
15599244	Superdex™ 75 Increase 3,2/300	1pcs	20079,00
15589244	Superdex™ 75 Increase 5/150 GL	1pcs	19058,00
15579244	Superdex™ 75 Increase 10/300 GL	1pcs	19058,00
15175665	Superdex™ 200 Increase 3,2/300	1pcs	20079,00
15142355	Superdex™ 200 Increase 5/150 GL	1pcs	19058,00
15182085	Superdex™ 200 Increase 10/300 GL	1pcs	19058,00



Distributor
GE Healthcare

HISCREEN CAPTO Q IMPRES PREPACKED COLUMNS

- 10cm bed height of HiScreen columns is designed to allow method optimisation and parameter screening
- Easily connected in series to achieve 20 cm bed height
- Small bed volume gives fast results and minimal sample/buffer consumption
- Reproducible results, scalable to BioProcess columns packed with the same chromatography media using the same linear fluid velocity



FEATURED PRODUCT

Cat No.	Description	Exchange	Pack qty	List price
11653074	HiScreen Capto Q ImpRes Prepacked Column	Anion exchange	1pcs	1099,00
11390212	HiScreen™ Capto™ SP ImpRes Prepacked IEX Columns	Cation exchange	1pcs	1099,00



Distributor
GE Healthcare

AKTA START CHROMATOGRAPHY SYSTEM

- Easily purify your tagged proteins, antibodies, etc
- Gain insights from real-time monitoring
- Evaluate and share your purification results
- Fraction collector and UNICORN start software are optional modules



FEATURED PRODUCT

Cat No.	Flow Rate	Dimensions, W x H, mm	Pack qty	List price
15269764	0,5 to 5 mL/min	340 x 360	1pcs	59295,00
Accessories				
15279764	Fraction Collector		1pcs	11628,00
15279784	GE Healthcare Unicorn™ Start 1.0 Control software designed for AKTA™ start		1pcs	7454,00



Distributor
GE Healthcare

HITRAP™ CAPTO™ IMPRES IEX COLUMNS

- Combined properties of high binding capacity, high flow rate and low back pressure reduces process cycle times and improves productivity at process scale
- Rigid agarose matrix enables high bed heights and purification of viscous samples at high flow rates
- Screening of selectivity, binding and elution conditions, as well as small scale purifications
- Strong quarternary amine anion exchanger
- High-throughput medium for intermediate and polishing protein purification



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
10734572	HiTrap™ Capto™ Q ImpRes IEX Columns	5pcs	1127,00
11673074	HiTrap™ Capto™ Q ImpRes IEX Columns	5pcs	2886,00
11623074	HiTrap™ Capto™ SP ImpRes IEX Columns	5pcs	1122,00
11633074	HiTrap™ Capto™ SP ImpRes IEX Columns	5pcs	2886,00



thermo scientific

SAVANT™ SPDI3IDDA SPEEDVAC™ CONCENTRATOR

- Adjustable vacuum level to prevent samples from freezing
- Chemical-resistant dual automatic bleeder valves
- Radiant glass cover for excellent chemical resistance and additional heating for faster drying
- Digital display of parameters for easy, reproducible selection
- Independent timers for both heat and run
- PTFE tubing
- Deep PTFE-coated chamber and vapour pathway for superior chemical resistance and easy cleaning
- Accepts all standard small SpeedVac™ rotors and also a two-place microtitre plate rotor



GREAT SAVINGS

Cat No.	Temperature Range, °C	Dimensions (L x W x H), mm	Pack qty	List price	Promo price
10031382	45 to 80	450 x 360 x 330	1pcs	†02145,00	70010,18

analytikjena

UV CROSSLINKERS

- Energy efficient
- Large LED display
- Easy to use control panel
- Utilising UV energy to xbond DNA
- The internal microprocessor constantly monitors irradiated energy while ensuring repeatable results. The large LED display continuously indicates the remaining time or energy setting. Each unit has a safety interlocking system to guard against accidental exposure to the user
- The CX-2000 (254nm) crosslinker features an interlocked pull-out drawer in contrast to the pull-down door of the CL-1000 series



FEATURED PRODUCT

Cat No.	Description	Wavelength, nm	Pack qty	List price
10001292	CL-1000	254	1pcs	16530,00
12062759	CX-2000	254	1pcs	14045,00

fisherbrand

POWER SUPPLY, MIDI 300V/4

- With nearly twice the current and power of an alternative market leading product, at 700mA and 150W, the Midi 300V/4 offers a specification comparable to any equivalent power supply presently available on the market. The Midi 300V/4 is suitable for use with all Fisherbrand horizontal sub-gel systems and verti-gel units. Microprocessor controlled, with four sets of power terminals that allow simultaneous operation of up to four electrophoresis units either at constant voltage or constant current. The timer function may be set to run continuously or up to a maximum 999 minutes with an alarm that sounds to signify termination of the run. A user friendly keypad houses a clear 3-digit LED to aid set up, as well as a convenient 'pause/resume' key, particularly useful during extended runs when it is necessary to access the gel tank to monitor buffer levels and sample migration.
- Stackable design
- In-built safety systems
- Automatic crossover between parameters
- Individual indication of control parameter
- Dual voltage compatibility
- Voltage, V: 2 to 300
- Current, mA: 1 to 700
- Power (max.), W: 150
- Dimensions, mm: 190 x 305 x 95



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12613546	Power supply Midi 300V/4	1pcs	5994,00

fisherbrand

VERTI-GEL ELECTROPHORESIS UNITS

- Incorporates a sealing system which is compatible with all major types of 80mm x 100mm and 100mm x 100mm precast gel. Fisherbrand™ Vertical Gel Tank's ultra soft silicone seals and pressure bars which surround the glass plates guarantee leak proof gel casting



FEATURED PRODUCT

Cat No.	Description	Plate dimensions, mm	Pack qty	List price
11843293	Verti-Gel Mini complete vertical system	100 x 100	1pcs	6163,00
12623546	Verti-Gel Maxi vertical gel tank	200 x 200	1pcs	11111,00
15126644	Verti-Gel Maxi complete system	200 x 200	1pcs	11632,00

Join the
conversation

Connect with us on Facebook, Twitter and YouTube for the latest news, prizes and offers

invitrogen

iBLOT™ 2 TRANSFER STACKS

- Increased shelf life
- Reduced and compact packaging
- Complete compatibility with iBLOT™ 2 Gel Transfer Device



SPECIAL PRICE

Cat No.	Membrane	Pore size	Size, mm	Pack qty	List price	Promo price
15296526	Nitrocellulose	0,2 µm	130 x 83	10pcs	1500,00	1125,00
15239296	Nitrocellulose	0,2 µm	80 x 80	10pcs	1234,00	925,50
15249296	PVDF	0,2 µm	130 x 83	10pcs	1674,00	1255,50
15286066	PVDF	0,2 µm	80 x 80	10pcs	1386,00	1039,50

Accessory

10505003	iBLOT 2 Blotting Roller	86	1pcs	225,00	168,75
----------	-------------------------	----	------	--------	--------

thermo scientific

SUPERSIGNAL™ WEST PICO PLUS CHEMILUMINESCENT SUBSTRATE

- Picogram to femtogram sensitivity
- High signal stability - incubated blots provide stable signal duration over the critical 4 hour time period with up to 24 hours of light output under optimal conditions
- Stable reagent 8-hour working solution stability; 1-year kit stability at room temperature
- Exceptional robustness - provides high performance outside of the recommended antibody dilution ranges, including the most common 1:5K to 1:10K secondary antibody dilutions of a 1 mg/mL stock solution



UP TO 30% OFF

Cat No.	Description	Pack qty	List price	Promo price
15659364	SuperSignal West Pico PLUS Chemiluminescent Substrate	20mL	1030,00	772,50
15669364	SuperSignal West Pico PLUS Chemiluminescent Substrate	200mL	2825,00	2118,75
15626144	SuperSignal West Pico PLUS Chemiluminescent Substrate	500mL	1830,00	1372,50
15679364	SuperSignal West Pico PLUS Chemiluminescent Substrate	1L	486,00	364,50

thermo scientific

THERMOFLEX CHILLER

- Temperature range from 5°C to 40°C
- Internal pressure relief valve
- Custom temperature and pressure alarms
- Integrated funnel for easy filling
- Full flow fluid filter ensures clean fluid to protect your application and maximise recirculation system life
- Innovative, patented packaging for rapid installation
- Quick start guide for seamless start-up in minutes
- Minimises space requirements, even when compared to chillers with smaller footprints
- Temperature stability to ±0,1°C accuracy



FEATURED PRODUCT

Cat No.	Cooling capacity (50 Hz at 20°C)	Dimensions H x W x D, cm	Pack qty	List price
13426983	750 W / 2651 BTU	69,2 x 36 x 62,4	1pcs	46645,00
13458126	2,2 kW / 7513 BTU	73,6 x 43,6 x 67,3	1pcs	62915,00

invitrogen

iBIND™ ACCESSORIES

- The iBind™ Flex Western System offers automated western blot processing and can accommodate various blot formats including midi, mini, and vertically cut strips. Based on Sequential Lateral Flow (SLF) technology, all blocking, antibody incubation and washing steps during immunodetection are hands-free.



The iBind Solution Kits serve as a combined blocking, washing, and antibody diluent solution. Each kit contains reagents sufficient for 10 midi blots or 20 mini blots

25% OFF

Cat No.	Title	Description	Pack qty	List price	Promo price
15481244	iBind Flex Midi Insert Replacement	To process one mini blot	1set	1590,00	1192,50
15401254	iBind Flex Mini Insert Replacement	To process up to two mini blots	1set	1570,00	1177,50
15411254	iBind Flex Multi-Strip Insert Replacement	To process up to six vertically cut strip blots	1set	1580,00	1185,00
15367238	iBind Flex Cards	Pre-marked plastic-backed glass-fibre matrix	10pcs	1942,00	1456,50
15317248	iBind Flex Fluorescent Detection (FD) Solution Kit	Compatible with downstream fluorescent detection	1kit	1112,00	834,00
15357248	iBind Flex Solution Kit	Compatible with downstream chemiluminescent or colorimetric detection	1kit	922,00	691,50

MERCK

LUMINATA™ WESTERN HRP CHEMILUMINESCENCE SUBSTRATES

- Horseradish peroxidase (HRP) catalyses the oxidation of luminol in the presence of peroxide, releasing photons (chemiluminescence) that are captured on x-ray films or by CCD cameras. Most HRP substrates are plagued by quickly decaying signal and short shelf life, but Luminata substrates are formulated differently to include stabilizing reagents
- Consistent results with less pipetting error
- The stabilisers added to Luminata substrates:
- Prevent spontaneous oxidation of luminol by peroxide in the absence of HRP
- Enhance the stability of both the luminol and peroxide, ensuring that Luminata substrates are stable for one year in their premixed format
- Stabilise, strengthen, and prolong the chemiluminescence signal once HRP starts to catalyse oxidation
- Three different sensitivities for a range of detection needs:
 - Luminata Classico premixed reagent for everyday detection needs
 - Luminata Crescendo premixed reagent for detection in the picogram range
 - Luminata Forte premixed reagent for detection in the femtogram range



SPECIAL PRICE

Cat No.	Description	Sensitivity	Pack qty	List price	Promo price
12634645	Luminata Forte Western HRP Substrate	~400fg	500mL	3730,71	3307,90
11790644	Luminata Crescendo Western HRP Substrate	~1 to 3pg	500mL	3027,86	2684,70
11770644	Luminata Classico Western HRP Substrate	~6pg	500mL	2490,00	2207,80

analytikjena

STAINLESS STEEL M-26V MODEL TRANSILLUMINATOR

UVP's range of benchtop transilluminators can be used with any UVP imaging system. They use 8W tubes and have a specially constructed filter that provides the most uniform light emission possible from a tube based unit. The units have a stainless steel top that is easy to clean and a hinged UV blocking lid to protect the user. Single wavelength units have high, medium and low intensity settings that make the units suitable for both visualisation and band cutting. The dual and 3UV wavelength units have a single intensity setting. Single and dual wavelength benchtop transilluminator



FEATURED PRODUCT

Cat No.	Description	Filter size, mm	Pack qty	List price
12823928	M-26V transilluminator, 302nm	210 x 260	1pcs	11850,00
12843928	LMS-26 transilluminator, 302/365/254nm	210 x 260	1pcs	13950,00
Accessories				
11778261	UV/white converter plate 210mm x 260mm		1pcs	2345,00
12813918	VisiBlue converter plate 210mm x 260mm		1pcs	3755,00
15482170	UV Blocking Face Shield		1pcs	573,60
12853908	Gel Trays - 270 x 230mm		1pcs	399,00
12853948	Gel Trays - 290 x 230mm		1pcs	515,00
11778221	UVGL-25 Compact PL lamp 1 x 365nm tube + 1 x 254nm tube		1pcs	2615,00
10712467	UVGL-58 Hand held PL lamp 1 x 365nm tube + 1 x 254nm tube (2UV)		1pcs	3175,00
11718261	UVP Ultra Violet Product UV light		1pcs	2415,00

analytikjena

BIODOC-IT™2 IMAGING SYSTEMS

- High resolution monochrome camera with 5.0 megapixels resolution; open lens is easily accessible for adjustment of aperture, zoom and focus
- Large 254mm touch screen clearly shows captured gel images
- New software interface with easy-to-use functions to simplify your gel imaging research
- Save images to the system, USB stick (included) or connect to a network for transfer of images to another PC
- Ethidium bromide filter standard and located in the new two-position filter tray; additional filters are available
- Unique viewer window allows visibility of gels without opening the door
- Wide front door with safety switch turns the UV off when the door is open
- Side access doors allow gels to be repositioned while viewing a live preview of gel on the screen
- Compact design reduces laboratory bench space needs



FEATURED PRODUCT

Cat No.	Description	Filter size, cm	UV length, nm	Pack qty	List price
15477204	M-26	21 x 26	302	1pcs	70000,00
15435843	LMS-26	21 x 26	254/302/365	1pcs	77030,00



Fisher BioReagents™

Discover. Synthesize. Analyze. Customize.

SPECIAL OFFER! Fisher BioReagents customers spending 2500 SEK+ in a single order can claim a free 500mL sample of FASTRun™ Tris SDS PAGE Running Buffer, 10X



FASTRun Tris SDS PAGE Running Buffer, 10X will significantly cut the time you spend waiting for your gels while providing comparable or better resolution and a broad MW separation range! It is designed to be used with readily available Tris-Glycine polyacrylamide gels, precast or homemade. Recommended voltage is 200 V. Even when run under your current conditions, FASTRun Tris SDS PAGE Running Buffer will reduce the time spent waiting for results.

Claim your free sample and view the full terms and conditions at eu.fishersci.com/go/freefastrun

thermo**scientific**

MICROPLATES, 96-WELL, POLYSTYRENE, MAXISORP™ / POLYSORP

- Choose from two binding surfaces for reliable and reproducible results in immunoassays
 - MaxiSorp™: MaxiSorp, Immulon 4 HBX, and Enhanced Binding (EB) are hydrophilic and ideal for antibody sandwich assay
 - PolySorp™, PolySorp, Immulon 1B, and Universal Binding (UB) for adsorption of hydrophobic molecules
- Available in clear, black and white
- Alphanumeric orientation system
- Fits standard equipment
- Dimensions 128mm x 86mm
- Well shape: flat
- Well volume (total): 400µl
- Well volume (working): 350µl
- White plates
 - Recommended for luminescence based immunoassays
 - Maximum reflection
 - Minimum autofluorescence and autoluminescence
- Black plates
 - Recommended for fluorescence based immunoassays
 - Minimum background in fluorescence reading
 - Minimum back scattered light
- Transparent plates
 - Intended for quantitative and qualitative solid phase immunoassays
 - Recommended for colorimetric methods



20% OFF

Cat No	Surface	Colour	Type	Pack qty	List price	Promo price
10547781	Maxisorp™	Clear	F96	60pcs	1395,00	1116,00
10498802	Polysorp™	Clear	F96	180pcs	3850,00	3080,00
10396181	Maxisorp™	White	F96	80pcs	2357,18	1885,74
10020581	Polysorp™	White	F96	80pcs	2338,40	1870,72
10030581	Maxisorp™	Black	F96	80pcs	3015,00	2412,00
10316191	Polysorp™	Black	F96	80pcs	3050,00	2440,00

thermo**scientific**

COVALENT BINDING OF BIOMOLECULES WITH IMMOBILISER AMINO IMMUNOASSAY MICROPLATES

- Reduce the amount of reagents used in the lab with Thermo Scientific™ plates and modules with covalent binding surfaces
- Increase the sensitivity of your immunoassay



20% OFF

Cat No.	Format	Colour	Volume total, µL	Pack qty	List price	Promo price
10316872	Solid 96-well plate	Clear	400	30pcs	5555,00	4444,00
10432612	Solid 96-well plate	White	400	30pcs	6188,78	4951,02
10052232	Solid 96-well plate	Black	400	30pcs	6188,78	4951,02
10432802	Solid module F8 (x12) in frame	Clear	400	30pcs	6405,00	5124,00

fisher**brand**

Focus on Fisherbrand

for all your **plasticware** needs



For full range refer to www.eu.fishersci.com/go/fisherbrand



AUTOMATION-COMPATIBLE PCR MICROPLATE

- Rigid frame with a small split prevents warping during and post-thermal cycling
- In addition, ultra-thin walls provide fast and consistent heat transfer during PCR
- Plate starts flat, stays flat!



SPECIAL PRICE

Cat No.	No. of wells	Sterile	Well colour	Well design	Pack qty	List price	Promo price
15531433	96	Yes	Clear	Conical	50pcs	1570,96	1238,00



DNA/RNA CELL STAINS



SPECIAL PRICE

Cat No.	Description	Sub-cellular localisation	Detection	Excitation / Emission, nm	Flow cytometer laser lines	Pack qty	List price	Promo price
10093764	FxCycle Violet Stain	Nucleic Acids	Fluorescent	358/461	405	1pcs	3730,00	2797,50
10003774	FxCycle Far Red Stain	Nucleic Acids	Fluorescent	640/658	633	1pcs	3705,00	2778,75
13433329	FxCycle PI/RNase Staining Solution	Nucleus	Fluorescent	535/617	488	100mL	866,00	649,50
12333553	NucBlue™ Fixed Cell ReadyProbes™ Reagent	Nucleus	Fluorescent	360/460	Blue	1pcs	1028,00	740,16
15139325	NucRed™ Dead 647 ReadyProbes™ Reagent	Nucleus	Fluorescent	642/661	Far-red	1pcs	1864,00	1342,08
14879680	ActinGreen™ 488 ReadyProbes™ Reagent	Cytoskeleton	Fluorescent	495/518	Green	1pcs	1884,00	1356,48
15119325	ActinRed™ 555 ReadyProbes™ Reagent	Cytoskeleton	Fluorescent	540/565	Red-orange	1pcs	1872,00	1347,84



CELL PROFILARATION ASSAY KITS

- Multiplexable - compatible with R-PE (and tandems) and fluorescent proteins
- Accurate - superior results compared to BrdU assays
- Fast - results in as little as 60 minutes



SPECIAL PRICE

Cat No.	Description	Label or Dye	Detection	Excitation / Emission, nm	Flow cytometer laser lines	Pack qty	List price	Promo price
15564394	Click-iT™ Plus EdU Flow Cytometry Assay Kit	Alexa Fluor™ 350	Edu	350/438	-	50assays	4985,00	4237,25
15267156	Click-iT™ Plus EdU Flow Cytometry Assay Kit	Alexa Fluor™ 488	Edu	495/519	488	50assays	4500,00	3825,00
15574394	Click-iT™ Plus EdU Flow Cytometry Assay Kit	Alexa Fluor™ 594	Edu	532 or 561/615	-	50assays	4985,00	4237,25
15234657	Click-iT™ Plus EdU Flow Cytometry Assay Kit	Alexa Fluor™ 647	Edu	650/670	633/635	50assays	4600,00	3910,00
15342202	Click-iT™ Plus EdU Flow Cytometry Assay Kit	Pacific Blue™	Edu	404/450	405	50assays	4505,00	3829,25
10220455	CellTrace™ Violet Cell Proliferation Kit, for flow cytometry	CellTrace™ Violet	Fluorescent	405/450	405	1kit	3520,00	2816,00
10013162	CellTrace CFSE Cell Proliferation Kit, for flow cytometry	CellTrace™ CFSE	Fluorescent	492/517	488	1kit	2415,00	1692,00
15421114	CellTrace Yellow Cell Proliferation Kit, for flow cytometry	CellTrace™ Yellow	Fluorescent	546/579	532, 561	1kit	3345,00	2676,00
15396613	CellTrace Far Red Cell Proliferation Kit, for flow cytometry	CellTrace™ Far Red	Fluorescent	630/661	633/635	1kit	3240,00	2568,00



SPECIALITY APPLICATION PARTICLES



FEATURED PRODUCT

Cat No.	Description	Diameter, µm	Application	Pack qty	List price
11878273	PHARM-TROL™ Count Precision Size Standards	15	Ready-to-use precision standards designed for regular use in liquid particle counters for USP <788> and <789> verification.	6pcs	2985,00
11868273		15		20pcs	6100,00
10584945	COUNT-CAL™ Count Precision Standards	1,999 ± 0,020	Packaged in single-use bottles and intended to be sampled directly from the bottle, standards were developed for manufacturers of parenteral drugs and ophthalmic solutions seeking interim verification of USP 788 and 789 (Particulate Matter in Injections and Particulate Matter in Ophthalmic Solutions).	6pcs	1660,00
10369384		5,010 ± 0,035		6pcs	1630,00
10156124		10,00 ± 0,08		6pcs	1630,00
10280035		15,02 ± 0,15		6pcs	1630,00
10702037		19,99 ± 0,25		6pcs	1660,00
10322265		24,61 ± 0,22		6pcs	1630,00
10046634		29,75 ± 0,56		6pcs	1660,00
10034294		50,2 ± 0,3		6pcs	1660,00
10116514		69,1 ± 0,9		6pcs	1660,00
11848273		Validex™ Count Precision Size Standards		5,00 ± 0,10	Excellent for instrument calibration and verification in the drinking water industry
11858273	10,0 ± 0,10		500mL	2910,00	
11829902	Ezy-Cal™ Count Precision Size Standards	1,999 ± 0,020	Designed for regular use in liquid particle counters to verify instrument calibration before each use, standards assure reproducibility of the particle counter by permitting a continuous record of its performance using a NIST™-traceable microsphere suspension with a precise concentration.	100mL	2065,00
11839902		5,010 ± 0,035		100mL	2030,00
11849902		10,00 ± 0,08		100mL	2050,00
11573283		15,02 ± 0,15		100mL	2050,00
11869902		19,99 ± 0,25		100mL	2095,00
11879902		24,61 ± 0,22		100mL	2050,00
11889902		29,75 ± 0,56		100mL	2095,00
11899902		50,00 ± 0,50		100mL	2095,00
11809912		69,1 ± 0,9		100mL	2095,00

thermo scientific

EVOLUTION™ 260 BIO UV-VISIBLE SPECTROPHOTOMETER

- Reliable, high quality performance and customer usability
- Built in bio assays including nucleic acid and protein concentration and colorimetric analysis, etc.
- Precise temperature control for kinetic measurements
- Long life Xenon flash lamp reduces maintenance and downtime, lowers cost of ownership
- Double beam design with 1nm spectral bandwidth for increased sample identification and quantification results
- Innovative CUE (customised user environment) scripting and INSIGHT software minimises user error and increases flexibility
- Software and accessories for performance verification and calibration enables verification in minutes instead of hours
- Wide range of accessories including temperature control options, sample changers, sippers, and solid sampling devices increases flexibility



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
10438278	Evolution™ 260 Bio UV-Vis spectrophotometer, PC control	1pcs	88820,00

JENWAY

SPECTROPHOTOMETER GENOVA PLUS

- 5nm spectral bandwidth
- Absorbance range -0,3 to 2,5A
- Pre-programmed for DNA/RNA analysis
- Fitted with 10mm x 10mm cuvette holder
- Pre-programmed methods for protein analysis
- Press to read xenon lamp
- Three year warranty including xenon lamp




FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12201012	Genova Plus life science spectrophotometer	1pcs	38758,98

thermo scientific

IMMUNOHISTO CHEMISTRY ANTIBODIES FOR PATHOLOGY



FEATURED PRODUCT

Cat No.	Description	Formulation	Pack qty	List price
12644527	CD3 (Early T-Cell Marker), Rabbit Monoclonal Antibody	-	500µL	3897,00
12643707	CD3 (Early T-Cell Marker), Rabbit Monoclonal Antibody	-	1mL	5458,50
12633707	CD3 (Early T-Cell Marker), Rabbit Monoclonal Antibody	-	12mL	4419,00
12612627	CD38 Ab-3, Mouse Monoclonal Antibody	-	100µL	1285,00
12632627	CD38 Ab-3, Mouse Monoclonal Antibody	-	500µL	2970,00
12693697	Ki-67, Rabbit Monoclonal Antibody	-	1mL	6045,00
12683697	Ki-67, Rabbit Monoclonal Antibody	-	100µL	1370,00
12632927	Myelin Basic Protein Ab-1, Rabbit Polyclonal Antibody	-	7mL	4077,00
12620247	Perforin Ab-2, Mouse Monoclonal Antibody	-	100µL	1175,00
12640247	Perforin Ab-2, Mouse Monoclonal Antibody	-	500µL	3267,00
12630247	Perforin Ab-2, Mouse Monoclonal Antibody	-	1mL	5290,00
12610247	Perforin Ab-2, Mouse Monoclonal Antibody	-	7mL	2820,00
12642197	Placental Alk. Phosphatase (PLAP) Ab-4, Mouse Monoclonal Antibody	-	100µL	1015,00
12652197	Placental Alk. Phosphatase (PLAP) Ab-4, Mouse Monoclonal Antibody	-	500µL	2555,00
12632197	Placental Alk. Phosphatase (PLAP) Ab-4, Mouse Monoclonal Antibody	-	1mL	3154,50
12622197	Placental Alk. Phosphatase (PLAP) Ab-4, Mouse Monoclonal Antibody	-	7mL	1730,00
12632687	Negative Control for Mouse IgG1 Ab-1, Mouse Monoclonal Antibody	Purified Ab with BSA and Azide	500µL	2196,00
12602067	TTF-1 (Thyroid Transcription Factor-1) Ab-1, Mouse Monoclonal Antibody	Purified Ab with BSA and Azide	1mL	4855,50
12670847	Troponin T, Cardiac Isoform Ab-1, Mouse Monoclonal Antibody	Purified Ab with BSA and Azide	500µL	2700,00
12650847	Troponin T, Cardiac Isoform Ab-1, Mouse Monoclonal Antibody	Purified Ab with BSA and Azide	1mL	4266,00
12620857	Troponin T, Cardiac Isoform Ab-1, Mouse Monoclonal Antibody	Ready-to-use for IHC (UVLP)	7mL	2440,00
12629356	Factor XIII A Ab-1, Mouse Monoclonal Antibody	Ready-to-use for IHC (UVLP)	7mL	2855,00
12651987	CD44 Std./H-CAM Ab-4, Mouse Monoclonal Antibody	Ready-to-use for IHC (UVLP)	7mL	2005,00
12630947	c-erbB-3/HER-3 Ab-6, Mouse Monoclonal Antibody	1 mg/mL, unlabeled	200µL	5125,50
12643667	Estrogen Receptor (SP1), Rabbit Monoclonal Antibody	Culture Supernatant with Azide	500µL	6579,00

thermo**scientific**

FINNPIPETTE™ F1 AND F2 GOOD LABORATORY PIPETTING (GLP) KITS

- Everything you need for safe, and reliable pipetting in one convenient kit!
- Each kit includes a selection of three or four FinnpiPETTE F1 or F2 pipettes, a stand, a rack of Finntip or Finntip Flex tips for each pipette and sample reagent reservoirs
- **Thermo Scientific FinnpiPETTE™ F1 pipettes:**
Upgrade to exceptional performance with comfort, redefined
- Volume adjustment is virtually effortless and silent, enabling also one-handed volume changes for convenience
- A textured adjustment knob and optional grip strips enable a firmer grip that makes a huge difference in your work
- Lowered pipetting forces now reduce Repetitive Stress Injury (RSI) even further without compromising results
- **Thermo Scientific FinnpiPETTE™ F2 pipettes:**
Robust and Durable
- Light, smooth pipetting action provides a low pipetting force for reduction of RSI (Repetitive Strain Injury)
- Advanced volume gearing mechanism significantly improves accuracy, precision and durability
- Patented super blow-out in volumes 50 µl and below for increased accuracy and precision
- Fully autoclavable offering high level sterility
- The F1 and F2 GLP kits offer a variety of single channels in a range of volumes that are packaged with compatible Finntips, and a convenient pipette stand to support good laboratory pipetting in your lab



15% OFF

Cat No.	Description	Pack qty	List price	Promo price
15541945	FinnpiPETTE GLP F1 Kit 1 (1-1000µL)	1pcs	6226,35	5292,40
15551945	FinnpiPETTE GLP F1 Kit 2 (0,2-1000µL)	1pcs	7522,33	6393,98
15561945	FinnpiPETTE GLP F1 Kit 3 (10-10000µL)	1pcs	6339,04	5388,18
15571945	FinnpiPETTE GLP F1 Kit 4 (2-1000µL)	1pcs	6226,35	5292,40
11835850	FinnpiPETTE GLP F2 Kit 1 (1-1000µL)	1set	5709,83	4853,36
11845850	FinnpiPETTE GLP F2 Kit 2 (0,2-1000µL)	1pcs	7015,20	5962,92
11855850	FinnpiPETTE GLP F2 Kit 3 (10-10000µL)	1set	5813,13	4941,16
11773027	FinnpiPETTE GLP F2 Kit 4 (2-1000µL)	1set	5709,83	4853,36

thermo**scientific**

FINNTIP™ PIPETTE TIPS

- Achieve excellence in a complete Thermo Scientific FinnpiPETTE™ system
- Finntip tips are designed and manufactured expressly for FinnpiPETTE models to deliver optimal performance, precision and accuracy. They meet ISO 9001, ISO 14001 and ISO 134851 standards for exceptional quality
- Available in standard, extended length, and wide orifice tips in several packaging options within volume range from 0.2 µL to 10 mL
- Easy Identification: colour-coding on each FinnpiPETTE and Finntip rack maximises efficiency



25% OFF

Cat No.	Description	Pack qty	List price	Promo price
11722544	Finntip 250 Universal, 10x96/rack	960pcs	740,00	532,50
11732544	Finntip 250 Universal, 10x96/rack, sterile	960pcs	880,00	660,00
11792504	Finntip 10, 10x96/rack	960pcs	785,00	588,75
11702514	Finntip 10, 10x96/rack, sterile	960pcs	960,00	720,00
11783913	Finntip 1000, 10x96/rack, sterile	960pcs	750,00	562,50
10399310	Finntip 1000, 10x96/rack	960pcs	960,00	720,00
11762544	Finntip 300, 10x96/rack	960pcs	700,00	525,00
10512481	Finntip 300, 10x96/rack, sterile	960pcs	890,00	667,50

thermo**scientific**

PIPETTOR TIPS, ART™ BARRIER

- Designed to aid in pipetting for a variety of applications
- Available in sterile, non-sterile and in various packaging options
- Maximum Sample Recovery: Using advanced low-retention technology, the tips hydrophobic inner surface significantly reduces sample loss to adhesion and maximises sample recovery
- Extended length to prevent contamination by keeping the pipettor shaft from entering the sample vessel
- Gel loading for loading acrylamide or agarose gels, or wide orifice for genomic DNA, fragile cell lines and other viscous materials



30% OFF

Cat No.	Description	Pack qty	List price	Promo price
12661480	ART™ 100E Hinged Racks, Sterile	960pcs	1062,70	743,89
10768735	ART™ 200 Barrier Tip, Hinged Racks, Sterile, Low Retention, 200µL	960pcs	1062,70	743,89
15193606	ART™ 300 Barrier Tip, Hinged Racked, Sterile, Low Retention, 300µL	960pcs	1062,70	743,89
15123636	ART™ 1000 REACH Barrier tip, extended length, low retention, universal, hinged rack 8 x 96 tips	768pcs	1293,28	905,30
11942689	ART™ 10, MicroPoint, reload insert, low retention, universal, hinged rack 10 x 96 tips	960pcs	1062,70	743,89
10421175	ART™ 10 REACH Low Retention, Hinged Racks, Sterile, 10µL	960pcs	1062,70	743,89
10032874	ART™ 20E MicroPoint, low retention, universal, hinged rack 10 x 96 tips	960pcs	1062,70	743,89
10276944	ART™ 20 Barrier Tip, Hinged racks, Sterile, Low Retention, 20µL	960pcs	1062,70	743,89
15163606	ART™ 200 REACH Barrier tip, extended length, low retention, universal, hinged rack 8 x 96 tips	768pcs	1412,62	988,83
10431185	ART™ 1000 Barrier Tip, Hinged Racks, Sterile, Low Retention, 1000µL	768pcs	1063,78	744,65

fisher**brand**

Focus on Fisherbrand

for all your liquid handling needs



For full range refer to www.eu.fishersci.com/go/fisherbrand

fisher**scientific**
by Thermo Fisher Scientific



TRANSFER PIPETTES

- Applications: subculturing, harvesting, sampling 15, 50, 100 mL tubes and flasks
- Remove supernatant and media from cell cultures
- Wash cells
- Post-centrifugation to remove all cell layers



GREAT SAVINGS

Cat No.	Description	Length, mm	Stem Diameter, mm	Overall Capacity, mL	Bulb Draw, mL	Drops per mL	Individually wrapped	Pack qty	List price	Promo price
15397853	Extra Long	300	8,6	23	7,3	20	Yes	100pcs	386,00	301,01
15367833	Extended Fine Tip Standard	157	4,6	5,8	3,5	50	No	500pcs	304,00	237,07
15387833	Extended Fine Tip Standard	157	4,6	5,8	3,5	50	Yes	400pcs	555,00	427,81
10374850	Extended Fine Tip Large Bulb	153	3,0	5,0	3,3	50	No	250pcs	257,00	200,42
15317833	Extended Fine Tip Large Bulb	153	3,0	5,0	3,3	50	Yes	400pcs	600,00	453,03
15357823	Extended Fine Tip Small Bulb	104	3,0	1,5	1,0	50	Yes	400pcs	458,00	353,04



HDPE VISCODISPO™ SAMPLER

- Disposable sampler for highly viscous media. No cross-contamination with other samples, and no complicated, time-consuming, expensive cleaning
- HDPE
- Diameter 21mm
- 20 pieces per pack, individually packed



GREAT SAVINGS

Cat No.	Capacity, mL	Length, mm	Sterile	Pack qty	List price	Promo price
10369212	100	500	No	20pcs	1836,81	1010,08
10674293	190	1000	No	20pcs	2254,26	1239,73
10286942	100	500	Yes	20pcs	2129,03	1170,77
10414812	190	1000	Yes	20pcs	2796,96	1537,98



LIQUIDISPO™ DISPOSABLE SAMPLER

- Avoid cross-contamination with other samples and complicated, time-consuming, expensive cleaning
- HDPE
- Diameter 21mm
- 20 pieces per pack, individually packed
- Ergonomic handle
- No seals
- Clean room manufactured
- Optionally sterilised by gamma rays



GREAT SAVINGS

Cat No.	Capacity, mL	Length, mm	Sterile	Pack qty	List price	Promo price
10013872	100	500	No	20pcs	1836,81	1010,08
10339022	190	1000	No	20pcs	2254,26	1239,73
10656003	100	500	Yes	20pcs	2129,03	1170,77
10475002	190	1000	Yes	20pcs	2796,96	1537,98



DISOPIPETTE DISPOSABLE SAMPLER

- Intended for sampling of low viscosity media
- All-layer samples can be obtained when the DisoPipette is inserted into a liquid while opened
- Target point samples are obtained by inserting while closed and opening the DisoPipette at the sampling point
- Two opening sizes (2,5 and 5,5mm) for different viscosities
- Clean room manufactured
- Sterile varieties are sterilised by gamma rays



FEATURED PRODUCT

Cat No.	Capacity, mL	Length, mm	Sterile	Pack qty	List price
15521741	100	500	No	20pcs	650,00
10767075	180	900	Yes	50pcs	1419,35
15482985	100	500	Yes	20pcs	760,00
11901501	200	1000	Yes	50pcs	1753,32



POLYSTYRENE GRADUATED PIPETTE

- For dairy, environment, industries, microbiology
- Graduated



SPECIAL PRICE

Cat No.	Capacity, mL	Length, mm	Sterile	Pack qty	List price	Promo price
15354734	1	270	Yes	4000pcs	4615,00	3616,46
15344734	10	294	Yes	1000pcs	2360,00	1849,84



POLYPROPYLENE FUNNEL STAND

- Funnel stand with reducing attachment for safe and effective use with different funnel sizes. Suitable for funnels of 40mm to 160mm diameter.
- Polypropylene, white, break-resistant
- For two funnels
- Each with two reducer inserts
- Steplessly height adjustable



UP TO 45% OFF

Cat No.	Dimensions, L x W x H, mm	Pack qty	List price	Promo price
10057980	300 x 140 x 480	1pcs	830,00	634,25
10393731	550 x 150 x 500	1pcs	1556,76	856,20



TITRATING BURETTE

- **Titration burette with shatter protection**
The titrating burette with splinter protection is a fast acting burette with Schellbach stripes based on Dr. Schilling's design with automatic zero point extraction and Bürkle dosing armature. The glass burette has a plastic sheathing which increases resistance to breaking and protects against splinters. The possibility of the burette breaking is much lower than with a conventional burette, and if it does break, there will be no dangerous splinters.
- **Titration burette automatic**
Burette with Schellbach stripes type Schilling / Bürkle with original Bürkle dosing armature. A universal aid in the laboratory, electroplating industry, water treatment plants, for field analyses, etc. Precise dosing using the press button, exact fine titration with the micro-screw and precise discharge jet.



FEATURED PRODUCT

Cat No.	Volume, mL	Pack qty	List price
With shatter protection			
15591751	15	1pcs	540,00
11981501	25	1pcs	540,00
11991501	50	1pcs	610,00
Automatic, borosil. Glass			
10503461	5	1pcs	570,00
10208230	10	1pcs	570,00
10701961	15	1pcs	628,70
10686331	25	1pcs	570,00
10513461	50	1pcs	665,00



CIMAREC+™ STIRRING HOTPLATE SERIES

- Stir-Trac features innovative engineering for slow-speed stirring, consistent speed control, and strong magnetic coupling
- Temperature adjustable in 1 degree increments
- Hot top warning system protects from accidental burns with prominent display when heating surface is above 50°C (122°F)
- Top plate choices in ceramic or aluminium
- Easy-to-read LED display for heating and stirring
- Raised display design protects electronics from spills
- Optional temperature probe and splash guard protection shield
- Variable sizes for different sample volumes
- Dimensions (L x W x H), mm: 330 x 208 x 97



FEATURED PRODUCT

Cat No.	Description	Max. Temperature	Speed range	Surface Area Heating, mm	Pack qty	List price
15403610	Stirrer	-	50 to 1500rpm	184 x 184	1pcs	2335,00
15493600	Hotplate	540°C	-	184 x 184	1pcs	2235,00
15413610	Stirring Hotplate EU	540°C	50 to 1500rpm	184 x 184	1pcs	3900,00



STIRRERS, OVERHEAD, DIGITAL

- Built for high volume/high viscosity lab or pilot plant applications
- High performance brushless DC motors
- Overload protection
- Hollow through shaft for easy impeller adjustment
- Clamp, stand and mixing impellers available separately
- Safety approved to EN 61326 CE / IEC 61010-1 and 61010-2-51, EMC EN 50081-1, EN 50082-1 ROHS, WEEE
- Three year warranty



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price	Promo price
15203807	Overhead stirrer, 40 to 2000rpm, 25L	1pcs	8829,85	
15293797	Overhead stirrer, high torque, 20 to 3000rpm, 60L	1pcs	14526,66	12778,38
Accessories				
11984985	Compact stand with base 457mm x 8mm (for 15203807 only)	1pcs	3090,00	2666,00
12937106	Pivoting blade impeller 400mm length x 60mm diameter	1pcs	870,84	788,08
12917106	Crossed blade impeller 400mm length x 50mm diameter	1pcs	934,43	875,91



SONIC DISMEMBRATOR

- Four models available:
- Model 50 Sonic Dismembrator: processes 0,2 to 50mL, 50 watts, continuous or remote operation, (footswitch optional) lightweight converter (340g) for hand-held operation
- Model 120 Sonic Dismembrator: processes 0,2 to 50mL, 120 watts, programmable from 1 second to 10 hours, overload protection, pulse mode (1 to 59 seconds), energy display (watts and Joules), lightweight converter (340g) for hand-held operation
- Model 505 Sonic Dismembrator: processes 0,5 to 500mL, 500 watts, programmable from 1 second to 10 hours, overload protection, pulse mode (1 to 59 seconds), energy display (watts and Joules)
- Model 705 Sonic Dismembrator: processes 0,2 to 1000mL, 700 watts, touch screen display, programmable from 1 second to 99 hours, save up to 10 different programs, overload protection, pulse mode (1 second to 24 hours), energy display (real-time watts and cumulative Joules), temperature monitoring (requires, optional temperature probe), full amplitude control from 1 to 100%



10% OFF

Cat No.	Description	Pack qty	List price	Promo price
12961151	Model 50: Includes generator, converter, 3mm probe, cables, wrench set, manual, 230V, 50/60Hz, 203mm x 190mm x 152mm (w x l x h)	1pcs	21025,00	18922,50
12337338	Model 120: Includes generator, converter, 3mm probe, cables, wrench set, manual, 230V, 50/60Hz, 203mm x 330mm x 152mm (w x l x h)	1pcs	34580,00	31122,00
12931151	Model 705: Includes generator, converter, 13mm probe, cables, wrench set, manual, 230V, 50/60Hz, 203mm x 387mm x 216mm (w x l x h)	1pcs	41744,83	37570,35

Accessories

12921181	Microtip probe 1,9mm diameter, 0,2mL to 5mL processing volume, for 50 and 120 watt sonicators	1pcs	2481,48	2233,33
12911181	Microtip probe 3,1mm diameter, 0,5mL to 15mL processing volume, for 50 and 120 watt sonicators	1pcs	2215,61	1994,05
12931181	Microtip probe 6,3mm diameter, 5mL to 50mL processing volume, for 120 watt sonicator	1pcs	2481,48	2233,33
12921171	Probe 12,7mm diameter, 10mL to 250mL with replaceable tip, for 500 and 700 watt sonicators	1pcs	4652,78	4187,50
12971161	Probe 19,05mm diameter with replaceable tip, 25mL to 500mL processing volume	1pcs	4741,40	4267,26
12901171	Probe, 25,4mm diameter with replaceable tip 50mL to 1000mL for 700 watt sonicator	1pcs	5007,27	4506,54
12951171	Replacement tip 12,7mm diameter, for 500 and 700 watt sonicators	1pcs	690,00	621,00
12961171	Replacement tip 19,05mm diameter, for 500 and 700 watt sonicators	1pcs	620,37	558,33
12971171	Replacement tip 25,4mm diameter, for 500 and 700 watt sonicators	1pcs	708,99	638,09



Distributor
GE Healthcare

FILTER UNITS, SYRINGE, 30MM, PURADISC™

- Exhibit very low protein binding, so they are suitable for protein recovery applications. Whatman™ 30mm Puradisc™ FP sterile cellulose acetate syringe filters are hydrophilic, making them functional with aqueous and alcoholic media despite having limited solvent resistance.



SPECIAL PRICE

Cat No.	Description	Pore size, µm	Colour	Sterile	Pack qty	List price	Promo price
11399304	Cellulose acetate, non-pyrogenic	0,2	Red	Yes	50pcs	1318,33	1213,00



RIGHT-TO-KNOW, SAFETY-VENTED WASH BOTTLES, 500ML

- Wide mouth 53mm for easy refilling
- Translucent LDPE bottles with colour-coded caps for easy identification
- Pre-printed with chemical name and formula, GHS pictogram and signal word, National Fire Protection Association (NFPA) Section 704 Four-colour Diamond, health hazards, Chemical Abstract System (CAS) number, and suggested protective clothing and equipment



FEATURED PRODUCT

Cat No.	Chemical label	Colour	Pack qty	List price
15509991	Ethanol	White	4pk	338,00
15579981	Sodium Hypochlorite (Bleach)	Yellow	4pk	330,00



SPATULA

- The spoon spatula is a double spatula with a 0,5 mL spoon and a 17mm spoon blade
- Double scoop spatula with rounded ends, content 1,0 and 0,2mL
- Available either in a large package with up to 100 pieces or packaged individually and sterilised by gamma rays



FEATURED PRODUCT

Cat No.	Description	Sterile	Pack qty	List price
15674104	Spoon spatula, LaboPlast, PS, 0,5 ml/17 mm	No	100pcs	162,48
15614114	Spoon spatula, SteriPlast, PS, 0,5ml/17mm	Yes	100pcs	395,93
15684104	Micro-spatula, LaboPlast, PS, 1,0/0,25mL	No	100pcs	162,48
15624114	Micro-spatula, SteriPlast, PS, 1,0/0,25mL	Yes	100pcs	395,93

thermo scientific

CELL CULTURE CENTRIFUGE

- Includes: Thermo Scientific SL8, 230V (Mfr PN: 15132559); TX-150 Rotor (Mfr PN: 15175233); TX-150 Roundbucket (set of 4) (Mfr PN: 15192633); Sealing lids (set of 4) (Mfr PN: 15152693); 50mL conical skirted and unskirted adapter (set of 4) (Mfr PN: 15184233); 15mL conical adapter (set of 4) (Mfr PN: 15183263)



10% OFF

Cat No.	Description	Pack qty	List price	Promo price
15287559	SL8 Cell culture package	1pcs	35256,84	31731,16

thermo scientific

CENTRIFUGE TUBES, 15ML AND 50ML, CONICAL, NUNC™

- Environmentally friendly recyclable plastic rack that reduces contamination from traditional foam packaging
- Leakproof to protect your samples



FEATURED PRODUCT

Cat No.	Volume, mL	Max. rcf, xg	Packing	Pack qty	List price
11397201	15	10500	Bulk	500pcs	1970,00
11307211	15	10500	Racked	500pcs	2405,00
11347201	50	17000	Bulk	500pcs	1634,06
11317211	50	17000	Racked	300pcs	1405,00

MERCK

FILTER DEVICE AMICON™ CENTRIFUGAL TUBES

- Concentrates 500µL down to 15µL
- Concentration factor 25x to 30x
- Fixed angle rotor: 14000xg for concentration spin; 1000xg for reverse spin high sample recovery (greater than 90% of dilute starting solution)
- Fast processing time - 10 to 30 minutes
- Reverse spin capability provides consistent recoveries



GREAT SAVINGS

Cat No.	MWCO, kDa	Pack qty	List price	Promo price
10795785	3	8pcs	825,00	437,95
10212584	3	24pcs	1900,00	1030,75
10028713	10	8pcs	526,86	437,95
10088753	10	24pcs	1950,00	1030,75
10505274	30	8pcs	523,43	435,10
10012584	30	24pcs	2010,00	1021,25
10378434	50	8pcs	523,43	435,10
10194304	50	24pcs	1775,00	1021,25
10527164	100	8pcs	523,43	435,10
10114694	100	24pcs	2035,00	1021,25

thermo scientific

MEDIA BOTTLES

- Reduce the risk of cross contamination with these sterile, disposable media bottles
- An inexpensive and shatterproof alternative to glass media bottles
- Nonpyrogenic and noncytotoxic to protect cells
- Moulded-in graduations won't rub off
- Made from clear, lightweight, polyethylene terephthalate copolyester (PETG) providing impact resistance and excellent gas-barrier properties
- Closure is high density polyethylene (HDPE)
- Bottles packed in a shrink-wrapped tray to stay clean
- Heat-shrink band around cap and bottle neck provides visual assurance of sterility for individual bottles until band is removed



GREAT SAVINGS

Cat No.	Description	Volume, mL	Pack qty	List price	Promo price
12074118	Media Bottles, Clear, PETG	125	24pcs	625,00	298,21
12084118	Media Bottles, Clear, PETG	250	24pcs	755,00	381,81
12064118	Media Bottles, Clear, PETG	500	12pcs	462,00	298,21
12094118	Media Bottles, Clear, PETG	1000	12pcs	745,00	416,48

burkle

ALL PURPOSE TRAYS

- These are practical all-purpose trays that can be nested when empty, also durable, shatterproof and resist flexing. The trays have a strong, round rim. Suitable as a spill trough for environmentally dangerous, aggressive substances, for experiments, as an etching tray or an instrument holder. Also available as a set of 4 and 6 trays.
- White polypropylene
- Trays can be stacked
- Temperature resistant up to +60°C (not suitable for temperature controlled storage or dishwashers)



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
11931651	Trays, polypropylene, four piece set: 0,5L/ 1,5L/ 3L and 10L	1set	515,00	390,69
11941651	Trays, polypropylene, six piece set: 0,5L/ 1,5L/ 3L, 10L, 21L and 39L	1pcs	1470,00	1124,36
10130501	Tray, all-purpose, polypropylene. 130mm x 180mm, 0,5L	1pcs	59,65	45,68
10028320	Tray, all-purpose, polypropylene, 180mm x 240mm, 1,5L	1pcs	134,00	81,27
10221771	Tray, all-purpose, polypropylene, 240mm x 300mm, 3L	1pcs	196,00	118,12
10129880	Tray, all-purpose, polypropylene, 300mm x 400mm, 10L	1pcs	251,00	191,85
10404121	Tray, all-purpose, polypropylene, 400mm x 500mm, 21L	1pcs	437,00	333,27
10363511	Tray, all-purpose, polypropylene, 500mm x 700mm, 39L	1pcs	640,00	486,42

thermo**scientific**

BIOLITE FLASKS, CELL CULTURE TREATED



- Polystyrene
- HPDE cap
- Sterile

GREAT SAVINGS

Cat No.	Culture area, cm ²	Cap Type	Suggested working volume, mL	Pack qty	List price	Promo price
11830765	25	Vented	7	200pcs	1380,50	555,00
11849594	25	Plug Seal	7	200pcs	1305,37	491,00
11884235	75	Vented	25	100pcs	1117,55	478,00
11840765	75	Plug Seal	25	100pcs	1108,16	426,00
11850845	175	Vented	50	40pcs	748,48	323,00
11845275	175	Plug Seal	50	40pcs	755,99	284,00

thermo**scientific**

BIOLITE MULTIDISHES, CELL CULTURE TREATED



- Polystyrene, flat bottom well
- With lid
- Sterile

GREAT SAVINGS

Cat No.	Number of Wells	Well Dimension, mm	External dimensions, mm	Culture area, cm ² /well	Sug-gested working vol/well, mL	Pack qty	List price	Promo price
11825275	6	35	85,4 x 127,6	9,6	3,0	50pcs	638,60	317,00
11889684	12	22,1	85,4 x 127,6	3,8	2,0	50pcs	684,62	317,00
11874235	24	15,6	85,4 x 127,6	1,9	1,0	50pcs	555,96	357,00
11820765	48	9,75	85,4 x 127,6	0,75	0,500	50pcs	664,89	317,00
11835275	96	6,4	85,4 x 127,6	0,32	0,200	50pcs	566,29	297,00



ECM MIMETIC AND DEFINED SURFACES



- Ready-to-use cultureware suitable as replacement for natural, self-coated Fibronectin for adult stem, progenitor, and primary cell types in defined media environments

SPECIAL PRICE

Cat No.	Coating	Description	Pack qty	List price	Promo price
15310001		6 well plate	10pcs	2673,20	2501,36
15350001	Corning™ Pure-Coat™ ECM Mimetic Cultureware Fibronectin Peptide	Flask, 75 cm ² , vented cap, canted neck	10pcs	4255,94	3279,58
15380001		Flask, 175 cm ² , vented cap, canted neck	10pcs	5049,10	3890,84
15319830		Multi-Flask 3-layer	8pcs	7173,11	6023,57
15329830		Multi-Flask 5-layer	6pcs	9380,17	7228,25
15300011		6 well plate	10pcs	2435,60	2279,04
15340011		Flask, 75 cm ² , vented cap, canted neck	10pcs	2910,79	2723,69
15360011	Corning™ PureCoat™ ECM Mimetic Collagen I Peptide	Flask, 175 cm ² , vented cap, canted neck	10pcs	3623,58	3390,66
15349830		Multi-Flask 3-layer	8pcs	6734,94	5804,50
15359830		Multi-Flask 5-layer	6pcs	8346,60	7009,09

thermo**scientific**

BIOLITE DISHES, CELL CULTURE TREATED



- Polystyrene
- With lid and airvent
- Sterile

GREAT SAVINGS

Cat No.	External Dimen-sions, mm	Internal Dimensions, mm	Culture area, cm ²	Suggested working volume, mL	Pack qty	List price	Promo price
11810765	35 x 10	34,3 x 9,3	9,4	3	500pcs	986,07	478,00
11844335	60 x 15	52,7 x 12,6	21,5	5	500pcs	1230,24	491,00
11815275	100 x 20	88 x 16	60,8	12,5	200pcs	761,62	310,00
11804125	150 x 20	138,5 x 15,4	148	35	120pcs	1389,89	581,00

memmert
Experts in Thermostatics

BATHS, WATER, HEATED, WNB SERIES



- Microprocessor PID temperature controller
- Digital display
- Over-temperature protection
- Integrated digital timer

5 % OFF

Cat No.	Capacity, L	Power, W	Dimensions external, W x D x H, mm	Pack qty	List price	Promo price
10738303	7	1200	356 x 468 x 337	1pcs	5446,88	5174,54
10329342	10	1200	356 x 578 x 337	1pcs	5916,44	5620,62
10573693	14	1800	436 x 578 x 347	1pcs	6367,21	6048,85
10185702	22	2000	436 x 578 x 405	1pcs	7146,68	6789,35
10433302	29	2400	516 x 818 x 343	1pcs	8217,27	7806,41
10003832	45	2800	516 x 818 x 401	1pcs	10198,81	9688,87



CELL CULTURE MEDIA & REAGENTS



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15373541	Dulbecco's Modification of Eagle's Medium/ Ham's F-12 50/50 Mix	6x500mL	600,00	455,21
15323609	Accutase™ Cell Detachment Solution, Liquid	100mL	362,00	326,23
15393571	Dulbecco's Phosphate-Buffered Saline, 1X	6x500mL	467,00	332,48
15333571	Dulbecco's Phosphate-Buffered Saline, 10X	6x500mL	1675,00	1132,02



Distributor
GE Healthcare

HYCLONE™ CELL CULTURE MEDIA



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
10397042	ADCF-Mab liquid medium with L-glutamine, HEPES; without phenol red	1000mL	526,00
10235312	ADCF-Mab liquid medium with HEPES; without L-glutamine, phenol red	1000mL	566,00
10433142	CDM4CHO liquid medium with L-glutamine	1000mL	618,00
10463522	CDM4CHO liquid medium without L-glutamine	1000mL	618,00
11506341	CDM4HEK293 liquid medium without phenol red, L-glutamine, HEPES	1000mL	579,00
10399663	CDM4MAb liquid medium with L-glutamine	1000mL	750,00
10771484	CDM4MAb liquid medium without L-glutamine	1000mL	665,00
10761634	CDM4NS0 liquid medium without L-glutamine	1000mL	838,00
15280496	HyCell TransFx-H liquid medium	1000mL	703,00
10443295	PF-CHO LS liquid medium with L-glutamine; without phenol red	1000mL	596,00
10483522	SFM4CHO liquid media with L-glutamine	1000mL	567,00
10376852	SFM4CHO liquid media without L-glutamine	1000mL	567,00
10403912	SFM4CHO-Utility liquid media with L-glutamine	1000mL	496,00
10500283	SFM4HEK293 liquid medium with L-glutamine, poloxamer 188, sodium bicarbonate; without protein, phenol red, HEPES	1000mL	624,00
10397092	SFM4MAb liquid medium with L-glutamine	1000mL	524,00
11536341	SFM4Transfx-293 liquid medium with poloxamer 188, sodium bicarbonate; without L-glutamine, phenol red, HEPES	1000mL	722,00
10673143	SFX-Insect liquid medium with L-glutamine, sodium bicarbonate	1000mL	526,00

fisherbrand ASPIRATION SYSTEM



- Chemical resistant oil-free aspiration station
- Diaphragm pump with long service life
- Rapid action coupling hose connectors with automatic isolation valves

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15567479	Aspiration Basic System, comfortable operation, free flow 15L/min, vacuum 150 mbar	1pcs	9661,88	8233,31
15348191	Aspiration Advantage System with 4L PP bottle, with vacuum regulator on front panel and vacuum gauge <ul style="list-style-type: none"> • Integrated chemical resistant diaphragm pump • Pumping speed 15L/min • Ultimate pressure mbar 100 	1pcs	18311,73	14078,17
10700854	HandVac basic pipettor kit, comprising a single channel plastic adapter for pipettor tip with ejector, silicone hose (inside diameter 6mm, wall thickness 3mm, length 2m) and rapid action coupling hose connector	1pcs	4658,29	3969,46
10412163	Eight-channel stainless steel adapter	1pcs	1095,54	852,51
13594340	Eight-channel plastic adapter with ejector	1pcs	1271,31	1002,60
15441934	Single channel plastic adapter for disposal tips with ejector	1pcs	437,14	336,96
15432055	Silicone hose, 2m length, inside diameter 6mm x wall thickness 3mm	2m	262,93	230,00
15451934	Table top stand for HandVac and accessories	1pcs	1165,36	906,84
12967914	In-line filter	1pcs	2469,26	2210,63
15461934	4L polypropylene receiving bottle includes lid, in-line filter, tubing and rapid action coupling hose connector	4L	3569,10	2743,95
15309804	Footswitch for Aspiration Advantage System, with on-off function	1pcs	1220,45	1093,87

fisherbrand PIPETTES, SEROLOGICAL, STRAIGHT, POLYSTYRENE, BULK WRAPPED



- DNase and RNase free
- Sterilisation using gamma irradiation
- A black magnifier strip runs lengthwise along the pipette making the meniscus obvious and reducing reading errors
- The 5mL, 10mL and 50mL sizes have both ascending and descending graduations
- Tip design on the 50mL size increases the effect of surface tension to minimise dripping

GREAT SAVINGS

Cat No.	Capacity, mL	Graduations, mL	Colour	Pack qty	List price	Promo price
11879660	1	0,01	Yellow	1000pcs	1014,75	617,88
11889660	2	0,01	Green	500pcs	689,20	419,81
11899660	5	0,1	Blue	500pcs	741,74	451,73
11809670	10	0,1	Orange	500pcs	896,27	510,25
11829670	50	0,5	Purple	100pcs	1308,35	822,89

fisherbrand PIPETTES, SEROLOGICAL, STRAIGHT, POLYSTYRENE, INDIVIDUALLY WRAPPED



- Non-pyrogenic and non-cytotoxic warranty
- DNase and RNase free
- Supplied in paper/plastic packaging, or plastic/plastic packaging
- Sterilisation using gamma irradiation

GREAT SAVINGS

Cat No.	Capacity, mL	Graduations, mL	Colour	Pack qty	List price	Promo price
Pipettes individually wrapped plastic/plastic						
11819660	1	0,01	Yellow	1000pcs	1025,05	639,54
11879650	2	0,01	Green	500pcs	690,23	388,17
11829660	5	0,1	Blue	200pcs	327,22	193,71
11839660	10	0,1	Orange	200pcs	346,15	221,92
11517752	25	0,2	Red	200pcs	757,20	438,14
11537752	50	0,5	Purple	100pcs	1256,84	787,65
Pipettes individually wrapped paper/plastic						
11849181	1	0,01	Yellow	1000pcs	1137,34	699,11
11859181	2	0,01	Green	500pcs	655,21	398,72
11869181	5	0,1	Blue	200pcs	334,82	229,18
11879181	10	0,1	Orange	200pcs	368,81	221,92
11839181	25	0,2	Red	200pcs	832,43	506,07
11889181	50	0,5	Purple	100pcs	1289,81	831,82

thermo scientific

HERACELL™ VIOS 160i AND 250i CO₂ INCUBATORS



- THRIVE active airflow system for optimal cell growth
- New sensor technology allowing shorter recovery times
- Exclusive condensation-free humidification system
- High temperature sterilisation with push button simplicity
- HEPA air filter for ISO 5 air quality conditions

FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15331035	VIOS 160i CO ₂ incubator single stainless steel chamber TC sensor	1pcs	86821,37
15363212	VIOS 160i CO ₂ incubator single stainless steel chamber IR sensor	1pcs	94531,52
15572710	VIOS 250i CO ₂ incubator single chamber, stainless steel inner chamber, TC180 CO ₂ sensor	1pcs	124639,61
15512720	VIOS 250i CO ₂ incubator single chamber, stainless steel inner chamber, IR180Si CO ₂ sensor	1pcs	137834,21

thermo scientific

SHAKER, CO₂ RESISTANT, STARTER KIT, MAXQ™



- Speed 0 to 300rpm, ±1rpm
- Shaking orbit 19mm
- Timer function 0 to 99h 50min or continuous mode
- Starter kit includes : shaker, universal platform, flask clamps (2 x 100mL, 4 x 250mL, 2 x 1L, 2 x 2L)

FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15351105	Shaker orbital, MaxQ 2000 CO ₂ Plus, carbon dioxide-resistant	1pcs	23261,93

fisherbrand

GENERAL PURPOSE WATER BATHS

- Temperature range from ambient to 99,9°C
- Temperature stability ±0,1°C and temperature uniformity ±0,05°C at 37°C
- Bright LCD controller with icon driven menu
- Stores three preset temperatures
- Automatic timer with option to set time to start procedure
- Global voltage - one part number works anywhere in the world
- Audible alarms help protect your work
- General purpose water baths include a polycarbonate gabled lid
- Shaking water baths with reciprocal oscillation from 30 to 200 oscillations per minute, include stainless steel gable cover, shaking tray



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
15375877	General purpose water bath 5L	1pcs	6060,00	4632,26
15385877	General purpose water bath 10L	1pcs	6850,00	5499,18
15395877	General purpose water bath 20L	1pcs	7745,00	6315,27
15305887	General purpose water bath 28L	1pcs	8375,00	6438,75
15325887	Shaking water bath 15L	1pcs	45265,00	12524,96
15335887	Shaking water bath 27L	1pcs	49250,00	16487,63

fisherbrand

Focus on Fisherbrand

for all your equipment needs



For full range refer to
www.eu.fishersci.com/go/fisherbrand

fisherscientific
by Thermo Fisher Scientific

Julabo
THE TEMPERATURE CONTROL COMPANY

PURA WATER BATH

- High-quality water baths with a temperature stability of ±0,2 °C. With a working temperature range of +25 °C to 99,9 °C
- Can be used for a wide variety of purposes, including routine laboratory tasks, holding samples at a constant temperature, incubation, materials testing, corrosion testing, temperature control of cultures, thermal testing of food, and dental applications
- Spray-water protected with a membrane keypad and have an LED temperature display for straightforward operation
- Highly durable and scratch-proof enamelled coating
- Integrated protection against running dry
- Integrated drain



FEATURED PRODUCT

Cat No.	Working temperature range (°C)	Bath opening / bath depth (W x L / D cm)	Filling volume liters	Ambient temperature	Dimensions W x L x H (cm)	Pack qty	List price
15690404	+18 to +99,9	30 x 27 / 17	1 to 14	5 to 40 °C	42 x 38 x 30	1pcs	7065,19
15699754	+18 to +99,10	76 x 27 / 17	2 to 30	5 to 40 °C	85 x 38 x 30	1pcs	9175,57

fisherbrand

MICROSCOPES, RESEARCH, WITH LED ILLUMINATION

- Ergonomically designed; professional frame with low position coaxial focusing allows one handed operation, making extended use comfortable
- Adjustable 3W LED illumination
- Built-in low drive mechanical stage with extended post for X and Y control (135mm x 140mm)
- 18mm eyepieces with rubber eyeshields
- 30° inclined ergonomic Siedentopf viewing head
- Five year warranty (excludes bulbs, cords and fuses)
- Interpupillary distance, mm: 55 to 75
- Inclination: 30° inclined
- Eyepieces: WF 10x/18mm
- Revolving nosepiece: four positions



15% OFF

Cat No.	Optical Heads	Objectives	Pack qty	List price	Promo price
15398824	Binocular head, Siedentopf	Achromat	1pcs	8388,41	7130,15
15328834		Super high contrast	1pcs	9108,03	7741,83
15358834		Plan	1pcs	10284,12	8741,50

invitrogen

MICROSCOPY REAGENTS



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15205739	ProLong™ Diamond Antifade Mountant	5x2mL	1960,00	1411,20
15247528	ProLong™ Diamond Antifade Mountant with DAPI	5x2mL	2115,00	1522,80
11549306	ProLong™ Gold Antifade Mountant with DAPI	10mL	2045,00	1472,40
11514147	BlockAid™ Blocking Solution	50mL	1454,00	1046,88
12363643	Image-iT™ Fixation/Permeabilization Kit	1pcs	2020,00	1454,40




GEN3™ MULTI-BARRIER PIPETTE TIPS

- Effectively separates liquid from the self-sealing barrier using an additional protective layer of conventional polyethylene filter material.
- Multi-barrier self-sealing filter prevents aerosol or liquid cross contamination
- Proprietary filter design prevents possible PCR inhibition
- Maxymum Recovery™ technology for ultra-low retention and high yield
- Hinged rack for convenience and ease of use
- Pre-sterilised; no need to autoclave
- Maxymum Recovery™ is included as a standard feature on all GEN3 multi-barrier tips
- Ultra-smooth surfaces
- Exclusive molding technology
- Reduced sample loss and saves costly reagents
- Maximises accuracy and precision



SPECIAL PRICE

Cat No.	Capacity, µL	Pack qty	List price	Promo price
12802283	10	4800pcs	4982,96	4172,00
12892273	20	4800pcs	4982,96	4172,00
12882273	200	4800pcs	4982,96	4172,00



Go to www.eu.fishersci.com to request literature specific to your industry

invitrogen

FLOID™ IMAGING SYSTEMS

- Routine fluorescent imaging teaching labs
- Simple installation, USB ports
- Intuitive software, high-resolution LCD display
- Manual mechanical stage
- Mechanical or fixed stage option
- DVI ports
- Networking capability
- 5-position objective turret
- 20x fixed objective
- Fluorescence channels: 3
- Monochrome camera
- Epifluorescence
- Transmitted light
- Teaching tool
- Fits in hood or on benchtop
- Associated printer
- Multilanguage user interface
- Integrated reagent selection guide



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15520390	Floid™ Imaging Systems	1pcs	126500,00	107525,00

invitrogen


STREPTAVIDIN CONJUGATES

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
10307962	Streptavidin, Alexa Fluor™ 488 conjugate	1mg	2830,00	2037,60
10164762	Streptavidin, Alexa Fluor™ 488 conjugate	0,5mL	3445,00	2480,40
10453442	Streptavidin, Alexa Fluor™ 555 conjugate	1mg	2930,00	2109,60
10728853	Streptavidin, Alexa Fluor™ 555 conjugate	0,5mL	3450,00	2484,00
10626153	Streptavidin, Alexa Fluor™ 594 conjugate	1mg	2930,00	2109,60
10237202	Streptavidin, Alexa Fluor™ 594 conjugate	0,5mL	3445,00	2480,40
10308062	Streptavidin, Alexa Fluor™ 647 conjugate	1mg	2920,00	2102,40
10349222	Streptavidin, Alexa Fluor™ 647 conjugate	0,5mL	3450,00	2484,00

invitrogen

SECONDARY DETECTION KITS & REAGENTS



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
10727464	Goat anti-Mouse IgG1 Cross-Adsorbed Secondary Antibody, Alexa Fluor 488	250µL	1858,00	1337,76
10082502	F(ab') ₂ -Goat anti-Rabbit IgG (H+L) Secondary Antibody, Alexa Fluor 488	250µL	1858,00	1337,76
10216302	Goat anti-Mouse IgG2a Secondary Antibody, Alexa Fluor 555	250µL	1584,00	1140,48
10256392	Goat anti-Mouse IgM Heavy Chain Secondary Antibody, Alexa Fluor 594	250µL	1568,00	1128,96
10789384	Goat anti-Mouse IgG1 Secondary Antibody, Alexa Fluor 647	250µL	1568,00	1128,96
10134883	Goat anti-Mouse IgG1 Secondary Antibody, Biotin-XX	250µL	1126,00	810,72
10206882	F(ab') ₂ -Goat anti-Rabbit IgG (H+L) Secondary Antibody, Alexa Fluor 555	250µL	1670,00	1202,40
10514963	F(ab') ₂ -Goat anti-Rabbit IgG (H+L) Secondary Antibody, Alexa Fluor 594	250µL	1670,00	1202,40
10739574	Goat anti-Rabbit IgG (H+L) Secondary Antibody, Alexa Fluor 647	250µL	1992,00	1434,24
10378262	Goat anti-Rabbit IgG (H+L) Secondary Antibody, Biotin-XX	250µL	1214,00	874,08

invitrogen

COUNTLESS™ II AUTOMATED CELL COUNTERS

- The Invitrogen™ Countess™ II and Countess™ II FL Automated Cell Counters contain advanced autofocusing and counting algorithms to allow you to quickly and accurately identify and count cells within a population while avoiding user variation associated with manual focusing mechanisms. The Countess II FL automated cell counter offers bright-field and two optional, user-interchangeable fluorescent channels. The fluorescent colours collected are determined by the insertion of individual light cubes.
- Counting mode: bright-field channel plus two optional, user-changeable fluorescent channels
- Slides: reusable and disposable slides
- Counting time: as little as 10 seconds
- Focus: auto focus with manual focus option
- User profiles: customise and save up to 10 user profiles
- User interface: intuitive touch screen with ability to use a mouse, if desired



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15307812	Countess II	1set	49900,00	41916,00
15311986	Invitrogen™ Countess™ II FL Reusable Slide	1pcs	2190,00	1839,60

invitrogen

SECONDARY DETECTION KITS & REAGENTS



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15642232	Alexa Fluor™ 488 Tyramide SuperBoost™ Kit, goat anti-mouse IgG	50pcs	2120,00	1526,40
15611892	Alexa Fluor™ 488 Tyramide SuperBoost™ Kit, goat anti-mouse IgG	150pcs	4185,00	3013,20
15631902	Alexa Fluor™ 488 Tyramide SuperBoost™ Kit, goat anti-rabbit IgG	50pcs	2120,00	1526,40
15661892	Alexa Fluor™ 488 Tyramide SuperBoost™ Kit, goat anti-rabbit IgG	150pcs	4105,00	2955,60
15621892	Alexa Fluor™ 555 Tyramide SuperBoost™ Kit, goat anti-mouse IgG	150pcs	4090,00	2944,80
15671892	Alexa Fluor™ 555 Tyramide SuperBoost™ Kit, goat anti-rabbit IgG	150pcs	4095,00	2948,40
15652232	Alexa Fluor™ 594 Tyramide SuperBoost™ Kit, goat anti-mouse IgG	50pcs	2110,00	1519,20
15631892	Alexa Fluor™ 594 Tyramide SuperBoost™ Kit, goat anti-mouse IgG	150pcs	4180,00	3009,60
15641902	Alexa Fluor™ 594 Tyramide SuperBoost™ Kit, goat anti-rabbit IgG	50pcs	2120,00	1526,40
15681892	Alexa Fluor™ 594 Tyramide SuperBoost™ Kit, goat anti-rabbit IgG	150pcs	4185,00	3013,20
15641892	Alexa Fluor™ 647 Tyramide SuperBoost™ Kit, goat anti-mouse IgG	150pcs	4090,00	2944,80
15691892	Alexa Fluor™ 647 Tyramide SuperBoost™ Kit, goat anti-rabbit IgG	150pcs	4090,00	2944,80
15597505	Biotin XX Tyramide SuperBoost™ Kit, goat anti-mouse IgG	150pcs	3655,00	2631,60
15651892	Biotin XX Tyramide SuperBoost™ Kit, goat anti-rabbit IgG	150pcs	4090,00	2944,80

fisherbrand

Focus on Fisherbrand



for all your
**sampling and
storage** needs

For full range refer to www.eu.fishersci.com/go/fisherbrand



fisherbrand

**TRACEABLELIVE™
WI-FI DATALOGGING
REFRIGERATOR/FREEZER
THERMOMETERS**

- The simple, efficient, and reliable way to stay connected to your critical environments
- Wireless Temperature Monitoring and more - connect to existing WiFi
- Alerts: temperature alarm, connectivity interruption, low battery
- Cloud-based data interface - no local software required; web, Android, Apple IOS data visibility
- Continuous data transmission and monitoring using SSL encrypted SECURE connection
- Receive MOBILE PUSH / E-MAIL / TEXT notifications
- Unlimited scalability - use TraceableLIVE™ in individual labs or an entire network of facilities
- Multiple point temperature calibration - traceable to NIST



GREAT SAVINGS

Cat No.		Pack qty	List price	Promo price
15621904	TraceableLIVE WiFi Therm, 1 bottle 1 bullet	1pcs	2812,61	2021,10
15611904	TraceableLIVE WiFi Therm, 2 bottle probes	1pcs	2042,00	2141,42

Note: Requires TraceableLIVE™ subscription, starting at \$8,99/month per device.

thermoscientific

**CRYOGENIC
STORAGE SYSTEMS, LN2,
LOCATOR PLUS™**

- Versatile: available in four sizes, suited for both manual and computerised inventory record-keeping methods; storage capacity up to 6000 vials!
- Convenient: unique hanging rack and gridded box design allows for vial inventory tracking and easy sample retrieval



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15568813	Locator Jr, Plus, 71L, includes 4 racks, dimensions (D x H): 558mm x 683mm	1pcs	34446,81
15578813	Locator 4 Plus, 121L, includes 4 racks, dimensions (D x H): 558mm x 953mm	1pcs	39593,17
15598813	Locator 6 Plus, 184L, includes 6 racks, dimensions (D x H): 660mm x 953mm	1pcs	44739,54
15588813	Locator 8 Plus, 121L, includes 8 racks, dimensions (D x H): 558mm x 953mm	1pcs	39593,17

thermo scientific

THERMO SCIENTIFIC™ HERAFREEZE™ SERIES FREEZERS -86°C

- For storage from -50°C to -86°C
- Protection of critical samples, with exceptional thermal performance, safety and security through state-of-the-art engineering
- Outstanding BTU (British Thermal Unit) reserve, and fast door recovery times to retain sample integrity
- Save up to 15% in energy usage with the energy savings mode (compared to high performance mode) or choose the high performance mode for applications requiring ultra tight temperature uniformity
- Easy-to-use, touch screen interface
- Sample storage capacity: available in five sizes, excellent capacity per footprint
- 24-hour freezer monitoring, a detailed event log, and a built-in USB port enables easy downloading of reports to a portable drive



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12893263	Thermo Scientific™HERAFreeze™ Series -86°C freezer, 421 liters, 230V/50Hz, 198 x 68.6 x 95.5 cm (H x W x D)	1pcs	105547,36
11650823	Thermo Scientific™HERAFreeze™ Series, 548 liters, 230V/50Hz, 198 x 82 x 95.5 cm (H x W x D)	1pcs	119324,21
12304393	Thermo Scientific™HERAFreeze™ Series, 682 liters, 230V/50Hz, 198 x 96.5 x 95.5 cm (H x W x D)	1pcs	128471,21
11670823	Thermo Scientific™HERAFreeze™ Series, 815 liters, 230V/50Hz, 198 x 110.8 x 95.5 cm (H x W x D)	1pcs	136672,00
11680823	Thermo Scientific™HERAFreeze™ Series, 949 liters, 230V/50Hz, 198 x 125.10 x 95.5 cm (H x W x D)	1pcs	147591,63

fisherbrand

ULTRA-LOW TEMPERATURE WIFI DATALOGGING THERMOMETER

- A simple, efficient, and reliable way to stay connected to your critical environments with smartphone, tablet or personal computer.
- Monitor remotely temperature and more - connect to existing WiFi
- Ideal for monitoring ultra-low temperature freezers and include a FEP cable
- Alerts: temperature alarm, connectivity interruption, low battery
- Cloud-based data interface - no local software required; web, Android, Apple IOS data visibility
- Continuous data transmission and monitoring using SSL encrypted SECURE connection
- Receive MOBILE PUSH / E-MAIL / TEXT notifications
- Unlimited scalability - use TraceableLIVE™ in individual labs or an entire network of facilities
- Multiple point temperature calibration - traceable to NIST



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
15651904	TraceableLIVE -90 WiFi Therm, 1 probe	1pcs	3528,00	2535,14
15661904	TraceableLIVE -90 WiFi Therm, 2 probes	1pcs	4116,00	2957,70

Note: Requires TraceableLIVE™ subscription, starting at \$8,99/month per device.

TENAK

ECO ALURACK FOR MICROTITRE PLATES FOR UPRIGHT FREEZERS

- A rack with a shelf height of 65 mm can be used with ex. 1 x 60 mm Matrix plate or 6 x 10 mm 96 microwell plate or 4 x 15mm 96 microwell plate or 2 x 96 deepwell plate. Always make sure that you have measured the exact height of the plate you want to store
- Our racks can be used with all freezer models. The racks are made of aluminum
- No. of shelves is height x depth = total no. of shelves (plates) per rack

FEATURED PRODUCT

Cat No.	Shelf height, mm	Dimensions, H x D x W, mm	No. of shelves	Pack qty	List price
12418250	16	153 x 540 x 135	9 x 6 = 54	1pcs	2050,00
12428250	16	221 x 540 x 135	13 x 6 = 78	1pcs	2225,00
12438250	16	272 x 540 x 135	16 x 6 = 96	1pcs	2375,00
12448250	16	323 x 540 x 135	19 x 6 = 114	1pcs	2575,00

TENAK

ECO ALURACK FOR CRYOBOXES FOR UPRIGHT FREEZERS

- The racks below are sold without cryoboxes
- They are also offered including standard white cardboard boxes with 9x9 cell dividers.
- Below No. of boxes is expressed by height x depth = No. of boxes



FEATURED PRODUCT

Cat No.	Box height, mm	Dimensions, H x D x W, mm	No. of boxes	Pack qty	List price
12442510	50	221 x 558 x 139	4 x 4 = 16	1pcs	1045,00
12462510	50	276 x 558 x 139	5 x 4 = 20	1pcs	1060,00
12492510	50	331 x 558 x 139	6 x 4 = 24	1pcs	1230,00
12442520	50	221 x 688 x 139	4 x 5 = 20	1pcs	1225,00
12462520	50	276x688x139	5x5 = 25	1pcs	1280,00



memmert
Experts in Thermostatics

Product Range



UNIVERSAL OVENS

VACUUM OVENS

STERILISERS

INCUBATORS

COOLED INCUBATORS

CO₂-INCUBATORS

CLIMATE CHAMBERS

HUMIDITY CHAMBERS

CONSTANT CLIMATE CHAMBERS

WATERBATHS

OILBATHS

100% ATMOSAFE. MADE IN GERMANY.

www.memmert.com | www.eu.fishersci.com





JERRYCANS

- Different liquids can be stored and decanted in a small space, it can also be stored lying on its side e.g. in a refrigerator. It is very stable due to ribbed structure and thick sides.
- PP, transparent
- With tamper-proof closure
- Canister without valve autoclavable
- With molded, indestructible litre graduation
- Canister support made of aluminium for secure stand as accessory



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
10464661	Flachmann, without thread. conn., 5L	1pcs	349,00	266,88
10757561	Flachmann, with thread. conn., 5L	1pcs	374,00	285,04
11971451	Jerry can and stopcock set	1pcs	482,00	368,30
10058241	Compact Jerry can and stopcock set	1pcs	735,00	444,79



SPRAY BOTTLES

- Media-contacting parts PP, PE, and stainless steel
- Bottle made of LDPE
- Nozzle diameter 0,6mm
- Spray volume per stroke: 1,2mL ±0,1
- Adjustable spray jet
- Temperature max. 30°C
- Spraying lance extension 30cm available as accessory



GREAT SAVINGS

Cat No.	Capacity, mL	Height, mm	Pack qty	List price	Promo price
11896840	250	220	1pcs	59,65	45,50
10190071	500	240	1pcs	67,90	51,55
11886840	1000	295	1pcs	68,85	53,19



WASH BOTTLES, ADJUSTABLE SPRAY

- For hot, cold, aqueous, and solvent based solutions
- Set to «fine mist» for TLC plates, «needle» for rinsing glassware, «big stream» for a powerful jet needed to wash filters
- Instant off-trigger release eliminates messy drips
- Chemical-resistant and virtually unbreakable bottles are made of high-density polyethylene (HDPE), spray made of polypropylene (PP)
- Can be gas or chemically sterilised and cleaned in the laboratory dishwasher
- Zero contamination
- The double headed 1000mL sized bottle has a side fill opening to help maintain purity of solutions



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
11795873	Double-headed adjustable spray wash bottle 1000mL	1pcs	475,00	380,02
11714356	Adjustable spray wash bottle 240mL	3pcs	205,00	154,26



BOTTLES, CARBOY, HDPE AND PP

- Space saving rectangular design
- Wide mouth closure
- Available in 9L or 20L, with or without spigot
- Graduated in litres and gallons
- Material:
 - HDPE, able to withstand low temperatures
 - Polypropylene, autoclavable



GREAT SAVINGS

Cat No.	Material	Capacity, L	Pack qty	List price	Promo price
13418039	HDPE with spigot	9	1pcs	725,00	408,79
13408039	HDPE with spigot	20	1pcs	835,00	408,79
13438039	PP with spigot	9	1pcs	630,00	370,79
13428039	PP with spigot	20	1pcs	905,00	531,62



SOLVENT PUMPS

- Solvent pump, Mini, for reliable filling of small quantities
- The Mini solvent pump is a high quality hand pump for small containers such as bottles, canisters and tanks up to approx. 10 litres. For reliable filling and dosing of solvents and readily flammable liquids.
- Ideal for 1 to 5L containers
- Connecting thread GL 45
- Parts in contact with medium manufactured from PTFE and stainless steel (1.4301, 1.4404)
- Including 2m grounding cable
- The following solvent pumps are also available. Please see our website for more information



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
15322265	Mini solvent pump	1pcs	3555,00	3489,74
10504772	Foot Operated Solvent Pump	1pcs	5275,00	4031,02
10238232	Manual Solvent Pump	1pcs	6080,00	3681,60
10482594	Hand Operated Solvent Pump for Tin-foil Canisters	1pcs	5065,00	3875,66
10223924	Foot Operated Solvent Pump for Burkle Discharge Hose/Stopcock	1pcs	8920,00	5402,26



NATURAL & FORCED CONVECTION INCUBATORS WITH SINGLE/TWIN DISPLAY AND TOUCHSCREEN

- Temperature range up to +80°C with natural convection or forced air circulation N/F
- Sizes from 30L to 750L with two model variants: SingleDISPLAY and TwinDISPLAY
- A wide range of options for programming and documentation using interfaces, integrated data logger and AtmoCONTROL software



5 % OFF

Cat No.	Model	Capacity, L	Temp Range	Convection type	Pack qty	List price	Promo price
12885873	IN30 SingleDISPLAY	32	5 to 80	Natural	1pcs	9199,96	8734,26
12895873	IN55 SingleDISPLAY	53	5 to 80	Natural	1pcs	11353,93	10786,23
12805883	IN75 SingleDISPLAY	74	5 to 80	Natural	1pcs	13448,16	12775,75
12815883	IN110 SingleDISPLAY	108	5 to 80	Natural	1pcs	15861,69	15068,61
12825883	IN160 SingleDISPLAY	161	5 to 80	Natural	1pcs	21186,48	20127,16
12835883	IN260 SingleDISPLAY	256	5 to 80	Natural	1pcs	28258,03	26845,13
12865883	IN30plus Twin-Display	32	5 to 80	Natural	1pcs	15354,56	14586,83
12875883	IN55plus Twin-Display	53	5 to 80	Natural	1pcs	17326,71	16460,37
12885883	IN75plus Twin-Display	74	5 to 80	Natural	1pcs	19543,03	18565,88
12895883	IN110plus TwinDisplay	108	5 to 80	Natural	1pcs	21965,95	20867,65
12805893	IN160plus TwinDisplay	161	5 to 80	Natural	1pcs	27469,17	26095,71
12815893	IN260plus TwinDisplay	256	5 to 80	Natural	1pcs	33855,17	32162,41
12845893	IF30 SingleDisplay	32	10 to 80	Forced	1pcs	11447,84	10875,45
12855893	IF55 SingleDisplay	53	10 to 80	Forced	1pcs	14781,70	14042,62
12865893	IF75 SingleDisplay	74	10 to 80	Forced	1pcs	16819,59	15978,61
12875893	IF110 SingleDisplay	108	10 to 80	Forced	1pcs	20078,32	19074,40
12885893	IF160 SingleDisplay	161	10 to 80	Forced	1pcs	24989,90	23740,41
12895893	IF260 SingleDisplay	256	10 to 80	Forced	1pcs	32775,18	31136,42



BIOLOGICAL SAFETY CABINETS SAFE 2020

- Featuring an advanced design, these cabinets comply with stringent regulatory standards and are certified to EN 12469
- Superior personnel and product protection with best-in-class air flow and filter technology
- Exceptionally comfortable ergonomic design prevents user fatigue and promotes safe working habits
- Maximised energy efficiency for environmental protection and significant cost savings
- SmartFlow™ digital technology maintains consistent airflow during normal filter loading or temporary airflow obstruction
- Sash opening: 20cm working height; can be fully opened to 77,3cm for thorough cleaning



FEATURED PRODUCT

Cat No.	Size	Pack qty	List price
10284804	0,9m	1pcs	89742,02
10462614	1,2m	1pcs	93000,76
10294804	1,5m	1pcs	117389,63
10422784	1,8m	1pcs	120845,58

Accessories

10603955	Adjustable stand 750-950mm for 0,9m	1pcs	6555,04
10171813	Adjustable stand 750-950mm for 1,20m	1pcs	7606,85
10294014	Adjustable stand 750-950mm for 1,50m	1pcs	8283,01
10690394	Adjustable stand 750-950mm for 1,80m	1pcs	8808,92



ACCUMET™ AE150 BENCHTOP PH METERS

- The essential functions needed for routine measurements without complicated features and processes at a budget-friendly price
- One to five point pH calibration
- Automatic buffer recognition with selection of two pH buffer sets
- pH slope and offset display
- Automatic or manual temperature compensation
- Ready icon to indicate measurement stability
- Electrode condition icon to indicate status
- 100 dataset storage with non-volatile memory
- pH range 0,00 to 14,00 pH
- pH resolution 0,01 pH
- Warranty: 3 years



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15514693	Benchtop pH Meter AE150 meter and stand	1pcs	3578,47
15534693	Benchtop pH Meter AE150 bio set with double junction	1pcs	4954,81

thermo scientific

ORION™ STAR™ A211 BENCHTOP PH METERS

Reliably measure pH, mV, ORP and temperature with the Thermo Scientific™ Orion™ Star A211 pH Benchtop Meter for advanced lab analysis and printer or computer interfacing

- pH Range, pH: -2,000 to 20,000
Resolution, pH: 0,1/ 0,01/ 0,001
Accuracy, pH: ±0,002
pH calibration: Up to 5 points, with editing
- mV Range, mV: ±2000,0 (Rel.mV)
Resolution, mV: 0,1/ 0,1 (Rel.mV)
Accuracy, mV: ±0,2/(Rel.mV) or ±0,05% of reading whichever is greater
- Temp. ORP: Yes (mV/RmV, E, ORP mode)
Range, °C: -5 to 105
Resolution, °C: 0,1
Accuracy, °C: ±0,1
Temperature: 1 point (offset calibration)
- General Memory: 2000 with time and date stamp
Data storage method data logging function: manual, ready (includes AUTO-READ™), timed
Inputs: pH electrode: BNC, reference pin, ATC probe: 8 pin mini-DIN



15% OFF

Cat No.	Description	Pack qty	List price	Promo price
11671449	Star A211 Meter, plus Triode 3-in-1 pH/ATC probe, ROSS™ solution kit (475mL each of pH 4, 7 and 10 buffers; storage solution, cleaning solution, and pH electrode storage bottle), electrode arm with holder and universal power adapter	1pcs	8400,00	7140,00
15585630	Star A211 pH meter, 8172BNWP ROSS Sure-Flow glass-body pH electrode, 927007MD stainless steel ATC probe, 810199 ROSS pH buffer and solution kit, electrode stand, universal power adapter, literature CD, printed quick start guide, computer interface cable and meter test certificate	1pcs	7762,54	7762,54
11621499	Meter with pH/temperature module, Triode pH/ATC probe, stirrer probe, ROSS™ pH buffers, electrode storage and cleaning solutions and bottle	1pcs	15495,00	13170,75
10495933	Combination pH Micro Electrode	1pcs	5865,00	4985,25
12493617	Orion ROSS Ultra Glass Triode pH/ATC Combination Electrodes with a BNC and 8-pin MiniDIN Connector and 1 meter cable	1pcs	4125,00	3506,25
10544973	Glass body combination electrode	1pcs	3545,00	3013,25



SEVENCOMPACT™ PH/ION METER, MICRO KIT

- pH meter Mettler S220-Micro benchtop, 5 calibration points, TRIS compatible probe included, PC connectivity
- The SevenCompact™ series combines precise electrochemical measurement technologies with flexibility, innovative design and ease-of-use.
- Optimal measuring specifications ensure great performance in nearly all pH or conductivity related lab measurements.
- IP54 water and dust resistance for a long lasting investment
- Protection for the connectors
- Transparent cover as an extra barrier



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15573423	SevenCompact pH/Ion meter S220-Micro Kit	1 unit	7529,31	6177,81

fisherbrand ISOTEMP™ CHILLER

- The Isotemp Recirculating Chiller platform offers a compact line of refrigerated recirculating chillers. Ideal for laboratory applications that require higher cooling and pumping capacities
- Intuitive digital controller for ease of use
- Customer controlled high and low temperature alarm limits
- Small footprint and quiet operation
- Stability: ±0,1°C
- Temperature Range: Setpoint Temperature Range: -10°C to +30°C



GREAT SAVINGS

Cat No.	Cooling capacity, W	Volume, L	Pack qty	List price	Promo price
15542541	1000	9,5	1pcs	25930,17	19935,27
11356641	1900	18,9	1pcs	34713,00	32879,24



FIVEEASY™ PLUS PH/MV METERS

- Compact, space saving design with ABS housing
- RS232 or USB ports to export data directly either to a printer or a PC for further processing
- ATC and MTC temperature compensation
- Five point (max.) linear/segmented calibration with four predefined buffer groups
- Ideal for wide range of laboratory applications

Specifications:

- pH Range: -2,00 to 16,00
- pH Resolution: 0,01/0,1
- pH Accuracy: ±0,01
- mV Range: -2000 to 2000mV
- mV Resolution: 1mV
- mV Accuracy: ±1mV
- Temperature Range: -5° to +105°C (23° to 221°F)
- Temperature Resolution: 0,1°C
- Temperature Accuracy: 0,3°C
- Memory Size: 200 measurements, current calibration



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
15583370	Mettler Toledo™ FiveEasy™ Plus FP20 pH/mV Meters	1pcs	7518,42	5780,21

fisherbrand HEATED AND REFRIGERATED BATH CIRCULATORS

- The robust refrigeration system ensures that your samples and applications can be cooled quickly
- Stainless steel tanks for durability and ease of cleaning
- Features energy savings mode
- Heater, kW: 2,0
- Cooling capacity, W (at 20°C): 900W



SPECIAL PRICE

Cat No.	Capacity, L	Temp Range, °C	Pack qty	List price	Promo price
11356651	8 to 15	-35 to +200	1pcs	40460,00	34205,37

fisherbrand

BATHS, CIRCULATORS, HEAT-REFRIGERATED, CLOSED SYSTEM

- Full range heating and cooling from -10°C to 80°C
- Small reservoir ensures fast temperature response time
- Small footprint
- Quiet operation
- Intuitive user interface allows user to choose between five temperature set-points
- Temperature range: -10 to 80°C
- Capacity: 2,8L
- Stability: ±0,1°C
- Heater capacity: 2000W
- Dimensions (w x d x h), mm: 227 x 436 x 619
- Mass: 29,9kg



SPECIAL PRICE

Cat No.	Pump pressure, mbar	Pump flow, max., L/min	Cooling capacity, W (at 20°C)	Pack qty	List price	Promo price
12019907	300	17	250	1pcs	24065,00	19745,38
11396651	300	17	500	1pcs	29280,00	24024,30

fisherbrand

EXTERNALLY THREADED CRYOGENIC STORAGE VIALS



50% OFF

Cat No.	Description	Volume, mL	Pack qty	List price	Promo price
11777939	Cryovial	1,2	100pcs	202,00	101,00
11787939	Cryovial	2	100pcs	202,00	101,00
11797939	Cryovial	5	50pcs	1104,28	800,00

thermoscientific

NUNC™ UNIVERSAL TUBES

- Use this versatile tube either inside a large storage system, or in smaller scale when used for traditional cryostorage in dewars



GREAT SAVINGS

Cat No.	Volume, mL	Thread	Code type	Packaging	Pack qty	List price	Promo price
15583724	1,8	External	2D	Barcoded rack	480pcs	6502,10	4400,00
13406113	1,8	External	uncoded	Latch rack	480pcs	3108,48	2055,69
13416113	1,8	External	2D	Bulk	480pcs	4000,64	2720,84
13426113	1,8	External	uncoded	Bulk	480pcs	2554,40	2177,17
15297368	2	Internal	uncoded	Latch rack	480pcs	3070,91	2618,00
15207378	2	Internal	2D	Bulk	480pcs	3944,29	3014,00
15217378	2	Internal	uncoded	Bulk	480pcs	2516,83	2145,00

fisherbrand

CRYOGENIC VIAL STORAGE BOXES



GREAT SAVINGS

Cat No.	Description	Dimensions (L x W x H), mm	Tube Array	For tubes, mL	Pack qty	List price	Promo price
11788029	Cryobox	133 x 133x 50	9 x 9	1,2 to 2	4pcs	248,00	202,00
11798029	Cryobox	133 x 133x 50	9 x 9	5	4pcs	268,00	268,00

fisherbrand

COLOUR CODERS FOR EXTERNALLY AND INTERNALLY THREADED CRYOGENIC STORAGE VIALS



GREAT SAVINGS

Cat No.	Colour	Pack qty	List price	Promo price
11927974	White	500pcs	935,00	601,94
11937974	Yellow	500pcs	935,00	601,94
11947974	Green	500pcs	935,00	601,94
11957974	Blue	500pcs	935,00	601,94
11967974	Red	500pcs	935,00	601,94

sartorius

PIPETTORS, ELECTRONIC, PICUS

- Ergonomic lightweight design to protect user from Repetitive Strain Injury (RSI)
- Electronic brake and piston control system guarantee accurate and precise pipetting
- Intuitive user interface
- Fast volume setting and menu navigation
- Hot key for quick access to favourite programs
- Calibration in 1, 2 or 3 points
- Autoclavable lower parts
- Li-Polymer battery allows charging in approx. 1 hour



5 % OFF

Cat No.	Volume range, µL	Increments, µL	Pack qty	List price	Promo price
Single channel					
13223019	0,2 to 10	0,02	1pcs	6526,86	6200,52
13233019	5 to 120	0,10	1pcs	6526,86	6200,52
13486358	10 to 300	0,20	1pcs	6526,86	6200,52
12645814	50 to 1000	1,00	1pcs	6526,86	6200,52
13263019	100 to 5000	5,00	1pcs	6526,86	6200,52
13273019	500 to 10000	10,00	1pcs	6526,86	6200,52
8 channel					
13283019	0,2 to 10	0,02	1pcs	10057,94	9555,04
13293019	5 to 120	0,10	1pcs	10057,94	9555,04
13203029	10 to 300	0,20	1pcs	10057,94	9555,04
13213029	50 to 1200	1,00	1pcs	10057,94	9555,04
12 channel					
13223029	0,2 to 10	0,02	1pcs	12030,09	11428,59
13233029	5 to 120	0,10	1pcs	12030,09	11428,59
13476368	10 to 300	0,20	1pcs	12030,09	11428,59
12617854	50 to 1200	1,00	1pcs	12030,09	11428,59

thermo scientific

STERILIN UNIVERSAL POLYSTYRENE CONTAINERS, 30ML

- 1/2-turn cap ensures improved end-user handling
- Available pre-filled with boric acid (0,4g) for the preservation of urine samples
- Suitable for centrifugation at 3800xg
- Aseptically manufactured
- In vitro use only
- OD x H, mm: 31,0 x 94,0
- Material container/cap: PS/PE
- 95kPa



GREAT SAVINGS

Cat No.	Description	Label	Pack qty	List price	Promo price
11339633	Quickstart universal	None	400pcs	995,00	776,10
11319153	Quickstart universal	Plain	400pcs	1000,00	778,16
11339133	Quickstart universal	Printed	400pcs	1010,00	632,55

fisherbrand

MICROCENTRIFUGES, MICRO 17/17R

- Capacity: 48mL (24 x 15/2mL tubes)
- Maximum speed/RCF: 13300rpm/17000xg
- Adjustable 99 minute timer and quick-spin feature for short runs
- Supplied with rotor and biocontainment lid



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
11526873	Microcentrifuge Micro 17R, refrigerated	1pcs	36770,00

fisherbrand

VORTEX MIXER

- Speed rpm 0 to 3000rpm
- Touch or continuous operation
- Technopolymer body for higher chemical resistance
- Orbit diameter: 4,5mm
- Dimensions, internal (w x d x h), mm: 153 x 170 x 125



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
13214789	Vortex mixer, speed 0 to 3000rpm	1pcs	1340,00	1017,86

thermo scientific

STERILIN POLYSTYRENE CONTAINERS, 7ML

- Aseptically manufactured
- Suitable for centrifugation at 7200xg
- Leak-tested in accordance with BS EN 14254 Annexe D
- Available pre-filled with boric acid (0,09g) for the preservation of urine samples
- Available with removable glass coverslip - ideal for work on the culture of intracellular organisms (e.g. Chlamydia trachomatis, Rickettsia, etc.)
- Material container/cap: PS/PE



GREAT SAVINGS

Cat No.	Description	Label	Pack qty	List price	Promo price
11399133	Bijou	Printed	700pcs	2800,00	2051,90

fisherbrand

MICROCENTRIFUGE GRADUATED TUBES

- Flat cap
- DNase, RNase free
- Autoclavable



GREAT SAVINGS

Cat No.	Capacity, mL	Colour	Type	Pack qty	List price	Promo price
11996955	0,6	Natural	Low retention	500pcs	358,00	272,17
11986955	1,5	Natural	Low retention	250pcs	84,65	62,49
11906965	2	Natural	Low retention	250pcs	77,30	60,76

fisherbrand

MICROCENTRIFUGE TUBES WITH LOCKING SNAP CAP

- Convenient, pierceable locking lids snap into place with one-handed operation
- Polypropylene (PP)
- Colour: Natural
- Sterile: No
- Graduated
- Flat and pierceable cap
- Locking lid to secure the sample and avoid undesired opening
- Marking area
- Autoclavable and can withstand temperatures to -86°C
- DNase/RNase and pyrogen free
- Easy one-hand use
- Applications include:
 - Sample preparation
 - Sample storage and transportation
 - Centrifugation



GREAT SAVINGS

Cat No.	Capacity, mL	Bottom shape	Max. Rcf (xg)	Diameter, mm	Length, mm	Pack qty	List price	Promo price
15422545	0,5	Conical	30000	10	30	500pcs	292,00	139,27
15432545	1,5	Conical	26000	13	39	500pcs	183,00	138,84
15442545	2	Round	25000	13	40	500pcs	251,00	160,08

fisherbrand

CONICAL CENTRIFUGE TUBES

- 15mL and 50mL Capacity
- Available in PP and PET
- Black graduation marks
- White marking area
- Flat cap and plug seal cap
- Gamma irradiation sterilisation
- Bulk or rack version



GREAT SAVINGS

Cat No.	Capacity, mL	Max. rcf, xg	Material	Dia. x H, mm	Sterile	Packaging	Pack qty	List price	Promo price
Cap flat top									
11755075	15	6000	PP	17 x 119	Yes	Bulk	500pcs	958,33	590,52
11512303	50	9400	PP	29 x 114	Yes	Bulk	500pcs	685,08	553,76
11849650	15	6000	PP	17 x 119	Yes	Rack	500pcs	834,46	613,32
11819650	50	9400	PP	29 x 114	Yes	Rack	500pcs	1215,64	713,64
Cap plug seal									
11765075	15	6000	PP	17 x 119	Yes	Bulk	500pcs	958,09	590,52
11809650	50	9400	PP	29 x 114	Yes	Bulk	500pcs	1210,49	699,11
11829650	50	9400	PP	29 x 114	No	Bulk	500pcs	968,39	566,96
11879640	15	1800	PET	17 x 119	Yes	Rack	500pcs	1045,65	613,32
11889640	15	6000	PP	17 x 119	Yes	Rack	500pcs	1061,11	609,52
11839650	50	1800	PET	29 x 114	Yes	Rack	500pcs	1421,28	858,04
11899640	50	9400	PP	29 x 114	Yes	Rack	500pcs	1205,33	703,76

fisherbrand

MINI TUBE ROTATOR

- LED indicator displaying speed and settings
- Speed range of 10 to 40rpm
- Includes assortment of nine paddles



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15534080	Mini tube rotator, 1,5mL to 50mL tubes	1pcs	2722,00	2310,98
Accessories				
15554140	Two paddles for 0,5mL/0,8mL tubes	2unit	102,00	83,71
15544140	Two paddles for 1,5mL/2mL tubes	2unit	102,00	83,71
15524140	Two paddles for 10mL/15mL tubes	2unit	102,00	83,71
15534140	Two paddles for 5mL/7mL tubes	2unit	102,00	83,71



CUTLERY, MULTIFUNCTION

- Traceability on each bag
- Sterile 10⁻³
- Material: polystyrene
- Colour: white



SPECIAL PRICE

Cat No.	Length, mm	Pack qty	List price	Promo price
11728623	180	500pcs	3620,00	2652,81

fisherbrand

Focus on Fisherbrand

for all your **weighing** needs



For full range refer to www.eu.fishersci.com/go/fisherbrand



fisherbrand

CELL STRAINERS

- Sterile, rapid, easy to use device for isolating primary cells to consistently obtain a uniform single-cell suspension from tissues
- Polypropylene frame with nylon mesh
- Sterile
- Individually packed
- Fits 50mL conical tubes
- Cell strainers are made of a strong nylon mesh available in three mesh sizes, 40µm, 70µm, and 100µm, for optimal performance in a variety of applications including stem cells and primary cells
- Fisherbrand cell strainers provide a fast, simple alternative to gauze filtration when dissociating cells from clumps or primary tissues



GREAT SAVINGS

Cat No.	Mesh size, µm	Pack qty	List price	Promo price
11587522	40	50pcs	1670,00	528,01
11597522	70	50pcs	590,00	528,01
11517532	100	50pcs	1055,00	528,01

fisherbrand

BALANCES, ANALYTICAL

- Weighing modes include weighing, parts counting and percentage weighing
- Large draft shield with sliding top and side doors
- Dimensions (w x d x h): 196mm x 320mm x 287mm
- Readability: 0,1mg
- Linearity: 0,2mg
- Repeatability: 0,1mg
- Two year warranty
- Communication interface: RS232



FEATURED PRODUCT

Cat No.	Capacity	Readability	Calibration	Pack qty	List price
15305113	210g	0,1mg	External	1pcs	9476,00
15395103	65g	0,1mg	Internal	1pcs	9753,00
15315113	210g	0,1mg	Internal	1pcs	11143,00



S-BLENDER-1 AND S-DILUTER-1



- S-Blender-1 is a high performance blender suitable for a range of sample dilutions between 50 and 400mL
- S-Diluter-1 increases your laboratory's productivity, helps to secure dilution operations and reduces the risk of contamination

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15589783	Gosselin™ S-Blender-1	1pcs	30550,00	27513,76
15599783	Gosselin™ S-Diluter-1	1pcs	42960,00	38690,38
11777919	Sterile bag PE Stomacher 400 full filter	1500pcs	8735,00	6132,76
11787919	Sterile bag PE Stomacher 400 lateral filter	1500pcs	8735,00	6132,76
12682255	Sterile bags PE blender	2500pcs	4255,00	3118,15
11522882	Blender Bag Clips	200pcs	4915,00	1294,22



MULTI PURPOSE TUBE ROTATOR



- Digital microprocessor control
- Six carousel options for maximum versatility
- Adjustable speed from 5-80rpm
- Adjustable angled head offers mixing angles from 95 to 180°
- Timer features both time up and time down counting modes
- Display will show last used settings, even after power has been turned off
- Two year warranty

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15524080	Multi purpose tube rotator	1pcs	3358,00	2711,72
Accessories				
15554090	50 place tube carousel for 1,5/2,0mL tubes	1pcs	4263,00	982,82
15564090	30 place tube carousel for 15mL tubes	1pcs	4263,00	982,82
15574090	15 place tube carousel for 50mL tubes	1pcs	4263,00	982,82
15584090	10 and 20 place tube carousel for 10-15mL and 5-7mL tubes	1pcs	4263,00	982,82
15594090	64 place drum tube carousel for 14mm dia. tubes	1pcs	4263,00	982,82
15504100	64 place drum tube carousel for 18,5mm dia. 15-20mL tubes	1pcs	4263,00	982,82



INCUBATORS, MICROBIOLOGICAL

Safe

- Safe view of samples through internal glass door, without impact on temperature
- Safe conditions with superior temperature uniformity
- Safe containment with automatic overtemperature alarm
- Alarm in case of temperature deviations

Easy

- Easy temperature setting with intuitive user interface
- Easy to remove 'One Click' shelf system
- Easy to clean interior with rounded corners
- Easy to read large display

Efficient

- Efficient small footprint to optimise bench top space
- Conveniently stackable without the need for tools or stacking kits
- Efficient utilisation of interior with flexible shelf system
- Incubators are supplied with two shelves



10% OFF

	10529070	10839810	10744262	12625306	12635306
Model No	IGS60	IGS100	IGS180	IGS400	IGS750
Capacity, L	75	117	194	400	750
Power rating, W	300	520	710	1200	1500
Current, A	3,7	2,3	3,1	5	6,5
Dimensions, external [w x d x h], mm	530 x 565 x 720	354 x 414 x 508	640 x 738 x 920	778 x 1545 x 770	1261 x 1545 x 770
Dimensions, internal [w x d x h], mm	354 x 414 x 508	464 x 414 x 608	464 x 589 x 708	544 x 1305 x 569	1044 x 1307 x 569
Number of shelves, supplied/max	2/13	2/16	2/19	2/39	2/39
List price	42432,34	17529,40	22444,90	38588,32	52806,55
Promo price	11189,11	15776,46	20200,41	34729,49	47525,90

thermo**scientific**

IMP REFRIGERATED INCUBATORS WITH PELTIER TECHNOLOGY

- Temperature range +5°C to +70°C and stability $\leq \pm 0,1^\circ\text{C}$
- Internal glass door for easy viewing of samples and access port as standard
- Programmable controller with timer function and over-temperature alarm



10% OFF

Cat No.	Model	Capacity, L	Pack qty	List price	Promo price
15608540	Heratherm IMP180	178	1pcs	51651,44	46486,30
15618540	Heratherm IMP400	381	1pcs	70433,78	63390,40

fisher**brand**

SHAKER, BOTTLE/TUBE ROLLER, VARIABLE SPEED, DIGITAL

- Adjustable speed to accommodate a variety of applications
- Stackable unit to save some bench space
- CO₂ resistant, optional remote control to operate while in incubator or biosafety cabinet
- Speed range 5 to 80rpm
- Temperature range 4°C to 60°C
- Six rollers supplied as standard, up to 11 can be used
- Continuous



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15376617	Bottle roller, global version	1pcs	8597,00	6609,43

Accessories

15386617	Additional rollers (pack of two)	2pk	275,00	226,90
15316627	Stacking kit	4set	619,00	477,14
15326627	Remote control	1pcs	756,00	582,74

stuart

INCUBATORS, SHAKING, SI SERIES

- Combined incubator shakers
- Digital display for temperature and speed
- Integrated timer
- Unique retractable platform for easy loading and unloading
- Angle adjustable accessory tube racks, with magnalock coupling system available
- Optimised for use with microtitre plates with an orbit of 1,5mm and speed range between 250rpm and 1250rpm
- With water trays to control the humidity levels in the chamber and minimise sample evaporation



GREAT SAVINGS

Cat No.	Temperature range, °C	Speed range, rpm	Pack qty	List price	Promo price
11986458	+5 to +60	30 to 300	1pcs	25530,00	24313,17

Accessories

10774337	Tube rack for 64 x 1,5mL microtubes	1pcs	4330,00	
10470295	Tube rack for 25 x 15mL microtubes	1pcs	4330,00	
10784337	Tube rack for 12 x 50mL microtubes	1pcs	4315,00	
10794337	Tube rack for 16 x 30mL microtubes	1pcs	4330,00	

thermo**scientific**

COMPACT MICROBIOLOGICAL INCUBATORS HERATHERM

- Minimal footprint for space restricted lab areas
- Temperatures at or below ambient
- High temperature accuracy
- Internal light facilitates sample observation
- Convection technology: mechanical convection
- Temperature range °C 17 - 40
- Spatial temperature deviation at 37 °C: $\pm 1,2^\circ\text{C}$
- Temperature deviation over time at 37 °C: $\pm 0,2^\circ\text{C}$
- Footprint m²: 0,12
- Dimensions
 - Interior, mm (W x H x D): 180 x 310 x 290
 - Exterior, mm (W x H x D): 260 x 415 x 4705
- Number of shelves supplied/max: 2/3
- Energy consumption at 37 °C: 14W



10% OFF

Cat No.	Capacity, L	Pack qty	List price	Promo price
10699640	18	1pcs	6796,08	6116,47



POLYPROPYLENE (PP) CENTRIFUGE BOTTLES: TRANSLUCENT

- High performance silicone O-ring caps are designed for high speed applications
- Standard screw caps are designed for low speed applications
- Bottles have extra-thick walls for durability and wide-mouth design for easy liquid handling
- Nonsterile
- Clear graduations and textured marking spot
- Autoclavable at 121°C, 15 psi for 20 minutes



GREAT SAVINGS

Cat No.	Capacity, mL	Max force, xg	Material	Pack qty	List price	Promo price
15158259	250	27500	Polypropylene	36pcs	4560,73	2998,87
15118219	250	27500	Polycarbonate	36pcs	4629,94	3044,37
15180835	250	13500	Polypropylene	36pcs	2286,56	1503,47

fisherbrand

**PETRI DISHES,
PLASTIC, ROUND**



- Inner packs of 20

GREAT SAVINGS

Cat No.	Ø x H, mm	Vent	Pack qty	List price	Promo price
Aseptic					
12694785	90 x 14,2	3	600pcs	725,00	408,79
12664785	90 x 16,2	3	600pcs	735,00	411,64
12644785	90 x 16,2	0	600pcs	780,00	478,80
12654785	90 x 16,2	1	600pcs	895,00	411,64
Sterile					
12604795	90 x 14,2	3	600pcs	810,00	532,57
12674785	90 x 16,2	3	600pcs	1275,00	561,93
12684785	90 x 16,2	1	600pcs	885,00	561,93

Scienceware
BET-BART PRODUCTS

PETRI DISH RACK

- Each rack has six columns and holds up to fifty-four petri dishes with lids
- Large centre divider has space for labelling
- Includes carrying handle that does not interfere with samples
- Up to four racks can be stacked and interlocked
- Finger access slots in base allow safe removal of stack of dishes
- Rack base is epoxy-coated steel, and all other parts are made of white polypropylene
- Racks are shipped flat and are easily assembled in minutes
- Steam autoclavable at 121°C



GREAT SAVINGS

Cat No.	Dimensions, L x W x H, mm	Pack qty	List price	Promo price
11745234	267 x 171 x 171	1pcs	945,00	555,90

Scienceware
BET-BART PRODUCTS

**MOLDED CONTACT
PLATE RACK**

- Compatible with RODAC™ and BD™ plates
- Rack is configured in 6 columns with bottom access slots for easy removal of plates from the rack
- Centre divider has a built-in handle and ample space for labelling
- Racks interlock for safe stacking up to four high
- Made from polypropylene, rack is autoclavable at 121°C and assembles easily



SPECIAL PRICE

Cat No.	Dimensions, L x W x H, mm	Positions	Pack qty	List price	Promo price
10566384	266 x 165 x 178	60	1pcs	870,00	683,92
13383599	330 x 210 x 178	42	1pcs	770,00	604,79

thermoscientific

**PETRI DISHES,
PLASTIC, ROUND**



- Inner packs of 20

20% OFF

Cat No.	Type	Dish diameter, mm	Pack qty	List price	Promo price
11769252	Clear Triple Vent Petri Dishes	90	500pcs	840,00	672,00
15370366	Aseptic standard Petri Dishes	90	500pcs	865,00	674,70
11339283	Triple Vent Petri Dishes	140	80pcs	1060,00	826,80
11339273	Square Petri Dish	100	120pcs	3065,00	2390,00

Scienceware
BET-BART PRODUCTS

**BACTERIAL CELL
SPREADERS**

- Autoclavable and may be flame sterilised
- Features triangular end with smooth rounded surfaces throughout
- Smaller 16mm width bar suitable for use with 6cm plates
- Larger 25mm width bar suitable for use with 9cm plates
- Both sizes 30cm long
- Autoclavable
- Sterile



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
11336954	Bacterial Cell Spreaders	1pcs	261,75	138,32

Scienceware
BET-BART PRODUCTS

**STACKABLE PETRI DISH
INCUBATION TRAY**

- Holds five 100mm or fifteen 60mm Petri dishes
- Made of polypropylene
- Includes two label areas for sample identification and holes at bottom to facilitate air flow
- Not recommended above 100°C; autoclaving may warp
- Disinfected using 70% ethanol or laboratory disinfectants



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
11734736	Stackable Petri Dish Incubation Tray	3pcs	780,00	547,40



SPECIFIC CULTURE MEDIA

- Laminin is a major component of basement membranes. It has numerous biological activities including promotion of cell adhesion, migration, chemotaxis, growth, and differentiation, including neurite outgrowth
- Corning™ basic Fibroblast Growth Factor (bFGF), human recombinant, 10µg, is a heparin-binding mitogenic protein that enhances the proliferation of a wide variety of cell types under serum-free or serum-reduced conditions
- The Corning Synthemax II-SC Substrate is an easy-to-use, self-coating material for creating a unique synthetic surface that mimics the natural cell environment
- Corning Dispase, 100mL, is a bacillus-derived neutral metalloprotease recommended for recovering cells cultured on Corning Matrigel Basement Membrane Matrix
- Corning Matrigel matrix is a solubilised basement membrane preparation extracted from the Engelbreth-Holm-Swarm (EHS) mouse sarcoma, a tumour rich in extracellular matrix proteins, including Laminin (a major component), Collagen IV, heparin sulfate proteoglycans, entactin/nidogen, and a number of growth factors



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
10256312	Laminin/Entactin High Concentration Extracellular Matrix Protein	10,50mg	5385,97	4150,42
11563530	Basic Fibroblast Growth Factors (bFGF), Human Recombinant	10ug	1718,87	1608,38
15319397	Synthemax™ II-SC Substrate	10mg	5291,42	3786,27
11553550	Dispase	100mL	4196,11	3233,54
11573560	Matrigel™ hESC-Qualified Matrix, *LDEV-Free	5mL	3569,99	2751,01



COOLRACK™ M MICROCENTRIFUGE TUBE MODULES

- Place the CoolRack module directly onto a temperature source between -196°C to >100°C and it will rapidly adapt to that temperature.
- Ensures ±0°C temperature uniformity across all tubes when cooling, snap freezing, heating or thawing
- All CoolRack modules may be autoclaved, high heat sterilised or decontaminated with bleach, alcohol or other disinfectants or lab detergents
- Suggested applications include cooling reagents such as restriction enzymes, dNTPs and antibodies, alcohol-free dry ice snap freezing of tissue, virus and bacteria samples and bench top cryogenic tube sorting in liquid nitrogen



SPECIAL PRICE

Cat No.	Compartment	Format	Pack qty	List price	Promo price
15592801	6	2 x 3	1pcs	674,74	529,53
15522811	15	3 x 5	1pcs	1629,64	1278,93
15562811	30	5 x 6	1pcs	2875,21	2256,45
15592811	90	6 x 15	1pcs	4328,38	3396,90



LB MEDIUM, DEHYDRATED

- Formulations for the E. coli growth media known as LB Medium typically contain 10g of tryptone and 5g of yeast extract per litre, and LB Medium contains 10g NaCl per litre

FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
11343229	Capsules	1kg	3250,00
11363229	Large Capsules	1kg	6282,00



DIGITAL DRY BATHS

- Digital control which displays time and temperature (+5°C to +130°C)
- Built-in temperature-sensing probe for improved temperature accuracy and control
- Interchangeable aluminium alloy blocks provide versatility and easy cleaning (blocks available separately*)



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
15397928	Dual block, overall dimensions 318 x 200 x 100mm	1pcs	3585,00	2756,17

* Blocks are available separately. For further information, please see our website: www.eu.fishersci.com



COOLBOX™ XT SYSTEMS

- Accommodates variety of sample formats and temperatures
- Modular design holds wide array of sample tubes and plates
- Enables cooling up to 16 hours or freezing up to 8 hours with one Corning XT Cooling Core
- Reduces risk of ice contamination
- Available in single- or double-capacity systems
- Internal cooling or freezing cores provide hours of 0,5°C to 4°C cooling or -20°C to 0°C sample freezing
- Can be used with dry ice



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15582791	Holder and one cooling core	1pcs	5470,15	4292,96
15502801	Holder and two cooling cores	1pcs	8428,33	6614,53

fisherbrand

ADVANCED IR VORTEX MIXER



- Speed rpm 0 to 3000rpm
- Infrared sensor, continuous
- Zinc alloy and technopolymer
- Orbit diameter: 4,5mm
- Dimensions, internal (w x d x h), mm: 150 x 165 x 130

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
13284769	Advanced IR Vortex Mixer	1pcs	1405,00	1212,21
Accessories				
12365107	Foam attachment for no. 19 micro vials and 1,5mL Eppendorf™ tubes	1pcs	510,00	374,01
11792583	Foam attachment for microtitre plate	1pcs	605,00	462,59

MERCK

TRANSFER MEMBRANES, PVDF, BLOTTING, IMMOBILON™

- Immobilon™-P
 - For general immunodetection and amino acid analysis
 - Rapid immunodetection method reduces detection time by two hours
 - High mechanical strength for easy handling
- Immobilon™-FL
 - The first transfer membrane optimised for fluorescence applications
 - Extremely low background
 - Compatible with all commonly used fluorescent probes
 - Can be used at all excitation and emission wavelengths
 - Ideal for multiplexing



GREAT SAVINGS

Cat No.	Description	Pore size, µm	Dimensions, mm	Pack qty	List price	Promo price
10344661	Immobilon™-P	0,45	265 x 3750	1roll	6600,00	3404,80
10452792	Immobilon™-FL	0,45	265 x 3750	1roll	4421,79	3920,65



Distributor
GE Healthcare

VACU-GUARD™ FILTER



- Protects pump systems from solvent vapour or gaseous contaminants and aqueous aerosols
- Features a hydrophobic PTFE membrane
- Designed for in-line use with stepped barb connections for 10 to 12mm ID hose

SPECIAL PRICE

Cat No.	Diameter, mm	Pack qty	List price	Promo price
11364634	47	100pcs	1455,63	1325,00



Distributor
GE Healthcare

WHATMAN™ ME MEMBRANE FILTERS



- Made from mixed cellulose esters
- Economical
- Specifically for aqueous solutions
- Hydrophilic
- For use up to 125°C
- Sterile, individually packed
- All membrane filters have a high-contrast grid (3,1mm)
- Excellent contrast for easier particle detection
- Grids are nontoxic and do not inhibit bacterial growth, ensuring sample integrity


9% OFF

Cat No.	Pore size, µm	Diameter, mm	Pack qty	List price	Promo price
11395364	0,45	50	10pcs	1278,04	1163,00



Distributor
GE Healthcare

WHATMAN™ NUCLEPORE™ POLYCARBONATE HYDROPHILIC MEMBRANES



- Offers a highly defined pore size cutoff in the sub-micron range, enabling particle deformity measurement and capture of extremely small particles.
- Track-etched for highly defined pore size and high pore density
- Treated with PVP to render it hydrophilic
- Low extractables, minimising sample contamination risk
- High chemical resistance and good thermal stability for a wide range of sample types
- Low, consistent ash and low tare weights
- Smooth flat surface for high particle visibility


SPECIAL PRICE

Cat No.	Pore size, µm	Diameter, mm	Pack qty	List price	Promo price
11342895	0,2	47	100pcs	1054,25	959,00

thermo scientific

CULTI-LOOPSTM

- Gain efficiencies in qualitative quality control of clinical, food, and pharmaceutical products. Thermo Scientific Culti-Loops are ready-to-use, disposable bacteriological loops containing stabilised, viable microorganisms in a unique gel matrix for the control of culture media and microbiology testing protocols
- Convenience: decrease time and unnecessary costs of laboratory resources with a direct-streak inoculating loop that requires no rehydration
- Performance: obtain enhanced viability and recovery with proprietary gel-preserved organisms
- Confidence: reach complete traceability with ATCC™ License Derivative strains



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12938858	Culti-Loops™ Bacillus subtilis ATCC® 6633™	5	736,58
12637639	Culti-Loops™ Candida albicans ATCC® 10231™	5	1065,12
11973272	Culti-Loops™ Enterococcus faecalis ATCC® 29212™	5	736,58
11984082	Culti-Loops™ Escherichia coli ATCC® 25922™	5	632,00
12637669	Culti-Loops™ Escherichia coli ATCC® 8739™	5	842,00
12998988	Culti-Loops™ Klebsiella pneumoniae subsp. Pneumoniae ATCC® 700603™	5	968,68
11983272	Culti-Loops™ Pseudomonas aeruginosa ATCC® 27853™	5	632,00
12905743	Culti-Loops™ Staphylococcus aureus subsp. aureus ATCC® 25923™	5	736,58
11934082	Culti-Loops™ Staphylococcus aureus subsp. aureus ATCC® 29213™	5	736,58
12985893	Culti-Loops™ Streptococcus pneumoniae ATCC® 49619	5	1429,70

thermo scientific

QUANTI-CULT PLUS™

- Ensure reproducible control of culture media for pharmaceutical materials with quantitative quality control test microorganisms. Thermo Scientific Quanti-Cult Plus provides specific and reproducible numbers of viable microorganisms, derived from authentic high-quality ATCC™ cultures, in a safe and ready-to-rehydrate vial. Each vial delivers up to 10 inocula of 0,1mL, each containing <100 CFU (colony forming units)
- Safety: reduce the risk of contamination and infection with a ready-to-rehydrate, two-vial system that requires zero organism handling.
- Performance: obtain enhanced viability and recovery with proprietary gel-preserved organisms that provide consistent counts for up to 8 hours after rehydration when stored at 2–8°C
- Confidence: reach complete traceability with ATCC™ License Derivative strains



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
12965993	Quanti-Cult Plus™ Aspergillus brasiliensis ATCC® 16404™	100 tests	5080,00
12975993	Quanti-Cult Plus™ Bacillus subtilis ATCC® 6633™	100 tests	5080,00
12995993	Quanti-Cult Plus™ Candida albicans ATCC® 10231™	100 tests	5080,00
12926003	Quanti-Cult Plus™ Clostridium sporogenes ATCC® 11437™	100 tests	5080,00
13137876	Quanti-Cult Plus™ Clostridium sporogenes ATCC® 19404™	100 tests	5413,13
12976037	Quanti-Cult Plus™ Escherichia coli ATCC® 8739™	100 tests	4148,10
12956027	Quanti-Cult Plus™ Pseudomonas aeruginosa ATCC® 9027™	100 tests	5080,00



AXYGEN™ MINI PLATE SPINNER CENTRIFUGE

- The Axysen Axyspin Plate Spinner is the first centrifuge designed specifically for quick and easy centrifugation of samples in PCR plates or microplates.
- Capacity for two standard plates
- Top loading design
- Exceptionally small footprint
- Starts and stops in seconds
- Use before and after thermal cycling to improve PCR yield
- Quickly spin down droplets and condensation
- Accepts skirted and non-skirted and all standard PCR plates
- Less than 1/4 the size of most plate centrifuges
- Maximum capacity: 2 standard PCR or microtitre microplates
- Ambient operating range 4° to 35°C



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15582995	Mini Plate Spinner Centrifuge, 230V, EU Plug	1pcs	4450,38	3691,00



MICROTITRE SHAKER, MINI

- High speed, small orbital action – ideal for microtitre plates
- Capacity for four or eight microtitre plates
- Digital selection of speed
- Built in digital timer
- Accessories available for mixing microcentrifuge tubes
- With BioCote™ antimicrobial protection
- Plate dimensions, mm: 220 x 220
- Plate positions: 4
- Orbit diameter: 1,5mm



SPECIAL PRICE

Cat No.	Temperature range, °C	Speed range, rpm	Pack qty	List price	Promo price
11486608	+4 to 40	250 to 1200	1pcs	9370,00	8922,94
Accessories					
11496608	Tube holder 1,5mL		1pcs	1225,00	1124,74
12944025	Tube holder 0,5mL		1pcs	1225,00	1124,74
12693335	Tube holder 0,2mL		1pcs	1840,00	1115,37



NOVASPEC PRO SPECTROPHOTOMETER

- Designed for general biotech laboratory use, this spectrophotometer is lightweight and easy to use. The Biochrom™ Novaspec Pro Spectrophotometer features a small footprint and rugged lightweight construction, ideal for laboratories where space is at a minimum
- Performs measurements in absorbance, % transmission and concentration
- Includes sorted methods for Bradford, BCA, Biuret and Lowry protein quantification
- Automatic calibration at start-up
- No lid or moving parts, this robust instrument is easy to use, clean and store
- Graphical display enables wavelength scans, kinetic assays and standard curves to be viewed
- Increase measurement flexibility with the optional test tube cell holder which allows measurement of 10 to 18mm test tubes
- The standard cell holder accepts 10mm glass and disposable cuvettes and is easily removed for cleaning
- Optional thermostatted cell holder is available for the added capability of measurements at 37°C
- Visible spectrophotometer covers 325 to 1100nm allowing you to see more of your sample



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15591166	Biochrom™ Novaspec Pro Spectrophotometer	1pcs	26939,49



WRITE-ON™ PAPER LABEL TAPES

- Use pencil, ball-point pen or solvent ink marker to create labels that remain legible under the most demanding laboratory conditions, paper label tapes resist moisture and most solvents and are great for use in low temperature applications



GREAT SAVINGS

Cat No.	Colour	Dimensions,m	Pack qty	List price	Promo price
15196324	White, Red, Yellow, Green, Blue, Orange	13,7	6pcs	890,00	699,78
13455226	White	13,7	4pcs	277,00	203,72
13405236	White	13,7	3pcs	274,00	189,50
13474166	White	36,5	6pcs	855,00	518,62
13475236	White	36,5	4pcs	600,00	459,11
13435246	White	36,5	3pcs	636,94	336,48
13485236	Red	36,5	4pcs	550,00	459,11
13435236	Red	36,5	6pcs	690,00	518,62
13445246	Red	36,5	3pcs	560,00	430,66
13434720	Yellow	36,5	6pcs	645,00	518,62
13455246	Yellow	36,5	3pcs	530,00	430,66
13495236	Yellow	36,5	4pcs	565,00	459,11
13415246	Green	36,5	4pcs	575,00	459,11
13475246	Blue	36,5	3pcs	775,00	430,66
13405246	Blue	36,5	4pcs	555,00	459,11
13465236	Blue	36,5	6pcs	605,00	518,62
13455236	Orange	36,5	6pcs	665,00	518,62
13425246	Orange	36,5	4pcs	660,00	459,11
13465246	Green	36,5	3pcs	636,94	336,48
13485246	Orange	36,5	3pcs	605,00	430,66
13495246	White	54,8	6pcs	895,00	703,87
10045180	White	54,8	4pcs	720,00	565,06

fisherbrand CUVETTES, POLYSTYRENE



- Material: PS
- Pathlength, mm: 10

GREAT SAVINGS

Cat No.	Capacity, µL	Description	Pack qty	List price	Promo price
11847832	1600	Cuvettes, visible	1000pcs	775,00	497,14
11837832	4000	Cuvettes, visible	1000pcs	740,00	503,41
11537692	1600	Cuvettes, disposable visible	100pcs	83,50	47,50
11547692	4000	Cuvettes, disposable visible	100pcs	83,50	47,50



BAG HOLDERS, BIOHAZARD, CLAVIES™

- Recommended for use with the clear, unprinted and standard biohazard bag ranges, these biohazard bag holders ensure the bag is securely supported on a stainless steel wire frame and polypropylene bottom. The bottom also serves as a containment tray in the event of a bag leakage. The tray is orange for easy identification as a biohazard containment system and can be autoclaved at 121°C to maintain sterility. For ease of transportation they are supplied unassembled, however assembly only takes a few minutes and once done becomes permanent



GREAT SAVINGS

Cat No.	Tray dimensions, mm	Height, mm	Bag size, mm	Pack qty	List price	Promo price
11322154	185 x 356	470	300 x 610	1pcs	675,07	357,08
11787864	356 x 356	546	610 x 760	1pcs	660,00	514,53
11797864	356 x 356	696	610 x 910	1pcs	720,00	562,44



MULTI-ROLL LABEL TAPE DISPENSER

- Provide smooth writing surface for multiple tapes of different colours and width, making labelling easy and fast
- Three plastic separating discs to assure smooth, snag free dispensing for up to four rolls of tape
- Includes one free roll of 12,7mm white Write-On™ labelling tape



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
10073213	Label Tape Dispenser	1pcs	417,00	341,88

thermo**scientific**

COLIFORM WATER BATHS PRECISION™

- Temperature controller provides $\pm 0,05^{\circ}\text{C}$ uniformity and $\pm 0,1^{\circ}\text{C}$ stability with stainless-steel gable cover
- Achieve outstanding temperature uniformity with perimeter-directed water flow
- Easily clean and maintain bath with coil-free design
- Optimise scheduling with auto-on and auto-off timers
- Help prevent bath overheating with low-fluid protection
- Easily operate and monitor with icon-based graphical display
- Accommodate taller labware with new extended-height, hinged lid
- Protect your work with audible alarms



FEATURED PRODUCT

Cat No.	Capacity, L	Heating capacity, W	Temperature Range, °C	Pack qty	List price
15305877	19	1200	5 to 100	1pcs	20040,00
15315877	35	1500	5 to 100	1pcs	22600,00

thermo**scientific**

CIRCULATING WATER BATHS PRECISION™

- Achieve enhanced temperature uniformity with perimeter-directed water flow
- Easily clean and maintain bath with coil-free internal design
- Optimise scheduling with auto-on and auto-off timers
- Accommodate taller labware with new hinged lid and extended height
- Help prevent bath overheating and damage with low-fluid protection
- Easily operate and monitor with icon-based graphical display
- Protect your work with audible alarms



FEATURED PRODUCT

Cat No.	Capacity, L	Heating capacity, W	Temperature Range, °C	Pack qty	List price
15345867	35	1500	5 to 100	1pcs	20015,00
15355867	89	1500	5 to 100	1pcs	25610,00

Scienceware
BEL-ART PRODUCTS

THERMOMETER, MAXIMUM REGISTERING, BI-METAL, DUAL SCALE

- Determine the highest temperature attained without dangerous mercury or inconvenient computer connections
- This durable bi-metal thermometer has a maximum registering pointer (red) that is moved up by the temperature measuring pointer (black) to indicate the highest temperature achieved, and remains there until manually reset.
- Threaded connection (1/4" NPT) screws directly into in-line processes; can be used to mount thermometer on included stand
- Waterproof and pressurised; withstands autoclave cycles and has a removable protective silicone sleeve to prevent moisture from making contact with the dial
- Constructed from durable 304 stainless steel with a glass lens; withstands pressure, high temperatures and rust
- Probe length: 50mm



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15519971	Thermometer, maximum registering, bi-metal, dual scale	1pcs	347,00

memmert
Experts in Thermostatics

BATHS, WATER, HEATED, WNE SERIES

- Fuzzy supported PID microprocessor controller
- Digital display
- 2 x PT100 sensors
- Calibration facility on controller
- Audio and visual alarm
- Over-temperature protection
- Shaking device option



5 % OFF

Cat No.	Capacity, L	Power, W	Dimensions, external, W x D x H, mm	Pack qty	List price	Promo price
10778873	7	1200	365 x 468 x 337	1pcs	6207,56	5897,18
10024092	10	1200	356 x 578 x 337	1pcs	6733,47	6396,80
10633913	14	1800	456 x 578 x 347	1pcs	7268,77	6905,33
10195892	22	2000	436 x 578 x 405	1pcs	8048,23	7645,82
10544463	29	2400	516 x 818 x 343	1pcs	9137,61	8680,73
10105282	45	2800	516 x 818 x 401	1pcs	11147,32	10589,95

fisherbrand

TABLETOP AUTOCLAVES

- Microprocessor control system provides nine fully automatic programs, including add water, sterilisation and dry
- Eight pre-installed sterilisation cycles; custom cycles are not an option for the SterilElite16 or 24
- Two separate water-level check mechanisms: water tank level and chamber water level
- Auto dry function
- Pressure door auto-lock device
- Chamber protected by pressure switch and steam pressure safety valve
- Vacuum/pressure release button cuts off the power supply, exhausts the chamber of pressure, and releases the chamber vacuum simultaneously



FEATURED PRODUCT

Cat No.	Capacity, L	Pack qty	List price
15515199	16	1 unit	32382,44
15525199	24	1 unit	41179,98



PLATE READER

- For examining cell monolayer, cell clones and for counting viral plaques or bacterial colonies
- Easy and safe reading of hemagglutination, hemagglutination inhibition, hemolysis and agglutination tests
- Clear deck for resting plates
- Large, adjustable mirror is perfect for viewing plate underside
- Mirror magnification is 3X actual plate size
- Made of polypropylene with a glass-clear acrylic deck for placement of plates and petri dishes



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
11328763	Microtitre Plate Reader	1 pcs	2745,13	1450,41



REFRIGERATORS AND FREEZERS, LABORATORY, GPS SERIES

- GPS series: medium and large capacity refrigerators and freezers for routine sample protection

Cat No.	Description	Capacity, L	Dimensions, mm	Pack qty	List price
12664527	GPS series freezer, auto defrost, 230V, 50Hz	400L	600 x 642 x 1822	1 pcs	16763,24
12629275	GPS series refrigerator, glass door, auto defrost, 230V, 50Hz	400L	600 x 642 x 1822	1 pcs	19824,76
12639275	GPS series refrigerator, solid door, auto defrost, 230V, 50Hz	400L	600 x 642 x 1822	1 pcs	15504,82
12684527	GPS series freezer, auto defrost, 230V, 50Hz	700L	700 x 817 x 1972	1 pcs	19618,16
12649275	GPS series refrigerator, glass door, auto defrost, 230V, 50Hz	700L	700 x 817 x 1972	1 pcs	23055,32
12659275	GPS series refrigerator, solid door, auto defrost, 230V, 50Hz	700L	700 x 817 x 1972	1 pcs	17176,45
12908794	GPS series freezer, auto defrost, 230V, 50Hz	1400L	1400 x 817 x 1972	1 pcs	27262,57
12647286	GPS series refrigerator, glass door, auto defrost, 230V, 50Hz	1400L	1400 x 817 x 1972	1 pcs	32333,80
12627286	GPS series refrigerator, solid door, auto defrost, 230V, 50Hz	1400L	1400 x 817 x 1972	1 pcs	26040,00

FEATURED PRODUCT



BIOHAZARD DISPOSAL BAGS

- Sterilisation patch darkens when subjected to steam sterilisation
- High temperature steam autoclavable to 135°C
- 0,3mm thick polypropylene, compatible with autoclave temperatures; less likely to melt/deform like PE based bags
- Pass ASTM 1922 Tear Resistance and ASTM 1709 Dart Impact Tests
- Include black biohazard symbol and precautionary procedures printed in four languages (English, Spanish, French, German)



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
11392154	Biohazard bag, red, 360mm x 480mm	200 pcs	1950,97	1030,70
11302164	Biohazard bag, red, 480mm x 580mm	200 pcs	2842,14	1501,58
11312164	Biohazard bag, red, 640mm x 890mm	200 pcs	4979,59	2630,72
11737927	Biohazard bag, red, 790mm x 970mm	200 pcs	4625,70	3595,69
10296674	Biohazard bag, red, 970mm x 1220mm	100 pcs	3715,00	2923,87



DEODORANT AUTOCLAVE PADS ODOCLAVE™

- Heat-release deodorant works automatically at autoclave temperatures
- Open pouch and place in bag prior to autoclaving
- Three-year shelf life or longer if stored in a sealed container
- 32mm square



SPECIAL PRICE

Cat No.	Fragrance	Pack qty	List price	Promo price
11966072	Pine	100 pcs	860,00	705,50



OUR CHEMICAL OFFERING GETS BIGGER



Whatever your project, we have a solution

Analysis



Research and development



A Thermo Fisher Scientific Brand



Life science applications



Production



fisherbrand

LABORATORY COATS, UNISEX



- Manufactured using a quality pre-shrunk polycotton fabric, these coats are suitable for industrial laundering at temperatures up to 85°C. Features include: two lower patch pockets, breast pocket with integral pen pocket, concealed stud fastened front, single rear vent for ease of movement. Not CE marked.
- Material: 65% Polyester/35% Cotton

UP TO 41% OFF

Cat No.	Size	Colour	Length, cm	Pack qty	List price	Promo price
11558102	S	White	132	1pcs	323,00	197,55
11568102	M	White	132	1pcs	297,84	213,31
11578102	L	White	132	1pcs	323,00	252,10
11588102	XL	White	132	1pcs	322,00	236,22
11598102	XXL	White	132	1pcs	318,00	187,62
11508112	XXXL	White	132	1pcs	319,00	221,45



COVERALLS, DUPONT™ EASYSAFE



- A certified Type 5/6 chemical protective coverall
- Based on a new optimised polyethylene non-woven fabric at low cost to meet budget requirements:
 - softer feel, touch the garment and feel the difference!
 - great balance of protection and comfort
- Optimised design
- Optimised packaging

SPECIAL PRICE

Cat No.	Size	Pack qty	List price	Promo price
15310541	S	1pcs	46,00	36,16
15320541	M	1pcs	45,85	45,85
15330541	L	1pcs	46,00	36,16
15340541	XL	1pcs	46,00	36,16
15350541	XXL	1pcs	46,00	36,20
15360541	XXXL	1pcs	46,00	36,16



A7 P+™ LABORATORY COATS



- Cat III PPE labcoat - warranted protection against chemical splashes
- Thumb loops to increase wrist protection and wrist coverage between gloves and labcoat
- Light SMS tissue to avoid heat and discomfort during work

SPECIAL PRICE

Cat No.	Size	Pack qty	List price	Promo price
13529020	S	15pcs	855,00	665,33
13549020	M	15pcs	855,00	665,33
13559020	L	15pcs	855,00	704,09
13569020	XL	15pcs	855,00	855,00
13589020	XXL	15pcs	855,00	665,33



SPECTACLES, SAFETY, 3M™ SECUREFIT™ 200 SERIES



- Design provides excellent coverage and good field of vision
- Offers excellent protection against ultraviolet (UV) radiation
- Lightweight and comfortable design (18g)
- Good impact resistance with anti-scratch coating for improved durability
- Special anti-fog coating which reduces fogging and ensures clearer vision
- 3M Pressure Diffusion Technology allows temples to naturally adjust to individual head sizes
- Self adjusting feature eliminates need to stock multiple frame sizes at worksites
- Secure, snug fit keeps eyewear in place without slippage even when the person is in motion
- Combines well with 3M™ 9300 Series respirators

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15282025	SecureFit™, clear, anti-scratch/anti-fog	1pcs	80,50	66,41
15292025	SecureFit™, grey, anti-scratch/anti-fog	1pcs	86,80	71,61
15202035	SecureFit™, amber, anti-scratch/anti-fog	1pcs	82,60	68,15



SECUREFIT SF400 SERIES SAFETY SPECTACLES



- Optical class 1 lens suitable for prolonged use
- Design provides excellent coverage and field of vision
- Protection against UV radiation
- Lightweight design (19g)
- Adjustable nose pads for secure and personalised fit
- Soft-touch dual moulded temples for increased comfort and grip

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15371994	SF401AF-EU, Clear PC	1pcs	72,27	59,62
15391994	SF402AF-EU, Grey PC	1pcs	72,52	59,83
15302004	SF403AF-EU, Amber PC	1pcs	72,52	59,83
15322004	SF408AS-EU, Mirror blue PC	1pcs	82,88	68,38
15312004	SF410AS-EU, Mirror I/O PC	1pcs	82,88	68,38



GOGGLEGEAR™ 500 SERIES SAFETY GOGGLES



- Can be disinfected with diluted bleach soaking or alcohol wipes without losing its anti-fog performance
- Offers excellent protection against UV radiation
- Adjustable elastic strap
- Indirect ventilation
- Low profile frame/shroud
- Polycarbonate lens

SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15603074	GoggleGear™ 500 Series Safety Goggles	1pcs	174,34	170,00
15693064	GoggleGear™ 500 Series Prescription Insert	1pcs	87,17	82,00



DISPOSABLE POLYPROPYLENE FACEMASKS

- Provides advanced protection
- Disposable, single use
- Choice of ear loops or tie closure
- Nonwoven polypropylene outer layer
- Thermally bonded nonwoven polypropylene inner layer
- Latex-free
- Choice of blue or white
- Universal size/fit



GREAT SAVINGS

Cat No.	Description	Colour	Thickness	Pack qty	List price	Promo price
15451932	Disposable face mask with ear loops	White	2-ply	500pcs	251,00	178,70
15461932	Disposable face mask with ear loops	Blue	3-ply	500pcs	319,00	287,20
15471932	Disposable face mask with ear loops	White	3-ply	500pcs	319,00	227,11



VERSAFLO™ POWERED AIR TURBO TR-600

- Comfortable airflow - three user selectable flow rates to maximise comfort
- Ease of use - no calibration, LED battery/filter indicator, audio vibration alarm
- Multiple applications - easy clean, suitable for shower decontamination and temporary submersion
- Safety - approved to EN 12941 Class TH2 or TH3 dependent on headtop
- The 3M™ Versaflo™ Powered Air Turbo TR-600 is a belt-mounted air purifying device. When combined with one of the 3M approved headtops, it forms a powerassisted system for respiratory protection against particles, nuisance odours (where applicable) and gases and vapours. The product complies with European EMC Directive 2004/108/EC



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15529305	Turbo unit, includes airflow indicator and spark arrestors	1pcs	5383,07	4441,03
15559255	Starter kit, includes everything required to get started	1pcs	12278,19	10129,51
15549305	Standard battery	1pcs	2395,74	1976,49
15589255	High capacity battery	1pcs	2950,74	2434,36
15599255	Battery charger cradle	1pcs	1169,89	965,16
15509265	Single station battery charger kit (Europe/UK)	1pcs	1784,87	1472,52
15519265	4-Station battery charger kit (Europe/UK)	1pcs	7899,53	6517,11
15559305	Particulate filter	5pcs	1423,98	1174,78
15539295	A1P filter - Organic gases and vapours	5pcs	2735,19	2256,53
15549295	ABE1P filter - Organic gases and vapours	5pcs	3464,27	2858,02
15559295	A2P filter - Organic gases and vapours	5pcs	3003,80	2478,14
15569295	ABE2K1HgP filter - Organic gases and vapours	5pcs	4119,40	3398,51
15558811	6051i indicator filter	8pcs	622,77	513,79
15568811	6055i indicator filter	8pcs	790,73	652,35



AURA™ PARTICULATE RESPIRATORS, SERIES 9300

- Unique three-panel design is lightweight and foldable
- Soft cover web on inner panel ensures greater comfort, evenly tensioned headbands relieve strap pressure
- Centrally positioned cool flow valve reduces heat build up and makes breathing easier and cooler



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
11717574	9310+ Premium dust respirator, provides protection against fine dusts. Individually packed, EN 149:2001 FFP1. Particulates Assigned Protection Factor 4 (NPF 4)	20pcs	615,00	507,38
12672545	9312+ Premium valved dust respirator, provides protection against fine dusts. Individually packed, EN 149:2001 FFP1. Particulates Assigned Protection Factor: 4 (NPF 4)	10pcs	307,00	253,28
10583051	9320 Premium dust/mist respirator, provides protection against fine dusts and water based mists. Individually packed, EN 149:2001 FFP2. Particulates Assigned Protection Factor: 10 (NPF 12)	20pcs	700,00	577,50
10116544	9322+ Premium valved dust/mist respirator, provides protection against fine dusts and water based mists. Individually packed, EN 149:2001 FFP2. Particulates Assigned Protection Factor: 10 (NPF 12)	10pcs	565,00	466,13
12602555	9330+ Particulate Respirator unvalved, Protection against dusts and mists, FFP3	20pcs	785,00	647,63
11346974	9332+ Premium valved dust/mist/metal fumes respirator, provides protection against fine dusts, water based mists and metal fumes. Individually packed, EN 149:2001 FFP3 Particulates Assigned Protection Factor: 20 (NPF 50)	10pcs	1042,44	860,01



RESPIRATORS, K100 SERIES

- Diamond-shaped front allows the respirator to maintain its shape in hot and humid conditions. Single loop strap allows easy adjustment on the face and head. For dust and/or non-volatile liquid particles



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
12922475	K101 unvalved, FFP1, AFP4	200pcs	1400,00	1155,00
12932475	K111 valved, FFP1, AFP4	100pcs	1086,30	896,20
12942475	K102 unvalved, FFP2, AFP10	200pcs	1846,20	1523,12
12952475	K122 valved, FFP2, AFP10	100pcs	1295,40	1068,71



RESPIRATORS, 8000 SERIES (CLASSIC)

- The 3M™ Respirators 8000 Series provide lightweight, comfortable and effective respiratory protection against dust and mist. The convex shape, twin strap design, nose foam and nose clip ensure comfortable wear over a range of face sizes
- The collapse resistant shell and proprietary 3M™ Cool Flow™ valve in the 8822 and 8812 respirators offer durable, comfortable protection particularly in hot and humid conditions



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
11703363	3M™ 8710 Cup-Shaped Unvalved Dust/Mist Respirator - small pack FFP1	20pcs	243,00	200,48
11723363	3M™ 8812 Cup-Shaped Valved Dust/Mist Respirator FFP1	10pcs	218,00	179,85
11713363	3M™ 8810 Cup-Shaped Unvalved Dust/Mist Respirator - small pack FFP2	20pcs	740,00	610,50
11733363	3M™ 8822 Cup-Shaped Valved Dust/Mist Respirator FFP2 3M™	10pcs	331,00	273,08



TOUCHNTUFF™ DISPOSABLE BLUE NITRILE GLOVES

- Provides three times the puncture resistance of natural latex and vinyl gloves
- Exceptionally soft feel
- No additives to contaminate product or irritate skin
- Rolled cuff, 0,1mm thickness at the palm
- Length: 240mm



UP TO 42% OFF

Cat No.	Size	Pack qty	List price	Promo price
11766803	6,5-7	100pcs	209,00	145,00
11776803	7,5-8	100pcs	250,00	145,00
11786803	8,5-9	100pcs	250,00	145,00
11726843	9,5-10	100pcs	249,00	145,00



REFRIGERATORS AND FREEZERS, LABORATORY, ES SERIES FMS

- ES series FMS: refrigerators, freezers and combination refrigerator/freezers for the safe storage of potentially flammable materials. Compatible with ATEX Directive 94/9/EC Annex VIII



FEATURED PRODUCT

Cat No.	Description	Capacity	Dimensions, mm	Pack qty	List price
15234737	Refrigerator, auto defrost ATEX	151L	595 x 595 x 845	1pcs	9701,08
15254737	Refrigerator, auto defrost ATEX	288L	595 x 634 x 1570	1pcs	13805,02
15284737	Combination, manual-F/ auto-R defrost - ATEX	263L	540 x 600 x 1680	1pcs	17965,31
15264737	Freezer, manual defrost	151L	595 x 595 x 845	1pcs	9701,08
15274737	Freezer, manual defrost	232L	595 x 634 x 1570	1pcs	13805,02



GLOVES, NITRILE, DISPOSABLE, KIMTECH SCIENCE™ STERLING™ NITRILE GLOVES

- Combines protection, comfort and sustainability for all-round care of the science and the environment
- Tested against over 40 chemicals and 10 cytotoxic drugs. Strength, quality and cleanliness tracked by certificate of analysis
- PPE Category III compliant, AQL 0.65, EN374-1 and -2, ISO 16604, food contact approved
- Available in two lengths: 240mm or 250mm



GREAT SAVINGS

Cat No.	Size	Length, mm	Pack qty	List price	Promo price
10634525	XS	240	150pcs	258,00	234,01
10765495	S	240	150pcs	237,00	118,57
10735685	M	240	150pcs	234,00	114,52
10682415	L	250	150pcs	236,00	133,03
10754175	XL	250	140pcs	238,00	115,50



TOUCHNTUFF™ DISPOSABLE NITRILE GLOVE

- Ambidextrous, nitrile gloves are powder free with textured fingers for enhanced grip throughout the palm. Beaded cuff, anthracite grey colour. Compliant with ASTM D 6319, Type II, ISO 9001
- Ansell Grip Technology™
- Enhanced grip minimises force required to grip dry, wet or oily objects
- Soft and durable for long wear periods
- Three times the puncture resistance of latex or PVC disposables
- Chemical splash resistance
- Protects against Type I skin allergies
- Silicone free
- 'Touchscreen-capable' glove



GREAT SAVINGS

Cat No.	Size	Pack qty	List price	Promo price
15569461	XS	100pcs	190,83	140,00
15579461	S	100pcs	190,83	140,00
15589461	M	100pcs	190,83	140,00
15599461	L	100pcs	190,83	140,00
15509471	XL	100pcs	190,83	140,00

GLOVES, CRYOGENIC -190°C, PPE CATEGORY II

- Mid-Arm gloves
- For use with materials stored with liquid nitrogen and in other very cold environments
- Superb insulation needs less thickness giving more dexterity
- Synthetic fibre outer for flexibility
- Water resistant
- The middle layer is of special polyolefin/polyester micro-fibre, which is almost twice as effective in thermal insulation as a similar thickness of down. This minimises heat transfer without being so thick as to be cumbersome. The inner layer is of cotton for comfortable wear
- They are made to EN 420 sizes and have been tested to EN388 (2, X, 2, X,) and EN 511-contact cold 2, convective cold 2, Flexibility PASS
- Avoid contact with alcohol-based products



FEATURED PRODUCT

Cat No.	Size	Length, mm	Pack qty	List price
11883801	S	300	1pr	1635,00
11883831	M	360	1pr	1690,00
11893831	L	380	1pr	1675,00



GLOVES, NITRILE, DISPOSABLE, XCEED™ INT



- Category III, EN 374-2; AQL .65, ISO 13485.
- Ansell Microflex™ non-sterile nitrile blue gloves with textured fingers. Ambidextrous with powder free, chlorinated, polymer coating. Standard examination
- Length 240mm, finger thickness 0,11mm

GREAT SAVINGS

Cat No.	Size	Pack qty	List price	Promo price
15289373	XS	250pcs	222,00	180,00
15269373	S	250pcs	227,00	180,00
15259373	M	250pcs	227,00	180,00
15249373	L	250pcs	227,00	180,00
15279373	XL	230pcs	222,00	180,00



HYFLEX™ ULTRA LIGHT LEVEL 2 CUT RESISTANT NYLON GLOVES

- For ultra fine mechanical applications in variety of industries where mid-level performance cut resistance is needed
- The thin liner uncoated (18 gauge knit design) will maintain the hand at a comfortable temperature during use
- Spandex reinforcement
- Snug/anatomical fit on the hand thanks to the spandex reinforcement in the tension areas of the hand (knuckles and palm)



GREAT SAVINGS

Cat No.	Size	Pack qty	List price	Promo price
15496160	6	12pr	1280,00	945,00
15416170	7	12pr	1280,00	945,00
15426170	8	12pr	1280,00	945,00
15436170	9	12pr	1280,00	945,00
15466170	10	12pr	1280,00	945,00
15476170	11	12pr	1280,00	945,00



HYFLEX™ PALM DIPPED MEDIUM DUTY, 10 GAUGE GLOVES W/INTERCEPT™ TECHNOLOGY

- Excellent flexibility and fit with high levels of cut protection
- The comfortable and flexible plaited liner also features an optimised lycra content that is soft and cool on the skin for all-day comfort.
- The palm dipped polyurethane (PU) coating provides enhanced grip in dry and light oil environments for secure handling and offers EN Abrasion Level 4 (ANSI Level 6) for extended glove life
- The PU coating is also silicone-free to protect the quality of industrial finishing processes



GREAT SAVINGS

Cat No.	Size	Pack qty	List price	Promo price
15569891	6	12pr	1580,00	1229,00
15579891	7	12pr	1580,00	1229,00
15589891	8	12pr	1580,00	1229,00
15599891	9	12pr	1580,00	1229,00
15549891	10	12pr	1580,00	1229,00
15559891	11	12pr	1580,00	1229,00

fisherbrand

DISPOSABLE HIGH-TACK SHOE COVERS

- High-tack, disposable shoe covers provide superior protection
- Single-use
- Polyurethane
- Blue
- Latex-free
- High-tack sole



GREAT SAVINGS

Cat No.	Description	Colour	Size	Pack qty	List price	Promo price
15421952	Shoe cover, high-tack	Blue	Universal	300pcs	625,00	471,91
15431952	Shoe cover, high-tack	Blue	X-Large	300pcs	655,00	494,56

fisherbrand

DISPOSABLE SHOE COVERS

- Spunbond polypropylene shoe covers protect against dirt, grime and certain dry particulates
- Single-use
- Blue
- Polypropylene
- Latex-free
- With or without anti-skid sole



SPECIAL PRICE

Cat No.	Description	Colour	Size	Pack qty	List price	Promo price
15481942	Shoe cover without anti-skid	Blue	Universal	300pcs	229,00	213,08
15491942	Shoe cover without anti-skid	Blue	X-Large	300pcs	240,00	178,20
15401952	Shoe cover with anti-skid	Blue	Universal	300pcs	284,00	263,38
15411952	Shoe cover with anti-skid	Blue	X-Large	300pcs	366,00	260,58

thermo

WASTE DISPOSAL SYSTEM

- Hinged cover seal and long funnel stem reduces evaporation
- Built-in vent facilitates disposal
- Constructed of chemical and break resistant plastic
- Each consists of a safety waste funnel with a standard Thermo Scientific closure screwed onto the waste container (bottle or carboy) for temporary disposal of most chemical and biological wastes. For temporary waste disposal only. Product is not a safety can as defined by OSHA guidelines



GREAT SAVINGS

Cat No.	Description	Pack qty	List price	Promo price
10084680	139,7mm safety waste funnel HDPE with 4L bottle	1pcs	1033,03	751,10
10098021	254mm safety waste funnel HDPE with 10L bottle	1pcs	663,02	511,07
11764434	139,7mm safety waste funnel HDPE PP closure	1pcs	481,30	371,00
11774434	254mm safety waste funnel HDPE PP closure	1pcs	511,82	394,52



Distributor
GE Healthcare

SURFACE PROTECTORS, BENCHKOTE, WHATMAN™

- Benchkote is an absorbent, impermeable material designed to protect laboratory surfaces against hazardous spills. The material features a high-quality, smooth, absorbent Whatman paper that quickly absorbs liquid spills and a laminated polyethylene layer that prevents flowthrough to the working surface. After use the sheet is incinerated or disposed of according to local regulations.
- Material is very strong, making it tear resistant when wet or dry.
- Smooth white surface can be written on with ink or pencil and lies flat.
- Suitable for saturation with disinfectant to protect benches where pathogens and other bacteria are present.
- Use polyethylene side up to collect deposits without absorption.
- Paper side quickly absorbs liquid spills, preventing liquids from going through to the work surface.
- Spillages are trapped in the absorbent paper.
- Benchkote can be incinerated after use; the polyethylene layer does not melt or drip but is rapidly consumed in the flames



SPECIAL PRICE

Cat No.	Dimensions, L x W	Pack qty	List price	Promo price
11300674	460mm x 570mm, sheets	50pcs	847,65	771,00
11353195	400mm x 570mm, sheets	100pcs	974,52	887,00

thermoscientific

CLEAN BENCH HERAGUARD™ ECO

- Ideal for laboratory applications where sample protection is required.
- Complies with EN 61010-1: 2011, EN 61326:2006
- Low energy consumption
- Transparent side windows made of safety glass with access port for tubing or a media tap.
- Standard melamine work top or optional stainless steel
- Intuitive control panel displays safety and performance data
- Easy maintenance: prefilter on the intake side can be changed without tools
- Optional UV disinfection light



FEATURED PRODUCT

Cat No.	Description	Pack qty	List price
15128902	Clean Bench Heraguard Eco 0,9m	1pcs	44485,97
15158902	Clean Bench Heraguard Eco 1,20m	1pcs	44298,15
15148902	Clean Bench Heraguard Eco 1,50m	1pcs	52862,90
15138902	Clean Bench Heraguard Eco 1,80m	1pcs	54008,62
15168902	Clean Bench Heraguard Eco 180m, increased internal depth 78,5cm and height 95cm	1pcs	78698,01



600MM WING DOOR SAFETY CABINET FOR ACIDS AND ALKALIS

- Safety elements assembled outside the storage compartment for increased protection against corrosion, easy to clean
- Two hermetically separated storage compartments
- Cylinder locking (integration in locking system possible)
- Adjusting aids to compensate for uneven floor
- Integrated metal-free air ducts ready for connection (NW 75) to a technical exhaust system with only one extraction air inlet
- Excellent venting - separate extraction air of storage compartments
- GS approval
- EN 14727 conformity (laboratory furniture)
- Weight, kg: 125
- Extraction air, NW: 75
- Air exchange 30 times, m³/h: 35
- Pressure drop inside the cabinet, Pa: 56



10% OFF

Cat No.	Dimensions exterior, W x D x H, mm	Dimensions interior, W x D x H, mm	Pack qty	List price	Promo price
12643515	600 x 605 x 1970	485 x 565 x 930	1pcs	18619,76	16757,78



SAFETY STORAGE, UNDERBENCH CABINET, TYPE 90, FLAMMABLE LIQUID

- The new underbench cabinets from Asecos with a special height of 800mm fit perfectly underneath workbenches and worktops. They are perfect for the storage of larger containers because of the maximised interior height
- Fire resistance of 90 minutes – type tested according to DIN EN 14470-1
- Three different models with doors or drawers available
- Cylinder lock with locking state indicator (red/green)
- Integrated air ducts ready for connection (NW 75) to a technical exhaust system
- Interior options include second level drawer, perforated insert, pull-out tray, collecting sump



SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15385569	Safety storage cabinet, type 90, underbench, with two drawers, two perforated inserts and two second-level drawers	1pcs	26499,18	23900,00



SAFETY STORAGE CABINET Q-CLASSIC-90

- Suitable for flammables, acids/alkalis, chemical and toxins
- Not suitable for compressed gases
- Cabinet prevents external fires igniting the contents for 90 minutes (EN 14470-1 and EN 14727)
- Requires venting or filtration. Supplied with 3 shelves and a sump



10% OFF

Cat No.	Colour	Dimensions, W x D x H, mm	Pack qty	List price	Promo price
15596845	Anthracite grey	1200 x 615 x 1955	1pcs	24984,62	22486,16
11954836	Yellow	893 x 615 x 1953	1pcs	25905,13	23314,62
11934796	Yellow	893 x 615 x 1953	1pcs	20843,86	18800,00



SAFETY CABINETS FOR ACIDS AND ALKALIS

- Cabinet for acids / alkalis with interior drawers
- Safety elements assembled outside the storage compartment for increased protection against corrosion, easy to clean
- Two hermetically separated storage compartments
- Cylinder locking (integration in locking system possible)
- Adjusting aids to compensate for uneven floor
- Integrated metal-free air ducts ready for connection (NW 75) to a technical exhaust system with only one extraction air inlet
- Excellent venting - separate extraction air of storage compartments
- GS approval
- EN 14727 conformity (laboratory furniture)



10% OFF

Cat No.	Description	Pack qty	List price	Promo price
15231035	Safety cabinets for acids and alkalis	1pcs	28951,99	26056,25



STORAGE CABINETS, SAFETY

- Professional storage and handling of hazardous materials
- 90 minute fire resistance (EN 14470-1 and EN 14470-2)
- 10 year warranty for Q-LINE cabinets
- Five year warranty for all other cabinets
- Delivery charge not included
- Dimensions, W x D x H, mm: 600 x 615 x 2050




SPECIAL PRICE

Cat No.	Description	Pack qty	List price	Promo price
15121866	G-ULTIMATE-90 Safety storage cabinet, light grey, with mounting rails, bottle holders and fold-away ramp, 250mm	1pcs	28447,65	25700,00
Accessory				
12173951	Recirculating air filter system	1pcs	12688,81	11500,00

Experts in stable incubation processes

BINDER Cooling incubators



EXTREMELY ROBUST AND LONG-LASTING QUALITY



FLEXIBLE, MIXED STACKING OPTIONS



KT Series | with thermoelectric cooling

Combines exceptional performance
with low energy consumption

Now available
through
Fisher Scientific
across Europe

BINDER
Best conditions for your success



HOMOGENEOUS TEMPERATURE DISTRIBUTION



EASY DECONTAMINATION

KB series | with compressor technology

Safe and reproducible incubation,
even at high ambient temperatures

SPOTLIGHT

ISSUE 2

Your **NEW** Fisher Scientific website

Faster Search | Improved Self-service!



YOU SPOKE



WE LISTENED

VISIT www.fishersci.se

© 2017 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at fishersci.com/trademarks.

In Sweden:
For customer service, call 031 352 32 00
To fax an order, use 031 352 32 50
To order online: fishersci.se

 **fisher**scientific
by Thermo Fisher Scientific