

# smart spend

Great promotional offers and new innovative products from February to March 2017

## Featured Product

Protect your hearing. Situational awareness.  
Enhanced attention

### 3M™ PELTOR™ ProTac™ III Headset



**3M**

See page 12 for more details

## Featured Product

Ideal for applications that require stable  
temperatures above or below ambient

### NEW Thermo Scientific™ RI150 and RI250 Refrigerated Incubators with Powerful Compressor Technology



thermo**scientific**

See page 5 for more details

## Featured Product

Muffle furnaces with flap door or lift door

### Muffle Furnaces with NEW Fibre Material Classified Non-carcinogenic



**Nabertherm**  
MADE IN GERMANY

See page 6 for more details

## Featured Product

The Fisherbrand microplate shaker is  
designed for convenience and flexibility

### Microplate shaker



**Fisherbrand**

See page 10 for more details

Try Aqualine™ product performance at 25% discount



## Fisher Chemical™ Aqualine™ Reagents for Karl Fischer Titration

- Your solution for water content determination by Karl Fischer titration
- Reliable – Fast and stable endpoints ensure reliable and accurate results
- Safe – Low toxicity and pyridine free
- Convenient – Available as a one or two component solution

25% discount on a selection of Aqualine volumetric products for high water content analysis



Cat No.	Description	Quantity	Price
<b>Single component reagents</b>			
10538810	Aqualine Complete 1 – Water equivalent 1mg H <sub>2</sub> O/mL	1L	111.75
10744761	Aqualine Complete 2 – Water equivalent 2mg H <sub>2</sub> O/mL	1L	108.75
10754761	Aqualine Complete 2 – Water equivalent 2mg H <sub>2</sub> O/mL	2.5L	237.75
10181570	Aqualine Complete 2 – Water equivalent 2mg H <sub>2</sub> O/mL	1L	93.00
10676262	Aqualine Complete 2 – Water equivalent 2mg H <sub>2</sub> O/mL	2.5L	204.75
<b>Reagents for aldehydes and ketones</b>			
11433813	Aqualine Complete 5K – Water equivalent 5mg H <sub>2</sub> O/mL	500mL	61.46
10092042	Aqualine Complete 5K – Water equivalent 5mg H <sub>2</sub> O/mL	1L	97.50
11443813	Aqualine Matrix K – Matrix K should be used in conjunction with Complete 5K	500mL	56.55
10205592	Aqualine Matrix K – Matrix K should be used in conjunction with Complete 5K	1L	84.00
<b>Two component reagents</b>			
10786291	Aqualine Titrant 2 – Water equivalent 2mg H <sub>2</sub> O/mL	1L	72.64
10284100	Aqualine Titrant 2 – Water equivalent 2mg H <sub>2</sub> O/mL	2.5L	114.75
10528810	Aqualine Titrant 5 – Water equivalent 5mg H <sub>2</sub> O/mL	1L	42.04
10172050	Aqualine Titrant 5 – Water equivalent 5mg H <sub>2</sub> O/mL	2.5L	105.75
10101580	Aqualine Solvent	1L	41.29
10264390	Aqualine Solvent	2.5L	117.00
10199511	Aqualine Solvent CM – Solvent for samples with high hydrocarbon content	1L	61.73
10510852	Aqualine Solvent CM – Solvent for samples with high hydrocarbon content	2.5L	122.25

30% off

invitrogen

## Invitrogen™ iBlot™2 Gel Transfer Device

- Complete protein transfer in seven minutes or less
- High detection sensitivity and even transfer
- Increased blotting reliability and reproducibility
- Flexible gel size formats and membrane types
- A simple, user-friendly system

Transfer proteins quickly and efficiently from gel to membrane in seven minutes or less, plus save 30% off list price.



Cat No.	Description	Price
15217995	Invitrogen iBlot 2 gel transfer device	961.80

35% off

## Applied Biosystems™ SimpliAmp™ Thermal Cycler

- Compact design - saves bench space
- Easy to operate - large, colour LCD touchscreen that simplifies operation
- Three independent temperature zones for PCR optimisation - features Applied Biosystems™ VeriFlex™ Block
- Cloud-enabled for secure remote access through Thermo Fisher Connect

Enjoy the consistent performance of an Applied Biosystems™ Thermal Cycler for less.



Cat No.	Description	Price
15224438	SimpliAmp thermal cycler	2593.50

Manual identification is now easier than ever! With Thermo Scientific™ RapID™ Systems you'll experience faster identification, reduced turnaround time, visible colour reactions and precise reporting for more organisms on more platforms

## Thermo Scientific™ RapID™ Systems

- One-Step inoculation; Decreased prep time allows for increased productivity
- Four hour incubation; Quicker time-to-result for a faster response
- Common procedure; Lessened materials usage saves time and money
- Visible colour reactions; Reduced subjectivity decreases repeat testing

The RapID Systems are panels for manual identification of microorganisms, they can be used in different contexts like food pathogen identifications and clinical pathogens as well.



thermo scientific

Cat No.	Description	Pack qty	Price
11464139	RapID ANA II	20 panels/kit	175.06
13245619	RapID CB Plus	20 panels/kit	324.37
11454139	RapID NH Plus	20 panels/kit	139.16
11484139	RapID NF Plus	20 panels/kit	147.29
11494139	RapID ONE	20 panels/kit	98.54
13255619	RapID Staph Plus	20 panels/kit	249.30
13225619	RapID SS/u	20 panels/kit	148.17
11474139	RapID STR	20 panels/kit	148.78
11404149	RapID Yeast Plus	20 panels/kit	175.80
10302213	RapID Inoculation Fluid 1mL	20 tubes/pack	40.11
13246239	McFarland Standard No. 3	1	54.02
12853693	RapID Spot Indole 15mL	1	18.24
12995973	RapID ANA II QC Set/3	5 loops/pack	309.80
12998264	ERIC RapID Compendium	1	17.85
12883693	RapID Inoculation Fluid 2mL	20 tubes/pack	16.51
10108582	McFarland Standard No. 4	1	54.02
13285559	RapID Nitrate A 15mL	1	18.54
10392223	RapID Nitrate B 15mL	1	18.54
12985973	RapID CB Plus QC Set/3	5 loops/pack	309.80
10302213	RapID Inoculation Fluid 1mL	20 tubes/pack	40.11
12853693	RapID Spot Indole 15mL	1	18.24
10399814	BactiDrop Oxidase 0.75mL	1	127.93
12905983	RapID NH Plus QC Set/4	5 loops/pack	309.80
10037212	McFarland Standard No. 1	1	54.02
12978244	RapID NF Plus QC Set/4	5 loops/pack	309.80
10413843	McFarland Standard No. 2	1	54.02
12945983	RapID ONE QC Set/4	5 loops/pack	309.80
12975973	RapID Staph Plus QC Set/4	5 loops/pack	309.80
12925983	RapID SS/u QC Set/5	5 loops/pack	309.80
12955983	RapID Yeast Plus QC Set/5	5 loops/pack	309.80

Buy the new laboratory stands and get 30% discount



## NEW LaboPlast™ Laboratory Stands Promotion

- Funnel stands/filtering stands; stand for Imhoff sedimentation funnel or pipettes
- With innovative magnet concept
- Easy and flexible height adjustment
- Adaptable with corresponding accessories
- Comfortable and smooth operation

The lab frame can be flexibly adapted with corresponding accessories, for example as a frame for Imhoff sedimentation funnels, a pipette stand or as a funnel frame.



Cat No.	Description	Price
10057980	Funnel stand/filtering stand, 2 funnel	50.93
10393731	Funnel stand/filtering stand, 4 funnel	85.50
15624144	Funnel stand/funnel holder with reducing insert	5.42
10740773	Funnel rack, carrying handle, 2 funnel	34.21
15614144	Sedimentation funnel holder for Imhoff funnel	6.60
15634144	Pipette stand, rotatable, with 40 openings	35.55
15644144	Replacement set for pipette stand, plates/drip pan	9.75

Buy a KBF series constant climate chamber and get IQ/OQ qualification folder for FREE



## KBF Series Constant Climate Chambers

- Long term stability test conditions in line with ICH Q1A and Q1B
- Independent from water supply
- Homogeneous light distribution over the entire usable space (KBF P, LQC series)
- Independent light dose control via spherical sensors (KBF LQC series)
- Support for GMP/GLP conformity with IQ/OQ qualification folders

Achieve the best test results with the BINDER KBF series.

Get a KBF series constant climate chamber with free IQ/OQ qualification folder.



Cat No.	Description	Price
15669105	KBF115-230V with Qualification folder IQ/OQ (8012-0867)	8,837.17
15679105	KBF240-230V with Qualification folder IQ/OQ (8012-0867)	10,017.28
15689105	KBF720-230V with Qualification folder IQ/OQ (8012-0867)	12,729.09
15659105	KBF1020-230V with Qualification folder IQ/OQ (8012-0867)	14,319.77
15619115	KBFLQC240-230V with Qualification folder IQ/OQ (8012-0869)	16,714.63
15629115	KBFLQC720-230V with Qualification folder IQ/OQ (8012-0869)	19,334.30
15699105	KBFP240-230V with Qualification folder IQ/OQ (8012-0868)	12,458.57
15609115	KBFP720-230V with Qualification folder IQ/OQ (8012-0868)	15,184.39

Buy a BINDER cooling incubator at a special price

## BINDER KB Series (Compressor) and KT Series (Peltier) Cooled Incubators

- KB series: Powerful cooling with well proven compressor cooling unit
- KT series: Energy efficient thermoelectric cooling with Peltier technology
- Reliable and reproducible incubation conditions
- Disinfection routine at 100°C
- High quality, corrosion resistant stainless steel inner chamber

Get a special price on BINDER Cooling Incubators. KB Series with compressor technology or KT Series with thermoelectric cooling.



Cat No.	Description	Price
<b>KB Series</b>		
15609105	KB023-230V cooled incubator, with powerful compressor, 23L	2,458.40
15649105	KB053-230V cooled incubator, with powerful compressor, 53L	2,669.50
15639105	KB115-230V cooled incubator, with powerful compressor, 115L	3,190.78
15629105	KB240-230V cooled incubator, with powerful compressor, 247L	4,580.09
15619105	KB400-230V cooled incubator, with powerful compressor, 400L	5,124.37
15699095	KB720-230V cooled incubator, with powerful compressor, 698L	6,967.72
<b>KT Series</b>		
15649115	KT053-230V cooled incubator with Peltier technology, 53L	2,250.83
15659115	KT115-230V cooled incubator with Peltier technology, 102L	2,816.93
15639115	KT170-230V cooled incubator with Peltier technology, 163L	3,477.96

## Featured Product

Ideal for applications that require stable temperatures above or below ambient

## NEW Thermo Scientific™ RI150 and RI250 Refrigerated Incubators with Powerful Compressor Technology

- Broad temperature range of +4°C to +60°C
- Adjustable over- and under-temperature alarm alerts operator of any unexpected temperature deviations
- Glass window in door for undisturbed inspection of samples
- Access port provided as standard
- Automatic defrost for uninterrupted long term sub-ambient applications

Choose a new incubator that provides a broad temperature range enabling it to be used across a wide range of applications.



Cat No.	Description	Price
15658175	RI150 refrigerated incubator with compressor, 220V EU, 150L	3,000.00
15668175	RI250 refrigerated incubator with compressor, 220V EU, 250L	3,260.00
15678175	Shelf accessory for RI150 and RI250 refrigerated incubators	99.00

## Featured Product



Muffle furnaces with flap door or lift door

### Muffle Furnaces with NEW Fibre Material Classified Non-carcinogenic

- T-max 1,100°C or 1,200°C
- Heating from two sides by ceramic heating plates (heating from three sides for muffle furnaces L 24/11 - LT 40/12)
- Insulation made of non-classified fibre material
- Housing made of sheets of textured stainless steel
- Dual shell housing for low external temperatures and high stability
- Optional flap door (L) which can be used as work platform or lift door (LT) with hot surface facing away from the operator
- Adjustable air inlet integrated in door
- Exhaust air outlet in rear wall of furnace
- Solid state relays provide for low noise operation

The muffle furnaces L 1/12 - LT 40/12 are the right choice for daily laboratory use. These models stand out for their excellent workmanship, advanced and attractive design, and high level of reliability. The muffle furnaces come equipped with either a flap door or lift door at no extra charge.

Note: other volume sizes and lift door versions are also available, please contact Customer Service.

Additional equipment:

- Chimney, chimney with fan or catalytic converter (not for L 1)
- Over-temperature limiter with adjustable cut-out temperature for thermal protection Class 2 in accordance with EN 60519-2 as temperature limiter to protect the furnace and load
- Protective gas connection to purge with non-flammable protective or reaction gases (not available in combination with chimney, chimney with fan or catalytic converter)
- Manual or automatic gas supply system



Cat No.	Description	Price
<b>T-max 1,100°C</b>		
15483392	L 3 / 11 muffle furnace, volume 5L, with flap door	1,329.43
15493392	L 9 / 11 muffle furnace, volume 9L, with flap door	1,793.19
15426732	L 40 / 11 muffle furnace, volume 40L, with flap door	2,968.03
<b>T-max 1,200°C</b>		
15446672	L 3 / 12 muffle furnace, volume 5L, with flap door	1,638.60
15426702	L 9 / 12 muffle furnace, volume 9L, with flap door	2,294.04
15446732	L 40 / 12 muffle furnace, volume 40L, with flap door	3,413.24

Data-logging thermometers with user-defined time intervals



### Traceable™ Excursion-Trac™ Data Logging Refrigerator/Freezer Thermometer

- Meets current CDC requirements for vaccine storage and monitoring
- Rolling memory structure, maintains most recent 525,600 temperature observations
- Stores up to 10 unique alarm events
- Memory can be cleared on unit after downloading data to USB stick
- Hassle-free retrieval of data - unit can remain in use while downloading and analysing data
- Status indicators - low battery, memory full, USB data transfer and active alarm state
- Supplied with stand, Velcro™, magnets, wall mount, batteries, Traceable™ Certificate



Cat No.	Range	Resolution	Accuracy	Probe	Price
15294016	-50 to 70°C (-58 to 158°F)	0.01°	±0.25°C	1 Bottle	98.80
15204026	-50 to 70°C (-58 to 158°F)	0.01°	±0.25°C	2 Bottles	105.95
15214026	-50 to 70°C (-58 to 158°F)	0.01°	±0.25°C	1 Bullet™	94.90
15224026	-50 to 70°C (-58 to 158°F)	0.01°	±0.25°C	2 Bullets™	102.70
15234026	-50 to 70°C (-58 to 158°F)	0.01°	±0.25°C	1 Bottle & Bullet™	105.95
15244026	-50 to 70°C (-58 to 158°F)	0.01°	±0.25°C	1 Stainless steel	105.95
15388754	-50 to 70°C (-58 to 158°F)	0.01°	±0.25°C	2 Stainless steel	110.50

The world's most popular refrigerator/freezer thermometers

## Traceable™ Refrigerator/Freezer Thermometer

- Triple display simultaneously shows MIN/MAX and current temperatures, °C/°F switchable
- Accurately monitor temperatures in freezers, water baths, heating blocks, incubators and refrigerators with this enclosed temperature-buffered sensor
- Bottle insulates sensor from rapid temperature changes when refrigerator door is opened
- Minimum/maximum monitors high and low readings overnight, on weekends, or for any time period
- Alarm is programmable in 1° increments, sounds continuously when temperature rises above/falls below set points and can be turned off manually
- Visual and audible signals continue even if temperature returns to non-alarm range
- High-impact, chemical-resistant ABS plastic case: 70mm x 108mm x 19mm, weight 113g
- Supplied with stand, 3m micro cable, wall mounting, Velcro™, magnetic strips, Traceable™ Certificate



11873460

11799735

11749745

Cat No.	Range	Resolution	Accuracy	Probe	Price
11873460	-50 to 70°C (-58 to 158°F)	1°	±1°C	Bottle (patented) specifically for refrigerator/freezer applications	35.43
11799735	-50 to 70°C (-58 to 158°F)	1°	±1°C	5mm diameter and 19mm long	32.89
11749745	-50 to 70°C (-58 to 158°F)	1°	±1°C	5mL vaccine bottle mimics vaccine temperature and fits in tray	41.37

The lab's most popular timer

## Traceable™ Four-Channel Alarm Timer

- Features: count up/down, time in/out, memories, alarm, 12/24 hour clock
- Large high display allows viewing from a distance
- Remarkable memory returns display to previously programmed countdown time at the touch of a button
- Extra-loud, high decibel alarms for 1 minute or can be silenced manually
- Timing capacity: 99 hours, 59 minutes, 59 seconds
- Resolution: 1 second
- Accuracy: 0.01%
- Supplied with clip, stand, magnet, battery, Traceable™ Certificate



Cat No.	Description	Price
11745863	Traceable™ Four-Channel Alarm Timer	20.41

Focus on Fisherbrand  
for all your  
**Traceable™** Products  
Keeping your research on track

- Lab tools
- Humidity monitors
- Spatulas
- Stop watches
- Timers
- Thermometers



For the full range visit [eu.fishersci.com/go/fisherbrand](http://eu.fishersci.com/go/fisherbrand)

Fisherbrand nutating mixers deliver gentle uniform mixing of fluids in tubes, plates, culture flasks and dishes

### Nutating mixer - fixed speed

- Budget friendly design
- Easy to use, plug and play operation
- Fixed-speed movement at 20° angle
- Operating temperature range: 5 to 40°C, usable in non-CO<sub>2</sub> incubators, refrigerators and chambers
- Fixed speed of 18 or 20rpm (voltage dependent)
- Two year warranty



Cat No.	Description	Price
15584070	Nutating mixer - fixed speed	195.70

Fisherbrand nutating mixers deliver gentle uniform mixing of fluids in tubes, plates, culture flasks and dishes

### Nutating mixer - variable speed

- Variable speed from 2 to 60rpm
- 306mm x 266mm platform
- Display will show last used settings
- Operating temperature range: 5 to 40°C, usable in non-CO<sub>2</sub> incubators, refrigerators and chambers
- Continues to run by last set speed
- Continuous or timed operation with automatic switch-off
- Audible and visual alarm when the set time reaches zero
- Two year warranty



Cat No.	Description	Price
15594070	Nutating mixer - variable speed	404.70

Fisherbrand 3D platform rotators provide smooth orbital motion for mixing in culture flasks, dishes, plates and tubes

### 3D Platform rotator

- Simple speed adjustment from 2 to 80rpm
- Gentle movement ideal for staining gels without destroying the edges or keeping fragile tissue intact during incubation
- Large platform accommodates most commonly used vessels with non-slip mat to prevent sample movement
- Low voltage power supply provides safe cold room operation and low energy consumption
- Operating temperature range of 5°C to 40°C allows for use in refrigerators or non-CO<sub>2</sub> incubators
- Continuous or timed operation with automatic switch-off
- Display will show last used settings, even after power has been turned off
- Two year warranty



Cat No.	Description	Price
15504080	3D Platform rotator	341.05

The Fisherbrand multi function 3D rotator provides flexibility of motion to include rotation, reciprocal rotation, and vibration

### Multi function 3D rotator

- Rotation speed: 5 to 100rpm for long wave motion
- Reciprocal rotation: 0 to 360° turning angle for long and medium wave motion
- Vibration: 0 to 5° turning angle, programmable in bursts of 0 to 5 seconds, for short wave motion
- Operating temperature range of 5 to 40°C allows for use in refrigerators or non-CO<sub>2</sub> incubators
- Continuous or timed operation
- Display will show last used settings
- Two year warranty



Cat No.	Description	Price
15514080	Multi function 3D Platform rotator	555.75



Fisherbrand mini tube rotators accommodate 1.5mL to 50mL tubes. Adjustable digital speed control and agitation mode allow for a wide variety of applications

### Mini tube rotator

- Small footprint fits easily into refrigerators and non-CO<sub>2</sub> incubators
- Digital RPM display
- Paddles are easily adjustable/interchanged
- Extra paddles are needed if using all same size tubes (except for 50mL size)
- Two year warranty



Cat No.	Description	Price
15534080	Mini tube rotator	173.93
<b>Optional accessories</b>		
15554140	Two paddles for 0.5mL/0.8mL tubes	9.05
15544140	Two paddles for 1.5mL/2mL tubes	9.05
15534140	Two paddles for 5mL/7mL tubes	9.05
15524140	Two paddles for 10mL/15mL tubes	9.05
15564140	Paddle for 50mL tubes	9.05

The Fisherbrand multi purpose tube rotator offers flexibility with multiple carousel options, adjustable speed, and an angled head for added versatility

### Multi purpose tube rotator

- Digital microprocessor control
- Six carousel options for maximum versatility
- Adjustable speed from 5 to 80rpm
- Adjustable angled head
- Timer features both time up and time down counting modes
- Display will show last used settings
- Two year warranty



Cat No.	Description	Price
15524080	Multi purpose tube rotator	281.20
<b>Optional accessories</b>		
15554090	50 place tube carousel for 1.5/2.0mL tubes	106.40
15564090	30 place tube carousel for 15mL tubes	106.40
15574090	15 place tube carousel for 50mL tubes	106.40
15584090	10 and 20 place tube carousel for 10 to 15mL and 5 to 7mL tubes	106.40
15594090	64 place drum tube carousel for 14mm dia. tubes	106.40
15504100	64 place drum tube carousel for 18.5mm dia. 15 to 20mL tubes	106.40

The Fisherbrand multi platform shaker has smooth orbital shaking and a wide variety of optional accessories for multiple applications

### Multi platform shaker

- Variable shaking speed: 50 to 450rpm
- Automatic load balancing system
- Continuous or timed operation with automatic switch-off
- Maximum flexibility with many options for interchangeable accessories
- Simple set up and easy operation
- Display will show last used settings
- Low voltage power supply provides for safe cold room operation and low energy consumption
- Operating temperature range of 5 to 40°C allows for use in cold rooms, refrigerators or non-CO<sub>2</sub> incubators
- Two year warranty



Cat No.	Description	Price
15594060	Multi platform shaker	487.82
<b>Required accessories: these items must also be purchased as appropriate</b>		
15524100	Flat platform with non-slip rubber mat	75.24
15534100	Platform with clamps for 6 x 250mL flasks/150mL beakers	139.65
15544100	Platform with clamps for 12 x 100mL flasks/150mL beakers	139.65
15554100	Universal platform with adjustable bars	161.50
15564100	25mL stainless steel flask clamp	17.72
15574100	50mL stainless steel flask clamp	18.81
15584100	100mL stainless steel flask clamp	18.81
15594100	250mL stainless steel flask clamp	21.14
15504110	500mL stainless steel flask clamp	24.99
15514110	1000mL stainless steel flask clamp	28.50
15524110	2000mL stainless steel flask clamp	39.19
15534110	3000mL stainless steel flask clamp	57.00
15544110	5000mL stainless steel flask clamp	71.25
15554110	Sticky mat 140mm x 140mm	28.50
15564110	500mL infusion bottle clamp	28.50
15574110	1000mL infusion bottle clamp	35.63
15584110	Funnel clamp	57.00
15594110	Fixed tube rack 40 x 14mm dia.	57.00
15504120	Fixed tube rack 40 x 16mm dia.	57.00
15514120	Adjustable angle tube rack 40 x 14mm dia.	64.13
15524120	Adjustable angle tube rack 40 x 16mm dia.	64.13

Note: Not for use in CO<sub>2</sub> incubators. Will void warranty.

These sleek and compact bottle/tube rollers are capable of rolling multiple sizes of tubes or bottles

### Digital CO<sub>2</sub> and humidity resistant bottle/tube roller

- Anti-corrosive and wear resistant construction
- Adjustable rollers
- Stacking kit allows you to stack units two high
- Optional remote control




Cat No.	Description	Price
15376617	Digital CO <sub>2</sub> and humidity resistant bottle/tube roller	720.01
<b>Optional accessories</b>		
15386617	Additional rollers (pack of two)	23.03
15316627	Stacking kit	51.83
15326627	Remote control	63.35

## Featured Product

The Fisherbrand microplate shaker is designed for convenience and flexibility

### Microplate shaker

- Variable shaking speed: 150 to 1000rpm
- Quick and easy fitting of any standard-depth microplates
- Continuous or timed operation, with automatic switch off
- Display will show last used settings
- Operating temperature of 5 to 40°C allows for use in refrigerators, cold rooms and non-CO<sub>2</sub> incubators
- Two year warranty



Tube rack

Cat No.	Description	Price
15504070	Microplate shaker	281.44
<b>Optional accessories</b>		
15594110	Fixed tube rack 14mm dia. – 40 tubes max	57.00
15504120	Fixed tube rack 16mm dia. – 40 tubes max	57.00
15514120	Adjustable tube rack 14mm dia. – 40 tubes max	64.13
15524120	Adjustable tube rack 16mm dia. – 40 tubes max	64.13

Fisherbrand open air rockers are ideal for cell culture work and can be used in a variety of environmental conditions

### Open air rocker

- The microprocessor control provides electronic tilt angle adjustment which allows the user to easily adjust rocking angle from 1 to 15° while unit is operating
- Precise speed control
- Display will show last user settings
- LED display
- Audible and visual signals will activate when system detects an obstruction or overload of the tray
- Slowly increases speed to desired set-point
- Operating conditions: usable in refrigerators, cold rooms or non-CO<sub>2</sub> incubators, maximum 80% relative humidity, non-condensing
- Two year warranty



Cat No.	Description	Price
15514070	Open air rocker	536.75
<b>Optional accessories</b>		
15584130	Stacking tray assembly for open air rocker	66.12
15594130	Dimpled rubber mat	27.36
15514140	Rubber mat	27.36

Precise stirring control with exceptional safety and temperature performance

### Thermo Scientific™ Cimarec+™ Digital Stirrers, Hotplates and Stirrer Hotplates

- Hotplates, stirrers and stirrer hotplates all available in three sizes: 10cm x 10cm, 18cm x 18cm and 25cm x 25cm
- Ceramic surfaces available for all sizes and aluminium available for 10cm x 10cm
- Easy to read digital LED display
- Automatic over-temperature protection and HOT TOP warning system
- Three year warranty to provide peace of mind

Receive special prices for the Cimarec+™ range for a limited period.



Cat No.	Description	Price
<b>Stirrers</b>		
15463600	Cimarec+ basic, stirrer, ceramic top, 10cm x 10cm	167.40
15403610	Cimarec+ basic, stirrer, ceramic top, 18cm x 18cm	174.60
15443610	Cimarec+ basic, stirrer, aluminum top, 18cm x 18cm	226.80
15453620	Cimarec+ basic, stirrer, ceramic top, 25cm x 25cm	292.50
<b>Hotplates</b>		
15443600	Cimarec+ basic, hotplate, ceramic top, 10cm x 10cm	153.90
15493600	Cimarec+ basic, hotplate, ceramic top, 18cm x 18cm	166.50
15443620	Cimarec+ basic, hotplate, ceramic top, 25cm x 25cm	284.40
<b>Stirrer hotplates</b>		
15483600	Cimarec+ basic, stirrer hotplate, ceramic top, 10cm x 10cm	271.80
15433610	Cimarec+ basic, stirrer hotplate, ceramic top, 18cm x 18cm	291.60
15453610	Cimarec+ basic, stirrer hotplate, aluminum top, 18cm x 18cm	350.10
15483620	Cimarec+ basic, stirrer hotplate, ceramic top, 25cm x 25cm	381.60

## Featured Product

**3M**

Protect your hearing. Situational awareness.  
Enhanced attention

### 3M™ PELTOR™ ProTac™ III Headset

- Hearing protector helps protect against harmful noise, whilst allowing ambient sounds to be heard
- Built-in level-dependent function allows the wearer to hear ambient sounds, such as conversations, machine sounds and warning signals
- 3.5mm jack to connect to your mobile phone or 2-way radio to hear the call right in your ear
- Auto power-off turns headset off after four hours of non-use to save energy, a low battery warning informs you when to change batteries
- Voice guided menu system so no need to take your headset off
- Level dependant microphone for ambient listening

You have two great reasons to protect your hearing – your ears. The 3M™ PELTOR™ ProTac™ III Headset with level-dependent function for ambient listening helps you focus on your work and communicate with colleagues in noisy environments, without removing the hearing protection. Stereo feedback will allow you to hear and better determine direction of alarms, warning cries, reversing vehicles and so on.



Cat No.	Description	Price
15638175	PELTOR ProTac III Headset, black, headband	69.95
15618175	PELTOR ProTac III Slim Headset, black, headband	64.95
15648175	PELTOR ProTac III Headset, black, helmet attachment	64.95
15628175	PELTOR ProTac III Slim Headset, black, helmet attachment	69.95

© 2017 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [www.fishersci.com/trademarks](http://www.fishersci.com/trademarks).

#### In the United Kingdom:

For customer service, call 01509 555000

To fax an order, use 01509 555111

To order online: [www.fishersci.co.uk](http://www.fishersci.co.uk)

12279\_UK



**Fisher  
Scientific**

A Thermo Fisher Scientific Brand