


# Your guide to chromatography media



Pure proteins today. Powerful results tomorrow.



GE format guide - For ordering information, visit [www.gelifesciences.com](http://www.gelifesciences.com)

Type of purification	Manual purification	Manual or system purification	System purification
<b>Format</b>	Spin columns	96-well plates	Gravity flow columns
<b>Format name</b>	SpinTrap™	MultiTrap™	GraviTrap™, MiniTrap™, MidiTrap™, PD10
<b>Use</b>	Screening and quick desalting of small sample quantities using a benchtop centrifuge	High-throughput screening and small scale purification using centrifuge or vacuum equipment	Simple one-step purification of proteins or sample desalting without the need for equipment
		Bottles of chromatography media (resins)	Small-column cartridges
		Lab pack	HiTrap™
		Batch purification and self-packing	Easy to use with a syringe, peristaltic pump, or a chromatography system
			Other columns
			HiScreen™, HiPrep, HiLoad, RESOURCE™, Tricorn™, and Precision
			Larger scale or high performance applications

Need further guidance?  
Download our free Purify App  
[www.gelifesciences.com/purify](http://www.gelifesciences.com/purify)



GE and GE monogram are trademarks of General Electric Company. Capto, GraviTrap, HiLoad, HiPrep, HiScreen, HiTrap, MabSelect SuRe, MidiTrap, MiniBeads, MiniTrap, MonoBeads, MultiTrap, RESOURCE, Sephacryl, Sephadex, Sepharose, SOURCE, SpinTrap, Superdex, Superose, and Tricorn are trademarks of General Electric Company or one of its subsidiaries. Strep-tag and Strep-tactin are trademarks of IBA GmbH. TALON is a trademark of Clontech Laboratories, Inc. Superflow is a trademark of Sterogene Bioseparations Inc. All other trademarks are the property of their respective owners. © 2015 General Electric Company—All rights reserved. First published July 2015. GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK. [www.gelifesciences.com](http://www.gelifesciences.com)