



Your guide to chromatography media

Pure proteins today. Powerful results tomorrow.

Click on the chromatography method you are interested in to jump to the corresponding selection guide.

AC
Affinity chromatography

SEC
Size exclusion chromatography

IEX
Ion exchange chromatography




GE format guide

Type of purification	Manual purification			Manual or system purification		System purification
Format	Spin columns	96-well plates	Gravity flow columns	Bottles of chromatography media (resins)	Small-column cartridges	Other columns
Format name	SpinTrap™	MultiTrap™	GraviTrap™, MiniTrap™, MidiTrap™, PD10	Lab pack	HiTrap™	HiScreen™, HiPrep™, HiLoad™, RESOURCE™, Tricorn™, and Precision
Use	Screening and quick desalting of small sample quantities using a benchtop centrifuge	High-throughput screening and small scale purification using centrifuge or vacuum equipment	Simple one-step purification of proteins or sample desalting without the need for equipment	Batch purification and self-packing	Easy to use with a syringe, peristaltic pump, or a chromatography system	Larger scale or high performance applications

Affinity chromatography

Guide to chromatography media


Click on the icons to be directed to the corresponding web page.



Size exclusion chromatography

Guide to chromatography media

Click on the icons to be directed to the corresponding web page.




Ion exchange chromatography

Home

Guide to chromatography media

Click on the icons to be directed to the corresponding web page.



Download the Purify app for help with choosing the right purification products.

www.gelifesciences.com/purify

Join ÄKTA™ today

www.gelifesciences.com/AKTAclub



For ordering information, visit www.gelifesciences.com

GE, GE monogram, Capto, GraviTrap, HiLoad, HiPrep, HiScreen, HiTrap, MabSelect SuRe, MidiTrap, MiniBeads, MiniTrap, MonoBeads, MultiTrap, RESOURCE, Sephacryl, Sephadex, Sepharose, SOURCE, SpinTrap, Superdex, Superose, and Tricorn are trademarks of General Electric Company. Strep-tag and Strep-tactin are trademarks of IBA GmbH. TALON is a trademark of Clontech Laboratories, Inc. Superflow is a trademark of Sterogene Bioseparations Inc. All other third-party trademarks are the property of their respective owners.

© 2015–2016 General Electric Company. First published Aug. 2015.

GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK. www.gelifesciences.com

29167217 AB 05/2016