



Fisherbrand
QUALITY. RELIABILITY. VALUE.

Dialysis

Dialysis Tubing

Dialysis Tubing

Diffuses low molecular weight compounds in aqueous solutions; retains proteins and bacteria in tubing

- Regenerated cellulose dialysis tubing
- Available range 3.5kDa to 14kDa
- Pack includes two locking membrane clamps and a pair of surgical scissors
- Tubing is supplied on a compact polystyrene spool that eliminates absorption of glycerin from membranes
- Dispenser box is packaged in a plastic bag with a desiccant for humidity control
- Supplied with individual lot data sheets from outside laboratory*



Cat. No	Description	Pack qty
15203857*	Dialysis tubing, 14.6mm diameter, 30m roll, 6 to 8kDa MWCO	1
15213857*	Dialysis tubing, 20.4mm diameter, 30m roll, 6 to 8kDa MWCO	1
15223857*	Dialysis tubing, 25.5mm diameter, 30m roll, 6 to 8kDa MWCO	1
15233857*	Dialysis tubing, 31.8mm diameter, 30m roll, 6 to 8kDa MWCO	1
15243857*	Dialysis tubing, 6.37mm diameter, 15m roll, 12 to 14kDa MWCO	1
15293847*	Dialysis tubing, 15.9mm diameter, 15m roll, 12 to 14kDa MWCO	1
15253857*	Dialysis tubing, 28.6mm diameter, 15m roll, 12 to 14kDa MWCO	1
15283847	Dialysis tubing, 6.37mm diameter, 30m roll, 12 to 14kDa MWCO	1
15273847	Dialysis tubing, 15.9mm diameter, 30m roll, 12 to 14kDa MWCO	1
15263847	Dialysis tubing, 21.0mm diameter, 30m roll, 12 to 14kDa MWCO	1
15253847	Dialysis tubing, 28.6mm diameter, 30m roll, 12 to 14kDa MWCO	1
15233847*	Dialysis tubing, 12.1mm diameter, 15m roll, 3,500Da MWCO	1
15263857*	Dialysis tubing, 29.3mm diameter, 15m roll, 3,500Da MWCO	1
15273857	Dialysis tubing clamp, sealing width up to 46mm	10
15243847	Dialysis tubing clamp, sealing width up to 65mm	10

2014_UK197_10599

© 2014 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at www.fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20

Germany: +49 2304 932-5 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478

Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00

Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303



**Fisher
Scientific**

A Thermo Fisher Scientific Brand