

Dropper Bottles and Assemblies



WHEATON® glass dropper bottles and assemblies are manufactured from **clear neutral tubular glass** and are the perfect primary packaging solution for a range of **diagnostic, pharmaceutical** and **healthcare** applications.

Made from chemically inert Type 1B Neutral Glass, the glass dropper bottles and dropper assemblies are suitable for key packaging applications such as blood grouping reagents and rare grouping sera. Glass dropper bottles are also an advantage for several in vitro diagnostic applications where false positives from plastic dropper bottles may be a concern.

The dropper assemblies are manufactured with pharmaceutical grade rubber bulbs and a choice of coloured caps to help differentiate products in kits, and are also manufactured with a controlled drop size for consistency.

The vials are also available to purchase either as trays of vials packed in a box, or trays of vials shrink wrapped onto a pallet to reduce regular shipping costs.

Key Features

- Manufactured from **Type 1B** Clear Neutral Glass
- Standard sizes available are **3ml, 4.5ml, 7ml** and **10ml**
- Controlled drop size dropper assembly
- Choice of **different coloured caps** on dropper assemblies allows for differentiation in diagnostic kits



Product Code	Product Description	Bottle Size (d x h mm)	Case Quantity	Packaging / Inner Packs
10ml Dropper Bottles and Assemblies				
W017922	10ml Clear Type 1B Tubular Glass Dropper Bottle	20.25 x 58.5	245	MOQ applies – 228 trays of 245 palletised
W017922/4	10ml Clear Type 1B Tubular Glass Dropper Bottle	20.25 x 58.5	980	4 trays of 245 in box
W017924	Black Cap / Black Bulb Dropper Assembly for 10ml	–	980	1 bag of 980 in box
W017925	White Cap / Black Bulb Dropper Assembly for 10ml	–	980	1 bag of 980 in box
W017926	Blue Cap / Black Bulb Dropper Assembly for 10ml	–	980	1 bag of 980 in box
W017927	Yellow Cap / Black Bulb Dropper Assembly for 10ml	–	980	1 bag of 980 in box
W017928	Red Cap / Black Bulb Dropper Assembly for 10ml	–	980	1 bag of 980 in box
W017929	Green Cap / Black Bulb Dropper Assembly for 10ml	–	980	1 bag of 980 in box
W017930	Grey Cap / Black Bulb Dropper Assembly for 10ml	–	980	1 bag of 980 in box
7ml Dropper Bottles and Assemblies				
W017931	7ml Clear Type 1B Tubular Glass Dropper Bottle	20.25 x 42.75	245	MOQ applies – 264 trays of 245 palletised
W017931/4	7ml Clear Type 1B Tubular Glass Dropper Bottle	20.25 x 42.75	980	4 trays of 245 in box
W017933	Black Cap / Black Bulb Dropper Assembly for 7ml	–	980	1 bag of 980 in box
W017934	White Cap / Black Bulb Dropper Assembly for 7ml	–	980	1 bag of 980 in box
W017935	Blue Cap / Black Bulb Dropper Assembly for 7ml	–	980	1 bag of 980 in box
W017936	Yellow Cap / Black Bulb Dropper Assembly for 7ml	–	980	1 bag of 980 in box
W017937	Red Cap / Black Bulb Dropper Assembly for 7ml	–	980	1 bag of 980 in box
W017938	Green Cap / Black Bulb Dropper Assembly for 7ml	–	980	1 bag of 980 in box
W017939	Grey Cap / Black Bulb Dropper Assembly for 7ml	–	980	1 bag of 980 in box
4.5ml Dropper Bottles and Assemblies				
W017940	4.5ml Clear Type 1B Tubular Glass Dropper Bottle	16.25 x 46.5	374	MOQ applies – 264 trays of 374 palletised
W017940/5	4.5ml Clear Type 1B Tubular Glass Dropper Bottle	16.25 x 46.5	1870	5 trays of 374 in box
W017942	Black Cap / White Bulb Dropper Assembly for 4.5ml	–	1870	1 bag of 1870 in box
W017943	White Cap / White Bulb Dropper Assembly for 4.5ml	–	1870	1 bag of 1870 in box
W017944	Blue Cap / White Bulb Dropper Assembly for 4.5ml	–	1870	1 bag of 1870 in box
W017945	Yellow Cap / White Bulb Dropper Assembly for 4.5ml	–	1870	1 bag of 1870 in box
W017946	Red Cap / White Bulb Dropper Assembly for 4.5ml	–	1870	1 bag of 1870 in box
3ml Dropper Bottles and Assemblies				
W017948	3ml Clear Type 1B Tubular Glass Dropper Bottle	16.25 x 35.5	374	MOQ applies – 324 trays of 374 palletised
W017948/5	3ml Clear Type 1B Tubular Glass Dropper Bottle	16.25 x 35.5	1870	5 trays of 374 in box
W017950	Black Cap / White Bulb Dropper Assembly for 3ml	–	1870	1 bag of 1870 in box
W017951	White Cap / White Bulb Dropper Assembly for 3ml	–	1870	1 bag of 1870 in box
W017952	Blue Cap / White Bulb Dropper Assembly for 3ml	–	1870	1 bag of 1870 in box
W017953	Yellow Cap / White Bulb Dropper Assembly for 3ml	–	1870	1 bag of 1870 in box
W017954	Red Cap / White Bulb Dropper Assembly for 3ml	–	1870	1 bag of 1870 in box
W017955	Green Cap / White Bulb Dropper Assembly for 3ml	–	1870	1 bag of 1870 in box
Lab Kits				
W212000	10ml Clear Glass Dropper Bottle with Separate Black Cap / Black Bulb Dropper Assembly	–	490	2 trays of 245 vials and 1 bag of 490 assemblies in box
W212001	7ml Clear Glass Dropper Bottle with Separate Black Cap / Black Bulb Dropper Assembly	–	490	2 trays of 245 vials and 1 bag of 490 assemblies in box
W212002	4.5ml Clear Glass Dropper Bottle with Separate White Cap / White Bulb Dropper Assembly	–	1122	4 trays of 374 vials and 1 bag of 1122 assemblies in box
W212003	3ml Clear Glass Dropper Bottle with Separate White Cap / White Bulb Dropper Assembly	–	1122	4 trays of 374 vials and 1 bag of 1122 assemblies in box



Scilabware is a company of **DWK Life Sciences**

For more information visit www.scilabware.com or email us with your questions enquiries@scilabware.com

Scilabware Limited, Unit 4, Riverside 2, Campbell Road, Stoke-on-Trent, ST4 4RJ, UK
Tel: +44 (0)1782 444406 | Fax: +44 (0)1782 940436